

Data Acquisition Toolbox

For Use with MATLAB®

Computation
└─

Visualization
└─

Programming
└─



User's Guide

Version 2

How to Contact The MathWorks:



508-647-7000

Phone



508-647-7001

Fax



The MathWorks, Inc.
3 Apple Hill Drive
Natick, MA 01760-2098

Mail



<http://www.mathworks.com>
<ftp.mathworks.com>
<comp.soft-sys.matlab>

Web
Anonymous FTP server
Newsgroup



support@mathworks.com
suggest@mathworks.com
bugs@mathworks.com
doc@mathworks.com
subscribe@mathworks.com
service@mathworks.com
info@mathworks.com

Technical support
Product enhancement suggestions
Bug reports
Documentation error reports
Subscribing user registration
Order status, license renewals, passcodes
Sales, pricing, and general information

Data Acquisition Toolbox User's Guide

© COPYRIGHT 1999 - 2000 by The MathWorks, Inc.

The software described in this document is furnished under a license agreement. The software may be used or copied only under the terms of the license agreement. No part of this manual may be photocopied or reproduced in any form without prior written consent from The MathWorks, Inc.

FEDERAL ACQUISITION: This provision applies to all acquisitions of the Program and Documentation by or for the federal government of the United States. By accepting delivery of the Program, the government hereby agrees that this software qualifies as "commercial" computer software within the meaning of FAR Part 12.212, DFARS Part 227.7202-1, DFARS Part 227.7202-3, DFARS Part 252.227-7013, and DFARS Part 252.227-7014. The terms and conditions of The MathWorks, Inc. Software License Agreement shall pertain to the government's use and disclosure of the Program and Documentation, and shall supersede any conflicting contractual terms or conditions. If this license fails to meet the government's minimum needs or is inconsistent in any respect with federal procurement law, the government agrees to return the Program and Documentation, unused, to MathWorks.

MATLAB, Simulink, Stateflow, Handle Graphics, and Real-Time Workshop are registered trademarks, and Target Language Compiler is a trademark of The MathWorks, Inc.

Other product or brand names are trademarks or registered trademarks of their respective holders.

Printing History: May 1999 First printing New for Version 1
November 2000 Second printing Revised for Version 2 (Release 12)

Preface

What Is the Data Acquisition Toolbox?	xvi
Exploring the Toolbox	xvi
Related Products	xvii
System Requirements	xvii
Associated Products	xvii
Using This Guide	xix
Expected Background	xix
Using the Documentation Examples	xx
How This Guide is Organized	xx
Installation Information	xxii
Toolbox Installation	xxii
Hardware and Driver Installation	xxii
Typographical Conventions	xxiii

Introduction to Data Acquisition

1

Overview	1-2
Anatomy of a Data Acquisition Experiment	1-3
The Data Acquisition System	1-4
Data Acquisition Hardware	1-6
Sensors	1-8
Signal Conditioning	1-11
The Computer	1-13

Software	1-13
The Analog Input Subsystem	1-16
Sampling	1-16
Quantization	1-20
Channel Configuration	1-24
Transferring Data from Hardware to System Memory	1-27
Making Quality Measurements	1-29
Accuracy and Precision	1-29
Noise	1-33
Matching the Sensor Range and A/D Converter Range	1-34
How Fast Should a Signal be Sampled?	1-35
Selected Bibliography	1-39

Getting Started with the Data Acquisition Toolbox

2

Toolbox Components	2-2
M-File Functions	2-3
The Data Acquisition Engine	2-4
The Hardware Driver Adaptor	2-7
Accessing Your Hardware	2-8
Acquiring Data	2-8
Outputting Data	2-9
Reading and Writing Digital Values	2-10
Understanding the Toolbox Capabilities	2-12
The Contents and Readme Files	2-12
Documentation Examples	2-12
The Quick Reference Guide	2-13
Demos	2-13
Examining Your Hardware Resources	2-17
General Toolbox Information	2-17

Adaptor-Specific Information	2-18
Device Object Information	2-19
Getting Help	2-20
The daqhelp Function	2-20
The propinfo Function	2-21

The Data Acquisition Session

3

Overview	3-2
Example: The Data Acquisition Session	3-3
Creating a Device Object	3-4
Creating an Array of Device Objects	3-5
Where Do Device Objects Exist?	3-6
Adding Channels or Lines	3-8
Mapping Hardware Channel IDs to MATLAB Indices	3-9
Configuring and Returning Properties	3-12
Property Types	3-12
Returning Property Names and Property Values	3-14
Configuring Property Values	3-18
Specifying Property Names	3-19
Default Property Values	3-19
daqpropedit: A Graphical Property Editor	3-20
Acquiring and Outputting Data	3-22
Starting the Device Object	3-23
Logging or Sending Data	3-23
Stopping the Device Object	3-24
Cleaning Up	3-25

Getting Started with Analog Input

4

Overview	4-2
Creating an Analog Input Object	4-3
Adding Channels to an Analog Input Object	4-4
Referencing Individual Hardware Channels	4-6
Example: Adding Channels for a Sound Card	4-7
Configuring Analog Input Properties	4-9
The Sampling Rate	4-9
Trigger Types	4-11
The Samples to Acquire per Trigger	4-12
Acquiring Data	4-13
Starting the Analog Input Object	4-13
Logging Data	4-14
Stopping the Analog Input Object	4-14
Analog Input Examples	4-15
Acquiring Data with a Sound Card	4-15
Acquiring Data with a National Instruments Board	4-19
Evaluating the Analog Input Object Status	4-22
Status Properties	4-22
The Display Summary	4-23

Doing More with Analog Input

5

Overview	5-2
Configuring and Sampling Input Channels	5-3
Input Channel Configuration	5-4
Sampling Rate	5-6

Channel Skew	5-7
Managing Acquired Data	5-9
Previewing Data	5-9
Extracting Data from the Engine	5-13
Returning Time Information	5-18
Configuring Analog Input Triggers	5-20
Defining a Trigger: Trigger Types and Conditions	5-21
Executing the Trigger	5-26
Trigger Delays	5-26
Repeating Triggers	5-30
How Many Triggers Occurred?	5-36
When Did the Trigger Occur?	5-37
Device-Specific Hardware Triggers	5-37
Configuring Events and Actions	5-46
Event Types	5-46
Recording and Retrieving Event Information	5-49
Creating and Executing Action Functions	5-52
Examples: Using Action Properties and Functions	5-53
Linearly Scaling the Data: Engineering Units	5-56
Example: Performing a Linear Conversion	5-57

Analog Output

6

Overview	6-2
Getting Started with Analog Output	6-3
Creating an Analog Output Object	6-3
Adding Channels to an Analog Output Object	6-4
Configuring Analog Output Properties	6-6
Outputting Data	6-9
Analog Output Examples	6-10
Evaluating the Analog Output Object Status	6-14

Managing Output Data	6-17
Queuing Data with putdata	6-17
Example: Queuing Data with putdata	6-19
Configuring Analog Output Triggers	6-21
Defining a Trigger: Trigger Types	6-22
Executing the Trigger	6-23
How Many Triggers Occurred?	6-23
When Did the Trigger Occur?	6-24
Device-Specific Hardware Triggers	6-25
Configuring Events and Actions	6-27
Event Types	6-27
Recording and Retrieving Event Information	6-29
Examples: Using Action Properties and Action Functions ...	6-31
Linearly Scaling the Data: Engineering Units	6-34
Example: Performing a Linear Conversion	6-35
Starting Multiple Device Objects	6-37

Digital Input/Output

7

Overview	7-2
Creating a Digital I/O Object	7-3
Adding Lines to a Digital I/O Object	7-4
Line and Port Characteristics	7-5
Referencing Individual Hardware Lines	7-8
Writing and Reading Digital I/O Line Values	7-10
Writing Digital Values	7-10
Reading Digital Values	7-12
Example: Writing and Reading Digital Values	7-13

Generating Timer Events	7-15
Timer Events	7-15
Starting and Stopping a Digital I/O Object	7-16
Example: Generating Timer Events	7-17
Evaluating the Digital I/O Object Status	7-18
The Display Summary	7-18

Saving and Loading the Session

8

Overview	8-2
Saving and Loading Device Objects	8-3
Saving Device Objects to an M-File	8-3
Saving Device Objects to a MAT-File	8-5
Logging Information to Disk	8-6
Specifying a Filename	8-7
Retrieving Logged Information	8-8
Example: Logging and Retrieving Information	8-10

Function Reference

9

Overview	9-2
Getting Command Line Function Help	9-2
Functions Grouped by Category	9-4
Functions Listed Alphabetically	9-8
addchannel	9-9
addline	9-14
addmuxchannel	9-17

analoginput	9-18
analogoutput	9-21
binvec2dec	9-24
clear	9-25
daqaction	9-27
daqfind	9-28
daqhelp	9-30
daqhwinfo	9-32
daqmem	9-34
daqpropedit	9-37
daqread	9-39
daqregister	9-43
daqreset	9-44
daqschool	9-45
dec2binvec	9-46
delete	9-48
digitalio	9-50
disp	9-52
flushdata	9-53
get	9-55
getdata	9-57
getsample	9-61
getvalue	9-62
ischannel	9-63
isdioline	9-64
isvalid	9-65
length	9-67
load	9-69
makenames	9-71
muxchanidx	9-72
obj2mfile	9-74
peekdata	9-76
propinfo	9-78
putdata	9-80
putsample	9-82
putvalue	9-83
save	9-84
set	9-85
setverify	9-88
showdaqevents	9-90

size	9-92
start	9-94
stop	9-95
trigger	9-97
waittilstop	9-98

Base Property Reference

10

Overview	10-2
Getting Command Line Property Help	10-2
Properties Grouped by Category	10-3
Analog Input Properties	10-3
Analog Output Properties	10-8
Digital I/O Properties	10-12
Properties Listed Alphabetically	10-14
BufferingConfig	10-15
BufferingMode	10-17
Channel	10-18
ChannelName	10-20
ChannelSkew	10-21
ChannelSkewMode	10-22
ClockSource	10-24
DataMissedAction	10-26
DefaultChannelValue	10-27
Direction	10-28
EventLog	10-29
HwChannel	10-31
HwLine	10-33
Index	10-34
InitialTriggerTime	10-36
InputOverRangeAction	10-38
InputRange	10-39
InputType	10-41
Line	10-43

LineName	10-44
LogFileName	10-45
Logging	10-46
LoggingMode	10-47
LogToDiskMode	10-48
ManualTriggerHwOn	10-49
MaxSamplesQueued	10-51
Name	10-52
NativeOffset	10-53
NativeScaling	10-55
OutputRange	10-56
Parent	10-58
Port	10-59
RepeatOutput	10-60
Running	10-61
RuntimeErrorAction	10-62
SampleRate	10-64
SamplesAcquired	10-66
SamplesAcquiredAction	10-67
SamplesAcquiredActionCount	10-68
SamplesAvailable	10-69
SamplesOutput	10-70
SamplesOutputAction	10-71
SamplesOutputActionCount	10-72
SamplesPerTrigger	10-73
Sending	10-74
SensorRange	10-75
StartAction	10-76
StopAction	10-77
Tag	10-79
Timeout	10-80
TimerAction	10-81
TimerPeriod	10-82
TriggerAction	10-83
TriggerChannel	10-84
TriggerCondition	10-85
TriggerConditionValue	10-89
TriggerDelay	10-90
TriggerDelayUnits	10-91
TriggerRepeat	10-92

TriggersExecuted	10-93
TriggerType	10-94
Type	10-96
Units	10-97
UnitsRange	10-98
UserData	10-99

11 | Device-Specific Property Reference

Overview	11-2
Getting Command Line Property Help	11-2
Properties Listed by Vendor	11-3
National Instruments Properties	11-3
ComputerBoards Properties	11-4
Agilent Technologies Properties	11-4
Sound Card Properties	11-5
Properties Listed Alphabetically	11-6
BitsPerSample	11-7
COLA	11-8
Coupling	11-9
DriveAISenseToGround	11-10
GroundingMode	11-11
InputMode	11-12
InputSource	11-14
NumMuxBoards	11-15
OutOfDataMode	11-16
RampRate	11-17
SourceMode	11-18
SourceOutput	11-19
Span	11-21
StandardSampleRates	11-23
Sum	11-24
TransferMode	11-25

Troubleshooting Your Hardware

A

Overview	A-2
Troubleshooting National Instruments Hardware	A-3
What Driver Are You Using?	A-3
Is Your Hardware Functioning Properly?	A-5
Troubleshooting ComputerBoards Hardware	A-6
What Driver Are You Using?	A-6
Is Your Hardware Functioning Properly?	A-8
Troubleshooting Agilent Technologies Hardware	A-9
What Driver Are You Using?	A-9
Is Your Hardware Functioning Properly?	A-10
Troubleshooting Sound Cards	A-12
Microphone and Sound Card Types	A-15
Testing with a Microphone	A-16
Testing with a CD Player	A-16
Running in Full Duplex Mode	A-18
Other Things to Try	A-19
Registering the Hardware Driver Adaptor	A-19
Contacting The MathWorks	A-20

Managing Your Memory Resources

B

Overview	B-2
How Much Memory Do You Need?	B-3
Example: Managing Memory Resources	B-4

C

Preface

What Is the Data Acquisition Toolbox?	xvi
Exploring the Toolbox	xvi
Related Products	xvii
System Requirements	xvii
Associated Products	xvii
Using This Guide	xix
Expected Background	xix
Using the Documentation Examples	xx
How This Guide is Organized	xx
Installation Information	xxii
Toolbox Installation	xxii
Hardware and Driver Installation	xxii
Typographical Conventions	xxiii

What Is the Data Acquisition Toolbox?

The Data Acquisition Toolbox is a collection of M-file functions and MEX-file dynamic link libraries (DLLs) built on the MATLAB® Technical Computing Environment. The toolbox provides you with these main features:

- A framework for bringing live, measured data into MATLAB using PC-compatible, plug-in data acquisition hardware
- Support for these popular hardware devices:
 - National Instruments E Series and 1200 Series boards
 - ComputerBoards (Measurement Computing Corporation) boards
 - Agilent Technologies E1432A/33A/34A VXI modules
 - Microsoft Windows sound cards

Additionally, you can use the Data Acquisition Toolbox Adaptor Kit to interface unsupported hardware devices to the toolbox.

- Support for analog input (AI), analog output (AO), and digital I/O (DIO) subsystems including simultaneous AI and AO conversions
- Data buffering for background acquisitions
- Data logging
- Event-driven acquisitions

Exploring the Toolbox

A list of the toolbox functions is available to you by typing

```
help daq
```

You can view the code for any function by typing

```
type function_name
```

You can view the help for any function by typing

```
daqhelp function_name
```

You can change the way any toolbox function works by copying and renaming the M-file, then modifying your copy. You can also extend the toolbox by adding your own M-files, or by using it in combination with other products such as the Signal Processing Toolbox or the Instrument Control Toolbox.

Related Products

System Requirements

The Data Acquisition Toolbox is a multiplatform product that you install on a host computer running Microsoft Windows 95, Windows 98, Windows 2000, or Windows NT 4.0.

The toolbox requires:

- MATLAB 6.0 (Release 12)
- A supported data acquisition device – for a complete listing of all supported devices, visit the Data Acquisition Toolbox section of the MathWorks Web site at <http://www.mathworks.com/products/daq/>
- Software such as drivers and support libraries, as required by your data acquisition device

Associated Products

The MathWorks provides several associated products that are especially relevant to the kinds of tasks you can perform with the Data Acquisition Toolbox. For more information about any of these products, see either:

- The online documentation for that product, if it is installed or if you are reading the documentation from the CD
- The MathWorks Web site, at <http://www.mathworks.com>; see the “products” section

Note The toolboxes listed below all include functions that extend MATLAB’s capabilities.

Products Associated with the Data Acquisition Toolbox

Product	Description
Control System Toolbox	Tool for modeling, analyzing, and designing control systems using classical and modern techniques
Database Toolbox	Tool for connecting to, and interacting with, most ODBC/JDBC databases from within MATLAB
Instrument Control Toolbox	Tool for communicating with instruments that support the GPIB (IEEE-488, HP-IB) interface, the VISA standard, or the serial port interface
MATLAB Report Generator	Tool for documenting information in MATLAB in multiple output formats
Neural Network Toolbox	A comprehensive environment for neural network research, design, and simulation within MATLAB
Signal Processing Toolbox	Tool for algorithm development, signal and linear system analysis, and time-series data modeling
Statistics Toolbox	Tool for analyzing historical data, modeling systems, developing statistical algorithms, and learning and teaching statistics
System Identification Toolbox	Tool for building accurate, simplified models of complex systems from noisy time-series data
Wavelet Toolbox	Tool for signal and image analysis, compression, and de-noising

Using This Guide

Expected Background

To use the Data Acquisition Toolbox, you should have some familiarity with:

- The basic features of MATLAB
- The capabilities of your hardware device
- The basic concepts associated with acquiring live, measured data

If You Are a New User

For a brief review of basic data acquisition concepts, you should start with Chapter 1, “Introduction to Data Acquisition.” Otherwise, start with Chapter 2, “Getting Started with the Data Acquisition Toolbox”, which provides simple examples that illustrate how to input and output data. Then read the appropriate chapter based on the hardware subsystem you are using. For example, if you are acquiring data with an analog input subsystem, you should read Chapter 4, “Getting Started with Analog Input.”

After you have successfully transferred data between your hardware device and MATLAB, you should read the appropriate reference material as needed. If you want detailed information about a specific function, refer to Chapter 9, “Function Reference.” If you want detailed information about a specific property, refer to Chapter 10, “Base Property Reference” or Chapter 11, “Device-Specific Property Reference.”

If You Are an Experienced User

Start with “Data Acquisition Toolbox 2.0” in the *Release Notes* for a description of new and modified toolbox features.

You should then read the appropriate reference material as needed. If you want detailed information about a specific function, refer to Chapter 9, “Function Reference.” If you want detailed information about a specific property, refer to Chapter 10, “Base Property Reference” or Chapter 11, “Device-Specific Property Reference.”

Using the Documentation Examples

When you encounter examples or code snippets in this book, you may want to try them for yourself. An easy way to do this is to cut the relevant text from the PDF or HTML versions of this guide and paste it into the MATLAB workspace. You access the PDF and HTML content with the Help browser.

Some examples are constructed as mini-applications that illustrate one or two important toolbox features, and serve as templates so you can see how to build applications that suit your specific needs. These examples are included as toolbox M-files and are treated as demos. You can list all Data Acquisition Toolbox demos by typing

```
hel p daqdem os
```

All documentation example M-files begin with daqdoc. To run an example, type the M-file name at the command line. Note that most examples are written for specific hardware devices. To use these examples with your hardware device, you should modify the creation function input arguments and the device object property values as needed.

How This Guide is Organized

The organization of this guide is described below.

Chapter	Description
Introduction to Data Acquisition	Provides you with general information about making measurements with data acquisition hardware. The topics covered should help you understand the specification sheet associated with your hardware.
Getting Started with the Data Acquisition Toolbox	Describes the toolbox components, and shows you how to access your hardware, examine your hardware resources, and get command line help.
The Data Acquisition Session	Describes all the steps you are likely to take when acquiring or outputting data.
Getting Started with Analog Input	Shows you how to perform basic data acquisition tasks using your analog input subsystem.

Chapter	Description
Doing More with Analog Input	Presents the complete analog input functionality.
Analog Output	Shows you how to perform data acquisition tasks using your analog output subsystem.
Digital Input/Output	Shows you how to perform data acquisition tasks using your digital I/O subsystem.
Saving and Loading the Session	Shows you how to save and load device objects, data, and event information using several disk file formats.
Function Reference	Presents a complete description of all toolbox functions.
Base Property Reference	Presents a complete description of all toolbox base properties.
Device-Specific Property Reference	Presents a complete description of all toolbox device-specific properties.
Troubleshooting Your Hardware	Presents tips to help you troubleshoot your hardware.
Managing Your Memory Resources	Describes how to allocate and return memory resources.
Glossary	Presents a glossary of data acquisition terms.

Installation Information

To acquire live, measured data into the MATLAB environment, or to output data from the MATLAB environment, you must install these three components:

- The Data Acquisition Toolbox
- Data acquisition hardware
- Software such as drivers and support libraries, as required by your data acquisition hardware

Toolbox Installation

To determine if the Data Acquisition Toolbox is installed on your system, type

```
ver
```

at the MATLAB prompt. MATLAB displays information about the versions of MATLAB you are running, including a list of installed add-on products and their version numbers. Check the list to see if the Data Acquisition Toolbox appears. For information about installing the toolbox, see the *MATLAB Installation Guide* for your platform.

If you experience installation difficulties and have Web access, look for the license manager and installation information at the MathWorks Web site (<http://www.mathworks.com>).

Hardware and Driver Installation

Installation of your hardware device, hardware drivers, and any other device-specific software is described in the documentation provided by your hardware vendor.

Note You need to install all necessary device-specific software provided by your hardware vendor in addition to the Data Acquisition Toolbox.

Typographical Conventions

This guide uses the following typographical conventions.

Item	Convention to Use	Example
Example code	Monospace font	To assign the value 5 to A, enter A = 5
Function names/syntax	Monospace font	The cos function finds the cosine of each array element. Syntax line example is MLGetVar ML_var_name
Keys	Boldface with an initial capital letter	Press the Return key.
Literal string (in syntax descriptions in Reference chapters)	Use monospace bold for literals.	f = freqspace(n, ' whole ')
Mathematical expressions	Variables in <i>italics</i> Functions, operators, and constants in standard text.	This vector represents the polynomial $p = x^2 + 2x + 3$
MATLAB output	Monospace font	MATLAB responds with A = 5
Menu names, menu items, and controls	Boldface with an initial capital letter	Choose the File menu.
New terms	<i>Italics</i>	An <i>array</i> is an ordered collection of information.
String variables (from a finite list)	<i>Monospace italics</i>	sysc = d2c(sysd, ' <i>method</i> ')

Introduction to Data Acquisition

Overview	1-2
Anatomy of a Data Acquisition Experiment	1-3
The Data Acquisition System	1-4
Data Acquisition Hardware	1-6
Sensors	1-8
Signal Conditioning	1-11
The Computer	1-13
Software	1-13
The Analog Input Subsystem	1-16
Sampling	1-16
Quantization	1-20
Channel Configuration	1-24
Transferring Data from Hardware to System Memory	1-27
Making Quality Measurements	1-29
Accuracy and Precision	1-29
Noise	1-33
Matching the Sensor Range and A/D Converter Range	1-34
How Fast Should a Signal be Sampled?	1-35
Selected Bibliography	1-39

Overview

Before you set up any data acquisition system, you should understand the physical quantities you want to measure, the characteristics of those physical quantities, the appropriate sensor to use, and the appropriate data acquisition hardware to use.

The purpose of this chapter is to provide you with some general information about making measurements with a data acquisition system. The topics covered should assist you in understanding the above considerations, and understanding the specification sheet associated with your hardware. Topics include:

- Anatomy of a data acquisition experiment
- The data acquisition system
- The analog input subsystem
- Making quality measurements
- Selected bibliography

Anatomy of a Data Acquisition Experiment

For each new data acquisition experiment, you need to perform these tasks:

- System setup
- Calibration
- Trials

System Setup

The first step in any data acquisition experiment is to install the hardware and software. Hardware installation consists of plugging a board into your computer or installing modules into an external chassis. Software installation consists of loading hardware drivers and application software onto your computer. After the hardware and software are installed, you can attach your sensors.

Calibration

After the hardware and software are installed and the sensors are connected, the data acquisition hardware should be *calibrated*. Calibration consists of providing a known input to the system and recording the output. For many data acquisition devices, calibration can be easily accomplished with software provided by the vendor.

Trials

After the hardware is set up and calibrated, you can begin to acquire data. You may think that if you completely understand the characteristics of the signal you are measuring, then you should be able to configure your data acquisition system and acquire the data.

In the real world however, your sensor may be picking up unacceptable noise levels and require shielding, or you may need to run the device at a higher rate, or perhaps you need to add an antialias filter to remove unwanted frequency components.

These real-world effects act as obstacles between you and a precise, accurate measurement. To overcome these obstacles, you need to experiment with different hardware and software configurations. In other words, you need to perform multiple data acquisition trials.

The Data Acquisition System

As a user of MATLAB and the Data Acquisition Toolbox, you are interested in measuring and analyzing physical phenomena. The purpose of any data acquisition system is to provide you with the tools and resources necessary to do so.

You can think of a data acquisition system as a collection of software and hardware that connects you to the physical world. A typical data acquisition system consists of these components:

- **Data acquisition hardware**

At the heart of any data acquisition system lies the data acquisition hardware. The main function of this hardware is to convert analog signals to digital signals, and to convert digital signals to analog signals.

- **Sensors and actuators (transducers)**

Sensors and actuators can both be *transducers*. A transducer is a device that converts input energy of one form into output energy of another form. For example, a microphone is a sensor that converts sound energy (in the form of pressure) into electrical energy, while a loudspeaker is an actuator that converts electrical energy into sound energy.

- **Signal conditioning hardware**

Sensor signals are often incompatible with data acquisition hardware. To overcome this incompatibility, the signal must be conditioned. For example, you may need to condition an input signal by amplifying it or by removing unwanted frequency components. Output signals may need conditioning as well. However, only input signal conditioning is discussed in this chapter.


- **The computer**

The computer provides a processor, a system clock, a bus to transfer data, and memory and disk space to store data.

- **Software**

Data acquisition software allows you to exchange information between the computer and the hardware. For example, typical software allows you to configure the sampling rate of your board, and acquire a predefined amount of data.

The data acquisition components, and their relationship to each other, are shown below.



The figure depicts the two important features of a data acquisition system:

- Signals are input to a sensor, conditioned, converted into bits that a computer can read, and analyzed to extract meaningful information.

For example, sound level data is acquired from a microphone, amplified, digitized by a sound card, and stored in MATLAB for subsequent analysis of frequency content.

- Data from a computer is converted into an analog signal and output to an actuator.

For example, a vector of data in MATLAB is converted to an analog signal by a sound card and output to a loudspeaker.


Data Acquisition Hardware

Data acquisition hardware is either internal and installed directly into an expansion slot inside your computer, or external and connected to your computer through an external cable. For example, VXI modules are installed in an external VXI chassis, and data is transferred between MATLAB and the module using an external link such as FireWire (IEEE 1394).

At the simplest level, data acquisition hardware is characterized by the *subsystems* it possesses. A subsystem is a component of your data acquisition hardware that performs a specialized task. Common subsystems include:

- Analog input
- Analog output
- Digital input/output
- Counter/timer

Hardware devices that consist of multiple subsystems, such as the one depicted below, are called *multifunction boards*.



Analog Input Subsystems

Analog input subsystems convert real-world analog input signals from a sensor into bits that can be read by your computer. Perhaps the most important of all the subsystems commonly available, they are typically multichannel devices offering 12 or 16 bits of resolution.

Analog input subsystems are also referred to as AI subsystems, A/D converters, or ADCs. Analog input subsystems are discussed in detail beginning in “The Analog Input Subsystem” on page 1-16.

Analog Output Subsystems

Analog output subsystems convert digital data stored on your computer to a real-world analog signal. These subsystems perform the inverse conversion of analog input subsystems. Typical acquisition boards offer two output channels with 12 bits of resolution, with special hardware available to support multiple channel analog output operations.

Analog output subsystems are also referred to as AO subsystems, D/A converters, or DACs.

Digital Input/Output Subsystems

Digital input/output (DIO) subsystems are designed to input and output digital values (logic levels) to and from hardware. These values are typically handled either as single bits or *lines*, or as a *port*, which typically consists of eight lines.

While most popular data acquisition cards include some digital I/O capability, it is usually limited to simple operations, and special dedicated hardware is often necessary for performing advanced digital I/O operations.

Counter/Timer Subsystems

Counter/timer (C/T) subsystems are used for event counting, frequency and period measurement, and pulse train generation.

Sensors

A sensor converts the physical phenomena of interest into a signal that is input into your data acquisition hardware. There are two main types of sensors based on the output they produce: digital sensors and analog sensors.

Digital sensors produce an output signal that is a digital representation of the input signal, and has discrete values of magnitude measured at discrete times. A digital sensor must output logic levels that are compatible with the digital receiver. Some standard logic levels include transistor-transistor logic (TTL) and emitter-coupled logic (ECL). Examples of digital sensors include switches and position encoders.

Analog sensors produce an output signal that is directly proportional to the input signal, and is continuous in both magnitude and in time. Most physical variables such as temperature, pressure, and acceleration are continuous in nature and are readily measured with an analog sensor. For example, the temperature of an automobile cooling system and the acceleration produced by a child on a swing all vary continuously.

The sensor you use depends on the phenomena you are measuring. Some common analog sensors and the physical variables they measure are listed below.

Table 1-1: Common Analog Sensors

Sensor	Physical Variable
Accelerometer	Acceleration
Microphone	Pressure
Pressure gauge	Pressure
Resistance temperature device (RTD)	Temperature
Strain gauge	Force
Thermocouple	Temperature

When choosing the best analog sensor to use, you must match the characteristics of the physical variable you are measuring with the characteristics of the sensor. The two most important sensor characteristics are:

- The sensor output
- The sensor bandwidth

Sensor Output

The output from a sensor can be an analog signal or a digital signal, and the output variable is usually a voltage although some sensors output current.

Current Signals. Current is often used to transmit signals in noisy environments because it is much less affected by environmental noise. The full scale range of the current signal is often either 4-20 mA or 0-20 mA. A 4-20 mA signal has the advantage that even at minimum signal value, there should be a detectable current flowing. The absence of this indicates a wiring problem.

Before conversion by the analog input subsystem, the current signals are usually turned into voltage signals by a current-sensing resistor. The resistor should be of high precision, perhaps 0.03% or 0.01% depending on the resolution of your hardware. Additionally, the voltage signal should match the signal to an input range of the analog input hardware. For 4-20 mA signals, a 50 ohm resistor will give a voltage of 1 V for a 20 mA signal by Ohm's law.

Voltage Signals. The most commonly interfaced signal is a voltage signal. For example, thermocouples, strain gauges, and accelerometers all produce voltage signals. There are three major aspects of a voltage signal that you need to consider:

- **Amplitude**

If the signal is smaller than a few millivolts, you may need to amplify it. If it is larger than the maximum range of your analog input hardware (typically ± 10 V), you will have to divide the signal down using a resistor network.

The amplitude is related to the sensitivity (resolution) of your hardware. Refer to "Accuracy and Precision" on page 1-29 for more information about hardware sensitivity.

- **Frequency**

Whenever you acquire data, you should decide the highest frequency you want to measure.

The highest frequency component of the signal determines how often you should sample the input. If you have more than one input, but only one analog input subsystem, then the overall sampling rate goes up in proportion to the number of inputs. Higher frequencies may be present as noise, which you can remove by filtering the signal before it is digitized.

If you sample the input signal at least twice as fast as the highest frequency component, then that signal will be uniquely characterized. However, this rate may not mimic the waveform very closely. For a rapidly varying signal, you may need a sampling rate of roughly 10 to 20 times the highest frequency to get an accurate picture of the waveform. For slowly varying signals, you need only consider the minimum time for a significant change in the signal.

The frequency is related to the bandwidth of your measurement. Bandwidth is discussed in the next section.

- **Duration**

How long do you want to sample the signal for? If you are storing data to memory or to a disk file, then the duration determines the storage resources required. The format of the stored data also affects the amount of storage space required. For example, data stored in ASCII format takes more space than data stored in binary format.

Sensor Bandwidth

In a real-world data acquisition experiment, the physical phenomena you are measuring has expected limits. For example, the temperature of your automobile's cooling system varies continuously between its low limit and high limit. The temperature limits, as well as how rapidly the temperature varies between the limits, depends on several factors including your driving habits, the weather, and the condition of the cooling system. The expected limits may be readily approximated, but there are an infinite number of possible temperatures that you can measure at a given time. As explained in "Quantization" on page 1-20, these unlimited possibilities are mapped to a finite set of values by your data acquisition hardware.

The *bandwidth* is given by the range of frequencies present in the signal being measured. You can also think of bandwidth as being related to the rate of change of the signal. A slowly varying signal has a low bandwidth, while a

rapidly varying signal has a high bandwidth. To properly measure the physical phenomena of interest, the sensor bandwidth must be compatible with the measurement bandwidth.

You may want to use sensors with the widest possible bandwidth when making any physical measurement. This is the one way to ensure that the basic measurement system is capable of responding linearly over the full range of interest. However, the wider the bandwidth of the sensor, the more you must be concerned with eliminating sensor response to unwanted frequency components.

Signal Conditioning

Sensor signals are often incompatible with data acquisition hardware. To overcome this incompatibility, the sensor signal must be conditioned. The type of signal conditioning required depends on the sensor you are using. For example, a signal may have a small amplitude and require amplification, or it may contain unwanted frequency components and require filtering. Common ways to condition signals include:

- Amplification
- Filtering
- Electrical isolation
- Multiplexing
- Excitation source

Amplification

Low-level signals – less than around 100 millivolts – usually need to be amplified. High-level signals may also require amplification depending on the input range of the analog input subsystem.

For example, the output signal from a thermocouple is small and must be amplified before it is digitized. Signal amplification allows you to reduce noise and to make use of the full range of your hardware thereby increasing the resolution of the measurement.

Filtering

Filtering removes unwanted noise from the signal of interest. A noise filter is used on slowly varying signals such as temperature to attenuate higher frequency signals that can reduce the accuracy of your measurement.

Rapidly varying signals such as vibration often require a different type of filter known as an antialiasing filter. An antialiasing filter removes undesirable higher frequencies that may lead to erroneous measurements.

Electrical Isolation

If the signal of interest may contain high-voltage transients that could damage the computer, then the sensor signals should be electrically isolated from the computer for safety purposes.

You can also use electrical isolation to make sure that the readings from the data acquisition hardware are not affected by differences in ground potentials. For example, when the hardware device and the sensor signal are each referenced to ground, problems occur if there is a potential difference between the two grounds. This difference can lead to a *ground loop*, which may lead to erroneous measurements. Using electrically isolated signal conditioning modules eliminates the ground loop and ensures that the signals are accurately represented.

Multiplexing

A common technique for measuring several signals with a single measuring device is multiplexing.

Signal conditioning devices for analog signals often provide multiplexing for use with slowly changing signals such as temperature. This is in addition to any built-in multiplexing on the DAQ board. The A/D converter samples one channel, switches to the next channel and samples it, switches to the next channel, and so on. Because the same A/D converter is sampling many channels, the effective sampling rate of each individual channel is inversely proportional to the number of channels sampled.

You must take care when using multiplexers so that the switched signal has sufficient time to settle. Refer to “Noise” on page 1-33 for more information about settling time.

Excitation Source

Some sensors require an excitation source to operate. For example, strain gauges, and resistive temperature devices (RTDs) require external voltage or current excitation. Signal conditioning modules for these sensors usually provide the necessary excitation. RTD measurements are usually made with a current source that converts the variation in resistance to a measurable voltage.

The Computer

The computer provides a processor, a system clock, a bus to transfer data, and memory and disk space to store data.

The processor controls how fast data is accepted by the converter. The system clock provides time information about the acquired data. Knowing that you recorded a sensor reading is generally not enough. You also need to know when that measurement occurred.

Data is transferred from the hardware to system memory via dynamic memory access (DMA) or interrupts. DMA is hardware controlled and therefore extremely fast. Interrupts may be slow due to the latency time between when a board requests interrupt servicing and when the computer responds. The maximum acquisition rate is also determined by the computer's bus architecture. Refer to "How Are Acquired Samples Clocked?" on page 1-23 for more information about DMA and interrupts.


Software

Regardless of the hardware you are using, you must send information to the hardware and receive information from the hardware. You send configuration information to the hardware such as the sampling rate, and receive information from the hardware such as data, status messages, and error messages. You may also need to supply the hardware with information so that you can integrate it with other hardware and with computer resources. This information exchange is accomplished with software.

There are two kinds of software:

- Driver software
- Application software

For example, suppose you are using the Data Acquisition Toolbox with a National Instruments AT-MIO-16E-1 board and its associated NI-DAQ driver. The relationship between you, the driver software, the application software, and the hardware is shown below.



The diagram illustrates that you supply information to the hardware, and you receive information from the hardware.

Driver Software

For data acquisition device, there is associated driver software that you must use. Driver software allows you to access and control the capabilities of your hardware. Among other things, basic driver software allows you to:

- Bring data on to and get data off of the board
- Control the rate at which data is acquired
- Integrate the data acquisition hardware with computer resources such as processor interrupts, DMA, and memory
- Integrate the data acquisition hardware with signal conditioning hardware
- Access multiple subsystems on a given data acquisition board
- Access multiple data acquisition boards

Application Software

Application software provides a convenient “front-end” to the driver software. Basic application software allows you to:

- Report relevant information such as the number of samples acquired
- Generate events
- Manage the data stored in computer memory
- Condition a signal
- Plot acquired data

With some application software, you can also perform analysis on the data. MATLAB and the Data Acquisition Toolbox provide you with these capabilities and more.

The Analog Input Subsystem

Many data acquisition hardware devices contain one or more subsystems that convert (digitize) real-world sensor signals into numbers your computer can read. Such devices are called analog input subsystems (AI subsystems, A/D converters, or ADCs). After the real-world signal is digitized, you can analyze it, store it in system memory, or store it to a disk file.

The function of the analog input subsystem is to *sample* and *quantize* the analog signal using one or more *channels*. You can think of a channel as a path through which the sensor signal travels. Typical analog input subsystems have eight or 16 input channels available to you. After data is sampled and quantized, it must be transferred to system memory.

Analog signals are continuous in time and in amplitude (within predefined limits). Sampling takes a “snapshot” of the signal at discrete times, while quantization divides the voltage (or current) value into discrete amplitudes. Sampling, quantization, channel configuration, and transferring data from hardware to system memory are discussed below.

Sampling

Sampling takes a “snapshot” of the sensor signal at discrete times. For most applications, the time interval between samples is kept constant (for example, sample every millisecond) unless externally clocked.

For most digital converters, sampling is performed by a sample and hold (S/H) circuit. An S/H circuit usually consists of a signal buffer followed by an electronic switch connected to a capacitor. The operation of an S/H circuit follows these steps:

- 1 At a given sampling instant, the switch connects the buffer and capacitor to an input.
- 2 The capacitor is charged to the input voltage.
- 3 The charge is held until the A/D converter digitizes the signal.
- 4 For multiple channels connected (multiplexed) to one A/D converter, the previous steps are repeated for each input channel.
- 5 The entire process is repeated for the next sampling instant.


A multiplexer, S/H circuit, and A/D converter are illustrated in the next section.

Hardware can be divided into two main categories based on how signals are sampled: *scanning* hardware, which samples input signals sequentially, and *simultaneous sample and hold* (SS/H) hardware, which samples all signals at the same time. These two types of hardware are discussed below.

Scanning Hardware

Scanning hardware samples a single input signal, converts that signal to a digital value, and then repeats the process for every input channel used. In other words, each input channel is sampled sequentially. A *scan* occurs when each input in a group is sampled once.

As shown below, most data acquisition devices have one A/D converter that is multiplexed to multiple input channels.




Therefore, if you use multiple channels, those channels cannot be sampled simultaneously and a time gap exists between consecutive sampled channels. This time gap is called the *channel skew*. You can think of the channel skew as the time it takes the analog input subsystem to sample a single channel.

Additionally, the maximum sampling rate your hardware is rated at typically applies for one channel. Therefore, the maximum sampling rate per channel is given by the formula

$$\text{Maximum sampling rate per channel} = \frac{\text{Maximum board rate}}{\text{Number of channels scanned}}$$

Typically, you can achieve this maximum rate only under ideal conditions. In practice, the sampling rate depends on several characteristics of the analog input subsystem including the settling time and the gain, as well as the channel skew. The sample period and channel skew for a multichannel configuration using scanning hardware is shown below.




If you cannot tolerate channel skew in your application, you must use hardware that allows simultaneous sampling of all channels. Simultaneous sample and hold hardware is discussed in the next section.

Simultaneous Sample and Hold Hardware

Simultaneous sample and hold (SS/H) hardware samples all input signals at the same time and holds the values until the A/D converter digitizes all the signals. For high-end systems, there can be a separate A/D converter for each input channel.

For example, suppose you need to simultaneously measure the acceleration of multiple accelerometers to determine the vibration of some device under test. To do this, you must use SS/H hardware since it does not have a channel skew. In general, you may need to use SS/H hardware if your sensor signal changes significantly in a time that is less than the channel skew, or if you need to use a transfer function or perform a frequency domain correlation.

The sample period for a multichannel configuration using SS/H hardware is shown below. Note that there is no channel skew.




Quantization

As discussed in the previous section, sampling takes a snapshot of the input signal at an instant of time. When the snapshot is taken, the sampled analog signal must be converted from a voltage value to a binary number that the computer can read. The conversion from an infinitely precise amplitude to a binary number is called *quantization*.

During quantization, the A/D converter uses a finite number of evenly spaced values to represent the analog signal. The number of different values is determined by the number of bits used for the conversion. Most modern converters use 12 or 16 bits. Typically, the converter selects the digital value that is closest to the actual sampled value.

The figure below shows a 1 Hz sine wave quantized by a 3-bit A/D converter.




The number of quantized values is given by $2^3 = 8$, the largest representable value is given by $111 = 2^2 + 2^1 + 2^0 = 7.0$, and the smallest representable value is given by $000 = 0.0$.

Quantization Error

There is always some error associated with the quantization of a continuous signal. Ideally, the maximum quantization error is ± 0.5 least significant bits (LSBs), and over the full input range, the average quantization error is zero.

As shown below, the quantization error for the previous sine wave is calculated by subtracting the actual signal from the quantized signal.




Input Range and Polarity

The *input range* of the analog input subsystem is the span of input values for which a conversion is valid. You can change the input range by selecting a different *gain* value. For example, National Instruments' AT-MIO-16E-1 board has eight gain values ranging from 0.5 to 100. Many boards include a programmable gain amplifier that allows you to change the device gain through software.

When an input signal exceeds the valid input range of the converter, an *overrange* condition occurs. In this case, most devices saturate to the largest representable value, and the converted data is almost definitely incorrect. The gain setting affects the precision of your measurement – the higher (lower) the gain value, the lower (higher) the precision. Refer to “How Are Range, Gain, and Measurement Precision Related?” on page 1-32 for more information about how input range, gain, and precision are related to each other.

An analog input subsystem can typically convert both *unipolar* signals and *bipolar* signals. A unipolar signal contains only positive values and zero, while a bipolar signal contains positive values, negative values, and zero.

Unipolar and bipolar signals are depicted below. Refer to the figure in “Quantization” on page 1-20 for an example of a unipolar signal.



In many cases, the signal polarity is a fixed characteristic of the sensor and you must configure the input range to match this polarity.

As you can see, it is crucial to understand the range of signals expected from your sensor so that you can configure the input range of the analog input subsystem to maximize resolution and minimize the chance of an overrange condition.

How Are Acquired Samples Clocked?

Samples are acquired from an analog input subsystem at a specific rate by a clock. Like any timing system, data acquisition clocks are characterized by their resolution and accuracy. Timing resolution is defined as the smallest time interval that you can accurately measure. The timing accuracy is affected by clock *jitter*. Jitter arises when a clock produces slightly different values for a given time interval.

For any data acquisition system, there are typically three clock sources that you can use: the onboard data acquisition clock, the computer clock, or an external clock. The Data Acquisition Toolbox supports all of these clock sources, depending on the requirements of your hardware.

The Onboard Clock. The onboard clock is typically a timer chip on the hardware board that is programmed to generate a pulse stream at the desired rate. The onboard clock generally has high accuracy and low jitter compared to the computer clock. You should always use the onboard clock when the sampling rate is high, and when you require a fixed time interval between samples. The onboard clock is referred to as the *internal clock* in this guide.

The Computer Clock. The computer (PC) clock is used for boards that do not possess an onboard clock. The computer clock is less accurate and has more jitter than the onboard clock, and is generally limited to sampling rates below 500 Hz. The computer clock is referred to as the *software clock* in this guide.

External Clock. An external clock is often used when the sampling rate is low and not constant. For example, an external clock source is often used in automotive applications where samples are acquired as a function of crank angle.

Channel Configuration

You can configure input channels in one of these two ways:

- Differential
- Single-ended

Your choice of input channel configuration may depend on whether the input signal is *floating* or *grounded*.

A floating signal uses an isolated ground reference and is not connected to the building ground. As a result, the input signal and hardware device are not connected to a common reference, which can cause the input signal to exceed the valid range of the hardware device. To circumvent this problem, you must connect the signal to the onboard ground of the device. Examples of floating signal sources include ungrounded thermocouples and battery devices.

A grounded signal is connected to the building ground. As a result, the input signal and hardware device are connected to a common reference. Examples of grounded signal sources include nonisolated instrument outputs and devices that are connected to the building power system.


Note For more information about channel configuration, refer to your hardware documentation.

Differential Inputs

When you configure your hardware for differential input, there are two signal wires associated with each input signal – one for the input signal and one for the reference (return) signal. The measurement is the difference in voltage between the two wires, which helps reduce noise and any voltage that is common to both wires.

As shown below, the input signal is connected to the positive amplifier socket (labeled +) and the return signal is connected to the negative amplifier socket

(labeled -). The amplifier has a third connector that allows these signals to be referenced to ground.




National Instruments recommends that you use differential inputs under any of these conditions:

- The input signal is low level (less than 1 volt).
- The leads connecting the signal are greater than 10 feet.
- The input signal requires a separate ground-reference point or return signal.
- The signal leads travel through a noisy environment.

Single-Ended Inputs

When you configure your hardware for single-ended input, there is one signal wire associated with each input signal, and each input signal is connected to the same ground. Single-ended measurements are more susceptible to noise than differential measurements due to differences in the signal paths.

As shown below, the input signal is connected to the positive amplifier socket (labeled +) and the ground is connected to the negative amplifier socket (labeled -).



National Instruments suggests that you can use single-ended inputs under any of these conditions:

- The input signal is high level (greater than 1 volt).
- The leads connecting the signal are less than 10 feet.
- The input signal can share a common reference point with other signals.

You should use differential input connectors for any input signal that does not meet the preceding conditions. You can configure many National Instruments boards for two different types of single-ended connections:

- Referenced single-ended (RSE) connection

The RSE configuration is used for floating signal sources. In this case, the hardware device itself provides the reference ground for the input signal.

- Nonreferenced single-ended (NRSE) connection

The NRSE input configuration is used for grounded signal sources. In this case, the input signal provides its own reference ground and the hardware device should not supply one.

Refer to your National Instruments hardware documentation for more information about RSE and NRSE connections.

Transferring Data from Hardware to System Memory

The transfer of acquired data from the hardware to system memory follows these steps:

- 1 Acquired data is stored in the hardware's first-in first-out (FIFO) buffer.
- 2 Data is transferred from the FIFO buffer to system memory using interrupts or DMA.

These steps happen automatically. Typically, all that's required from you is some initial configuration of the hardware device when it is installed.

The FIFO Buffer

The FIFO buffer is used to temporarily store acquired data. The data is temporarily stored until it can be transferred to system memory. The process of transferring data into and out of an analog input FIFO buffer is given below:

- 1 The FIFO buffer stores newly acquired samples at a constant sampling rate.
- 2 Before the FIFO buffer is filled, the software starts removing the samples. For example, an interrupt is generated when the FIFO is half full, and signals the software to extract the samples as quickly as possible.
- 3 Since servicing interrupts or programming the DMA controller can take up to a few milliseconds, additional data is stored in the FIFO for future retrieval. For a larger FIFO buffer, longer latencies can be tolerated.
- 4 The samples are transferred to system memory via the system bus (for example, PCI bus or AT bus). After the samples are transferred, the software is free to perform other tasks until the next interrupt occurs. For example, the data can be processed or saved to a disk file. As long as the average rates of storing and extracting data are equal, acquired data will not be missed and your application should run smoothly.

Interrupts

The slowest but most common method to move acquired data to system memory is for the board to generate an interrupt request (IRQ) signal. This signal can be generated when one sample is acquired or when multiple samples

are acquired. The process of transferring data to system memory via interrupts is given below:

- 1 When data is ready for transfer, the CPU stops whatever it is doing and runs a special interrupt handler routine that saves the current machine registers, and then sets them to access the board.
- 2 The data is extracted from the board and placed into system memory.
- 3 The saved machine registers are restored, and the CPU returns to the original interrupted process.

The actual data move is fairly quick, but there is a lot of overhead time spent saving, setting up, and restoring the register information. Therefore, depending on your specific system, transferring data by interrupts may not be a good choice when the sampling rate is greater than around 5 kHz.

DMA

Direct memory access (DMA) is a system whereby samples are automatically stored in system memory while the processor does something else. The process of transferring data via DMA is given below:

- 1 When data is ready for transfer, the board directs the system DMA controller to put it into in system memory as soon as possible.
- 2 As soon as the CPU is able (which is usually very quickly), it stops interacting with the data acquisition hardware and the DMA controller moves the data directly into memory.
- 3 The DMA controller gets ready for the next sample by pointing to the next open memory location.
- 4 The previous steps are repeated indefinitely, with data going to each open memory location in a continuously circulating buffer. No interaction between the CPU and the board is needed.

Your computer supports several different DMA channels. Depending on your application, you can use one or more of these channels. For example, simultaneous input and output with a sound card requires one DMA channel for the input and another DMA channel for the output.

Making Quality Measurements

For most data acquisition applications, you need to measure the signal produced by a sensor at a specific rate.

In many cases, the sensor signal is a voltage level that is proportional to the physical phenomena of interest (for example, temperature, pressure, or acceleration). If you are measuring slowly changing (quasi-static) phenomena like temperature, a slow sampling rate usually suffices. If you are measuring rapidly changing (dynamic) phenomena like vibration or acoustic measurements, a fast sampling rate is required.

To make high-quality measurements, you should follow these rules:

- Maximize the precision and accuracy
- Minimize the noise
- Match the sensor range to the A/D range


Accuracy and Precision

Whenever you acquire measured data, you should make every effort to maximize its accuracy and precision. The quality of your measurement depends on the accuracy and precision of the entire data acquisition system, and can be limited by such factors as board resolution or environmental noise.


In general terms, the *accuracy* of a measurement determines how close the measurement comes to the true value. Therefore, it indicates the correctness of the result. The *precision* of a measurement reflects how exactly the result is determined without reference to what the result means. The *relative precision* indicates the uncertainty in a measurement as a fraction of the result.

For example, suppose you measure a table top with a meter stick and find its length to be 1.502 meters. This number indicates that the meter stick (and your eyes) can resolve distances down to at least a millimeter. Under most circumstances, this is considered to be a fairly precise measurement with a relative precision of around 1/1500. However, suppose you perform the measurement again and obtain a result of 1.510 meters. After careful consideration, you discover that your initial technique for reading the meter stick was faulty because you did not read it from directly above. Therefore, the first measurement was not accurate.


Precision and accuracy are illustrated below.




Not precise
Not accurate



Precise
Not accurate



Not precise
Accurate



Precise
Accurate

For analog input subsystems, accuracy is usually limited by calibration errors while precision is usually limited by the A/D converter. Accuracy and precision are discussed in more detail below.

Accuracy

Accuracy is defined as the agreement between a measured quantity and the true value of that quantity. Every component that appears in the analog signal path affects system accuracy and performance. The overall system accuracy is given by the component with the worst accuracy.

For data acquisition hardware, accuracy is often expressed as a percent or a fraction of the least significant bit (LSB). Under ideal circumstances, board accuracy is typically ± 0.5 LSB. Therefore, a 12-bit converter has only 11 usable bits.

Many boards include a programmable gain amplifier, which is located just before the converter input. To prevent system accuracy from being degraded, the accuracy and linearity of the gain must be better than that of the A/D converter. The specified accuracy of a board is also affected by the sampling rate and the *settling time* of the amplifier. The settling time is defined as the time required for the instrumentation amplifier to settle to a specified accuracy. To maintain full accuracy, the amplifier output must settle to a level given by the magnitude of 0.5 LSB before the next conversion, and is on the order of several tenths of a millisecond for most boards.

Settling time is a function of sampling rate and gain value. High rate, high gain configurations require longer settling times while low rate, low gain configurations require shorter settling times.

Precision

The number of bits used to represent an analog signal determines the precision (resolution) of the device. The more bits provided by your board, the more precise your measurement will be. A high precision, high resolution device divides the input range into more divisions thereby allowing a smaller detectable voltage value. A low precision, low resolution device divides the input range into fewer divisions thereby increasing the detectable voltage value.

The overall precision of your data acquisition system is usually determined by the A/D converter, and is specified by the number of bits used to represent the analog signal. Most boards use 12 or 16 bits. The precision of your measurement is given by

$$\text{Precision} = \text{one part in } 2^{(\text{number of bits})}$$

The precision in volts is given by

$$\text{Precision} = \text{voltage range} / 2^{(\text{number of bits})}$$

For example, if you are using a 12-bit A/D converter configured for a 10 volt range, then

$$\text{Precision} = 10 \text{ volts} / 2^{(12)}$$

This means that the converter can detect voltage differences at the level of 0.00244 volts (2.44 mV).

How Are Range, Gain, and Measurement Precision Related?

When you configure the input range and gain of your analog input subsystem, the end result should maximize the measurement resolution and minimize the chance of an overrange condition. The actual input range is given by the formula

$$\text{Actual input range} = \text{Input range} / \text{Gain}$$

The relationship between gain, actual input range, and precision for a unipolar and bipolar signal having an input range of 10 V is shown below.

Table 1-2: Relationship Between Input Range, Gain, and Precision

Input Range	Gain	Actual Input Range	Precision (12 Bit A/D)
0 to 10 V	1.0	0 to 10 V	2.44 mV
	2.0	0 to 5 V	1.22 mV
	5.0	0 to 2 V	0.488 mV
	10.0	0 to 1 V	0.244 mV
-5 to 5 V	0.5	-10 to 10 V	4.88 mV
	1.0	-5 to 5 V	2.44 mV
	2.0	-2.5 to 2.5 V	1.22 mV
	5.0	-1.0 to 1.0 V	0.488 mV
	10.0	-0.5 to 0.5 V	0.244 mV

As shown in the table, the gain affects the precision of your measurement. If you select a gain that decreases the actual input range, then the precision increases. Conversely, if you select a gain that increases the actual input range, then the precision decreases. This is because the actual input range varies but the number of bits used by the A/D converter remains fixed.

Note With the Data Acquisition Toolbox, you do not have to specify the range and gain. Instead, you simply specify the actual input range desired.

Noise

Noise is considered to be any measurement that is not part of the phenomena of interest. Noise can be generated within the electrical components of the input amplifier (internal noise), or it can be added to the signal as it travels down the input wires to the amplifier (external noise). Techniques that you can use to reduce the effects of noise are described below.

Removing Internal Noise

Internal noise arises from thermal effects in the amplifier. Amplifiers typically generate a few microvolts of internal noise, which limits the resolution of the signal to this level. The amount of noise added to the signal depends on the bandwidth of the input amplifier.

To reduce internal noise, you should select an amplifier with a bandwidth that closely matches the bandwidth of the input signal.

Removing External Noise

External noise arises from many sources. For example, many data acquisition experiments are subject to 60 Hz noise generated by a.c. power circuits. This type of noise is referred to as *pick-up* or *hum*, and appears as a sinusoidal interference signal in the measurement circuit. Another common interference source is fluorescent lighting. These lights generate an arc at twice the power line frequency (120 Hz).

Noise is added to the acquisition circuit from these external sources because the signal leads act as aerials picking up environmental electrical activity.

Much of this noise is common to both signal wires. To remove most of this common-mode voltage, you should:

- Configure the input channels in differential mode. Refer to “Channel Configuration” on page 1-24 for more information about channel configuration.
- Use signal wires that are twisted together rather than separate.
- Keep the signal wires as short as possible.
- Keep the signal wires as far away as possible from environmental electrical activity.

Filtering

Filtering also reduces signal noise. For many data acquisition applications, a low-pass filter is beneficial. As the name suggests, a low-pass filter passes the lower frequency components but attenuates the higher frequency components. The cut-off frequency of the filter must be compatible with the frequencies present in the signal of interest and the sampling rate used for the A/D conversion.

A low-pass filter that's used to prevent higher frequencies from introducing distortion into the digitized signal is known as an antialiasing filter if the cut-off occurs at the Nyquist frequency. That is, the filter removes frequencies greater than one-half the sampling frequency. These filters generally have a sharper cut-off than the normal low-pass filter used to condition a signal. Antialiasing filters are specified according to the sampling rate of the system and there must be one filter per input signal.

Matching the Sensor Range and A/D Converter Range

When sensor data is digitized by an A/D converter, you must be aware of these two issues:

- The expected range of the data produced by your sensor. This range depends on the physical phenomena you are measuring and the output range of the sensor.
- The range of your A/D converter. For many devices, the hardware range is specified by the gain and polarity.

You should select the sensor and hardware ranges such that the maximum precision is obtained, and the full dynamic range of the input signal is covered.

For example, suppose you are using a microphone with a dynamic range of 20 dB to 140 dB and an output sensitivity of 50 mV/Pa. If you are measuring street noise in your application, then you might expect that the sound level never exceeds 80 dB, which corresponds to a sound pressure magnitude of 200 mPa and a voltage output from the microphone of 10 mV. Under these conditions, you should set the input range of your data acquisition card for a maximum signal amplitude of 10 mV, or a little more.

How Fast Should a Signal be Sampled?

Whenever a continuous signal is sampled, some information is lost. The key objective is to sample at a rate such that the signal of interest is well characterized and the amount of information lost is minimized.


If you sample at a rate that is too slow, then the signal is *undersampled*, and *aliasing* can occur. Aliasing can occur for both rapidly varying signals and slowly varying signals. For example, suppose you are measuring temperature once a minute. If your acquisition system is picking up a 60 Hz hum from an a.c power supply, then that hum will appear as constant noise level if you are sampling at 30 Hz.

Aliasing occurs when the sampled signal contains frequency components greater than one-half the sampling rate. The frequency components could originate from the signal of interest in which case you are undersampling and should increase the sampling rate. The frequency components could also originate from noise in which case you may need to condition the signal using a filter. The rule used to prevent aliasing is given by the *Nyquist theorem*, which states that:

- An analog signal can be uniquely reconstructed, without error, from samples taken at equal time intervals.
- The sampling rate must be equal to or greater than twice the highest frequency component in the analog signal. A frequency of one-half the sampling rate is called the Nyquist frequency.

However, if your input signal is corrupted by noise, then aliasing can still occur.

For example, suppose you configure your A/D converter to sample at a rate of 4 samples per second (4 S/s or 4 Hz), and the signal of interest is a 1 Hz sine wave. Since the signal frequency is one-fourth the sampling rate, then according to the Nyquist theorem, it should be completely characterized. However, if a 5 Hz sine wave is also present, then these two signals cannot be distinguished. In other words, the 1 Hz sine wave produces the same samples as the 5 Hz sine wave when the sampling rate is 4 S/s. This situation is shown below.





In a real-world data acquisition environment, you may need to condition the signal by filtering out the high frequency components.

Even though the samples appear to represent a sine wave with a frequency of one-fourth the sampling rate, the actual signal could be any sine wave with a frequency of

$$(n \pm 0.25) \times (\text{Sampling rate})$$

where n is zero or any positive integer. For this example, the actual signal may be at a frequency of 3 Hz, 5 Hz, 7 Hz, 9 Hz, and so on. The relationship $0.25 \times (\text{Sampling rate})$ is called the *alias* of a signal that may be at another frequency. In other words, aliasing occurs when one frequency assumes the identity of another frequency.

If you sample the input signal at least twice as fast as the highest frequency component, then that signal may be uniquely characterized but this rate would not mimic the waveform very closely. As shown below, to get an accurate picture of the waveform, you need a sampling rate of roughly 10 to 20 times the highest frequency.



As shown in the top figure, the low sampling rate produces a sampled signal that appears to be a triangular waveform. As shown in the bottom figure, a higher fidelity sampled signal is produced when the sampling rate is higher. In the latter case, the sampled signal actually looks like a sine wave.

How Can Aliasing be Eliminated?

The primary considerations involved in antialiasing are the sampling rate of the A/D converter and the frequencies present in the sampled data. To eliminate aliasing, you must:

- Establish the useful bandwidth of the measurement.
- Select a sensor with sufficient bandwidth.
- Select a low-pass anti-aliasing analog filter that can eliminate all frequencies exceeding this bandwidth.
- Sample the data at a rate at least twice that of the filter's upper cutoff frequency.

Selected Bibliography


- 1 *Transducer Interfacing Handbook – A Guide to Analog Signal Conditioning*, edited by Daniel H. Sheingold; Analog Devices Inc., Norwood, Massachusetts, 1980.
- 2 Bentley, John P., *Principles of Measurement Systems, Second Edition*; Longman Scientific and Technical, Harlow, Essex, UK, 1988.
- 3 Bevington, Philip R., *Data Reduction and Error Analysis for the Physical Sciences*, McGraw-Hill, New York, NY, 1969.
- 4 Carr, Joseph J., *Sensors*; Prompt Publications, Indianapolis, Indiana, 1997.
- 5 *The Measurement, Instrumentation, and Sensors Handbook*, edited by John G. Webster; CRC Press, Boca Raton, FL, 1999.
- 6 *PCI-MIO E Series User Manual, January 1997 Edition*; Part Number 320945B-01, National Instruments, Austin, TX, 1997.

Getting Started with the Data Acquisition Toolbox

Toolbox Components	2-2
M-File Functions	2-3
The Data Acquisition Engine	2-4
The Hardware Driver Adaptor	2-7
 Accessing Your Hardware	 2-8
Acquiring Data	2-8
Outputting Data	2-9
Reading and Writing Digital Values	2-10
 Understanding the Toolbox Capabilities	 2-12
The Contents and Readme Files	2-12
Documentation Examples	2-12
The Quick Reference Guide	2-13
Demos	2-13
 Examining Your Hardware Resources	 2-17
General Toolbox Information	2-17
Adaptor-Specific Information	2-18
Device Object Information	2-19
 Getting Help	 2-20
The daqhelp Function	2-20
The propinfo Function	2-21

Toolbox Components

The Data Acquisition Toolbox consists of three distinct components: M-file functions, the data acquisition engine, and hardware driver adaptors. As shown below, these components allow you to pass information between MATLAB and your data acquisition hardware.



The preceding diagram illustrates how information flows from component to component. Information consists of:

- **Property values**

You can control the behavior of your data acquisition application by configuring property values. In general, you can think of a property as a characteristic of the toolbox or of the hardware driver that can be manipulated to suit your needs.

- **Data**

You can acquire data from a sensor connected to an analog input subsystem and store it in MATLAB, or output data from MATLAB to an actuator connected to an analog output subsystem. Additionally you can transfer values (1's and 0's) between MATLAB and a digital I/O subsystem.

- **Events**

An event occurs at a particular time after a condition is met and may result in one or more actions that you specify. Events can be generated only after you configure the associated properties. Some of the ways you can use events include initiating analysis after a predetermined amount of data is acquired, or displaying a message to the MATLAB workspace after an error occurs.

M-File Functions

To perform any task with your data acquisition application, you must call M-file functions from the MATLAB environment. Among other things, these functions allow you to:

- Create device objects, which provide a gateway to your hardware's capabilities and allow you to control the behavior of your application.
- Acquire or output data.
- Configure property values.
- Evaluate your acquisition status and hardware resources.

For a listing of all Data Acquisition Toolbox functions, refer to Chapter 9, "Function Reference." You can also display all the toolbox functions by typing

```
hel p daq
```

The Data Acquisition Engine

The data acquisition engine (or just *engine*) is a MEX-file dynamic link library (DLL) file that:

- Stores the device objects and associated property values that control your data acquisition application
- Controls the synchronization of events
- Controls the storage of acquired or queued data

While the engine performs these tasks, you can use MATLAB for other tasks such as analyzing acquired data. In other words, the engine and MATLAB are *asynchronous*. The relationship between acquiring data, outputting data, and data flow is described below.

The Flow of Acquired Data


Acquiring data means that data is flowing from your hardware device into the data acquisition engine where it is temporarily stored in memory. The data is stored temporarily because it can be overwritten. The rate at which the data is overwritten depends on several factors including the available memory, the rate at which data is acquired, and the number of hardware channels from which data is acquired.

The stored data is not automatically available in the MATLAB workspace. Instead, you must explicitly extract data from the engine using the `getdata` function.

The flow of acquired data consists of these two independent steps:

- 1 Data acquired from the hardware is stored in the engine.
- 2 Data is extracted from the engine and stored in MATLAB, or output to a disk file.

These two steps are illustrated below.




The Flow of Output Data

Outputting data means that data is flowing from the data acquisition engine to the hardware device. However, before data is output, you must queue it in the engine with the `putdata` function. The amount of data that you can queue depends on several factors including the available memory, the number of hardware channels to which data is output, and the size of each data sample.

The flow of output data consists of these two independent steps:

- 1 Data from MATLAB is queued in the engine.
- 2 Data queued in the engine is output to the hardware.

These two steps are illustrated below.



The Hardware Driver Adaptor

The hardware driver adaptor (or just *adaptor*) is the interface between the data acquisition engine and the hardware driver. The adaptor's main purpose is to pass information between MATLAB and your hardware device via the driver.

Hardware drivers are provided by your device vendor. For example, if you are acquiring data using a National Instruments E Series board, then the appropriate version of the NI-DAQ driver must be installed on your platform. Note that hardware drivers are not installed as part of the Data Acquisition Toolbox. The supported vendors and the adaptor names used by the toolbox are listed below.

Table 2-1: Supported Vendors and Adaptor Names

Vendor	Adaptor Name
National Instruments	ni daq
ComputerBoards	cbi
Agilent Technologies	hpe1432
Windows sound cards	wi nsound

Note You can use the Data Acquisition Toolbox Adaptor Kit to interface unsupported hardware devices to the toolbox. Refer to the MathWorks Web site at <http://www.mathworks.com/products/daq/> for more information about the adaptor kit.

As described in “Examining Your Hardware Resources” on page 2-17, you can list the supported adaptor names with the `daqhwinfo` function. For a list of vendor driver requirements and limitations, refer to “Data Acquisition Toolbox 2.0” in the *Release Notes*.

Accessing Your Hardware

Perhaps the most effective way to get started with the Data Acquisition Toolbox is to connect to your hardware, and input or output data. This section provides simple examples that show you how to:

- Acquire data from analog input channels
- Output data to analog output channels
- Read values from and write values to digital I/O lines

Each example illustrates a typical *data acquisition session*. The data acquisition session comprises all the steps you are likely to take when acquiring or outputting data using a supported hardware device. You should keep these steps in mind when constructing your own data acquisition applications.

Note that the analog input and analog output examples use a sound card, while the digital I/O example uses a National Instruments PCI-6024E board. If you are using a different supported hardware device, you should modify the adaptor name and the device ID supplied to the creation function as needed.

If you want detailed information about any functions that are used, refer to Chapter 9, “Function Reference.” If you want detailed information about any properties that are used, refer to Chapter 10, “Base Property Reference.”

Acquiring Data

If you have a sound card installed, you can run the following example, which acquires one second of data from two analog input hardware channels, and then plots the acquired data.

You should modify this example to suit your specific application needs. If you want detailed information about acquiring data, refer to Chapter 5, “Doing More with Analog Input.”

1. Create a device object – Create the analog input object `ai` for a sound card.

```
ai = analoginput('winsound');
```

2. Add channels – Add two hardware channels to `ai`.

```
addchannel(ai, 1:2);
```

3. Configure property values – Configure the sampling rate to 44.1 kHz and collect 1 second of data (44,100 samples) for each channel.

```
set(ai, 'SampleRate', 44100)
set(ai, 'SamplesPerTrigger', 44100)
```

4. Acquire data – Start the acquisition. When all the data is acquired, ai automatically stops executing.

```
start(ai)
data = getdata(ai);
plot(data)
```

5. Clean up – When you no longer need ai, you should remove it from memory and from the MATLAB workspace.

```
delete(ai)
clear ai
```

Outputting Data

If you have a sound card installed, you can run the following example, which outputs 1 second of data to two analog output hardware channels.

You should modify this example to suit your specific application needs. If you want detailed information about outputting data, refer to Chapter 6, “Analog Output.”

1. Create a device object – Create the analog output object ao for a sound card.

```
ao = analogoutput('winsound');
```

2. Add channels – Add two hardware channels to ao.

```
addchannel(ao, 1:2);
```

3. Configure property values – Configure the sampling rate to 44.1 kHz for each channel.

```
set(ao, 'SampleRate', 44100)
```

4. Output data – Create 1 second of output data, and queue the data in the engine for eventual output to the analog output subsystem. You must queue one column of data for each hardware channel added.

```
data = sin(linspace(0, 2*pi, 44100)');  
putdata(ao, [data data])
```

Start the output. When all the data is output, ao automatically stops executing.

```
start(ao)
```

Block the MATLAB command line until ao stops running.

```
while strcmp(ao.Running, 'On')  
end
```

5. Clean up – When you no longer need ao, you should remove it from memory and from the MATLAB workspace.

```
delete(ao)  
clear ao
```

Reading and Writing Digital Values

If you have a supported National Instruments board with at least eight digital I/O lines, you can run the following example, which outputs digital values, and then reads back those values.

You should modify this example to suit your specific application needs. If you want detailed information about reading and writing digital values, refer to Chapter 7, “Digital Input/Output.”

1. Create a device object – Create the digital I/O object di o for a National Instruments PCI-6024E board with hardware ID 1.

```
di o = digital i o('ni daq', 1);
```

2. Add lines – Add eight hardware lines to di o, and configure them for output.

```
addline(di o, 0: 7, 'out');
```

3. Read and write values – Create an array of output values, and write the values to the digital I/O subsystem. Note that reading and writing digital I/O line values typically does not require that you configure specific property values.

```
pval = [1 1 1 1 0 1 0 1];  
putvalue(dio, pval)  
gval = getvalue(dio);
```

4. Clean up – When you no longer need `dio`, you should remove it from memory and from the MATLAB workspace.

```
delete(dio)  
clear dio
```

Note Digital line values are usually not transferred at a specific rate. Although some specialized boards support clocked I/O, the Data Acquisition Toolbox does not support this functionality.

Understanding the Toolbox Capabilities

In addition to the printed and online documentation, the Data Acquisition Toolbox provides these resources to help you understand the product capabilities:

- The Contents and Readme files
- Documentation examples
- The Quick Reference Guide
- Demos

The Contents and Readme Files

The Contents M-file lists the toolbox functions and demos. You can display this information by typing

```
hel p daq
```

The Readme HTML file contains toolbox modifications and bug fixes. You can display this information by typing

```
whatsnew daq
```

Documentation Examples

This guide provides detailed examples that show you how to acquire or output data. These examples are collected in the example index, which is available through the Help browser

Some examples are constructed as mini-applications that illustrate one or two important features of the toolbox and serve as templates so you can see how to build applications that suit your specific needs. These examples are included as toolbox M-files and are treated as demos. You can list all Data Acquisition Toolbox demos by typing

```
hel p daqdemos
```

All documentation example M-files begin with daqdoc. To run an example, type the M-file name at the command line. Note that most analog input (AI) and analog output (AO) examples are written for sound cards. To use these examples with your hardware device, you should modify the adaptor name and the device ID supplied to the creation function as needed.

Additionally, most documentation examples are written for clocked subsystems. However, some supported hardware devices – particularly ComputerBoards devices – do not possess onboard clocks. If the AI or AO subsystem of your hardware device does not have an onboard clock, then these examples will not work. To use the documentation examples, you can:

- Input single values using the `getsample` function, or output single values using the `putsample` function.
- Configure the `ClockSource` property to `Software`.

The Quick Reference Guide

The Quick Reference Guide gives you a brief but complete overview of the toolbox capabilities, functions, and properties. You may find it useful to print this guide and keep it handy when using the toolbox. You can access this guide through the Help browser.

Demos

The toolbox includes a large collection of demos, which are available through the `daqschool` interface. To launch this interface, type `daqschool` at the command line.

```
daqschool
```


You can also launch individual demos by typing the appropriate M-file name at the command line. Alternatively, you can launch a toolbox demo through the MATLAB Demo window by typing

```
demo toolbox 'Data Acquisition'
```

The data acquisition demos, organized by device object, are listed below.

Note Most analog input and analog output demos require that you have a sound card installed. The digital I/O demos require that you have a supported National Instruments board with digital I/O capabilities.

daqschool is shown below.



Common Demos

The common demos illustrate features that are common to all supported device objects. These demos are listed below.

Demo Name	Description
demodaq_action	Introduction to action functions.
demodaq_intro	Introduction to the Data Acquisition Toolbox.
demodaq_save	Save and load device objects.
daqtimerplot	Example action function that plots acquired data.

Analog Input Demos

The analog input demos are listed below.

Demo Name	Description
daqrecord	Record data from the specified adaptor.
daqscope	Example oscilloscope.
demoai_channel	Introduction to analog input channels.
demoai_fft	Fast Fourier transform (FFT) of an incoming signal.
demoai_intro	Introduction to analog input objects.
demoai_logging	Demonstrate logging data to a disk file.
demoai_trigger	Demonstrate the use of immediate, manual, and software triggers.

Analog Output Demos

The analog output demos are listed below.

Demo Name	Description
daqfngen	Example function generator.
daqplay	Output data to the specified adaptor.
daqsong	Output data from HANDEL.MAT to a sound card.
demoao_channel	Introduction to analog output channels.
demoao_intro	Introduction to analog output objects.
demoao_trigger	Demonstrate the use of immediate and manual triggers.

Digital I/O Demos

The digital I/O demos are listed below.

Demo Name	Description
demodio_intro	Introduction to digital I/O objects.
demodio_line	Introduction to digital I/O lines.
diopanel	Example panel for transferring values between MATLAB and a digital I/O subsystem.

Examining Your Hardware Resources

You can examine the data acquisition hardware resources visible to the toolbox with the `daqhwinfo` function. Hardware resources include installed boards, hardware drivers, and adaptors. The information returned by `daqhwinfo` depends on the supplied arguments, and is divided into these three categories:

- General toolbox information
- Adaptor-specific information
- Device object information

If you configure hardware parameters using a vendor tool such as National Instruments' Measurement and Automation Explorer or ComputerBoards' InstaCal, then `daqhwinfo` will return this configuration information. For example, if you configure your ComputerBoards device for 16 single-ended channels using InstaCal, then `daqhwinfo` returns this hardware configuration. However, the toolbox does not preserve configuration information that is not directly associated with your hardware. For example, channel name information is not preserved. Refer to Appendix A, "Troubleshooting Your Hardware" for more information about using vendor tools.

General Toolbox Information

To display general information about the Data Acquisition Toolbox

```
out = daqhwinfo
out =
    ToolboxName: 'Data Acquisition Toolbox'
    ToolboxVersion: [1x39 char]
    MatlabVersion: '6.0.0.74976 (R12)'
    InstalledAdaptors: {4x1 cell}
```

The `InstalledAdaptors` field lists the hardware driver adaptors installed on your system. To display the installed adaptors

```
out.InstalledAdaptors
ans =
    'cbi'
    'hpe1432'
    'ni daq'
    'winsound'
```

This information tells you that an adaptor is available for ComputerBoards, Agilent Technologies, National Instruments, and sound card devices.

Note Toolbox adaptors are available to you only if the associated hardware driver is installed on your platform.

Adaptor-Specific Information

To display hardware information for a particular vendor, you must supply the adaptor name as an argument to `daqhwinfo`. The supported vendors and adaptor names are given in “The Hardware Driver Adaptor” on page 2-7. For example, to display hardware information for the `winSound` adaptor

```
out = daqhwinfo('winSound')
out =
    AdaptorDllName: 'd:\v6\toolbox\daq\daq\private\mmwinSound.dll'
    AdaptorDllVersion: 'Version 2.0 (R12) 05-Oct-2000'
    AdaptorName: 'winSound'
    BoardNames: {'AudioPCI Record'}
    InstalledBoardIDs: {'0'}
    ObjectConstructorName: {'analoginput('winSound',0) [1x26 char]}
```

The `ObjectConstructorName` field lists the subsystems supported by the installed sound cards, and the syntax for creating a device object associated with a given subsystem. To display the device object constructor names available for the `AudioPCI Record` board

```
out.ObjectConstructorName(:)
ans =
    'analoginput('winSound',0)'
    'analogoutput('winSound',0)'
```

This information tells you that the sound card supports analog input and analog output objects. To create an analog input object for the sound card

```
ai = analoginput('winSound');
```

To create an analog output object for the sound card

```
ao = analogoutput('winSound');
```

Device Object Information

To display hardware information for a specific device object, you must supply the device object as an argument to `daqhwinfo`. The hardware information for the analog input object `ai` created in the preceding section is given below.

```
out = daqhwinfo(ai)
out =
    AdaptorName: 'windsound'
    Bits: 16
    Coupling: {'AC Coupled'}
    DeviceName: 'AudioPCI Record'
    DifferentialIDs: []
    Gains: []
    ID: '0'
    InputRanges: [-1 1]
    MaxSampleRate: 44100
    MinSampleRate: 8000
    NativeDataType: 'int16'
    Polarity: {'Bipolar'}
    SampleType: 'SimultaneousSample'
    SingleEndedIDs: [1 2]
    SubsystemType: 'AnalogInput'
    TotalChannels: 2
    VendorDriverDescription: 'Windows Multimedia Driver'
    VendorDriverVersion: '1.0'
```

Among other things, this information tells you that the minimum sampling rate is 8 kHz, the maximum sampling rate is 44.1 kHz, and there are two hardware channels that you can add to the analog input object.

Getting Help

The Data Acquisition Toolbox provides you with these help resources:

- The HTML and PDF versions of this guide, and the Quick Reference Guide, which are available through the Help browser
- M-file function help, which you can display with the `help` command (since some toolbox functions are overloaded, you may need to specify the appropriate pathname as well)
- The `daqhelp` function
- The `property` function

The `daqhelp` Function

You can use the `daqhelp` function to:

- Display command line help for functions and properties
- List all the functions and properties associated with a specific device object

A device object need not exist for you to obtain this information. For example, to display all the functions and properties associated with an analog input object, as well as the constructor M-file help

```
daqhelp analoginput
```

To display help for the `SampleRate` property

```
daqhelp SampleRate
```

You can also display help for an existing device object. For example, to display help for the `BitsPerSample` property for an analog input object associated with a sound card

```
ai = analoginput('winsound');  
out = daqhelp(ai, 'BitsPerSample');
```

The propinfo Function

You can use the `propinfo` function to return the characteristics of Data Acquisition Toolbox properties. For example, you can find the default value for any property using this function. `propinfo` returns a structure containing the fields shown below.

Table 2-2: `propinfo` Fields

Field Name	Description
Type	The property data type. Possible values are <code>acti on functi on</code> , <code>any</code> , <code>doubl e</code> , and <code>stri ng</code> .
Constrai nt	The type of constraint on the property value. Possible values are <code>acti on functi on</code> , <code>bounded</code> , <code>enum</code> , and <code>none</code> .
Constrai ntVal ue	The property value constraint. The constraint can be a range of valid values or a list of valid string values.
Defaul tVal ue	The property default value.
ReadOnl y	If the property is read-only, a 1 is returned. Otherwise, a 0 is returned.
ReadOnl yRunni ng	If the property is read-only while the device object is running, a 1 is returned. Otherwise, a 0 is returned.
Devi ceSpeci fi c	If the property is device-specific, a 1 is returned. If a 0 is returned, the property is supported for all device objects of a given type.

For example, to return the characteristics for all the properties associated with the analog input object `ai` created in the preceding section

```
AIinfo = propinfo(ai);
```

The characteristics for the `TriggerType` property are displayed below.

```
AIinfo.TriggerType
ans =
           Type: 'string'
      Constraint: 'Enum'
ConstraintValue: {3x1 cell}
      DefaultValue: 'Immediate'
          ReadOnly: 0
      ReadOnlyRunning: 1
      DeviceSpecific: 0
```

This information tells you that:

- The property value data type is a string.
- The property value is constrained as an enumerated list of values.
- There are three possible property values.
- The default value is `Immediate`.
- The property is not read-only.
- You cannot configure the property while the device object is running.
- The property is supported for all analog input objects.

To display the property value constraints

```
AIinfo.TriggerType.ConstraintValue
ans =
      'Manual'
      'Immediate'
      'Software'
```


The Data Acquisition Session

Overview	3-2
Example: The Data Acquisition Session	3-3
 Creating a Device Object	 3-4
Creating an Array of Device Objects	3-5
Where Do Device Objects Exist?	3-6
 Adding Channels or Lines	 3-8
Mapping Hardware Channel IDs to MATLAB Indices	3-9
 Configuring and Returning Properties	 3-12
Property Types	3-12
Returning Property Names and Property Values	3-14
Configuring Property Values	3-18
Specifying Property Names	3-19
Default Property Values	3-19
daqpropedit: A Graphical Property Editor	3-20
 Acquiring and Outputting Data	 3-22
Starting the Device Object	3-23
Logging or Sending Data	3-23
Stopping the Device Object	3-24
 Cleaning Up	 3-25

Overview

The data acquisition session comprises all the steps you are likely to take when acquiring or outputting data. These steps are:

- 1 **Create a device object** – You create a device object using the `analoginput`, `analogoutput`, or `digitalio` creation function. Device objects are the basic toolbox elements you use to access your hardware device.
- 2 **Add channels or lines** – After a device object is created, you must add channels or lines to it. Channels are added to analog input and analog output objects, while lines are added to digital I/O objects. Channels and lines are the basic hardware device elements with which you acquire or output data.
- 3 **Configure properties** – To establish the device object behavior, you assign values to properties using the `set` function or dot notation.

You can configure many of the properties at any time. However, some properties are configurable only when the device object is not running. Conversely, depending on your hardware settings and the requirements of your application, you may be able to accept the default property values and skip this step.

- 4 **Acquire or output data** – To acquire or output data, you must execute the device object with the `start` function. While the device object is running, it behaves according to the previously configured or default property values.

After data is acquired, you must extract it from the engine with the `getdata` function. Before you can output data, you must queue it in the engine with the `putdata` function.

- 5 **Clean up** – When you no longer need the device object, you should remove it from memory using the `delete` function, and remove it from the MATLAB workspace using the `clear` command.

The data acquisition session is used in many of the documentation examples included in this guide. Note that the fourth step is treated differently for digital I/O objects since they do not store data in the engine. Therefore, only analog input and analog output objects are discussed in this section.

Example: The Data Acquisition Session

This example illustrates the basic steps you take during a data acquisition session using an analog input object. You can run this example by typing `daqdoc3_1` at the MATLAB command line.

1. Create a device object – Create the analog input object `AI` for a sound card. The installed adaptors and hardware IDs are found with `daqhwinfo`.

```
AI = analoginput('winsound');
%AI = analoginput('ni daq', 1);
%AI = analoginput('cbi', 1);
```

2. Add channels – Add two channels to `AI`.

```
addchannel(AI, 1:2);
%addchannel(AI, 0:1); % For NI and CBI
```

3. Configure property values – Configure the sampling rate to 11.025 kHz and define a 2 second acquisition.

```
set(AI, 'SampleRate', 11025)
set(AI, 'SamplesPerTrigger', 22050)
```

4. Acquire data – Start `AI`, wait until the requested samples are acquired, and extract all the data from the engine. Before `start` is issued, you may want to begin inputting data from a microphone or a CD player.

```
start(AI)
while strcmp(AI.Running, 'On')
end
data = getdata(AI);
```

Plot the data and label the figure axes.

```
plot(data)
xlabel('Samples')
ylabel('Signal (Volts)')
```

5. Clean up – When you no longer need `AI`, you should remove it from memory and from the MATLAB workspace.

```
delete(AI)
clear AI
```

Creating a Device Object

Device objects are the toolbox components you use to access your hardware device. They provide a gateway to the functionality of your hardware, and allow you to control the behavior of your data acquisition application. Each device object is associated with a specific hardware subsystem.

To create a device object, you call M-file functions called *object creation functions* (or *object constructors*). These M-files are implemented using MATLAB's object-oriented programming capabilities, which are described in "MATLAB Classes and Objects" in the Help browser. The device object creation functions are listed below.

Table 3-1: Device Object Creation Functions

Function	Description
anal ogi nput	Create an analog input object.
anal ogout put	Create an analog output object.
di gi tal i o	Create a digital I/O object.

Before you can create a device object, the associated hardware driver adaptor must be registered. Adaptor registration occurs automatically. However, if for some reason an adaptor is not automatically registered, then you must do so manually with the `daqregister` function. Refer to "Registering the Hardware Driver Adaptor" on page A-19 for more information.


You can find out how to create device objects for a particular vendor and subsystem with the `ObjectConstructorName` field of the `daqhwinfo` function. For example, to find out how to create an analog input object for an installed National Instruments board, you supply the appropriate adaptor name to `daqhwinfo`.

```
out = daqhwinfo('ni daq');
out.ObjectConstructorName(:)
ans =
    'anal ogi nput (' ni daq', 1) '
    'anal ogout put (' ni daq', 1) '
    'di gi tal i o(' ni daq', 1) '
```

The constructor syntax tells you that you must supply the adaptor name and the hardware ID to the `analoginput` function

```
ai = analoginput('ni daq', 1);
```

The association between device objects and hardware subsystems is shown below.



Creating an Array of Device Objects

In MATLAB, you can create an array from existing variables by concatenating those variables together. The same is true for device objects. For example, suppose you create the analog input object `ai` and the analog output object `ao` for a sound card

```
ai = analoginput('winsound');
ao = analogoutput('winsound');
```

You can now create a device object array consisting of `ai` and `ao` using the usual MATLAB syntax. To create the row array `x`

```
x = [ai ao]
```

Index:	Subsystem:	Name:
1	Analog Input	winsound0-AI
2	Analog Output	winsound0-AO

To create the column array `y`

```
y = [ai; ao];
```

Note that you cannot create a matrix of device objects. For example, you cannot create the matrix

```
z = [ai ao; ai ao];
```

```
??? Error using ==> analoginput/vertcat
```

Only a row or column vector of device objects can be created.

Depending on your application, you may want to pass an array of device objects to a function. For example, using one call to the `set` function, you can configure both `ai` and `ao` to the same property value.

```
set(x, 'SampleRate', 44100)
```

Refer to Chapter 9, “Function Reference” to see which functions accept a device object array as an input argument.

Where Do Device Objects Exist?


When you create a device object, it exists in both the MATLAB workspace and the data acquisition engine. For example, suppose you create the analog input object `ai` for a sound card and then make a copy of `ai`.

```
ai = analoginput('winsound');  
newai = ai;
```

The copied device object `newai` is identical to the original device object `ai`. You can verify this by setting a property value for `ai` and returning the value of the same property from `newai`.

```
set(ai, 'SampleRate', 22050);  
get(newai, 'SampleRate')  
ans =  
    22050
```

As shown below, `ai` and `newai` return the same property value because they both reference the same device object in the data acquisition engine.



If you delete either the original device object or a copy, then the engine device object is also deleted. In this case, you cannot use any copies of the device object that remain in the workspace since they are no longer associated with any hardware. Device objects that are no longer associated with hardware are called *invalid objects*. The example below illustrates this situation.

```
delete(ai);
newai
newai =
Invalid Data Acquisition object.
This object is not associated with any hardware and
should be removed from your workspace using CLEAR.
```

You should remove invalid device objects from the workspace with the `clear` command.

Adding Channels or Lines

Channels and *lines* are the basic hardware device elements with which you acquire or output data

After you create a device object, you must add channels or lines to it. Channels are added to analog input and analog output objects, while lines are added to digital I/O objects. The channels added to a device object constitute a *channel group*, while the lines added to a device object constitute a *line group*.

The functions associated with adding channels or lines to a device object are listed below.

Table 3-2: Functions Associated with Adding Channels or Lines

Functions	Description
addchannel	Add hardware channels to an analog input or analog output object.
addline	Add hardware lines to a digital I/O object.
addmuxchannel	Add channels when using a National Instruments AMUX-64T multiplexer.

For example, to add two channels to an analog input object associated with a sound card, you must supply the appropriate hardware channel identifiers (IDs) to addchannel .


```
ai = analoginput('winsound');  
addchannel(ai, 1:2)
```

Note You cannot acquire or output data with a device object that does not contain channels or lines. Similarly, you cannot acquire or output data with channels or lines that are not contained by a device object.

You can think of a device object as a channel or line container that reflects the common functionality of a particular device. The common functionality of a device applies to all channels or lines that it contains. For example, the sampling rate of an analog input object applies to all channels contained by

that object. In contrast, the channels and lines contained by the device object reflect the functionality of a particular channel or line. For example, you can configure the input range (gain and polarity) on a per-channel basis.

The relationship between an analog input object and the channels it contains is shown below.



For digital I/O objects, the diagram would look the same except that lines would be substituted for channels.

Mapping Hardware Channel IDs to MATLAB Indices

When you add channels to a device object, the resulting channel group consists of a mapping between hardware channel IDs and MATLAB indices.

Hardware channel IDs are numeric values defined by the hardware vendor that uniquely identify a channel. For National Instruments and ComputerBoards hardware, the channel IDs are “zero-based” (begin at zero). For Agilent Technologies hardware and sound cards, the channel IDs are “one-based” (begin at one). However, when you reference channels, you use the MATLAB indices and not the hardware IDs. Given this, you should keep in mind that MATLAB is one-based. You can return the vendor’s hardware IDs with the `daqhwinfo` function.

For example, suppose you create the analog input object `ai` for a National Instruments board and you want to add the first three differential channels.

```
ai = analoginput('ni daq', 1);
```

To return the hardware IDs, supply the device object to `daqhwinfo`, and examine the `DifferentialIDs` field.

```
out = daqhwinfo(ai)
out.DifferentialIDs
ans =
      0      1      2      3      4      5      6      7
```


The first three differential channels have IDs 0, 1, and 2, respectively.

```
addchannel(ai, 0:2);
```

The index assigned to a hardware channel depends on the order in which you add it to the device object. In the above example, the channels are automatically assigned the MATLAB indices 1, 2, and 3, respectively. You can change the hardware channels associated with the MATLAB indices using the `HwChannel` property. For example, to swap the order of the second and third hardware channels

```
ai.Channel(2).HwChannel = 2;
ai.Channel(3).HwChannel = 1;
```

The original and modified index assignments are shown below.



Note If you are using scanning hardware, then the MATLAB indices define the scan order; index 1 is sampled first, index 2 is sampled second, and so on.

For digital I/O objects, the diagram would look the same except that lines would be substituted for channels.

Configuring and Returning Properties


You define and evaluate the behavior of your data acquisition application with device object properties. You define your application behavior by assigning values to properties with the set function or the dot notation. You evaluate your application configuration and status by displaying property values with the get function or the dot notation.

Property Types

Data Acquisition Toolbox properties are divided into two main types:

- **Common properties** – Common properties apply to every channel or line contained by a device object.
- **Channel/Line properties** – Channel/line properties are configured for individual channels or lines.

The relationship between an analog input object, the channels it contains, and their properties is shown below.




For digital I/O objects, the diagram would look the same except that lines would be substituted for channels.

Common properties and channel/line properties are subdivided into these two categories:

- **Base properties** – Base properties apply to all supported hardware subsystems of a given type, such as analog input. For example, the `SampleRate` property is supported for all analog input subsystems regardless of the vendor.
- **Device-specific properties** – Device-specific properties apply only to specific hardware devices. For example, the `BitsPerSample` property is supported only for sound cards. Note that base properties can have device-specific values. For example, the `InputType` property has a different set of values for each supported hardware vendor.

The relationship between common properties, channel/line properties, base properties, and device-specific properties is shown below.



For a complete description of all properties, refer to Chapter 10, “Base Property Reference” or Chapter 11, “Device-Specific Property Reference.”

Returning Property Names and Property Values

Once the device object is created, you can use the `set` function to return all configurable properties to a variable or to the command line. Additionally, if a property has a finite set of string values, then `set` also returns these values. You can use the `get` function to return one or more properties and their current values to a variable or to the command line.

The syntax used to return common and channel/line properties is described below. The examples are based on the analog input object `ai` created for a sound card and containing two channels.

```
ai = analoginput('winsound');  
addchannel(ai, 1:2);
```

Common Properties

To return all configurable common property names and their possible values for a device object, you must supply the device object to `set`. For example, all configurable common properties for `ai` are shown below. The base properties are listed first, followed by the device-specific properties.

```
set(ai)  
BufferingConfig  
BufferingMode: [ {Auto} | Manual ]  
Channel  
ChannelSkew  
ChannelSkewMode: [ {None} ]  
ClockSource: [ {Internal} ]  
DataMissedAction  
InputOverRangeAction  
InputType: [ {AC-Coupled} ]  
LogFileName  
LoggingMode: [ Disk | {Memory} | Disk&Memory ]  
LogToDiskMode: [ {Overwrite} | Index ]  
ManualTriggerHwOn: [ {Start} | Trigger ]  
Name  
RuntimeErrorAction  
SampleRate  
SamplesAcquiredAction  
SamplesAcquiredActionCount  
SamplesPerTrigger
```

```

StartAction
StopAction
Tag
Timeout
TimerAction
TimerPeriod
TriggerAction
TriggerChannel
TriggerCondition: [ {None} ]
TriggerConditionValue
TriggerDelay
TriggerDelayUnits: [ {Seconds} | Samples ]
TriggerRepeat
TriggerType: [ Manual | {Immediate} | Software ]
UserData

```

WINSOUND specific properties:

```

BitsPerSample
StandardSampleRates: [ Off | {On} ]

```

To return all common properties and their current values for a device object, you must supply the device object to get. For example, all common properties for ai are shown below. The base properties are listed first, followed by the device-specific properties.

```

get(ai)
    BufferingConfig = [512 30]
    BufferingMode = Auto
    Channel = [2x1 ai channel]
    ChannelSkew = 0
    ChannelSkewMode = None
    ClockSource = Internal
    DataMissedAction = daqaction
    EventLog = []
    InitialTriggerTime = [0 0 0 0 0 0]
    InputOverRangeAction =
    InputType = AC-Coupled
    LogFileName = logfile.daq
    Logging = Off
    LoggingMode = Memory
    LogToDiskMode = Overwrite

```

```
ManualTriggerHwOn = Start
Name = winsound0-AI
Running = Off
RuntimeErrorAction = daqaction
SampleRate = 8000
SamplesAcquired = 0
SamplesAcquiredAction =
SamplesAcquiredActionCount = 1024
SamplesAvailable = 0
SamplesPerTrigger = 8000
StartAction =
StopAction =
Tag =
Timeout = 1
TimerAction =
TimerPeriod = 0.1
TriggerAction =
TriggerChannel = [1x0 ai channel]
TriggerCondition = None
TriggerConditionValue = 0
TriggerDelay = 0
TriggerDelayUnits = Seconds
TriggerRepeat = 0
TriggersExecuted = 0
TriggerType = Immediate
Type = Analog Input
UserData = []

WINSOUND specific properties:
BitsPerSample = 16
StandardSampleRates = On
```

To display the current value for one property, you supply the property name to `get`.

```
get(ai, 'SampleRate')
ans =
    8000
```

To display the current values for multiple properties, you include the property names as elements of a cell array.


```
get(ai, {'StandardSampleRates', 'Running'})
ans =
    'On'      'Off'
```

You can also use the dot notation to display a single property value.

```
ai.TriggerType
ans =
Immediate
```

Channel and Line Properties

To return all configurable channel (line) property names and their possible values for a single channel (line) contained by a device object, you must use the `Channel (Line)` property. For example, to display the configurable channel properties for the first channel contained by `ai`

```
set(ai.Channel(1))
Channel Name
HwChannel
InputRange
SensorRange
Units
UnitsRange
```

All channel properties and their current values for the first channel contained by `ai` are shown below.

```
get(ai.Channel(1))
Channel Name = Left
HwChannel = 1
Index = 1
InputRange = [-1 1]
NativeOffset = 1.5259e-005
NativeScaling = 3.0518e-005
Parent = [1x1 analog input]
SensorRange = [-1 1]
Type = Channel
Units = Volts
UnitsRange = [-1 1]
```

As described in the preceding section, you can also return values for a specified number of channel properties with the `get` function or the dot notation.

Configuring Property Values

You configure property values with the `set` function or the dot notation. In practice, you can configure many of the properties at any time while the device object exists. However, some properties are not configurable while the object is running. Use the `propinfo` function, or refer to Chapter 10, “Base Property Reference” for information about when a property is configurable.

The syntax used to configure common and channel/line properties is described below. The examples are based on the analog input object `ai` created in “Returning Property Names and Property Values” on page 3-14.

Common Properties

You can configure a single property value using the `set` function

```
set(ai, 'TriggerType', 'Manual')
```

or the dot notation

```
ai.TriggerType = 'Manual';
```

To configure values for multiple properties, you can supply multiple property name/property value pairs to `set`.

```
set(ai, 'SampleRate', 44100, 'Name', 'Test 1- winsound')
```

Note that you can configure only one property value at a time using the dot notation.

Channel and Line Properties

To configure channel (line) properties for one or more channels (lines) contained by a device object, you must use the `Channel (Line)` property. For example, to configure the `SensorRange` property for the first channel contained by `ai`, you can use the `set` function

```
set(ai.Channel(1), 'SensorRange', [-2 2])
```

or the dot notation

```
ai.Channel(1).SensorRange = [-2 2];
```

To configure values for multiple channel or line properties, you supply multiple property name/property value pairs to set.

```
set(ai.Channel(1), 'SensorRange', [-2 2], 'Channel Name', 'Chan1')
```

To configure multiple property values for multiple channels

```
chs = ai.Channel(1:2);  
set(chs, {'SensorRange', 'Channel Name'}, {[-2 2], 'Chan1'; [0 4],  
'Chan2'});
```

Specifying Property Names

Device object property names are presented in this guide using mixed case. While this makes the names easier to read, you can use any case you want when specifying property names. Additionally, you need use only enough letters to identify the property name uniquely, so you can abbreviate most property names. For example, you can configure the `SampleRate` property any of these ways.

```
set(ai, 'SampleRate', 44100);  
set(ai, 'sampleRate', 44100);  
set(ai, 'sampler', 44100);
```

However, when you include property names in an M-file, you should use the full property name. This practice can prevent problems with future releases of the Data Acquisition Toolbox if a shortened name is no longer unique because of the addition of new properties.

Default Property Values

If you do not explicitly define a value for a property, then the default value is used. All configurable properties have default values. However, the default value for a given property may vary based on the hardware you are using. Additionally, some default values are calculated by the engine and depend on the values set for other properties. If the hardware driver adaptor specifies a default value for a property, then that value takes precedence over the value defined by the toolbox.

If a property has a finite set of string values, then the default value is enclosed by {} (curly braces). For example, the default value for the `LoggingMode` property is `Memory`.

```
set(ai, 'LoggingMode')  
[ Disk | {Memory} | Disk&Memory ]
```

You can also use the `propinfo` function, or refer to Chapter 10, “Base Property Reference” or Chapter 11, “Device-Specific Property Reference” to find the default value for any property.

daqpropedit: A Graphical Property Editor


The toolbox provides the Data Acquisition Property Editor, which is a graphical user interface (GUI) for accessing properties. You launch the property editor with the `daqpropedit` function. The GUI is designed so that you can:

- List all existing device objects as well as the channels or lines they contain.
- Configure property values.
- Display the names and current values for configurable properties.
- Display property characteristics.
- Display property help.

For example, to configure property values for the analog input object `ai` using `daqpropedit`

```
ai = analoginput('winsound');  
addchannel(ai, 1:2);  
daqpropedit(ai)
```

The Data Acquisition Property Editor is shown below.



Data Acquisition Property Editor

File Options Help

- winsound0-AI (1)
Channel 1: Left
Channel 2: Right

BitsPerSample: 16

BitsPerSample: 16
BufferingConfig: [512 30]
BufferingMode: 'Auto'
Channel: [2x1 aichannel]
ChannelSkew: 0
ChannelSkewMode: 'None'
ClockSource: 'Internal'
DataMissedAction: 'daquestion'
EventLog: []

Type:	double	Default Value:	16
Constraint:	Bounded	Constraint Value:	[8 16]
Read Only:	0 (False)	Read Only when Running:	1 (True)
Device Specific:	1 (True)		

BITSPERSAMPLE [8 | 16]

BitsPerSample specifies the number of bits the sound card uses to each sample.

The value of BitsPerSample cannot be modified while the object is

— List of device objects, channels, and lines

— Property to configure and its current value

— List of configurable properties and their current values

— Property characteristics

— Property help

Acquiring and Outputting Data

After you configure the device object, you can acquire or output data. Acquiring and outputting data involves these three steps:

- 1 Starting the device object
- 2 Logging data or sending data
- 3 Stopping the device object

As data is being transferred between MATLAB and your hardware, you can think of the device object as being in a particular *state*. Two types of states are defined for the Data Acquisition Toolbox:

- **Running** – For analog input objects, *running* means that data is being acquired from an analog input subsystem. However, the acquired data is not necessarily saved to memory or a disk file. For analog output objects, running means that data queued in the engine is ready to be output to an analog output subsystem.

The running state is indicated by the `Running` property for both analog input and analog output objects. `Running` can be `On` or `Off`.

- **Logging or Sending** – For analog input objects, *logging* means that data acquired from an analog input subsystem is being stored in the engine or saved to a disk file. The logging state is indicated by the `Logging` property. `Logging` can be `On` or `Off`.

For analog output objects, *sending* means the data queued in the engine is being output to an analog output subsystem. The sending state is indicated by the `Sending` property. `Sending` can be `On` or `Off`.

`Running`, `Logging`, and `Sending` are read-only properties that are automatically set to `On` or `Off` by the engine. When `Running` is `Off`, `Logging` and `Sending` must be `Off`. When `Running` is `On`, `Logging` and `Sending` are set to `On` only when a trigger occurs.

Note Digital I/O objects also possess a running state. However, since they do not store data in the engine, the logging and sending states do not exist.

Starting the Device Object

You start a device object with the `start` function. For example, to start the analog input object `ai`

```
ai = analoginput('winsound')
addchannel(ai, 1:2)
start(ai)
```

After `start` is issued, the `Running` property is automatically set to `On`, and both the device object and hardware device execute according to the configured and default property values.

While you are acquiring data with an analog input object, you can preview the data with the `peekdata` function. `peekdata` takes a “snapshot” of the most recent data but does not remove data from the engine. For example, to preview the most recent 500 samples acquired by each channel contained by `ai`

```
data = peekdata(ai, 500);
```

Since previewing data is usually a low priority task, `peekdata` does not guarantee that all requested data is returned. You can preview data at any time while the device object is running.

Logging or Sending Data

While the device object is running, you can:

- Log data acquired from an analog input subsystem to the engine (memory) or to a disk file.
- Output data queued in the engine to an analog output subsystem.

However, before you can log or send data, a trigger must occur. You configure an analog input or analog output trigger with the `TriggerType` property. All the examples presented in this section use the default `TriggerType` value of `Immediate`, which executes the trigger immediately after the `start` function is issued. For a detailed description of triggers, refer to “Configuring Analog Input Triggers” on page 5-20 or “Configuring Analog Output Triggers” on page 6-21.

Extracting Logged Data

When a trigger occurs for an analog input object, the Logging property is automatically set to On and data acquired from the hardware is logged to the engine or a disk file. You extract logged data from the engine with the `getdata` function. For example, to extract 500 samples for each channel contained by `ai`

```
data = getdata(ai, 500);
```

`getdata` blocks the MATLAB command line until all the requested data is returned to the workspace. You can extract data any time after the trigger occurs.

Sending Queued Data

For analog output objects, you must queue data in the engine with the `putdata` function before it can be output to the hardware. For example, to queue 8000 samples in the engine for each channel contained by the analog output object `ao`

```
ao = analogoutput('winsound');  
addchannel(ao, 1:2);  
data = sin(linspace(0, 2*pi, 8000))';  
putdata(ao, [data data])
```

Before the queued data can be output, you must start the analog output object.

```
start(ao)
```

When a trigger occurs, the Sending property is automatically set to On and the queued data is sent to the hardware.

Stopping the Device Object

An analog input (AI) or analog output (AO) object can stop under one of these conditions:

- You issue the `stop` function.
- The requested number of samples is acquired (AI) or sent (AO).
- A run-time hardware error occurs.
- A timeout occurs.

When the device object stops, the Running, Logging, and Sending properties are automatically set to Off. At this point, you can reconfigure the device object or immediately issue another `start` command using the current configuration.

Cleaning Up

When you no longer need a device object, you should clean up the MATLAB environment by removing the object from memory (the engine) and from the workspace. These are the steps you take to end a data acquisition session.

You remove device objects from memory with the `delete` function. For example, to delete the analog input object `ai` created in the preceding section

```
delete(ai)
```

A deleted device object is invalid, which means that you cannot connect it to the hardware. In this case, you should remove the object from the MATLAB workspace. To remove device objects and other variables from the MATLAB workspace, use the `clear` command.

```
clear ai
```

If you use `clear` on a device object that is connected to hardware, the object is removed from the workspace but remains connected to the hardware. You can restore cleared device objects to MATLAB with the `daqfind` function.

Getting Started with Analog Input

Overview	4-2
Creating an Analog Input Object	4-3
Adding Channels to an Analog Input Object	4-4
Referencing Individual Hardware Channels	4-6
Example: Adding Channels for a Sound Card	4-7
Configuring Analog Input Properties	4-9
The Sampling Rate	4-9
Trigger Types	4-11
The Samples to Acquire per Trigger	4-12
Acquiring Data	4-13
Starting the Analog Input Object	4-13
Logging Data	4-14
Stopping the Analog Input Object	4-14
Analog Input Examples	4-15
Acquiring Data with a Sound Card	4-15
Acquiring Data with a National Instruments Board	4-19
Evaluating the Analog Input Object Status	4-22
Status Properties	4-22
The Display Summary	4-23

Overview

Analog input (AI) subsystems convert real-world analog input signals from a sensor into bits that can be read by your computer. Perhaps the most important of all the subsystems commonly available, AI subsystems are typically multi-channel devices offering 12 or 16 bits of resolution. The Data Acquisition Toolbox provides access to analog input devices through an analog input object.

The purpose of this chapter is to show you how to perform simple analog input tasks using just a few functions and properties. After reading this chapter, you should be able to use the Data Acquisition Toolbox to configure your own analog input session. Topics include:

- Creating an analog input object
- Adding channels to an analog input object
- Configuring analog input properties
- Acquiring data
- Analog input examples
- Evaluating the analog input object status

Creating an Analog Input Object

You create an analog input object with the `analoginput` function. `analoginput` accepts the adaptor name and the hardware device ID as input arguments. For a list of supported adaptors, refer to “The Hardware Driver Adaptor” on page 2-7. The device ID refers to the number associated with your board when it is installed. Some vendors refer to the device ID as the device number or the board number. The device ID is optional for sound cards with an ID of 0. Use the `daqhwinfo` function to determine the available adaptors and device IDs.

Each analog input object is associated with one board and one analog input subsystem. For example, to create an analog input object associated with a National Instruments board with device ID 1

```
ai = analoginput('ni daq', 1);
```

The analog input object `ai` now exists in the MATLAB workspace. You can display the class of `ai` with the `whos` command.

```
whos ai
      Name      Size      Bytes  Class

      ai         1x1         1332  analoginput object
```

```
Grand total is 52 elements using 1332 bytes
```

Once the analog input object is created, the properties listed below are automatically assigned values. These general purpose properties provide descriptive information about the object based on its class type and adaptor.

Table 4-1: Descriptive Analog Input Properties

Property Name	Description
Name	Specify a descriptive name for the device object.
Type	Indicate the device object type.

You can display the values of these properties for `ai` with the `get` function.

```
get(ai, {'Name', 'Type'})
ans =
      'ni daq1- AI'      'Analog Input'
```

Adding Channels to an Analog Input Object

After creating the analog input object, you must add hardware channels to it. As shown by the figure in “Adding Channels or Lines” on page 3-8, you can think of a device object as a container for channels. The collection of channels contained by the device object is referred to as a *channel group*. As described in “Mapping Hardware Channel IDs to MATLAB Indices” on page 3-9, a channel group consists of a mapping between hardware channel IDs and MATLAB indices (see below).

When adding channels to an analog input object, you must follow these rules:

- The channels must reside on the same hardware device. You cannot add channels from different devices, or from different subsystems on the same device.
- The channels must be sampled at the same rate.

You add channels to an analog input object with the `addchannel` function. `addchannel` requires the device object and at least one hardware channel ID as input arguments. You can optionally specify MATLAB indices, descriptive channel names, and an output argument. For example, to add two hardware channels to the device object `ai` created in the preceding section

```
chans = addchannel(ai, 0:1);
```

The output argument `chans` is a *channel object* that reflects the channel array contained by `ai`. You can display the class of `chans` with the `whos` command.

```
whos chans
  Name      Size      Bytes  Class
  ----      -
  chans     2x1        512    ai channel object
```

```
Grand total is 7 elements using 512 bytes
```

You can use `chans` to easily access channels. For example, you can easily configure or return property values for one or more channels. As described in “Referencing Individual Hardware Channels” on page 4-6, you can also access channels with the `Channel` property.

Once you add channels to an analog input object, the properties listed below are automatically assigned values. These properties provide descriptive information about the channels based on their class type and ID.

Table 4-2: Descriptive Analog Input Channel Properties

Property Name	Description
HwChannel	Specify the hardware channel ID.
Index	Indicate the MATLAB index of a hardware channel.
Parent	Indicate the parent (device object) of a channel.
Type	Indicate a channel.

You can display the values of these properties for chans with the `get` function.

```
get(chans, {'HwChannel', 'Index', 'Parent', 'Type'})
ans =
    [0]    [1]    [1x1 analog input]    'Channel'
    [1]    [2]    [1x1 analog input]    'Channel'
```

If you are using scanning hardware, then the MATLAB indices define the scan order; index 1 is sampled first, index 2 is sampled second, and so on.

Note The number of channels you can add to a device object depends on the specific board you are using. Some boards support adding channels in any order and adding the same channel multiple times, while other boards do not. Additionally, each channel may have its own input range, which is verified with each acquired sample. The collection of channels you add to a device object is sometimes referred to as a *channel gain list* or a *channel gain queue*. For scanning hardware, these channels define the scan order.

Referencing Individual Hardware Channels

As described in the preceding section, you can access channels with the `Channel` property or with a channel object. To reference individual channels, you must specify either MATLAB indices or descriptive channel names.

MATLAB Indices

Every hardware channel contained by an analog input object has an associated MATLAB index that is used to reference the channel. When adding channels with the `addchannel` function, index assignments can be made automatically or manually. In either case, the channel indices start at 1 and increase monotonically up to the number of channel group members.

For example, the analog input object `ai` created in the preceding section had the MATLAB indices 1 and 2 automatically assigned to the hardware channels 0 and 1, respectively. To manually swap the hardware channel order, you supply the appropriate index to `chans` and use the `HwChannel` property.

```
chans(1).HwChannel = 1;  
chans(2).HwChannel = 0;
```

Alternatively, you can use the `Channel` property.

```
ai.Channel(1).HwChannel = 1;  
ai.Channel(2).HwChannel = 0;
```

Note that you can also use `addchannel` to specify the required channel order.

```
chans = addchannel(ai, [1 0]);
```

Descriptive Channel Names

Choosing a unique, descriptive name can be a useful way to identify and reference channels – particularly for large channel groups. You can associate descriptive names with hardware channels using the `addchannel` function. For example, suppose you want to add 16 single-ended channels to `ai`, and you want to associate the name `Tri gChan` with the first channel in the group.

```
ai.InputType = 'SingleEnded';  
addchannel(ai, 0, 'Tri gChan');  
addchannel(ai, 1:15);
```


Alternatively, you can use the Channel Name property.

```
ai . InputType = 'SingleEnded';
addchannel (ai , 0: 15);
ai . Channel (1) . Channel Name = 'TrigChan';
```

You can now use the channel name to reference the channel.

```
ai . TrigChan . InputRange = [- 10 10];
```

Example: Adding Channels for a Sound Card

Suppose you create the analog input object `ai` for a sound card.

```
ai = analoginput('winsound');
```

Most sound cards have just two hardware channels that you can add. If one channel is added, the sound card is said to be in *mono* mode. If two channels are added, the sound card is said to be in *stereo* mode. However, the rules for adding these two channels differ from those of other data acquisition devices. These rules are described below.

Mono Mode

If you add one channel to `ai`, the sound card is said to be in mono mode and the channel added must have a hardware ID of 1.

```
addchannel (ai , 1);
```

At the software level, mono mode means that data is acquired from channel 1. At the hardware level, you generally cannot determine the actual channel configuration and data can be acquired from channel 1, channel 2, or both depending on your sound card. Channel 1 is automatically assigned the descriptive channel name `Mono`.

```
ai . Channel . Channel Name
ans =
Mono
```

Stereo Mode

If you add two channels to `ai`, the sound card is said to be in stereo mode. You can add two channels using two calls to `addchannel` provided channel 1 is added first.

```
addchannel ( ai , 1 );  
addchannel ( ai , 2 );
```

Alternatively, you can use one call to `addchannel` provided channel 1 is specified as the first element of the hardware ID vector.

```
addchannel ( ai , 1: 2 );
```

Stereo mode means that data is acquired from both hardware channels. Channel 1 is automatically assigned the descriptive name `Left` and channel 2 is automatically assigned the descriptive name `Right`.

```
ai . Channel . Channel Name  
ans =  
    ' Left '  
    ' Right '
```

While in stereo mode, if you want to delete one channel, then that channel must be channel 2. If you try to delete channel 1, an error is returned.

```
delete(ai . Channel (2))
```

The sound card is now in mono mode.

Configuring Analog Input Properties

After hardware channels are added to the analog input object, you should configure property values. As described in “Configuring and Returning Properties” on page 3-12, the Data Acquisition Toolbox supports two basic types of properties for analog input objects: common properties and channel properties. Common properties apply to all channels contained by the device object while channel properties apply to individual channels.

The properties you configure depend on your particular analog input application. For many common applications, there is a small group of properties related to the basic setup that you will typically use. These basic setup properties control the sampling rate, define the trigger type, and define the samples to be acquired per trigger. Analog input properties related to the basic setup are given below.

Table 4-3: Analog Input Basic Setup Properties

Property Name	Description
<code>SampleRate</code>	Specify the per-channel rate at which analog data is converted to digital data.
<code>SamplesPerTrigger</code>	Specify the number of samples to acquire for each channel group member for each trigger that occurs.
<code>TriggerType</code>	Specify the type of trigger to execute.

The Sampling Rate

You control the rate at which an analog input subsystem converts analog data to digital data with the `SampleRate` property. `SampleRate` must be specified as samples per second. For example, to set the sampling rate for each channel of your National Instruments board to 100,000 samples per second (100 kHz)

```
ai = analoginput('ni daq', 1);
addchannel(ai, 0:1);
set(ai, 'SampleRate', 100000)
```

Data acquisition boards typically have predefined sampling rates that you can set. If you specify a sampling rate that does not match one of these predefined values, there are two possibilities:

- If the rate is within the range of valid values, then the engine automatically selects a valid sampling rate. The rules governing this selection process are described in the `SampleRate` reference pages in Chapter 10, “Base Property Reference.”
- If the rate is outside the range of valid values, then an error is returned.

Note For some sound cards, you can set the sampling rate to any value between the minimum and maximum values defined by the hardware. You can enable this feature with the `StandardSampleRates` property. Refer to Chapter 11, “Device-Specific Property Reference” for more information.

For hardware that supports simultaneous sampling of channels (sound cards and Agilent Technologies devices), the maximum sampling rate for each channel is given by the maximum board rate. For scanning hardware (most National Instruments and ComputerBoards devices), the per-channel sampling rate is given by the maximum hardware rate divided by the number of channels contained by the device object.

After setting a value for `SampleRate`, you should find out the actual rate set by the engine.

```
ActualRate = get(ai, 'SampleRate');
```

Alternatively, you can use the `setverify` function, which sets a property value and returns the actual value set.

```
ActualRate = setverify(ai, 'SampleRate', 100000);
```

You can find the range of valid sampling rates for your hardware with the `propinfo` function.

```
ValidRates = propinfo(ai, 'SampleRate');  
ValidRates.ConstraintValue  
ans =  
    1.0e+005 *  
    0.0000    2.0000
```

Trigger Types

For analog input objects, a trigger is defined as an event that initiates data logging to memory or to a disk file. Defining an analog input trigger involves specifying the trigger type with the `TriggerType` property. The `TriggerType` values that are supported for all hardware are given below.

Table 4-4: Analog Input `TriggerType` Property Values

TriggerType Value	Description
{Immediate}	The trigger occurs just after the start function is issued.
Manual	The trigger occurs just after you manually issue the <code>trigger</code> function.
Software	The trigger occurs when the associated trigger condition is satisfied. Trigger conditions are given by the <code>TriggerCondition</code> property.

Many devices have additional hardware trigger types, which are available to you through the `TriggerType` property. For example, to return all the trigger types for the analog input object `ai` created in the preceding section

```
set(ai, 'TriggerType')
[ Manual | {Immediate} | Software | HwDigital ]
```

This information tells you that the National Instruments board also supports a hardware digital trigger. For a description of device-specific trigger types, refer to “Device-Specific Hardware Triggers” on page 5-37, or the `TriggerType` reference pages in Chapter 10, “Base Property Reference.”

Note Triggering can be a complicated issue and it has many associated properties. For detailed information about triggering, refer to “Configuring Analog Input Triggers” on page 5-20.

The Samples to Acquire per Trigger

When a trigger executes, a predefined number of samples are acquired for each channel group member and logged to the engine or a disk file. You specify the number of samples to acquire per trigger with the `SamplesPerTrigger` property.

The default value of `SamplesPerTrigger` is calculated by the engine such that 1 second of data is collected, and is based on the default value of `SampleRate`. In general, to calculate the acquisition time for each trigger, you apply the formula

$$\text{acquisition time (seconds)} = \text{samples per trigger} / \text{sampling rate (in Hz)}$$

For example, to acquire 5 seconds of data per trigger for each channel contained by `ai`

```
set(ai, 'SamplesPerTrigger', 500000)
```

To continually acquire data, you set `SamplesPerTrigger` to `inf`.

```
set(ai, 'SamplesPerTrigger', inf)
```

A continuous acquisition will stop only if you issue the `stop` function, or an error occurs.

Acquiring Data

After you configure the analog input object, you can acquire data. Acquiring data involves these three steps:

- 1 Starting the analog input object
- 2 Logging data
- 3 Stopping the analog input object

Starting the Analog Input Object

You start an analog input object with the `start` function. For example, to start the analog input object `ai`

```
ai = analoginput('winsound')
addchannel(ai, 1:2)
start(ai)
```

After `start` is issued, the `Running` property is automatically set to `On`, and both the device object and hardware device execute according to the configured and default property values.

While you are acquiring data with an analog input object, you can preview the data with the `peekdata` function. `peekdata` takes a “snapshot” of the most recent data but does not remove data from the engine. For example, to preview the most recent 500 samples acquired by each channel contained by `ai`

```
data = peekdata(ai, 500);
```

Since previewing data is usually a low priority task, `peekdata` does not guarantee that all requested data is returned. You can preview data at any time while the device object is running. However, you cannot use `peekdata` in conjunction with hardware triggers since the device is idle until the hardware trigger is received.

Logging Data

While the analog input object is running, you can log acquired data to the engine (memory) or to a disk file. However, before you can log data a trigger must occur. You configure an analog input trigger with the `TriggerType` property. For a detailed description of triggers, refer to “Configuring Analog Input Triggers” on page 5-20.

When the trigger occurs, the `Logging` property is automatically set to `On` and data acquired from the hardware is logged to the engine or a disk file. You extract logged data from the engine with the `getdata` function. For example, to extract all logged samples for each channel contained by `ai`

```
data = getdata(ai);
```

`getdata` blocks the MATLAB command line until all the requested data is returned to the workspace. You can extract data any time after the trigger occurs. You can also return sample-time pairs with `getdata`. For example, to extract 500 sample-time pairs for each channel contained by `ai`

```
[data, time] = getdata(ai, 500);
```

`time` is an `m-by-1` array containing relative time values for all `m` samples. Time is measured relative to the time the first sample is logged, and is measured continuously until the acquisition stops. `getdata` is described in more detail in Chapter 9, “Function Reference.”

You can log data to disk with the `LoggingMode` property. You can replay data saved to disk with the `daqread` function. Refer to “Logging Information to Disk” on page 8-6 for more information about `LoggingMode` and `daqread`.

Stopping the Analog Input Object

An analog input object can stop under one of these conditions:

- You issue the `stop` function.
- The requested number of samples is acquired.
- A run-time hardware error occurs.
- A timeout occurs.

When the device object stops, the `Running` and `Logging` properties are automatically set to `Off`. At this point, you can reconfigure the device object or immediately issue another `start` command using the current configuration.

Analog Input Examples


This section illustrates how to perform basic data acquisition tasks using analog input subsystems and the Data Acquisition Toolbox. For most data acquisition applications, you must follow these basic steps:

- 1 Install and connect the components of your data acquisition hardware. At a minimum, this involves connecting a sensor to a plug-in or external data acquisition device.
- 2 Configure your data acquisition session. This involves creating a device object, adding channels, setting property values, and using specific functions to acquire data.
- 3 Analyze the acquired data using MATLAB.

Simple data acquisition applications using a sound card and a National Instruments board are given below.

Acquiring Data with a Sound Card

Suppose you must verify that the fundamental (lowest) frequency of a tuning fork is 440 Hz. To perform this task, you will use a microphone and a sound card to collect sound level data. You will then perform a fast Fourier transform (FFT) on the acquired data to find the frequency components of the tuning fork. The setup for this task is shown below.



Configuring the Data Acquisition Session

For this example, you will acquire 1 second of sound level data on one sound card channel. Since the tuning fork vibrates at a nominal frequency of 440 Hz, you can configure the sound card to its lowest sampling rate of 8000 Hz. Even at this lowest rate, you should not experience any aliasing effects since the tuning fork will not have significant spectral content above 4000 Hz, which is the Nyquist frequency. After you set the tuning fork vibrating and place it near the microphone, you will trigger the acquisition one time using a manual trigger.

You can run this example by typing `daqdoc4_1` at the MATLAB command line.

1. Create a device object – Create the analog input object `AI` for a sound card. The installed adaptors and hardware IDs are found with `daqhwinfo`.

```
AI = analoginput('winsound');
```

2. Add channels – Add one channel to `AI`.

```
chan = addchannel(AI, 1);
```

3. Configure property values – Assign values to the basic setup properties, and create the variables `blocksize` and `Fs`, which are used for subsequent analysis. The actual sampling rate is retrieved since it may be set by the engine to a value that differs from the specified value.

```
duration = 1; %1 second acquisition
set(AI, 'SampleRate', 8000)
ActualRate = get(AI, 'SampleRate');
set(AI, 'SamplesPerTrigger', duration*ActualRate)
set(AI, 'TriggerType', 'Manual')
blocksize = get(AI, 'SamplesPerTrigger');
Fs = ActualRate;
```

4. Acquire data – Start `AI`, issue a manual trigger, and extract all data from the engine. Before start is issued, you should begin inputting data from the tuning fork into the sound card.

```
start(AI)
trigger(AI)
data = getdata(AI);
```

5. Clean up – When you no longer need AI, you should remove it from memory and from the MATLAB workspace.

```
delete(AI)
clear AI
```

Analyzing the Data

For this example, analysis consists of finding the frequency components of the tuning fork and plotting the results. To do so, the function `daqdocfft` was created. This function calculates the FFT of data, and requires the values of `SampleRate` and `SamplesPerTrigger` as well as data as inputs.

```
[f, mag] = daqdocfft(data, Fs, blocksize);
```


`daqdocfft` outputs the frequency and magnitude of data, which you can then plot. `daqdocfft` is shown below.

```
function [f, mag] = daqdocfft(data, Fs, blocksize)
% [F, MAG]=DAQDOCFFT(X, FS, BLOCKSIZE) calculates the FFT of X
% using sampling frequency FS and the SamplesPerTrigger
% provided in BLOCKSIZE

xfft = abs(fft(data));
mag = 20*log10(xfft);
mag = mag(1:blocksize/2);
f = (0:length(mag)-1)*Fs/blocksize;
f = f(:);
```

The results are given below.

```
plot(f, mag)
grid on
ylabel('Magnitude (dB)')
xlabel('Frequency (Hz)')
```



The plot shows the fundamental frequency around 440 Hz and the first overtone around 880 Hz. A simple way to find actual fundamental frequency is


```
[ymax, maxi ndex] = max(mag);
maxi ndex
maxi ndex =
    441
```

The answer is 441 Hz.

Note The fundamental frequency is not always the frequency component with the largest amplitude. A more sophisticated approach involves fitting the observed frequencies to a harmonic series to find the fundamental frequency.

Acquiring Data with a National Instruments Board

Suppose you must verify that the nominal frequency of a sine wave generated by a function generator is 1.00 kHz. To perform this task, you will input the function generator signal into a National Instruments board. You will then perform a fast Fourier transform (FFT) on the acquired data to find the nominal frequency of the generated sine wave. The setup for this task is shown below.



Configuring the Data Acquisition Session

For this example, you will acquire 1 second of data on one input channel. The board is set to a sampling rate of 10 kHz, which is well above the frequency of interest. After you connect the input signal to the board, you will trigger the acquisition one time using a manual trigger.

You can run this example by typing `daqdoc4_2` at the MATLAB command line.

1. Create a device object – Create the analog input object `AI` for a National Instruments board. The installed adaptors and hardware IDs are found with `daqhwinfo`.

```
AI = analoginput('ni daq', 1);
```

2. Add channels – Add one channel to `AI`.

```
chan = addchannel(AI, 0);
```

3. Configure property values – Assign values to the basic setup properties, and create the variables `blocksize` and `Fs`, which are used for subsequent analysis. The actual sampling rate is retrieved since it may be set by the engine to a value that differs from the specified value.

```
duration = 1; %1 second acquisition
set(AI, 'SampleRate', 10000)
ActualRate = get(AI, 'SampleRate');
set(AI, 'SamplesPerTrigger', duration*ActualRate)
set(AI, 'TriggerType', 'Manual')
blocksize = get(AI, 'SamplesPerTrigger');
Fs = ActualRate;
```

4. Acquire data – Start AI, issue a manual trigger, and extract all data from the engine. Before start is issued, you should begin inputting data from the function generator into the data acquisition board.

```
start(AI)
trigger(AI)
data = getdata(AI);
```

5. Clean up – When you no longer need AI, you should remove it from memory and from the MATLAB workspace.

```
delete(AI)
clear AI
```

Analyzing the Data


For this experiment, analysis consists of finding the frequency of the input signal and plotting the results. You can find the signal frequency with `daqdocfft`.

```
[f, mag] = daqdocfft(data, Fs, blocksize);
```

This function, which is shown in “Analyzing the Data” on page 4-17, calculates the FFT of data, and requires the values of `SampleRate` and `SamplesPerTrigger` as well as data as inputs. `daqdocfft` outputs the frequency and magnitude of data, which you can then plot.

The results are given below.

```
plot(f, mag)
grid on
ylabel('Magnitude (dB)')
xlabel('Frequency (Hz)')
title('Frequency Output by Function Generator')
```



This plot shows the nominal frequency around 1000 Hz. A simple way to find actual frequency is shown below.

```
[ymax, maxi ndex] = max(mag);
maxi ndex
maxi ndex =
    994
```

The answer is 994 Hz.

Evaluating the Analog Input Object Status

You can evaluate the status of an analog input object by:

- Returning the values of certain properties
- Invoking the display summary

Status Properties

The properties associated with the status of your analog input object allow you to evaluate:

- If the device object is running
- If data is being logged to the engine or to a disk file
- How much data has been acquired
- How much data is available to be extracted from the engine

The analog input status properties are given below.

Table 4-5: Analog Input Status Properties

Property Name	Description
Loggi ng	Indicate if data is being logged to memory or to a disk file.
Runni ng	Indicate if the device object is running.
Sampl esAcqui red	Indicate the number of samples acquired per channel.
Sampl esAvai l abl e	Indicate the number of samples available per channel in the data acquisition engine.

When you issue the start function, Runni ng is automatically set to 0n. When the trigger executes, Loggi ng is automatically set to 0n and Sampl esAcqui red keeps a running count of the total number of samples per channel that have been logged to the engine or a disk file. Sampl esAvai l abl e tells you how many samples per channel are available to be extracted from the engine with the getdata function.

When the requested number of samples are acquired, `SamplesAcquired` reflects this number, and both `Running` and `Logging` are automatically set to `Off`. When you extract all the samples from the engine, `SamplesAvailable` is 0.

The Display Summary

You can invoke the display summary by typing the analog input object at the command line. The displayed information reflects many of the basic setup properties described in “Configuring Analog Input Properties” on page 4-9, and is designed so you can quickly evaluate the status of your data acquisition session. The display is divided into two main sections: general summary information and channel summary information.

General Summary Information

The general display summary includes the device object type and the hardware device name, followed by this information:

- Acquisition parameters
 - The sampling rate
 - The number of samples to acquire per trigger
 - The acquisition duration for each trigger
 - The destination for logged data
- Trigger parameters
 - The trigger type
 - The number of triggers, including the number of triggers already executed
- The engine status
 - Whether the engine is logging data, waiting to start, or waiting to trigger
 - The number of samples acquired since starting
 - The number of samples available to be extracted with `getdata`

Channel Summary Information

The channel display summary includes property values associated with:

- The hardware channel mapping
- The channel name
- The engineering units

The display summary for the example given in “Acquiring Data with a Sound Card” on page 4-15 before start is issued is shown below.

General display summary

Display Summary of Analog Input (AI) Object Using 'AudioPCI Record'.

Acquisition Parameters: 8000 samples per second on each channel.
8000 samples per trigger on each channel.
1 sec. of data to be logged per trigger.
Log data to 'Memory' on trigger.

Trigger Parameters: 1 'Manual' trigger(s) on TRIGGER.

Engine status: Waiting for START.
0 samples acquired since starting.
0 samples available for GETDATA.

Channel display summary

AI object contains channel(s):

Index:	Channel Name:	HwChannel:	InputRange:	SensorRange:	UnitsRange:	Units:
1	'Mono'	1	[- 1 1]	[- 1 1]	[- 1 1]	'Volts'

You can use the Channel property to display only the channel summary information.

AI.Channel

Doing More with Analog Input

Overview	5-2
Configuring and Sampling Input Channels	5-3
Input Channel Configuration	5-4
Sampling Rate	5-6
Channel Skew	5-7
Managing Acquired Data	5-9
Previewing Data	5-9
Extracting Data from the Engine	5-13
Returning Time Information	5-18
Configuring Analog Input Triggers	5-20
Defining a Trigger: Trigger Types and Conditions	5-21
Executing the Trigger	5-26
Trigger Delays	5-26
Repeating Triggers	5-30
How Many Triggers Occurred?	5-36
When Did the Trigger Occur?	5-37
Device-Specific Hardware Triggers	5-37
Configuring Events and Actions	5-46
Event Types	5-46
Recording and Retrieving Event Information	5-49
Creating and Executing Action Functions	5-52
Examples: Using Action Properties and Functions	5-53
Linearly Scaling the Data: Engineering Units	5-56
Example: Performing a Linear Conversion	5-57

Overview

This chapter presents the complete analog input functionality available to you with the Data Acquisition Toolbox. Properties and functions are presented in a way that reflects the typical procedures you will use to configure your analog input data acquisition session. Topics include:

- Configuring and sampling input channels
- Managing acquired data
- Configuring analog input triggers
- Configuring events and actions
- Linearly scaling the data

Configuring and Sampling Input Channels

The hardware you are using has characteristics that satisfy your specific application needs. Some of the most important hardware characteristics are related to configuring:

- The input channel type
- The sampling rate
- The channel skew (scanning hardware only)

Properties associated with configuring and sampling input channels are given below.

Table 5-1: Analog Input Properties Related to Sampling Channels

Property Name	Description
Channel Skew	Specify the time between consecutive scanned hardware channels.
Channel SkewMode	Specify how the channel skew is determined.
InputType	Specify the analog input hardware channel configuration.
SampleRate	Specify the per-channel rate at which analog data is converted to digital data.

Input Channel Configuration

You can configure your hardware input channels with the `InputType` property. The device-specific values for this property are given below.

Table 5-2: `InputType` Property Values

Vendor	InputType Value
National Instruments	{ Differential } SingleEnded NonReferencedSingleEnded
ComputerBoards	{ Differential } SingleEnded
Agilent Technologies	Differential
Sound Cards	AC-Coupled

The `InputType` value determines the number of hardware channels you can add to a device object. You can return the channel IDs with the `daqhwinfo` function. For example, suppose you create the analog input object `ai` for a National Instruments board. To display the differential channel IDs

```
ai = analoginput('ni daq', 1);
hwinfo = daqhwinfo(ai);
hwinfo.DifferentialIDs
ans =
    0     1     2     3     4     5     6     7
```

In contrast, the single-ended channel IDs would be numbered 0 through 15.

Note If the `InputType` value is changed, and that change decreases the number of channels contained by the analog input object, then a warning is returned and all channels are deleted.

National Instruments Devices

For National Instruments devices, `InputType` can be `Differential`, `SingleEnded`, or `NonReferencedSingleEnded`. Channels configured for differential input are not connected to a fixed reference such as earth, and input signals are measured as the difference between two terminals. Channels configured for single-ended input are connected to a common ground, and input signals are measured with respect to this ground. Channels configured for nonreferenced single-ended input are connected to their own ground reference, and input signals are measured with respect to this reference. The ground reference is tied to the negative input of the instrumentation amplifier.

The number of channels that you can add to a device object depends on the `InputType` property value. Most National Instruments boards have 16 or 64 single-ended inputs and 8 or 32 differential inputs, which are interleaved in banks of 8. This means that for a 64 channel board with single-ended inputs, you can add all 64 channels. However, if the channels are configured for differential input, you can only add channels 0-7, 16-23, 32-39, and 48-55.

ComputerBoards Devices

For `ComputerBoards` devices, `InputType` can be `Differential` or `SingleEnded`. Channels configured for differential input are not connected to a fixed reference such as earth, and the input signals are measured as the difference between two terminals. Channels configured for single-ended input are connected to a common ground, and input signals are measured with respect to this ground.

Agilent Technologies Devices

For Agilent Technologies devices, the only valid `InputType` value is `Differential`. Channels configured for differential input are not connected to a fixed reference such as earth, and the input signals are measured as the difference between two terminals.

Sound Cards

For sound cards, the only valid `InputType` value is `AC-Coupled`. When input channels are AC-coupled, they are connected so that constant (DC) signal levels are suppressed, and only nonzero AC signals are measured.

Sampling Rate

The sampling rate is defined as the per-channel rate (in samples/second) that an analog input subsystem converts analog data to digital data. You specify the sampling rate with the `SampleRate` property.

The maximum rate at which channels are sampled depends on the type of hardware you are using. If you are using simultaneous sample and hold (SS/H) hardware such as a sound card, then the maximum sampling rate for each channel is given by the maximum board rate. For example, suppose you create the analog input object `ai` for a sound card and configure it for stereo operation. If the device has a maximum rate of 48.0 kHz, then the maximum sampling rate per channel is 48.0 kHz.

```
ai = analoginput('winsound');
addchannel(ai, 1:2);
set(ai, 'SampleRate', 48000)
```

If you are using scanning hardware such as a National Instruments board, then the maximum sampling rate your hardware is rated at typically applies for one channel. Therefore, the maximum sampling rate per channel is given by the formula

$$\text{Maximum sampling rate per channel} = \frac{\text{Maximum board rate}}{\text{Number of channels scanned}}$$

For example, suppose you create the analog input object `ai` for a National Instruments board and add 10 channels to it. If the device has a maximum rate of 100 kHz, then the maximum sampling rate per channel is 10 kHz.

```
ai = analoginput('ni daq', 1);
set(ai, 'InputType', 'SingleEnded');
addchannel(ai, 0:9);
set(ai, 'SampleRate', 10000)
```

Typically, you can achieve this maximum rate only under ideal conditions. In practice, the sampling rate depends on several characteristics of the analog input subsystem including the settling time, the gain, and the channel skew. Channel skew is discussed in the next section.


Note Whenever the SampleRate value is changed, the BufferingConfig property value is recalculated by the engine if the BufferingMode property is set to Auto. Since BufferingConfig indicates the memory used by the engine, you should monitor this property closely.

Channel Skew

Many data acquisition devices have one A/D converter that is multiplexed to all input channels. If you sample multiple input channels from scanning hardware, then each channel is sampled sequentially following this procedure:

- 1 A single input channel is sampled.
- 2 The analog signal is converted to a digital value.
- 3 The process is repeated for every input channel being used.

Since these channels cannot be sampled simultaneously, a time gap exists between consecutively sampled channels. This time gap is called the *channel skew*. The channel skew and the sample period is illustrated below.



As shown in the preceding figure, a scan occurs when all channels in a group are sampled once and the scan rate is defined as the rate at which every channel in the group is sampled. The properties associated with configuring the channel skew are given below.

Table 5-3: Channel Skew Properties

Property Name	Description
Channel Skew	Specify the time between consecutive scanned hardware channels.
Channel SkewMode	Specify how the channel skew is determined.

Channel Skew and Channel SkewMode are configurable only for scanning hardware and not for simultaneous sample and hold (SS/H) hardware. For SS/H hardware, Channel SkewMode can only be None, and Channel Skew can only be 0. The values for Channel SkewMode are given below.

Table 5-4: ChannelSkewMode Property Values

ChannelSkewMode Value	Description
None	No channel skew is defined. This is the only valid value for simultaneous sample and hold (SS/H) hardware.
Equi sampl e	The channel skew is automatically calculated as $[(\text{sampling rate})(\text{number of channels})]^{-1}$.
Manual	The channel skew must be set with the Channel Skew property.
Mi ni mum	The channel skew is given by the smallest value supported by the hardware.

If Channel SkewMode is Mi ni mum or Equi sampl e, then Channel Skew indicates the appropriate read-only value. If Channel SkewMode is set to Manual , you must specify the channel skew with Channel Skew.

Managing Acquired Data

At the core of any analog input application lies the data you acquire from a sensor and input into your computer for subsequent analysis. The role of the analog input subsystem is to convert analog data to digitized data that can be read by the computer. There are two ways to manage acquired data:

- Preview the data with the `peekdata` function.
- Extract the data from the engine with the `getdata` function.

After data is extracted from the engine, you can analyze it, save it to disk, etc. In addition to these two functions, there are several properties associated with managing acquired data. These properties are listed below.

Table 5-5: Analog Input Data Management Properties

Property Name	Description
<code>SamplesAcquired</code>	Indicate the number of samples acquired per channel.
<code>SamplesAvailable</code>	Indicate the number of samples available per channel in the data acquisition engine.
<code>SamplesPerTrigger</code>	Specify the number of samples to acquire for each channel group member for each trigger that occurs.

Previewing Data


Before you extract and analyze acquired data, you may want to examine (preview) the data as it is being acquired. Previewing the data allows you to determine if the hardware is performing as expected and if your acquisition process is configured correctly. Once you are convinced that your system is in order, you may still want to monitor the data even as it is being analyzed or saved to disk.

Previewing data is managed with the `peekdata` function. For example, to preview the most recent 1000 samples acquired for the analog input object `ai`

```
data = peekdata(ai, 1000);
```


After start is issued, you can call peekdata. peekdata is a *nonblocking* function since it immediately returns control to MATLAB. Therefore, samples may be missed or repeated.

When a peekdata call is processed, the most recent samples requested are immediately returned, but the data is not extracted from the engine. In other words, peekdata provides a “snapshot” of the most recent requested samples. This situation is illustrated below.



■ Data stored in engine

If another peekdata call is issued, then once again, only the most recent requested samples are returned. This situation is illustrated below.



■ Data stored in engine

Rules for Using peekdata

Using peekdata to preview data follows these rules:

- You can call peekdata before a trigger executes. Therefore, peekdata is useful for previewing data before it is logged to the engine or a disk file.
- In most cases, you will call peekdata while the device object is running. However, you can call peekdata once after the device object stops running.
- If the specified number of preview samples is greater than the number of samples currently acquired, all available samples are returned with a warning message stating that the requested number of samples were not available.

For more information about peekdata, refer to its reference pages in Chapter 9, “Function Reference.”

Example: Polling the Data Block

Under certain circumstances, you may want to poll the data block. Polling the data block is useful when calling peekdata since this function does not block execution control. For example, you can issue peekdata calls based on the number of samples acquired by polling the SamplesAcquired property.

You can run this example by typing daqdoc5_1 at the MATLAB command line.

1. Create a device object – Create the analog input object AI for a sound card. The available adaptors and hardware IDs are found with daqhwinfo.

```
AI = analoginput('winsound');
%AI = analoginput('ni daq', 1);
%AI = analoginput('cbi', 1);
```

2. Add channels – Add one hardware channel to AI.

```
addchannel(AI, 1);
%addchannel(AI, 0); % For NI and CBI
```

3. Configure property values – Define a 10 second acquisition, set up a plot, and store the plot handle and title handle in the variables P and T, respectively.

```
duration = 10; % Ten second acquisition
ActualRate = get(AI, 'SampleRate');
set(AI, 'SamplesPerTrigger', duration*ActualRate)
figure
```

```
set(gcf, 'doublebuffer', 'on') %Reduce plot flicker
P = plot(zeros(1000, 1));
T = title(sprintf('Number of peekdata calls: '), num2str(0));
xlabel('Samples'), axis([0 1000 -1 1]), grid on
```

4. Acquire data – Start AI and update the display for each 1000 samples acquired by polling `Sampl esAcqui red`. The `drawnow` command forces MATLAB to update the plot. Since `peekdata` is used, all acquired data may not be displayed.

```
start(AI)
i = 1;
while AI.Sampl esAcqui red < AI.Sampl esPerTri gger
    while AI.Sampl esAcqui red < 1000*i
        end
        data = peekdata(AI, 1000);
        set(P, 'ydata', data);
        set(T, 'String', [sprintf('Number of peekdata calls: '),
            num2str(i)]);
        drawnow
        i = i + 1;
    end
while strcmp(AI.Runni ng, 'On')
end
```

5. Clean up – When you no longer need AI, you should remove it from memory and from the MATLAB workspace.

```
delete(AI)
clear AI
```

As you run this example, you may not preview all 80,000 samples stored in the engine. This is because the engine may store data faster than it can be displayed, and `peekdata` does not guarantee that all requested samples are processed.

Extracting Data from the Engine

Many data acquisition applications require that data is acquired at a fixed (often high) rate, and that the data is processed in some way immediately after it is collected. For example, you may want to perform an FFT on the acquired data and then save it to disk. When processing data, you must extract it from the engine. If acquired data is not extracted in a timely fashion, it may be overwritten.

Data is extracted from the engine with the `getdata` function. For example, to extract 1000 samples for the analog input object `ai`

```
data = getdata(ai, 1000);
```


In addition to returning acquired data, `getdata` can return relative time, absolute time, and event information. As shown below, `data` is an m -by- n array containing acquired data where m is the number of samples and n is the number of channels.

$$\begin{bmatrix} d_{11} & d_{12} & & d_{1n} \\ d_{21} & d_{22} & & d_{2n} \\ d_{31} & d_{32} & \dots & d_{3n} \\ \vdots & \vdots & & \vdots \\ \vdots & \vdots & & \vdots \\ \vdots & \vdots & & \vdots \\ d_{m1} & d_{m2} & & d_{mn} \end{bmatrix}$$

Extracted data. Each column represents a separate input channel.


`getdata` is considered a *blocking* function since it returns control to MATLAB only when the requested data is available. Therefore, samples are not missed or repeated. When a trigger executes, acquired data fills the engine. When a `getdata` call is processed, the requested samples are returned when the data is available, and then extracted from the engine.

As shown below, if a fraction of the data stored in the engine is extracted, then `getdata` always extracts the oldest data.



■ Data stored in engine

If another `getdata` call is issued, then once again, the oldest samples are extracted.



■ Data stored in engine

■ Data has been extracted from the engine

Rules for Using `getdata`

Using `getdata` to extract data stored in the engine follows these rules:

- If the requested number of samples is greater than the samples to be acquired, then an error is returned.
- If the requested data is not returned in the expected amount of time, an error is returned. The expected time to return data is given by the time it takes the engine to fill one data block plus the time specified by the `Timeout` property.
- You can issue `^C` (**Control+C**) while `getdata` is blocking. This will not stop the acquisition but will return control to MATLAB.
- The `SamplesAcquired` property keeps a running count of the total number of samples per channel that have been acquired.
- The `SamplesAvailable` property tells you how many samples you can extract from the engine per channel.
- MATLAB supports math operations only for the double data type. Therefore, if you extract data using the native data type of your hardware (typically `int16`), you must convert the data to doubles before performing math operations.

For more information about `getdata`, refer to its reference pages in Chapter 9, “Function Reference.”

Example: Previewing and Extracting Data

Suppose you have a data acquisition application that is particularly time consuming. By previewing the data, you can ascertain whether the acquisition is proceeding as expected without acquiring all the data. If it is not, then you can abort the session and diagnose the problem. This example illustrates how you might use `peekdata` and `getdata` together in such an application.

You can run this example by typing `daqdoc5_2` at the MATLAB command line.

1. Create a device object – Create the analog input object `AI` for a sound card. The installed adaptors and hardware IDs are found with `daqhwinfo`.

```
AI = analoginput('winsound');
%AI = analoginput('ni daq', 1);
%AI = analoginput('cbit', 1);
```

2. Add channels – Add one hardware channel to AI.

```
chan = addchannel (AI, 1);
%chan = addchannel (AI, 0); % For NI and CBI
```

3. Configure property values – Define a 10-second acquisition, set up the plot, and store the plot handle in the variable P. The amount of data to display is given by preview.

```
duration = 10; % Ten second acquisition
set(AI, 'SampleRate', 8000)
ActualRate = get(AI, 'SampleRate');
set(AI, 'SamplesPerTrigger', duration*ActualRate)
preview = duration*ActualRate/100;
subplot(211)
set(gcf, 'doublebuffer', 'on')
P = plot(zeros(preview, 1)); grid on
xlabel('Samples')
ylabel('Signal Level (Volts)')
```

4. Acquire data – Start AI and update the display using peekdata every time an amount of data specified by preview is stored in the engine by polling SamplesAcquired. The drawnow command forces MATLAB to update the plot. After all data is acquired, it is extracted from the engine. Note that whenever peekdata is used, all acquired data may not be displayed.

```
start(AI)
while AI.SamplesAcquired < preview
end
while AI.SamplesAcquired < duration*ActualRate
    data = peekdata(AI, preview);
    set(P, 'ydata', data)
    drawnow
end
```



Extract all the acquired data from the engine, and plot the data.

```
data = getdata(AI);
subplot(212), plot(data), grid on
title('All Acquired Data')
xlabel('Samples')
ylabel('Signal level (volts)')
```

5. Clean up – When you no longer need AI, you should remove it from memory and from the MATLAB workspace.

```
delete(AI)  
clear AI
```

The data is shown below.



Returning Time Information

You can return relative time and absolute time information with the `getdata` function. Relative time is associated with the extracted data. Absolute time is associated with the first trigger executed.

Relative Time

To return data and relative time information for the analog input object `ai`

```
[data, time] = getdata(ai);
```

`time` is an `m`-by-1 array of relative time values where `m` is the number of samples returned. `time = 0` corresponds to the first sample logged by the data acquisition engine, and time is measured continuously until the acquisition is stopped.

The relationship between the samples acquired and the relative time for each sample is shown below for `m` samples and `n` channels.

Data array. Each column represents one channel

Relative time array

$$\begin{bmatrix} d_{11} & d_{12} & \dots & d_{1n} \\ d_{21} & d_{22} & \dots & d_{2n} \\ d_{31} & d_{32} & \dots & d_{3n} \\ \vdots & \vdots & \vdots & \vdots \\ d_{m1} & d_{m2} & \dots & d_{mn} \end{bmatrix}$$

$$\begin{bmatrix} t_1 \\ t_2 \\ t_3 \\ \vdots \\ t_m \end{bmatrix}$$

Absolute Time

To return data, relative time information, and the absolute time of the first trigger for the analog input object `ai`

```
[data, time, abstime] = getdata(ai);
```

The absolute time is returned using MATLAB's clock format.

```
[year month day hour minute seconds]
```

The absolute time from the `getdata` call is

```
abstime
abstime =
1.0e+003 *
    1.9990    0.0020    0.0190    0.0130    0.0260    0.0208
```

To convert the clock vector to a more convenient form


```
t = fix(abstime);
sprintf('%d: %d: %d', t(4), t(5), t(6))
ans =
13: 26: 20
```

The absolute time of the first trigger is also recorded by the `InitialTriggerTime` property.

Note that absolute times are recorded by the `EventLog` property for each trigger executed. You can always find the absolute time associated with a data sample by adding its relative time to the absolute time of the associated trigger. Refer to “Recording and Retrieving Event Information” on page 5-49 for more information about returning absolute time information with the `EventLog` property.

Configuring Analog Input Triggers

An analog input trigger is defined as an event that initiates data logging. You can log data to the engine (memory) and to a disk file. As shown in the figure below, when a trigger occurs, the Logging property is automatically set On and data is stored in the specified target.



When defining a trigger, you must specify the trigger type. Additionally, you may need to specify one or more of these parameters:

- A trigger condition and trigger condition value
- The number of times to repeat the trigger
- A trigger delay
- An action function to execute when the trigger event occurs

Properties associated with analog input triggers are given below.

Table 5-6: Analog Input Trigger Properties

Property Name	Description
InitialTriggerTime	Indicate the absolute time of the first trigger.
ManualTriggerHwOn	Specify that the hardware device starts when a manual trigger is issued.
TriggerAction	Specify the M-file action function to execute when a trigger occurs.

Table 5-6: Analog Input Trigger Properties (Continued)

Property Name	Description
Tri ggerChannel	Specify the channels serving as trigger sources.
Tri ggerCondi ti on	Specify the condition that must be satisfied before a trigger executes.
Tri ggerCondi ti on Val ue	Specify one or more voltage values that must be satisfied before a trigger executes.
Tri ggerDel ay	Specify the delay value for data logging.
Tri ggerDel ayUni ts	Specify the units in which trigger delay data is measured.
Tri ggerRepeat	Specify the number of additional times the trigger executes.
Tri ggersExecuted	Indicate the number of triggers that execute.
Tri ggerType	Specify the type of trigger to execute.

Except for Tri ggerAct i on, these trigger-related properties are discussed in the following sections. Tri ggerAct i on is discussed in “Configuring Events and Actions” on page 5-46.

Defining a Trigger: Trigger Types and Conditions

Defining a trigger for an analog input object involves specifying the trigger type with the Tri ggerType property. You can think of the trigger type as the source of the trigger. For some trigger types, you may need to specify a trigger condition and a trigger condition value. Trigger conditions are specified with the Tri ggerCondi ti on property, while trigger condition values are specified with the Tri ggerCondi ti onVal ue property.

The analog input `TriggerType` and `TriggerCondition` values are given below.

Table 5-7: Analog Input `TriggerType` and `TriggerCondition` Values

TriggerType Value	TriggerCondition Value	Description
{Immediate}	None	The trigger occurs just after you issue the start function.
Manual	None	The trigger occurs just after you manually issue the trigger function.
Software	{Rising}	The trigger occurs when the signal has a positive slope when passing through the specified value.
	Falling	The trigger occurs when the signal has a negative slope when passing through the specified value.
	Leaving	The trigger occurs when the signal leaves the specified range of values.
	Entering	The trigger occurs when the signal enters the specified range of values.

For some devices, additional trigger types and trigger conditions are available. Refer to the `TriggerType` and `TriggerCondition` reference pages in Chapter 10, “Base Property Reference” for these device-specific values.

Trigger types are grouped into two main categories:

- Device-independent triggers
- Device-specific hardware triggers

The trigger types shown above are device-independent triggers since they are available for all supported hardware. For these trigger types, the action that initiates the trigger event involves satisfying a trigger condition in the engine (software trigger type), or issuing a toolbox function (start or trigger).

Conversely, device-specific hardware triggers depend on the specific hardware device you are using. For these trigger types, the action that initiates the trigger event involves an external analog or digital signal.

Device-specific hardware triggers for National Instruments, ComputerBoards, and Agilent Technologies devices are discussed in “Device-Specific Hardware Triggers” on page 5-37. Device-independent triggers are discussed below.

Immediate Trigger. If `TriggerType` is `Immediate` (the default value), the trigger occurs immediately after the `start` function is issued. You can configure an analog input object for continuous acquisition by use an immediate trigger and setting `SamplesPerTrigger` or `TriggerRepeat` to `inf`. Trigger repeats are discussed in “Repeating Triggers” on page 5-30.

Manual Trigger. If `TriggerType` is `Manual`, the trigger occurs just after you issue the `trigger` function. A manual trigger may provide you with more control over the data that is logged. For example, if the acquired data is noisy, you can preview the data using `peekdata`, and then manually execute the trigger after you observe that the signal is well-behaved.

Software Trigger. If `TriggerType` is `Software`, the trigger occurs when a signal satisfying the specified condition is detected on one of the hardware channels specified by the `TriggerChannel` property. The trigger condition is specified as either a voltage value and slope, or a range of voltage values using the `TriggerCondition` and `TriggerConditionValue` properties.

Example: Voice Activation Using a Software Trigger

This example demonstrates how to configure an acquisition with a sound card based on voice activation. The sample rate is set to 44.1 kHz and data is logged after an acquired sample has a value greater than or equal to 0.2 volt and a rising (positive) slope. A portion of the acquired data is then extracted from the engine and plotted.

You can run this example by typing `daqdoc5_3` at the MATLAB command line.

1. Create a device object – Create the analog input object `AIVoice` for a sound card. The installed adaptors and hardware IDs are found with `daqhwinfo`.

```
AIVoice = analoginput('winsound');
%AIVoice = analoginput('ni daq', 1);
%AIVoice = analoginput('cbit', 1);
```

2. Add channels – Add one hardware channel to AI Voice.

```
chan = addchannel (AI Voice, 1);
%chan = addchannel (AI Voice, 0); % For NI and CBI
```

3. Configure property values – Define a 2-second acquisition and configure a software trigger. The source of the trigger is chan, and the trigger executes when a rising voltage level has a value of at least 0.2 volt.

```
duration = 2; % two second acquisition
set (AI Voice, ' SampleRate', 44100)
Actual Rate = get (AI Voice, ' SampleRate' );
set (AI Voice, ' SamplesPerTrigger', Actual Rate*duration)
set (AI Voice, ' TriggerChannel', chan)
set (AI Voice, ' TriggerType', ' Software' )
set (AI Voice, ' TriggerCondition', ' Rising' )
set (AI Voice, ' TriggerConditionValue', 0.2)
```

4. Acquire data – Start AI Voice, acquire the specified number of samples, and extract the first 1000 samples from the engine as sample-time pairs. Display the number of samples remaining in the engine.

```
start (AI Voice)
while strcmp (AI Voice. Running, ' On' )
end
[data, time] = getdata (AI Voice, 1000);
remsamp = num2str (AI Voice. SamplesAvailable);
disp([' Number of samples remaining in engine: ', remsamp])
```


Plot all extracted data.

```
plot (time, data)
drawnow
xlabel (' Time (sec.) ' )
ylabel (' Signal Level (Volts) ' )
grid on
```

5. Clean up – When you no longer need AI Voice, you should remove it from memory and from the MATLAB workspace.

```
delete (AI Voice)
clear AI Voice
```

Note that when using software triggers, you must specify the `TriggerType` value before the `TriggerCondition` value. The output from this example is shown below.



The first logged sample has a signal level value of at least 0.2 volt, and this value corresponds to time = 0. Note that after you issue the `getdata` function, 87,200 samples remain in the engine.

```
AI Voice. SamplesAvailable
ans =
    87200
```

Executing the Trigger

For an analog input trigger to occur, you must follow these steps:

- 1 Configure the appropriate trigger properties.
- 2 Issue the start function.
- 3 Issue the trigger function if `TriggerType` value is `Manual`.

Once the trigger occurs, data logging is initiated. The device object and hardware device stop executing when the requested samples are acquired, a run-time error occurs, or you issue the stop function.

Note After a trigger occurs, the number of samples specified by `SamplesPerTrigger` is acquired for each channel group member before the next trigger can occur.


Trigger Delays

Trigger delays allow you to control exactly when data is logged after a trigger occurs. You can log data either before the trigger or after the trigger. Logging data before the trigger occurs is called *pretriggering*, while logging data after a trigger occurs is called *postriggering*.

You configure trigger delays with the `TriggerDelay` property. Pretriggers are specified by a negative `TriggerDelay` value, while postriggers are specified by a positive `TriggerDelay` value. You can delay data logging in time or in samples using the `TriggerDelayUnits` property. When `TriggerDelayUnits` is set to `Samples`, data logging is delayed by the specified number of samples. When the `TriggerDelayUnits` property is set to `Seconds`, data logging is delayed by the specified number of seconds.

Capturing Pretrigger Data

In some circumstances, you may want to capture data before the trigger occurs. Such data is called pretrigger data. When capturing pretrigger data, the `SamplesPerTrigger` property value includes the data captured before and after the trigger occurs. Capturing pretrigger data is illustrated in below.



■ Data stored in engine


You can capture pretrigger data for manual triggers and software triggers. If `TriggerType` is `Manual`, and the `trigger` function is issued before the trigger delay passes, then a warning is returned and the trigger is ignored (the trigger event does not occur).

You cannot capture pretrigger data for immediate triggers or device-specific hardware triggers.

Note Pretrigger data has negative relative time values associated with it. This is because time = 0 corresponds to the time the trigger event occurs and data logging is initiated.

Capturing Postrigger Data

In some circumstances, you may want to capture data after the trigger occurs. Such data is called postrigger data. When capturing postrigger data, the `Sampl esPerTri gger` property value and the number of postrigger samples are equal. Capturing postrigger data is illustrated below.



You can capture postrigger using any supported trigger type.

Example: Voice Activation and Pretriggers

This example modifies `daqdoc5_3` such that 500 pretrigger samples are acquired. You can run this example by typing `daqdoc5_4` at the MATLAB command line.

1. Create a device object – Create the analog input object `AI Voi ce` for a sound card. The installed adaptors and hardware IDs are found with `daqhwi nfo`.

```
AI Voi ce = anal ogi nput ( ' wi nsound' );
%AI Voi ce = anal ogi nput ( ' ni daq' , 1 );
%AI Voi ce = anal ogi nput ( ' cbi ' , 1 );
```

2. Add channels – Add one hardware channel to `AI Voi ce`.

```
chan = addchannel (AI Voi ce, 1);
%chan = addchannel (AI Voi ce, 0); % For NI and CBI
```

3. Configure property values – Define a 2-second acquisition, and configure a software trigger. The source of the trigger is chan, and the trigger executes when a rising voltage level has a value of at least 0.2 volt. Additionally, 500 pretrigger samples are collected.

```
duration = 2; % two second acquisition
set(AIVoice, 'SampleRate', 44100)
ActualRate = get(AIVoice, 'SampleRate');
set(AIVoice, 'SamplesPerTrigger', ActualRate*duration)
set(AIVoice, 'TriggerChannel', chan)
set(AIVoice, 'TriggerType', 'Software')
set(AIVoice, 'TriggerCondition', 'Rising')
set(AIVoice, 'TriggerConditionValue', 0.2)
set(AIVoice, 'TriggerDelayUnits', 'Samples')
set(AIVoice, 'TriggerDelay', -500)
```

4. Acquire data – Start AIVoice, acquire the specified number of samples, and extract the first 1000 samples from the engine as sample-time pairs.

```
start(AIVoice);
while strcmp(AIVoice.Running, 'On')
end
[data, time] = getdata(AIVoice, 1000);
```


Plot all the extracted data.

```
plot(time, data)
xlabel('Time (sec.)')
ylabel('Signal Level (Volts)')
grid on
```

5. Clean up – When you no longer need AIVoice, you should remove it from memory and from the MATLAB workspace.

```
delete(AIVoice)
clear AIVoice
```

The output from this example is shown below. Note that the pretrigger data constitutes half of the 1000 samples extracted from the engine. Additionally, pretrigger data has negative time associated with it since time = 0 corresponds to the time the trigger event occurs and data logging is initiated.



Repeating Triggers

You can configure triggers to occur once (one-shot acquisition) or multiple times. You control trigger repeats with the `TriggerRepeat` property. If `TriggerRepeat` is set to its default value of 0, then the trigger occurs once. If `TriggerRepeat` is set to a positive integer value, then the trigger is repeated the specified number of times. If `TriggerRepeat` is set to `inf`, then the trigger repeats continuously and you can stop the device object only by issuing the `stop` function.

Example: Voice Activation and Repeating Triggers

This example modifies `daqdoc5_3` such that two triggers are issued. The specified amount of data is acquired for each trigger and stored in separate variables. The `Timeout` value is set to five seconds. Therefore, if `getdata` does not return the specified number of samples in the time given by the `Timeout` property plus the time required to acquire the data, the acquisition will be aborted.

You can run this example by typing `daqdoc5_5` at the MATLAB command line.

1. Create a device object – Create the analog input object `AIVoice` for a sound card. The installed adaptors and hardware IDs are found with `daqhwinfo`.

```
AIVoice = analoginput('winsound');
%AIVoice = analoginput('ni daq', 1);
%AIVoice = analoginput('cbi', 1);
```

2. Add channels – Add one hardware channel to `AIVoice`.

```
chan = addchannel(AIVoice, 1);
%chan = addchannel(AIVoice, 0); % For NI and CBI
```

3. Configure property values – Define a 1-second total acquisition time and configure a software trigger. The source of the trigger is `chan`, and the trigger executes when a rising voltage level has a value of at least 0.2 volt. Additionally, the trigger is repeated once when the trigger condition is met.

```
duration = 0.5; % One-half second acquisition for each trigger
set(AIVoice, 'SampleRate', 44100)
ActualRate = get(AIVoice, 'SampleRate');
set(AIVoice, 'Timeout', 5)
set(AIVoice, 'SamplesPerTrigger', ActualRate*duration)
set(AIVoice, 'TriggerChannel', chan)
set(AIVoice, 'TriggerType', 'Software')
set(AIVoice, 'TriggerCondition', 'Rising')
set(AIVoice, 'TriggerConditionValue', 0.2)
set(AIVoice, 'TriggerRepeat', 1)
```

4. Acquire data – Start `AI Voice`, acquire the specified number of samples, extract all the data from the first trigger as sample-time pairs, and extract all the data from the second trigger as sample-time pairs. Note that you can extract the data acquired from both triggers with the command `getdata(AI Voice, 44100)`.

```
start(AI Voice)
while strcmp(AI Voice.Running, 'On')
end
[d1, t1] = getdata(AI Voice);
[d2, t2] = getdata(AI Voice);
```


Plot the data for both triggers.

```
subplot(211), plot(t1, d1), grid on, hold on
axis([t1(1)-0.05 t1(end)+0.05 -0.8 0.8])
xlabel('Time (sec)'), ylabel('Signal level (Volts)'),
title('Voice Activation First Trigger')
subplot(212), plot(t2, d2), grid on
axis([t2(1)-0.05 t2(end)+0.05 -0.8 0.8])
xlabel('Time (sec)'), ylabel('Signal level (Volts)'),
title('Voice Activation Second Trigger')
```

5. Clean up – When you no longer need `AI Voice`, you should remove it from memory and from the MATLAB workspace.


```
delete(AI Voice)
clear AI Voice
```

The data acquired for both triggers is shown below.



As described in “Extracting Data from the Engine” on page 5-13, if you do not specify the amount of data to extract from the engine with `getdata`, then the amount of data returned is given by the `SamplesPerTrigger` property. You can return data from multiple triggers with one call to `getdata` by specifying the appropriate number of samples. When you return data that spans multiple triggers, a NaN is inserted in the data stream between trigger events. Therefore, an extra “sample” (the NaN) is stored in the engine and returned by `getdata`. Identifying these NaNs allows you to locate where and when each trigger was issued in the data stream.

The figure below illustrates the data stored by the engine during a multiple-trigger acquisition. The data acquired for each trigger is given by the `SamplesPerTrigger` property value. The relative trigger times are shown on the `Time` axis where first trigger time corresponds to t_1 (0 seconds by definition), the second trigger time corresponds to t_2 , and so on.




The following code modifies `daqdoc5_5` so that multiple-trigger data is extracted from the engine with one call to `getdata`.

```
returndata = ActualRate*duration*(AIVoice.TriggerRepeat + 1);
start(AIVoice)
[d, t] = getdata(AIVoice, returndata);
```

Plot the data.

```
plot(t, d)
xlabel('Time (sec.)')
ylabel('Signal Level (Volts)')
title('Voice Activation for Both Triggers')
grid on
```

The multiple-trigger data is shown below.



You can find the relative trigger times by searching for NaNs in the returned data. You can find the index location of the NaN in `d` or `t` using the `isnan` function.

```
index = find(isnan(d))
index =
    22051
```

With this information, you can find the relative time for the second trigger.

```
t2time = t(index+1)
t2time =
    0.5980
```

How Many Triggers Occurred?

You can find out how many triggers occurred with the `TriggersExecuted` property value. The trigger number for each trigger executed is also recorded by the `EventLog` property. A convenient way to access event log information is with the `showdaevents` function.

For example, suppose you create the analog input object `ai` for a sound card and add one channel to it. `ai` is configured to acquire 40,000 samples with five triggers using the default sampling rate of 8000 Hz.

```
ai = analoginput('winsound');
ch = addchannel(ai, 1);
set(ai, 'TriggerRepeat', 4);
start(ai)
```

`TriggersExecuted` returns the number of triggers executed.

```
ai.TriggersExecuted
ans =
     5
```

`showdaevents` returns information for all the events that occurred while `ai` was executing.

```
showdaevents(ai)
 1 Start                ( 10: 22: 04,  0 )
 2 Trigger#1            ( 10: 22: 04,  0 )      Channel: N/A
 3 Trigger#2            ( 10: 22: 05, 8000 )      Channel: N/A
 4 Trigger#3            ( 10: 22: 06, 16000 )     Channel: N/A
 5 Trigger#4            ( 10: 22: 07, 24000 )     Channel: N/A
 6 Trigger#5            ( 10: 22: 08, 32000 )     Channel: N/A
 7 Stop                 ( 10: 22: 09, 40000 )
```

For more information about recording and retrieving events, refer to “Recording and Retrieving Event Information” on page 5-49.

When Did the Trigger Occur?

You can find the absolute time of the first trigger event with the `InitialTriggerTime` property value. The absolute time is returned using MATLAB's clock format.

[year month day hour minute seconds]

For example, the absolute time of the first trigger event for the preceding example is

```
abstime = ai.InitialTriggerTime
abstime =
1.0e+003 *
    1.9990    0.0040    0.0170    0.0100    0.0220    0.0041
```

To convert the clock vector to a more convenient form, you can use the `sprintf` function.

```
t = fix(abstime);
sprintf('%d: %d: %d', t(4), t(5), t(6))
ans =
10: 22: 4
```

You can also use the `showdaqevents` function to return the absolute time of each trigger event. For more information about trigger events, refer to “Recording and Retrieving Event Information” on page 5-49.

Device-Specific Hardware Triggers

Many data acquisition devices possess the ability to accept a hardware trigger. Hardware triggers are processed directly by the hardware and can be either a digital signal or an analog signal. Hardware triggers are often used when speed is required since a hardware device can process an input signal much faster than software.

The device-specific hardware triggers are presented to you as additional property values. Hardware triggers for National Instruments, ComputerBoards, and Agilent Technologies devices are discussed below and in Chapter 10, “Base Property Reference.”

Note that the available hardware trigger support depends on the board you are using. Refer to your hardware documentation for detailed information about its triggering capabilities.

National Instruments

When using National Instruments (NI) hardware, there are additional trigger types and trigger conditions available to you. These device-specific property values fall into two categories: hardware digital triggering and hardware analog triggering.

The device-specific trigger types and trigger conditions are described below and in Chapter 10, “Base Property Reference.”


Table 5-8: Analog Input TriggerType and TriggerCondition Property Values for NI Hardware

TriggerType Value	TriggerCondition Value	Description
HwDi gi tal	None	The trigger occurs when the falling edge of a digital (TTL) signal is detected.
HwAnal ogChannel or HwAnal ogPi n	{ AboveHi ghLevel }	The trigger occurs when the analog signal is above the specified value.
	Bel owLowLevel	The trigger occurs when the analog signal is below the specified value.
	Hi ghHysteresi s	The trigger occurs when the analog signal is greater than the specified high value with hysteresis given by the specified low value.
	I nsi deRegi on	The trigger occurs when the analog signal is inside the specified region.
	LowHysteresi s	The trigger occurs when the analog signal is less than the specified low value with hysteresis given by the specified high value.

Hardware Digital Triggering. If `TriggerType` is `HwDigital`, the trigger occurs when the falling edge of a digital (TTL) signal is detected. The following example illustrates how to configure a hardware digital trigger.

```
ai = analoginput('ni daq', 1);
addchannel(ai, 0:7);
set(ai, 'TriggerType', 'HwDigital')
```

The diagram below illustrates how you connect a digital trigger signal to an MIO-16E Series board. PFI0/TRIG1 corresponds to pin 11.



MIO-16E Series board

Hardware Analog Triggering. If `TriggerType` is `HwAnalogPin`, the trigger is given by a low-range analog signal (typically between -10 to 10 volts). For example, to trigger your acquisition when the trigger signal is between -4 volts and 4 volts


```
ai = analoginput('ni daq', 1);
addchannel(ai, 0:7);
set(ai, 'TriggerType', 'HwAnalogPin')
set(ai, 'TriggerCondition', 'InsideRegion')
set(ai, 'TriggerConditionValue', [-4.0 4.0])
```

If `TriggerType` is `HwAnalogChannel`, the trigger is given by an analog signal and the trigger channel is defined as the first channel in the channel group (MATLAB index of one). The valid range of the analog trigger signal is given by the full-scale range of the trigger channel. The following example illustrates how to configure such a trigger where the trigger channel is assigned the

descriptive name `Tri gChan` and the default `Tri ggerCondi ti on` property value is used.

```
ai = anal ogi nput ( ' ni daq' , 1 );
addchannel ( ai , 0: 7 );
set ( ai . Channel ( 1 ) , ' Channel Name' , ' Tri gChan' )
set ( ai , ' Tri ggerType' , ' HwAnal ogChannel ' )
set ( ai , ' Tri ggerCondi ti onVal ue' , 0. 2 )
```

The diagram below illustrates how you can connect analog trigger signals to an MIO-16E Series board.



ComputerBoards

When using ComputerBoards (Measurement Computing Corporation) hardware, there are additional trigger types and trigger conditions available to you. These device-specific property values fall into two categories: hardware digital triggering and hardware analog triggering.

The device-specific trigger types and trigger conditions are described below and in Chapter 10, “Base Property Reference.”

Table 5-9: Analog Input TriggerType and TriggerCondition Values for CBI Hardware

TriggerType Value	TriggerCondition Value	Description
HwDi gi tal	GateHi gh	The trigger occurs as long as the digital signal is high.
	GateLow	The trigger occurs as long as the digital signal is low.
	Tri gHi gh	The trigger occurs when the digital signal is high.
	Tri gLow	The trigger occurs when the digital signal is low.
	Tri gPosEdge	The trigger occurs when the positive (rising) edge of the digital signal is detected.
	{ Tri gNegEdge}	The trigger occurs when the negative (falling) edge of the digital signal is detected.


Table 5-9: Analog Input TriggerType and TriggerCondition Values for CBI Hardware (Continued)

TriggerType Value	TriggerCondition Value	Description
HwAnalog	{Tri gAbove}	The trigger occurs when the analog signal makes a transition from below the specified value to above.
	Tri gBel ow	The trigger occurs when the analog signal makes a transition from above the specified value to below.
	GateNegHys	The trigger occurs when the analog signal is more than the specified high value. The acquisition stops if the analog signal is less than the specified low value.
	GatePosHys	The trigger occurs when the analog signal is less than the specified low value. The acquisition stops if the analog signal is more than the specified high value.
	GateAbove	The trigger occurs as long as the analog signal is more than the specified value.
	GateBel ow	The trigger occurs as long as the analog signal is less than the specified value.
	GateInWi ndow	The trigger occurs as long as the analog signal is within the specified range of values.
	GateOutWi ndow	The trigger occurs as long as the analog signal is outside the specified range of values.

Hardware Digital Triggering. If Tri ggerType is HwDi gi tal , the trigger is given by a digital (TTL) signal. For example, to trigger your acquisition when the positive edge of a digital signal is detected

```
ai = analoginput(' cbi ' , 1);
addchannel(ai, 0: 7);
set(ai, ' Tri ggerType' , ' HwDi gi tal ' )
set(ai, ' Tri ggerCondi ti on' , ' Tri gPosEdge' )
```

The diagram below illustrates how you connect a digital trigger signal to a PCI-DAS1602/16 board. A/D External Trigger corresponds to pin 45.




PCI-DAS1602/16 board

Hardware Analog Triggering. If `TriggerType` is `HwAnalog`, the trigger is given by an analog signal. For example, to trigger your acquisition when the trigger signal is between -4 volts and 4 volts

```
ai = analoginput('cbi', 1);
addchannel(ai, 0:7);
set(ai, 'TriggerType', 'HwAnalog');
set(ai, 'TriggerCondition', 'GateInWindow');
set(ai, 'TriggerConditionValue', [-4.0 4.0]);
```

The diagram below illustrates how you connect an analog trigger signal to a PCI-DAS1602/16 board. AI Ch 0-7 corresponds to pins 2-17, while Analog Trigger In corresponds to pin 43.



PCI-DAS1602/16 board

Agilent Technologies

When using Agilent Technologies hardware, there are additional trigger types and trigger conditions available to you. These device-specific property values fall into two categories: hardware digital triggering and hardware analog triggering.

The device-specific trigger types and trigger conditions are described below and in Chapter 10, “Base Property Reference.”

Table 5-10: Analog Input TriggerType and TriggerCondition Property Values for Agilent Hardware

TriggerType Value	TriggerCondition Value	Description
HwDi gi tal	{Posi ti veEdge}	The trigger occurs when the positive (rising) edge of a digital signal is detected.
	Negat i veEdge	The trigger occurs when the negative (falling) edge of a digital signal is detected.
HwAnal og	{Ri si ng}	The trigger occurs when the analog signal has a positive slope when passing through the specified range of values.
	Fal l i ng	The trigger occurs when the analog signal has a negative slope when passing through the specified range of values.
	Leavi ng	The trigger occurs when the analog signal leaves the specified range of values.
	Enter i ng	The trigger occurs when the analog signal enters the specified range of values.

Note that when TriggerType is HwAnal og, the trigger condition values are all specified as two-element vectors. Setting two trigger levels prevents the module from triggering repeatedly due to a noisy signal.

Hardware Digital Triggering. If `TriggerType` is `HwDigital`, the trigger is given by a digital (TTL) signal. For example, to trigger your acquisition when the negative edge of a digital signal is detected

```
ai = analoginput('hpe1432', 8);  
addchannel(ai, 1:16);  
set(ai, 'TriggerType', 'HwDigital');  
set(ai, 'TriggerCondition', 'NegativeEdge');
```

Hardware Analog Triggering. If `TriggerType` is `HwAnalog`, the trigger is given by an analog signal. For example, to trigger your acquisition when the trigger signal is between -4 volts and 4 volts

```
ai = analoginput('hpe1432', 8);  
addchannel(ai, 1:16);  
set(ai, 'TriggerType', 'HwAnalog');  
set(ai, 'TriggerCondition', 'Entering');  
set(ai, 'TriggerConditionValue', [-4.0 4.0]);  
set(ai, 'TriggerChannel', ai.Channel(1:4));
```

Configuring Events and Actions

You can enhance the power and flexibility of your analog input application by utilizing *events*. An event occurs at a particular time after a condition is met and may result in one or more actions.

While the analog input object is running, you can use events to display a message, display data, analyze data, or perform just about any other action. Actions are controlled through *action properties* and *action functions*. All event types have an associated action property. Action functions are M-file functions that you construct to suit your specific data acquisition needs.

You execute an action when a particular event occurs by specifying the name of the M-file action function as the value for the associated action property. Note that `daqacti on` is the default value for some action properties.

Event Types

The analog input event types and associated action properties are described below.

Table 5-11: Analog Input Action Properties

Event Type	Property Name
Data missed	DataMi ssedActi on
Input overrange	I nputOverRangeActi on
Run-time error	Runt i meErrorActi on
Samples acquired	Sampl esAcqui redActi on
	Sampl esAcqui redActi onCount
Start	StartActi on
Stop	StopActi on
Timer	Ti merActi on
	Ti merPeri od
Trigger	Tri ggerActi on

Data Missed Event. A data missed event is generated immediately after acquired data is missed. In most cases, data is missed because:

- The engine cannot keep up with the rate of acquisition.
- The driver wrote new data into the hardware's FIFO buffer before the previously acquired data was read. You can usually avoid this problem by increasing the size of the memory block with the `BufferingConfig` property.

This event executes the action function specified for the `DataMissedAction` property. The default value for `DataMissedAction` is `daqaction`, which displays the event type and the device object name. When a data missed event occurs, the analog input object is automatically stopped.

Input Overrange Event. An input overrange event is generated immediately after an overrange condition is detected for any channel group member. An overrange condition occurs when an input signal exceeds the range specified by the `InputRange` property.

This event executes the action function specified for the `InputOverRangeAction` property. Overrange detection is enabled only when an action function is specified for `InputOverRangeAction`, and the analog input object is running.

Run-time Error Event. A run-time error event is generated immediately after a run-time error occurs. Additionally, a toolbox error message is automatically displayed to the MATLAB workspace. If an error occurs that is not explicitly handled by the toolbox, then the hardware-specific error message is displayed.

This event executes the action function specified for `RuntimeErrorAction`. The default value for `RuntimeErrorAction` is `daqaction`, which displays the event type, the time the event occurred, the device object name, and the error message.

Run-time errors include hardware errors and timeouts. Run-time errors do not include configuration errors such as setting an invalid property value.

Samples Acquired Event. A samples acquired event is generated immediately after a predetermined number of samples is acquired.

This event executes the action function specified for the `SamplesAcquiredAction` property every time the number of samples specified by `SamplesAcquiredActionCount` is acquired for each channel group member.

You should use `SamplesAcquiredAction` if you must access each sample that is acquired. However, if you are performing a CPU-intensive task with the data, then system performance may be adversely affected. If you do not have this requirement, you may want to use the `TimerAction` property.

Start Event. A start event is generated immediately after the start function is issued. This event executes the action function specified for `StartAction`. When the `StartActionM`-file has finished executing, `Running` is automatically set to `On` and the device object and hardware device begin executing.

Stop Event. A stop event is generated immediately after the device object and hardware device stop running. This occurs when:

- The stop function is issued.
- The requested number of samples is acquired.
- A run-time error occurs.

A stop event executes the action function specified for `StopAction`. Under most circumstances, the action function is not guaranteed to complete execution until sometime after the device object and hardware device stop running, and the `Running` property is set to `Off`.

Timer Event. A timer event is generated whenever the time specified by the `TimerPeriod` property passes. This event executes the action function specified for `TimerAction`. Time is measured relative to when the device object starts running.

Some timer events may not be processed if your system is significantly slowed or if the `TimerPeriod` value is too small. For example, a common application for timer events is to display data. However, since displaying data is a CPU-intensive task, some of these events may be dropped. To guarantee that events are not dropped, use the `SamplesAcquiredAction` property.

Trigger Event. A trigger event is generated immediately after a trigger occurs. This event executes the action function specified for the `TriggerAction` property. Under most circumstances, the action function is not guaranteed to complete execution until sometime after `Logging` is set to `On`.

Recording and Retrieving Event Information

While the analog input object is running, certain information is automatically recorded in the EventLog property for some of the event types listed in the preceding section. EventLog is a structure that contains two fields: Type and Data. The Type field contains the event type. The Data field contains event-specific information. Events are recorded in the order in which they occur. The first EventLog entry reflects the first event recorded, the second EventLog entry reflects the second event recorded, and so on.

The event types recorded in EventLog for analog input objects, as well as the values for the Type and Data fields, are given below.

Table 5-12: Analog Input Event Information Stored in EventLog

Event Type	Type Field Value	Data Field Value
Data missed	DataMi ssed	Rel Sampl e
Input overrange	OverRange	Rel Sampl e
		Channel
		OverRange
Run-time error	Error	AbsTi me
		Rel Sampl e
		Stri ng
Start	Start	AbsTi me
		Rel Sampl e
Stop	Stop	AbsTi me
		Rel Sampl e
Trigger	Tri gger	AbsTi me
		Rel Sampl e
		Channel
		Tri gger

Samples acquired events and timer events are not stored in EventLog.

Note Unless a run-time error occurs, EventLog records a start event, trigger event, and stop event for each data acquisition session.

The Data field values are described below.

The AbsTime Field. AbsTime is used by the run-time error, start, stop, and trigger events to indicate the absolute time the event occurred. The absolute time is returned using MATLAB's clock format.

day-month-year hour:minute:second

The Channel Field. Channel is used by the input overrange event and the trigger event. For the input overrange event, Channel indicates the index number of the input channel that experienced an overrange signal. For the trigger event, Channel indicates the index number for each input channel serving as a trigger source.

The OverRange Field. OverRange is used by the input overrange event, and can be On or Off. If OverRange is On, then the input channel experienced an overrange signal. If OverRange is Off, then the input channel no longer experienced an overrange signal.

The RelSample Field. RelSample is used by all events stored in EventLog to indicate the sample number that was acquired when the event occurred. RelSample is 0 for the start event and for the first trigger event regardless of the trigger type. RelSample is NaN for any event that occurs before the first trigger executes.

The String Field. String is used by the run-time error event to store the descriptive message that is generated when a run-time error occurs. This message is also displayed at the MATLAB command line.

The Trigger Field. Trigger is used by the trigger event to indicate the trigger number. For example, if three trigger events occur, then Trigger is 3 for the third trigger event. The total number of triggers executed is given by the TriggersExecuted property.

Example: Retrieving Event Information

Suppose you want to examine the events logged for the example given by “Example: Voice Activation Using a Software Trigger” on page 5-23. You can do this by accessing the EventLog property.

```
events = AIVoice.EventLog
events =
3x1 struct array with fields:
    Type
    Data
```

By examining the contents of the Type field, you can list the events that occurred while AI Voice was running.

```
{events.Type}
ans =
    'Start'    'Trigger'    'Stop'
```

To display information about the trigger event, you must access the Data field, which stores the absolute time the trigger occurred, the number of samples acquired when the trigger occurred, the index of the trigger channel, and the trigger number.

```
trigdata = events(2).Data
trigdata =
    AbsTime: [1999 4 15 18 12 5.8615]
    RelSample: 0
    Channel: 1
    Trigger: 1
```

Alternatively, you can use the showdaqevents function to display event information.

```
showdaqevents(events, 2)
    2 Trigger#1          ( 18:12:05, 0 )      Channel: 1
```

Creating and Executing Action Functions

When using action functions, you should be aware of these execution rules:

- Action functions execute in the order in which they are issued.
- All action functions except those associated with timer events are guaranteed to execute.
- Action function execution may be delayed if the action involves a CPU-intensive task such as updating a figure.

You can specify the action function to be executed when a specific event type occurs by including the name of the M-file as the value for the associated action property. For example, to execute the action function `myaction` for the analog input object `ai` every time 1000 samples are acquired

```
ai.SamplesAcquiredActionCount = 1000;
ai.SamplesAcquiredAction = 'myaction';
```

M-file action functions require at least two input arguments. The first argument is the device object. The second argument is a variable that captures the event information given in Table 5-12, Analog Input Event Information Stored in EventLog, on page 5-49. This event information pertains only to the event that caused the action function to execute. The function header for `myaction` is shown below.

```
function myaction(obj, event)
```

You can pass additional parameters to the action function by including them as elements of a cell array. For example, to pass the MATLAB variable `time` to `myaction`

```
time = datestr(now, 0);
ai.SamplesAcquiredActionCount = 1000;
ai.SamplesAcquiredAction = {'myaction', time};
```

The corresponding function header is

```
function myaction(obj, event, time)
```

If you pass additional parameters to the action function, then they must be included in the function header after the two required arguments.

Specifying a Toolbox Function as an Action

In addition to specifying your own action function, you can specify the start, stop, or trigger toolbox functions as actions. For example, to configure ai to stop running when an overrange condition occurs

```
ai.InputOverRangeAction = 'stop';
```

Examples: Using Action Properties and Functions

This section provides examples that show you how to create action functions and configure action properties.

Displaying Event Information with an Action Function

This example illustrates how action functions allow you to easily display event information. The example uses daqacti on to display information for trigger, run-time error, and stop events. The default SampleRate and SamplesPerTrigger values are used, which results in a 1-second acquisition for each trigger executed.

You can run this example by typing daqdoc5_6 at the MATLAB command line.

1. Create a device object – Create the analog input object AI for a sound card. The installed adaptors and hardware IDs are found with daqhwinfo.

```
AI = analoginput('winsound');
%AI = analoginput('ni daq', 1);
%AI = analoginput('cbi', 1);
```

2. Add channels – Add one hardware channel to AI.

```
chan = addchannel(AI, 1);
%chan = addchannel(AI, 0); % For NI and CBI
```

3. Configure property values – Repeat the trigger three times, find the time for the acquisition to complete, and define daqacti on as the M-file to execute when a trigger, run-time error, or stop event occurs.

```
set(AI, 'TriggerRepeat', 3)
time = (AI.SamplesPerTrigger/AI.SampleRate) * (AI.TriggerRepeat+1);
set(AI, 'TriggerAction', 'daqacti on')
set(AI, 'RuntimeErrorAction', 'daqacti on')
set(AI, 'StopAction', 'daqacti on')
```

4. Acquire data – Start AI and wait for it to stop running. The pause command allows the output from daqacti on to be displayed.

```
start(AI)
pause(time+1)
```

5. Clean up – When you no longer need AI, you should remove it from memory and from the MATLAB workspace.

```
delete(AI)
clear AI
```

Passing Additional Parameters to an Action Function

This example illustrates how additional arguments are passed to the action function. Timer events are generated every 0.5 second to display data using the action function daqdoc5_7plot.

You can run this example by typing daqdoc5_7 at the MATLAB command line.

1. Create a device object – Create the analog input object AI for a sound card. The installed adaptors and hardware IDs are found with daqhwinfo.

```
AI = analoginput('winsound');
%AI = analoginput('ni daq', 1);
%AI = analoginput('cbi', 1);
```

2. Add channels – Add one hardware channel to AI.

```
chan = addchannel(AI, 1);
%chan = addchannel(AI, 0); % For NI and CBI
```

3. Configure property values – Define a 10-second acquisition and execute the M-file daqdoc5_7plot every 0.5 seconds. Note that the variables bsize, P, and T are passed to the action function.

```
duration = 10; % Ten second duration
set(AI, 'SampleRate', 22050)
ActualRate = get(AI, 'SampleRate');
set(AI, 'SamplesPerTrigger', duration*ActualRate)
set(AI, 'TimerPeriod', 0.5)
bsize = (AI.SampleRate)*(AI.TimerPeriod);
figure
P = plot(zeros(bsize, 1));
T = title(['Number of action functions calls: ', num2str(0)]);
```



```
xlabel('Samples'), ylabel('Signal (Volts)')  
grid on  
set(gcf, 'doublebuffer', 'on')  
set(AI, 'TimerAction', {'daqdoc5_7plot', bsize, P, T})
```

4. Acquire data – Start AI. The `drawnow` command in `daqdoc5_7plot` forces MATLAB to update the display.

```
start(AI)  
pause(duration+0.5)
```

5. Clean up – When you no longer need AI, you should remove it from memory and from the MATLAB workspace.

```
delete(AI)  
clear AI
```

Linearly Scaling the Data: Engineering Units

The Data Acquisition Toolbox provides you with a way to linearly scale analog input signals from your sensor. You can associate this scaling with specific engineering units such as volts or Newtons that you may want to apply to your data. When specifying engineering units, there are three important considerations:

- The expected data range produced by your sensor. This range depends on the physical phenomena you are measuring and the maximum output range of the sensor.
- The range of your analog input hardware. For many devices, this range is specified by the gain and polarity. You can return valid input ranges with the `daqwinfo` function.
- The engineering units associated with your acquisition. By default, most analog input hardware converts data to voltage values. However, after the data is digitized, you may want to define a linear scaling that represents specific engineering units when data is returned to MATLAB.

The properties associated with engineering units and linearly scaling acquired data are given below.

Table 5-13: Analog Input Engineering Units Properties

Property Name	Description
SensorRange	Specify the range of data you expect from your sensor.
InputRange	Specify the range of the analog input subsystem.
Units	Specify the engineering units label.
UnitsRange	Specify the range of data as engineering units.

Note If supported by the hardware, you can set the engineering units properties on a per-channel basis. Therefore, you can configure different engineering unit conversions for each hardware channel.

Linearly scaled acquired data is given by the formula

$$\text{scaled value} = (\text{A/D value})(\text{units range})/(\text{sensor range})$$

Note The above formula assumes you are using symmetric units range and sensor range values, and represents the simplest scenario. If your units range or sensor range values are asymmetric, then you must also include the appropriate offset in the formula.

The A/D value is constrained by the `InputRange` property, which reflects the gain and polarity of your hardware channels, and is usually returned as a voltage value. You should choose an input range that utilizes the maximum dynamic range of your A/D subsystem. The best input range is the one that most closely encompasses the expected sensor range. If the sensor signal is larger than the input range, then the hardware will usually clip (saturate) the signal.

The units range is given by the `UnitsRange` property, while the sensor range is given by the `SensorRange` property. `UnitsRange` is specified as a voltage value, while `SensorRange` is specified as an engineering unit such as Newtons or g's (1 g = 9.80 m/s/s). These property values control the scaling of data when it is extracted from the engine with the `getdata` function. You can find the appropriate units range and sensor range from your sensor's specification sheet.

For example, suppose `SensorRange` is [- 1 1] and `UnitsRange` is [- 10 10]. If an A/D value is 2.5, then the scaled value is (2.5)(20/2) or 25, in the appropriate units.

Example: Performing a Linear Conversion

This example illustrates how to configure the engineering units properties for an analog input object connected to a National Instruments PCI-6024E board.

An accelerometer is connected to a device which is undergoing a vibration test. Your job is to measure the acceleration and the frequency components of the device while it is vibrating. The accelerometer has a range of ± 50 g's, a voltage sensitivity of 99.7 mV/g, and a resolution of 0.00016 g.

The accelerometer signal is input to a Tektronix TDS 210 digital oscilloscope, and to channel 0 of the data acquisition board. By observing the signal on the scope, the maximum expected range of data from the sensor is ± 200 mV, which corresponds to approximately ± 2 g's. Given this constraint, you should configure the board's input range to the ± 500 mV, which is the closest input range that encompasses the expected data range.

You can run this example by typing `daqdoc5_8` at the MATLAB command line.

1. Create a device object – Create the analog input object `AI` for a National Instruments board. The installed adaptors and hardware IDs are found with `daqhwinfo`.

```
AI = analoginput('ni daq', 1);
```

2. Add channels – Add one hardware channel to `AI`.

```
chan = addchannel(AI, 0);
```

3. Configure property values – Configure the sampling rate to 200 kHz and define a two second acquisition.

```
duration = 2;
ActualRate = setverify(AI, 'SampleRate', 200000);
set(AI, 'SamplesPerTrigger', duration*ActualRate)
```

Configure the engineering units properties. This example assumes you are using a National Instruments PCI-6024E board or an equivalent hardware device. `SensorRange` is set to the maximum accelerometer range in volts, and `UnitsRange` is set to the corresponding range in g's.

```
set(chan, 'SensorRange', [-5 5])
set(chan, 'InputRange', [-0.5 0.5])
set(chan, 'UnitsRange', [-50 50])
set(chan, 'Units', 'g' 's (1g = 9.80 m/s/s)')
```

4. Acquire data – Start the acquisition and wait until all the data is acquired.

```
start(AI)
while strcmp(AI.Running, 'On')
end
```

Extract and plot all the acquired data.

```
data = getdata(AI);  
subplot(2, 1, 1), plot(data)
```

Calculate and display the frequency information.

```
Fs = ActualRate;  
blocksize = duration*ActualRate;  
[f, mag] = daqdocfft(data, Fs, blocksize);  
subplot(2, 1, 2), plot(f, mag)
```

5. Clean up – When you no longer need AI, you should remove it from memory and from the MATLAB workspace.

```
delete(AI)  
clear AI
```


Analog Output

Overview	6-2
Getting Started with Analog Output	6-3
Creating an Analog Output Object	6-3
Adding Channels to an Analog Output Object	6-4
Configuring Analog Output Properties	6-6
Outputting Data	6-9
Analog Output Examples	6-10
Evaluating the Analog Output Object Status	6-14
Managing Output Data	6-17
Queuing Data with putdata	6-17
Example: Queuing Data with putdata	6-19
Configuring Analog Output Triggers	6-21
Defining a Trigger: Trigger Types	6-22
Executing the Trigger	6-23
How Many Triggers Occurred?	6-23
When Did the Trigger Occur?	6-24
Device-Specific Hardware Triggers	6-25
Configuring Events and Actions	6-27
Event Types	6-27
Recording and Retrieving Event Information	6-29
Examples: Using Action Properties and Action Functions	6-31
Linearly Scaling the Data: Engineering Units	6-34
Example: Performing a Linear Conversion	6-35
Starting Multiple Device Objects	6-37

Overview

Analog output subsystems convert digital data stored on your computer to a real-world analog signal. These devices perform the inverse conversion of analog input subsystems. Typical plug-in acquisition boards offer two output channels with 12 bits of resolution, with special hardware available to support multiple channel analog output operations. The Data Acquisition Toolbox provides access to analog output subsystems through an analog output object.

The purpose of this chapter is to show you how to perform data acquisition tasks using your analog output hardware. Topics include:

- Getting started with analog output
- Managing output data
- Configuring analog output triggers
- Configuring events and actions
- Linearly scaling the data
- Starting multiple device objects

Getting Started with Analog Output

The purpose of this section is to show you how to use the Data Acquisition Toolbox to perform basic tasks with your analog output (AO) hardware. This is accomplished by describing the most important properties and functions required for an analog output data acquisition session. In addition, several device-specific examples are provided as well as ways to evaluate the status of the analog output object.

After reading this section, you will be able to perform basic analog output tasks suited to your own data acquisition applications.

Creating an Analog Output Object

You create an analog output object with the `analogoutput` function. `analogoutput` accepts the adaptor name and the hardware device ID as input arguments. For a list of supported adaptors, refer to “The Hardware Driver Adaptor” on page 2-7. The device ID refers to the number associated with your board when it is installed. Some vendors refer to the device ID as the device number or the board number. The device ID is optional for sound cards with an ID of 0. Use the `daqhwinfo` function to determine the available adaptors and device IDs.

Each analog output object is associated with one board and one analog output subsystem. For example, to create an analog output object associated with a National Instruments board with device ID 1

```
ao = analogoutput('ni daq', 1);
```

The analog output object `ao` now exists in the MATLAB workspace. You can display the class of `ao` with the `whos` command.

```
whos ao
  Name      Size      Bytes  Class
  ao        1x1        1334   analogoutput object

Grand total is 53 elements using 1334 bytes
```

Once the analog output object is created, the properties listed below are automatically assigned values. These general purpose properties provide descriptive information about the object based on its class type and adaptor.

Table 6-1: Descriptive Analog Output Properties

Property Name	Description
Name	Specify a descriptive name for the device object.
Type	Indicate the device object type.

You can display the values of these properties for `ao` with the `get` function.

```
get ( ao, { ' Name' , ' Type' } )
ans =
    ' ni daq1- A0'      ' Anal og  Out put '
```

Adding Channels to an Analog Output Object

After creating the analog output object, you must add hardware channels to it. As shown by the figure in “Adding Channels or Lines” on page 3-8, you can think of a device object as a container for channels. The collection of channels contained by the device object is referred to as a *channel group*. As described in “Mapping Hardware Channel IDs to MATLAB Indices” on page 3-9, a channel group consists of a mapping between hardware channel IDs and MATLAB indices (see below).

When adding channels to an analog output object, you must follow these rules:

- The channels must reside on the same hardware device. You cannot add channels from different devices, or from different subsystems on the same device.
- The channels must be sampled at the same rate.

You add channels to an analog output object with the `addchannel` function. `addchannel` requires the device object and at least one hardware channel ID as input arguments. You can optionally specify MATLAB indices, descriptive channel names, and an output argument. For example, to add two hardware channels to the device object `ao` created in the preceding section

```
chans = addchannel ( ao, 0: 1 ) ;
```

The output argument `chans` is a *channel object* that reflects the channel array contained by `ao`. You can display the class of `chans` with the `whos` command.

```
whos chans
      Name      Size      Bytes  Class
      ----      -
chans      2x1      512  aochannel object
```

```
Grand total is 7 elements using 512 bytes
```

You can use `chans` to easily access channels. For example, you can easily configure or return property values for one or more channels. As described in “Referencing Individual Hardware Channels” on page 4-6, you can also access channels with the `Channel` property.

Once you add channels to an analog output object, the properties listed below are automatically assigned values. These properties provide descriptive information about the channels based on their class type and ID.

Table 6-2: Descriptive Analog Output Channel Properties

Property Name	Description
<code>HwChannel</code>	Specify the hardware channel ID.
<code>Index</code>	Indicate the MATLAB index of a hardware channel.
<code>Parent</code>	Indicate the parent (device object) of a channel.
<code>Type</code>	Indicate a channel.

You can display the values of these properties for `chans` with the `get` function.

```
get(chans, {'HwChannel', 'Index', 'Parent', 'Type'})
ans =
      [0]      [1]      [1x1 analogoutput]      'Channel'
      [1]      [2]      [1x1 analogoutput]      'Channel'
```

To reference individual channels, you must specify either MATLAB indices or descriptive channel names. Refer to “Referencing Individual Hardware Channels” on page 4-6 for more information.

Configuring Analog Output Properties

After hardware channels are added to the analog output object, you should configure property values. As described in “Configuring and Returning Properties” on page 3-12, the Data Acquisition Toolbox supports two basic types of properties for analog output objects: common properties and channel properties. Common properties apply to all channels contained by the device object while channel properties apply to individual channels.

The properties you configure depend on your particular analog output application. For many common applications, there is a small group of properties related to the basic setup that you will typically use. These basic setup properties control the sampling rate and define the trigger type. Analog output properties related to the basic setup are given below.

Table 6-3: Analog Output Basic Setup Properties

Property Name	Description
Sampl eRate	Specify the per-channel rate at which digital data is converted to analog data.
Tri ggerType	Specify the type of trigger to execute.

Setting the Sampling Rate

You control the rate at which an analog output subsystem converts digital data to analog data is controlled with the `Sampl eRate` property. `Sampl eRate` must be specified as samples per second. For example, to set the sampling rate for each channel of your National Instruments board to 100,000 samples per second (100 kHz)

```
ao = anal ogout put ( ' ni daq' , 1 ) ;  
addchannel ( ao, 0: 1 ) ;  
set ( ao, ' Sampl eRate' , 100000)
```

Data acquisition boards typically have predefined sampling rates that you can set. If you specify a sampling rate that does not match one of these predefined values, there are two possibilities:

- If the rate is within the range of valid values, then the engine automatically selects a valid sampling rate. The rules governing this selection process are described in the `SampleRate` reference pages in Chapter 10, “Base Property Reference.”
- If the rate is outside the range of valid values, then an error is returned.

Note For some sound cards, you can set the sampling rate to any value between the minimum and maximum values defined by the hardware. You can enable this feature with the `StandardSampleRates` property. Refer to Chapter 11, “Device-Specific Property Reference” for more information.

Most analog output subsystems allow simultaneous sampling of channels. Therefore, the maximum sampling rate for each channel is given by the maximum board rate.

After setting a value for `SampleRate`, you should find out the actual rate set by the engine.

```
ActualRate = get(ao, 'SampleRate');
```

Alternatively, you can use the `setverify` function, which sets a property value and returns the actual value set.

```
ActualRate = setverify(ao, 'SampleRate', 100000);
```

You can find the range of valid sampling rates for your hardware with the `propinfo` function.

```
ValidRates = propinfo(ao, 'SampleRate');
ValidRates.ConstraintValue
ans =
    1.0e+005 *
    0.0000    2.0000
```

Defining a Trigger

For analog output objects, a trigger is defined as an event that initiates the output of data from the engine to the analog output hardware.

Defining a trigger for an analog output object involves specifying the trigger type. Trigger types are specified with the `TriggerType` property. The valid `TriggerType` values that are supported for all hardware are given below.

Table 6-4: Analog Output `TriggerType` Property Values

TriggerType Values	Description
{Immediate}	The trigger occurs just after you issue the start function.
Manual	The trigger occurs just after you manually issue the trigger function.

Most devices have hardware-specific trigger types, which are available to you through the `TriggerType` property. For example, to see all the trigger types (including hardware-specific trigger types) for the analog output object `ao` created in the preceding section

```
set (ao, 'TriggerType')
[ Manual | {Immediate} | HwDigital ]
```

This information tells you that the National Instruments board also supports a hardware digital trigger. For a description of device-specific trigger types, refer to “Device-Specific Hardware Triggers” on page 6-25, or the `TriggerType` reference pages in Chapter 10, “Base Property Reference.”

Outputting Data

After you configure the analog output object, you can output data. Outputting data involves these three steps:

- 1 Queuing data
- 2 Starting the analog output object
- 3 Stopping the analog output object

Queuing Data in the Engine

Before you can start the device object, data must be queued in the engine. Data is queued in the engine with the `putdata` function. For example, to queue one second of data for each channel contained by the analog output object `ao`

```
ao = analogoutput('winsound');
addchannel(ao, 1:2);
data = sin(linspace(0, 2*pi, 8000))';
putdata(ao, [data data])
```

`putdata` is a blocking function, and will not return execution control to MATLAB until the specified data is queued. `putdata` is described in detail in “Managing Output Data” on page 6-17 and in Chapter 9, “Function Reference.”

Starting the Analog Output Object

You start an analog output object with the `start` function. For example, to start the analog output object `ao`

```
start(ao)
```

After `start` is issued, the `Running` property is automatically set to `On`, and both the device object and hardware device execute according to the configured and default property values. While the device object is running, you can continue to queue data.

However, running does not necessarily mean that data is being output from the engine to the analog output hardware. For that to occur, a trigger must execute. When the trigger executes, the `Sending` property is automatically set to `On`. Analog output triggers are described on “Defining a Trigger” on page 6-8 and “Configuring Analog Output Triggers” on page 6-21.

Stopping the Analog Output Object

An analog output object can stop under one of these conditions:

- You issue the stop function.
- The queued data is output.
- A run-time hardware error occurs.
- A timeout occurs.

When the device object stops, the `Running` and `Sending` properties are automatically set to `Off`. At this point, you can reconfigure the device object or immediately queue more data, and issue another `start` command using the current configuration.

Analog Output Examples


This section illustrates how to perform basic data acquisition tasks using analog output subsystems and the Data Acquisition Toolbox. For most data acquisition applications using analog output subsystems, you must follow these basic steps:

- 1 Install and connect the components of your data acquisition hardware. At a minimum, this involves connecting an actuator to a plug-in or external data acquisition device.
- 2 Configure your data acquisition session. This involves creating a device object, adding channels, setting property values, and using specific functions to output data.

Simple data acquisition applications using a sound card and a National Instruments board are given below.

Outputting Data with a Sound Card

In this example, sine wave data is generated in MATLAB, output to the D/A converter on the sound card, and sent to a speaker. The setup is shown below.



You can run this example by typing `daqdoc6_1` at the MATLAB command line.

1. Create a device object – Create the analog output object `A0` for a sound card. The installed adaptors and hardware IDs are found with `daqhwinfo`.

```
A0 = analogoutput('winsound');
```

2. Add channels – Add one channel to `A0`.

```
chan = addchannel(A0, 1);
```

3. Configure property values – Define an output time of four seconds, assign values to the basic setup properties, generate data to be queued, and queue the data with one call to `putdata`.

```
duration = 4;
set(A0, 'SampleRate', 8000)
set(A0, 'TriggerType', 'Manual')
ActualRate = get(A0, 'SampleRate');
len = ActualRate*duration;
data = sin(linspace(0, 2*pi, len))';
putdata(A0, data)
```

4. Output data – Start A0, issue a manual trigger, and wait for the device object to stop running.


```
start(A0)
trigger(A0)
while strcmp(A0.Running, 'On')
end
```

5. Clean up – When you no longer need A0, you should remove it from memory and from the MATLAB workspace.

```
delete(A0)
clear A0
```

Outputting Data with a National Instruments Board

In this example, sine wave data is generated in MATLAB, output to the D/A converter on a National Instruments board, and displayed with an oscilloscope. The setup is shown below.



You can run this example by typing `daqdoc6_2` at the MATLAB command line.

1. Create a device object – Create the analog output object A0 for a National Instruments board. The installed adaptors and hardware IDs are found with `daqhwinfo`.

```
A0 = analogoutput('ni daq', 1);
```

2. Add channels – Add one channel to A0.

```
chan = addchannel(A0, 0);
```

3. Configure property values – Define an output time of four seconds, assign values to the basic setup properties, generate data to be queued, and queue the data with one call to putdata.

```
duration = 4;
set(A0, 'SampleRate', 10000)
set(A0, 'TriggerType', 'Manual')
ActualRate = get(A0, 'SampleRate');
len = ActualRate*duration;
data = sin(linspace(0, 2*pi, len))';
putdata(A0, data)
```

4. Output data – Start A0, issue a manual trigger, and wait for the device object to stop running.

```
start(A0)
trigger(A0)
while strcmp(A0.Running, 'On')
end
```

5. Clean up – When you no longer need A0, you should remove it from memory and from the MATLAB workspace.

```
delete(A0)
clear A0
```

Evaluating the Analog Output Object Status

You can evaluate the status of an analog output object by:

- Returning the values of certain properties
- Invoking the display summary

Status Properties

The properties associated with the status of your analog output object allow you to evaluate:

- If the device object is running
- If data is being output from the engine
- How much data is queued in the engine
- How much data has been output from the engine

These properties are given below.

Table 6-5: Analog Output Status Properties

Property Name	Description
Runni ng	Indicate if the device object is running.
Sampl esAvai l abl e	Indicate the number of samples available per channel in the engine.
Sampl esOut put	Indicate the number of samples output per channel from the engine.
Sendi ng	Indicate if data is being sent (output) to the hardware device.

When data is queued in the engine, `Sampl esAvai l abl e` is updated to reflect the total number of samples per channel that was queued. When `start` is issued, `Runni ng` is automatically set to `0n`.

When the trigger executes, `Sendi ng` is automatically set to `0n` and `Sampl esOut put` keeps a running count of the total number of samples per channel output from the engine to the hardware. Additionally,

`SamplesAvailable` tells you how many samples per channel are still queued in the engine and ready to be output to the hardware.

When all the queued data is output from the engine, both `Running` and `Sending` are automatically set to `Off`, `SamplesAvailable` is 0, and `SamplesOutput` reflects the total number of samples per channel that was output.

The Display Summary

You can invoke the display summary by typing the analog output object at the command line. The information displayed reflects many of the basic setup properties described in “Configuring Analog Output Properties” on page 6-6, and is designed so you can quickly evaluate the status of your data acquisition session. The display is divided into two main sections: general summary information and channel summary information.

General Summary Information

The general display summary includes the device object type and the hardware device name, followed by this information:

- Output parameters – The sampling rate
- Trigger parameters – The trigger type
- The engine status
 - Whether the engine is sending data, waiting to start, or waiting to trigger
 - The total time required to output the queued data
 - The number of samples queued by putdata
 - The number of samples sent to the hardware device

Channel Summary Information

The channel display summary includes property values associated with:

- The hardware channel mapping
- The channel name
- The engineering units

The display summary shown below is for the example given in “Outputting Data with a Sound Card” on page 6-11 prior to issuing the start function.

General display summary

Display Summary of Analog Output (AO) Object Using 'AudioPCI Playback'.

Output Parameters: 8000 samples per second on each channel.

Trigger Parameters: 1 'Immediate' trigger on START.

Engine status: Waiting for START.

0 total sec. of data currently queued for START

0 samples currently queued by PUTDATA.

0 samples sent to output device since START.

Channel display summary

A0 object contains channel(s):

Index:	Channel Name:	HwChannel:	OutputRange:	UnitsRange:	Units:
1	'Mono'	1	[- 1 1]	[- 1 1]	'Volts'

You can use the Channel property to display only the channel summary information.

A0.Channel

Managing Output Data

At the core of any analog output application lies the data you want to send from a computer to an output device such as an actuator. The role of the analog output subsystem is to convert digitized data to analog data for subsequent output.

Before you can output data to the analog output subsystem, it must be queued in the engine. Queuing data is managed with the `putdata` function. In addition to this function, there are several properties associated with managing output data. These properties are given below.

Table 6-6: Analog Output Data Management Properties

Property Name	Description
<code>MaxSamplesQueued</code>	Indicate the maximum number of samples that can be queued in the engine.
<code>RepeatOutput</code>	Specify the number of additional times queued data is output.
<code>Timeout</code>	Specify an additional waiting time to queue data.

Queuing Data with `putdata`

Before data can be sent to the analog output hardware, you must queue it in the engine. Queuing data is managed with the `putdata` function. One column of data is required for each channel contained by the analog output object. For example, to queue one second of data for each channel contained by the analog output object `ao`

```
ao = analogoutput('winsound');  
addchannel(ao, 1:2);  
data = sin(linspace(0, 2*pi, 8000))';  
putdata(ao, [data data])
```

A data source consisting of m samples and n channels is illustrated below.

$$\begin{bmatrix} d_{11} & d_{12} & \dots & d_{1n} \\ d_{21} & d_{22} & & d_{2n} \\ d_{31} & d_{32} & \dots & d_{3n} \\ \vdots & \vdots & & \vdots \\ \vdots & \vdots & & \vdots \\ \vdots & \vdots & & \vdots \\ d_{m1} & d_{m2} & & d_{mn} \end{bmatrix}$$

Data source. Each column represents a separate output channel.

Rules for Using putdata

Using `putdata` to queue data in the engine follows these rules:

- You must queue data in the engine before starting the analog output object.
- If the value of the `RepeatOutput` property is greater than 0, then all queued data is automatically requeued until the `RepeatOutput` value is reached. You must configure `RepeatOutput` before `start` is issued.
- While the analog output object is running, you can continue to queue data unless `RepeatOutput` is greater than 0.
- You can queue data in the engine until the value specified by the `MaxSamplesQueued` property is reached, or the limitations of your hardware or computer are reached.

For more information about `putdata`, refer to its reference pages in Chapter 9, “Function Reference.”

Rules for Queuing Data

Data to be queued in the engine follows these rules:

- Data is output as soon as a trigger occurs.
- An error is returned if a NaN is included in the data stream.

- You can use the native data type of the hardware. Note that MATLAB supports math operations only for the double data type. Therefore, to use math functions on native data, you must convert it to doubles.
- If the data is not within the range of the `UnitsRange` property, then it is clipped to the maximum or minimum value specified by `UnitsRange`. Refer to “Linearly Scaling the Data: Engineering Units” on page 6-34 for more information about clipping.

Example: Queuing Data with `putdata`

This example illustrates how you can use `putdata` to queue 8000 samples, and then output the data a total of five times using the `RepeatOutput` property.

You can run this example by typing `daqdoc6_3` at the MATLAB command line.

1. Create a device object – Create the analog output object `A0` for a sound card. The installed adaptors and hardware IDs are found with `daqhwinfo`.

```
A0 = analogoutput('winsound');
%A0 = analogoutput('ni daq', 1);
%A0 = analogoutput('cbi', 1);
```

2. Add channels – Add one channel to `A0`.

```
chans = addchannel(A0, 1);
%chans = addchannel(A0, 0); % For NI and CBI
```

3. Configure property values – Define an output time of one second, assign values to the basic setup properties, generate data to be queued, and issue two `putdata` calls. Since the queued data is repeated four times and two `putdata` calls are issued, a total of 10 seconds of data is output.

```
duration = 1;
set(A0, 'SampleRate', 8000)
ActualRate = get(A0, 'SampleRate');
len = ActualRate*duration;
set(A0, 'RepeatOutput', 4)
data = sin(linspace(0, 2*pi, len));
putdata(A0, data)
putdata(A0, data)
```

4. Output data – Start A0 and wait for the device object to stop running.


```
start(A0)
while strcmp(A0. Running, ' On' )
end
```

5. Clean up – When you no longer need A0, you should remove it from memory and from the MATLAB workspace.

```
delete(A0)
clear A0
```

Configuring Analog Output Triggers

An analog output trigger is defined as an event that initiates the output of data. As shown in the figure below, when a trigger occurs, the `Sendi ng` property is automatically set to `On` and queued data is output from the engine to the hardware subsystem.



Properties associated with analog output triggers are given below.

Table 6-7: Analog Output Trigger Properties

Property Name	Description
<code>Ini ti al Tri ggerTi me</code>	Indicate the absolute time of the first trigger.
<code>Tri ggerActi on</code>	Specify the M-file action function to execute when a trigger occurs.
<code>Tri ggersExecuted</code>	Indicate the number of triggers that execute.
<code>Tri ggerType</code>	Specify the type of trigger to execute.

Except for `Tri ggerActi on`, these trigger-related properties are discussed in the following sections. `Tri ggerActi on` is discussed in “Configuring Events and Actions” on page 6-27.

Defining a Trigger: Trigger Types

Defining a trigger for an analog output object involves specifying the trigger type with the `TriggerType` property. You can think of the trigger type as the source of the trigger. The analog output `TriggerType` values are given below.

Table 6-8: Analog Output `TriggerType` Property Values

TriggerType Value	Description
<code>{Immediate}</code>	The trigger occurs just after you issue the <code>start</code> function.
<code>Manual</code>	The trigger occurs just after you manually issue the <code>trigger</code> function.

Trigger types can be grouped into two main categories:

- Device-independent triggers
- Device-specific hardware triggers

The trigger types shown above are device-independent triggers since they are available for all supported hardware. For these trigger types, the action that initiates the trigger event involves issuing a toolbox function (`start` or `trigger`). Conversely, device-specific hardware triggers depend on the specific hardware device you are using. For these trigger types, the action that initiates the trigger event involves an external digital signal.

Device-specific hardware triggers for National Instruments and Agilent Technologies devices are discussed in “Device-Specific Hardware Triggers” on page 6-25. Device-independent triggers are discussed below.

Immediate Trigger. If `TriggerType` is `Immediate` (the default value), the trigger occurs immediately after the `start` function is issued. You can configure an analog output object for continuous output, by using an immediate trigger and setting `RepeatOutput` to `inf`.

Manual Trigger. If `TriggerType` is `Manual`, the trigger occurs immediately after the `trigger` function is issued.

Executing the Trigger

For an analog output trigger to occur, you must follow these steps:

- 1 Queue data in the engine.
- 2 Configure the appropriate trigger properties.
- 3 Issue the start function.
- 4 Issue the `trigger` function if `TriggerType` is `Manual`.

Once the trigger occurs, queued data is output to the hardware, and the device object stops executing when all the queued data is output.

Note Only one trigger event can occur for analog output objects.

How Many Triggers Occurred?

For analog output objects, only one trigger can occur. You can determine if the trigger event occurred by returning the value of the `TriggersExecuted` property. If `TriggersExecuted` is 0, then the trigger event did not occur. If `TriggersExecuted` is 1, then the trigger event occurred. Event information is also recorded by the `EventLog` property. A convenient way to access event log information is with the `showdaqevents` function.

For example, suppose you create the analog output object `ao` for a sound card and add one channel to it. `ao` is configured to output 8,000 samples using the default sampling rate of 8000 Hz.

```
ao = analogoutput('winsound');  
addchannel(ao, 1);  
data = sin(linspace(0, 1, 8000))';  
putdata(ao, data)  
start(ao)
```

`TriggersExecuted` returns the number of triggers executed.

```
ao.TriggersExecuted
ans =
    1
```

You can use `showdaqevents` to return information for all events that occurred while `ao` was executing.

```
showdaqevents(ao)
    1 Start          ( 10: 43: 25, 0 )
    2 Trigger        ( 10: 43: 25, 0 )
    3 Stop           ( 10: 43: 26, 8000 )
```

For more information about recording and retrieving event information, refer to “Recording and Retrieving Event Information” on page 6-29.

When Did the Trigger Occur?

You can return the absolute time of the trigger with the `InitialTriggerTime` property. Absolute time is returned as a clock vector in the form

[year month day hour minute seconds]

For example, the absolute time of the trigger event for the preceding example is

```
abstime = ao.InitialTriggerTime
abstime =
1.0e+003 *
    1.9990    0.0040    0.0170    0.0100    0.0430    0.0252
```

To convert the clock vector to a more convenient form, you can use the `sprintf` function.

```
t = fix(abstime);
sprintf('%d: %d: %d', t(4), t(5), t(6))
ans =
10: 43: 25
```

As shown in the preceding section, you can also evaluate the absolute time of the trigger event with the `showdaqevents` function.

Device-Specific Hardware Triggers

Most data acquisition devices possess the ability to accept a hardware trigger. Hardware triggers are processed directly by the hardware and are typically transistor-transistor logic (TTL) signals. Hardware triggers are used when speed is required since a hardware device can process an input signal much faster than software.

The device-specific hardware triggers are presented to you as additional property values. Hardware triggers for National Instruments and Agilent Technologies devices are discussed below and in Chapter 10, “Base Property Reference.”

Note that the available hardware trigger support depends on the board you are using. Refer to your hardware documentation for detailed information about its triggering capabilities.

National Instruments

When using National Instruments hardware, there is an additional analog output trigger type available to you: digital triggering.

If `TriggerType` is set to `HwDigital`, the trigger is given by an external TTL signal that is input directly into the hardware device. The following example illustrates how to configure a hardware digital trigger.

```
ao = analogoutput('ni daq', 1);  
addchannel(ao, 0:1);  
set(ao, 'TriggerType', 'HwDigital')
```

With this trigger configuration, `ao` will not start outputting data until the TTL signal is detected by the hardware on the appropriate pin.

The diagram below illustrates how you can connect a digital trigger signal to an MIO-16E Series board. PFI6/WFTRIG corresponds to pin 5.



MIO-16E Series board

Agilent Technologies

When using Agilent Technologies hardware, there are additional analog output trigger types that you must be aware of: digital triggering on a positive edge and digital triggering on a negative edge.

If `TriggerType` is `HwDigitalPos`, the trigger source is the positive edge of a digital signal. If `TriggerType` is `HwDigitalNeg`, the trigger source is the negative edge of a digital signal.

In both cases, the digital signal is an external TTL signal that is input directly into the hardware device. The example below illustrates how to configure such a trigger.

```
ao = analogoutput('hpe1432', 8);
addchannel(ao, 1);
set(ao, 'TriggerType', 'HwDigitalPos')
```

With this trigger configuration, `ao` will not start outputting data until the TTL signal is detected by the hardware.

Configuring Events and Actions

You can enhance the power and flexibility of your analog output application by utilizing *events*. An event occurs at a particular time after a condition is met and may result in one or more actions.

While the analog output object is running, you can use events to display a message, display data, analyze data, or perform just about any other action. Actions are controlled through *action properties* and *action functions*. All event types have an associated action property. Action functions are M-file functions that you construct to suit your specific data acquisition needs.

You execute an action when a particular event occurs by specifying the name of the M-file action function as the value for the associated action property. Refer to “Creating and Executing Action Functions” on page 5-52 to learn how to create action functions. Note that `daqacti on` is the default value for some action properties.

Event Types

The analog output event types and associated action properties are described below.

Table 6-9: Analog Output Action Properties

Event Type	Property Name
Run-time error	<code>Runti meErrorActi on</code>
Samples output	<code>Sampl esOutputActi on</code>
	<code>Sampl esOutputActi onCount</code>
Start	<code>StartActi on</code>
Stop	<code>StopActi on</code>
Timer	<code>Ti merActi on</code>
	<code>Ti merPeri od</code>
Trigger	<code>Tri ggerActi on</code>

Run-time Error Event. A run-time error event is generated immediately after a run-time error occurs. This event executes the action function specified for `RunTimeErrorAction`. Additionally, a toolbox error message is automatically displayed to the MATLAB workspace. If an error occurs that is not explicitly handled by the toolbox, then the hardware-specific error message is displayed.

The default value for `RunTimeErrorAction` is `daqaction`, which displays the event type, the time the event occurred, the device object name, and the error message.

Run-time errors include hardware errors and timeouts. Run-time errors do not include configuration errors such as setting an invalid property value.

Samples Output Event. A samples output event is generated immediately after the number of samples specified by the `SamplesOutputActionCount` property is output for each channel group member. This event executes the action function specified for `SamplesOutputAction`.

Start Event. A start event is generated immediately after the start function is issued. This event executes the action function specified for `StartAction`. When the action function has finished executing, `Runni ng` is automatically set to `On` and the device object and hardware device begin executing.

Stop Event. A stop event is generated immediately after the device object and hardware device stop running. This occurs when:

- The stop function is issued.
- The requested number of samples is output.
- A run-time error occurs.

A stop event executes the action function specified for `StopAction`. Under most circumstances, the action function is not guaranteed to complete execution until sometime after the device object and hardware device stop running, and the `Runni ng` property is set to `Off`.

Timer Event. A timer event is generated whenever the time specified by the `TimerPeriod` property passes. This event executes the action function specified for `TimerAction`. Time is measured relative to when the device object starts running.

Some timer events may not be processed if your system is significantly slowed or if the `TimerPeriod` value is too small. For example, a common application for

timer events is to display data. However, since displaying data is a CPU-intensive task, some of these events may be dropped. To guarantee that events are not dropped, you may want to use the `SamplesOutputAction` property.

Trigger Event. A trigger event is generated immediately after a trigger occurs. This event executes the action function specified for `TriggerAction`. Under most circumstances, the action function is not guaranteed to complete execution until sometime after `Sending` is set to `On`.

Recording and Retrieving Event Information

While the analog output object is running, certain information is automatically recorded in the `EventLog` property for some of the event types listed in the preceding section. `EventLog` is a structure that contains two fields: `Type` and `Data`. The `Type` field contains the event type. The `Data` field contains event-specific information. Events are recorded in the order in which they occur. The first `EventLog` entry reflects the first event recorded, the second `EventLog` entry reflects the second event recorded, and so on.

The event types recorded in `EventLog` for analog output objects, as well as the values for the `Type` and `Data` fields, are given below.

Table 6-10: Analog Output Event Information Stored in EventLog

Event Type	Type Field Value	Data Field Value
Run-time error	Error	AbsTime
		Rel Sample
		String
Start	Start	AbsTime
		Rel Sample
Stop	Stop	AbsTime
		Rel Sample
Trigger	Trigger	AbsTime
		Rel Sample

Samples output events and timer events are not stored in EventLog.

Note Unless a run-time error occurs, EventLog records a start event, a trigger event, and stop event for each data acquisition session.

The Data field values are described below.

The AbsTime Field. AbsTime is used by all analog output events stored in EventLog to indicate the absolute time the event occurred. The absolute time is returned using MATLAB's clock format.

day-month-year hour:minute:second

The RelSample Field. RelSample is used by all events stored in EventLog to indicate the sample number that was output when the event occurred. RelSample is 0 for the start event and for the first trigger event regardless of the trigger type. RelSample is NaN for any event that occurs before the trigger executes.

The String Field. String is used by the run-time error event to store the descriptive message that is generated when a run-time error occurs. This message is also displayed at the MATLAB command line.

Example: Retrieving Event Information

Suppose you want to examine the events logged for the example given by “Example: Queuing Data with putdata” on page 6-19. You can do this by accessing the EventLog property.

```
events = A0.EventLog
events =
3x1 struct array with fields:
    Type
    Data
```

By examining the contents of the Type field, you can list the events that were recorded while A0 was running.

```
{events.Type}
ans =
    'Start'    'Trigger'    'Stop'
```

To display information about the trigger event, you must access the `Data` field which stores the absolute time the trigger occurred and the number of samples output when the trigger occurred.

```
trigdata = events(2).Data
trigdata =
    AbsTime: [1999 4 16 9 53 19.9508]
    Rel Sample: 0
```

Alternatively, you can use the `showdaqevents` function to display event information.

```
showdaqevents(events, 2)
2 Trigger ( 09:53:19, 0 )
```

Examples: Using Action Properties and Action Functions

Examples showing how to create action functions and configure action properties are given below.

Displaying the Number of Samples Output

This example illustrates how to generate samples output events. You can run this example by typing `daqdoc6_4` at the MATLAB command line. The `daqdoc6_4disp` action function displays the number of events that were output from the engine whenever the samples output event occurred.

1. Create a device object – Create the analog output object `A0` for a sound card. The installed adaptors and hardware IDs are found with `daqhwinfo`.

```
A0 = analogoutput('winsound');
%A0 = analogoutput('ni daq', 1);
%A0 = analogoutput('cbi', 1);
```

2. Add channels – Add two channels to `A0`.

```
chans = addchannel(A0, 1:2);
%chans = addchannel(A0, 0:1); % For NI and CBI
```

3. Configure property values – Configure the trigger to repeat four times, specify daqdoc6_4di sp as the M-file action function to execute whenever 8000 samples are output, generate data to be queued, and queue the data with one call to putdata.

```
set(A0, 'SampleRate', 8000)
ActualRate = get(A0, 'SampleRate');
set(A0, 'RepeatOutput', 4)
set(A0, 'SamplesOutputActionCount', 8000)
freq = get(A0, 'SamplesOutputActionCount');
set(A0, 'SamplesOutputAction', 'daqdoc6_4di sp')
data = sin(linspace(0, 2*pi, 3*freq));
putdata(A0, [data data])
```

4. Output data – Start A0.

```
start(A0)
```

5. Clean up – When you no longer need A0, you should remove it from memory and from the MATLAB workspace.

```
delete(A0)
clear A0
```

Displaying EventLog Information

This example illustrates how action functions allow you to easily display information stored in the EventLog property. You can run this example by typing daqdoc6_5 at the MATLAB command line. The daqdoc6_5di sp action function displays the absolute time and relative sample associated with the start, trigger, and stop events.

1. Create a device object – Create the analog output object A0 for a sound card. The installed adaptors and hardware IDs are found with daqhwi nfo.

```
A0 = analogoutput('winsound');
%A0 = analogoutput('ni daq', 1);
%A0 = analogoutput('cbi', 1);
```

2. Add channels – Add one channel to A0.

```
chan = addchannel(A0, 1);
%chan = addchannel(A0, 0); % For NI and CBI
```

3. Configure property values – Specify daqdoc6_5di.sp as the M-file action function to execute when the start, trigger, and stop events occur, generate data to be queued, and queue the data with one call to putdata.

```
set(A0, 'SampleRate', 8000)
ActualRate = get(A0, 'SampleRate');
set(A0, 'StartAction', 'daqdoc6_5di.sp')
set(A0, 'TriggerAction', 'daqdoc6_5di.sp')
set(A0, 'StopAction', 'daqdoc6_5di.sp')
data = sin(linspace(0, 2*pi, ActualRate));
data = [data data data];
time = (length(data)/A0.SampleRate);
putdata(A0, data)
```

4. Output data – Start A0.

```
start(A0)
```

5. Clean up – When you no longer need A0, you should remove it from memory and from the MATLAB workspace.

```
delete(A0)
clear A0
```

Linearly Scaling the Data: Engineering Units

The Data Acquisition Toolbox provides you with a way to linearly scale data as it is being queued in the engine. You can associate this scaling with specific engineering units such as volts or Newtons that you may want to apply to your data.

The properties associated with engineering units and linearly scaling output data are given below.

Table 6-11: Analog Output Engineering Units Properties

Property Name	Description
OutputRange	Specify the range of the analog output hardware subsystem.
Units	Specify the engineering units label.
UnitsRange	Specify the range of data as engineering units.

For many devices, the output range is expressed in terms of the gain and polarity.

Note You can set the engineering units properties on a per-channel basis. Therefore, you can configure different engineering unit conversions for each hardware channel.

Linearly scaled output data is given by the formula

$$\text{scaled value} = (\text{original value})(\text{output range})/(\text{units range})$$

The units range is given by the UnitsRange property, while the output range is given by the OutputRange property. UnitsRange controls the scaling of data when it is queued in the engine with the putdata function. OutputRange specifies the gain and polarity of your D/A subsystem. You should choose an output range that encompasses the output signal, and that utilizes the maximum dynamic range of your hardware.

For sound cards, you may have to adjust the volume control to obtain the full-scale output range of the device. Refer to “Troubleshooting Sound Cards” on page A-12 to learn how to access the volume control for your sound card.

For example, suppose `OutputRange` is `[- 10 10]`, and `UnitsRange` is `[- 5 5]`. If a queued value is 2.5, then the scaled value is $(2.5)(20/10)$ or 5, in the appropriate units.

Note The data acquisition engine always *clips* out-of-range values. Clipping means that an out-of-range value is fixed to either the minimum or maximum value that is representable by the hardware. Clipping is equivalent to saturation.

Example: Performing a Linear Conversion

This example illustrates how to configure the engineering units properties for an analog output object connected to a National Instruments PCI-6024E board.

The queued data consists of a 4 volt peak-to-peak sine wave. The `UnitsRange` property is configured so that queued data is scaled to the `OutputRange` property value, which is fixed at ± 10 volts. This scaling utilizes the maximum dynamic range of the analog output hardware.

You can run this example by typing `daqdoc6_6` at the MATLAB command line.

1. Create a device object – Create the analog output object `A0` for a National Instruments board. The installed adaptors and hardware IDs are found with `daqhwinfo`.

```
A0 = analogoutput('ni daq', 1);
```

2. Add channels – Add one hardware channel to `A0`.

```
chan = addchannel(A0, 0);
```

3. Configure property values – Create the data to be queued.

```
freq = 2;
w = 2*pi*freq;
t = linspace(0, 10, 10000);
data = 2*sin(w*t)';
```

Configure the sampling rate to 2 kHz, configure the trigger to repeat two times, and scale the data to cover the full output range of the D/A converter. Since the peak-to-peak amplitude of the queued data is 4, `UnitsRange` is set to `[-2 2]`, which scales the output data to 20 volts peak-to-peak.

```
set(A0, 'SampleRate', 2000)
set(A0, 'RepeatOutput', 2)
set(chan, 'UnitsRange', [-2 2])
```

Queue the data with one call to `putdata`.

```
putdata(A0, data)
```

4. Output data – Start A0 and wait until all the data is output.

```
start(A0)
while strcmp(A0.Running, 'On')
end
```

5. Clean up – When you no longer need A0, you should remove it from memory and from the MATLAB workspace.

```
delete(A0)
clear A0
```

Starting Multiple Device Objects

With the Data Acquisition Toolbox, you can start multiple device objects. You may find this feature useful when simultaneously using your hardware's analog output (AO) and analog input (AI) subsystems. For example, suppose you create the analog input object `ai` and the analog output object `ao` for a sound card, and add one channel to each device object.

```
ai = analoginput('winsound');
addchannel(ai, 1);
ao = analogoutput('winsound');
addchannel(ao, 1);
```

You should use manual triggers when starting multiple device objects since this trigger type executes faster than other trigger types with the exception of hardware triggers. Additionally, to synchronize the input and output of data, you should configure the `ManualTriggerHwOn` property to `Trigger` for `ai`.

```
set([ai ao], 'TriggerType', 'Manual')
set(ai, 'ManualTriggerHwOn', 'Trigger')
```

Configure `ai` for continuous acquisition, call the action function `qmoredata` whenever 1000 samples are output, and call `daqaction` when `ai` and `ao` stop running.

```
set(ai, 'SamplesPerTrigger', inf)
set(ao, 'SamplesOutputAction', {'qmoredata', ai})
set(ao, 'SamplesOutputActionCount', 1000)
set([ai ao], 'StopAction', 'daqaction')
```

As shown below, the action function `qmoredata` extracts data from the engine and then queues it for output.

```
function qmoredata(obj, event, ai)
data = getdata(ai, 1000);
putdata(obj, data)
```

Queue data in the engine, start the device objects, and execute the manual triggers.

```
data = zeros(4000, 1);  
putdata(ao, data)  
start([ai ao])  
trigger([ai ao])
```

Note Device objects cannot execute simultaneously unless you use an external hardware trigger.

You can determine the starting time for each device object with the `InitialTriggerTime` property. The difference, in seconds, between the starting times for `ai` and `ao` is

```
ai time = ai.InitialTriggerTime  
aotime = ao.InitialTriggerTime  
delta = abs(aotime - ai time);  
sprintf('%d', delta(6))  
ans =  
2.288818e-005
```

Note that this number depends on the specific platform you are using. To stop both device objects

```
stop([ai ao])
```

The output from `daqacti on` is shown below.

```
Stop event occurred at 13:00:25 for the object: winsound0-A0.  
Stop event occurred at 13:00:25 for the object: winsound0-AI.
```

Digital Input/Output

Overview	7-2
Creating a Digital I/O Object	7-3
Adding Lines to a Digital I/O Object	7-4
Line and Port Characteristics	7-5
Referencing Individual Hardware Lines	7-8
Writing and Reading Digital I/O Line Values	7-10
Writing Digital Values	7-10
Reading Digital Values	7-12
Example: Writing and Reading Digital Values	7-13
Generating Timer Events	7-15
Timer Events	7-15
Starting and Stopping a Digital I/O Object	7-16
Example: Generating Timer Events	7-17
Evaluating the Digital I/O Object Status	7-18
The Display Summary	7-18

Overview

Digital I/O (DIO) subsystems are designed to input and output digital values to and from hardware. These values are typically handled either as single bits or *lines*, or as a *port*, which typically consists of eight lines. While most popular data acquisition boards include some digital I/O capability, it is usually limited to simple operations and special dedicated hardware is required for performing advanced digital I/O operations. The Data Acquisition Toolbox provides access to digital I/O subsystems through a digital I/O object.

Note The Data Acquisition Toolbox does not directly support buffered digital I/O or handshaking (latching). However, you can write your own M-code to support this functionality. Buffered digital I/O means that the data associated with a digital I/O subsystem is stored in the engine. Handshaking allows your digital I/O subsystem to input or output data after receiving a digital pulse.

The purpose of this chapter is to show you how to perform data acquisition tasks using your digital I/O hardware. Topics include:

- Creating a digital I/O object
- Adding lines to a digital I/O object
- Writing and reading digital I/O line values
- Generating timer events
- Evaluating the digital I/O object status

Creating a Digital I/O Object

You create a digital I/O object with the `digitalio` function. `digitalio` accepts the adaptor name and the hardware device ID as input arguments. For a list of supported adaptors, refer to “The Hardware Driver Adaptor” on page 2-7. The device ID refers to the number associated with your board when it is installed. Some vendors refer to the device ID as the device number or the board number. Use the `daqhwinfo` function to determine the available adaptors and device IDs.

Each digital I/O object is associated with one board and one digital I/O subsystem. For example, to create a digital I/O object associated with a National Instruments board with device ID 1

```
di o = digitalio('ni daq', 1);
```

The digital I/O object `di o` now exists in the MATLAB workspace. You can display the class of `di o` with the `whos` command.

```
whos di o
      Name      Size      Bytes  Class

di o          1x1          1308  digitalio object
```

```
Grand total is 40 elements using 1308 bytes
```

Once the digital I/O object is created, the properties listed below are automatically assigned values. These general purpose properties provide descriptive information about the object based on its class type and adaptor.

Table 7-1: Descriptive Digital I/O Properties

Property Name	Description
Name	Specify a descriptive name for the device object.
Type	Indicate the device object type.

You can display the values of these properties for `di o` with the `get` function.

```
get(di o, {'Name', 'Type'})
ans =
      'ni daq1- DI0'      'Digital IO'
```

Adding Lines to a Digital I/O Object

After creating the digital I/O object, you must add lines to it. As shown by the figure in “Adding Channels or Lines” on page 3-8, you can think of a device object as a container for lines. The collection of lines contained by the device object is referred to as a *line group*. A line group consists of a mapping between hardware line IDs and MATLAB indices (see below).

When adding lines to a digital I/O object, you must follow these rules:

- The lines must reside on the same hardware device. You cannot add lines from different devices, or from different subsystems on the same device.
- You can add a line only once to a given digital I/O object. However, a line can be added to as many different digital I/O objects as you desire.
- You can add lines that reside on different ports to a given digital I/O object.

You add lines to a digital I/O object with the `addline` function. `addline` requires the device object, at least one hardware line ID, and the direction (input or output) of each added line as input arguments. You can optionally specify port IDs, descriptive line names, and an output argument. For example, to add eight output lines from port 0 to the device object `dio` created in the preceding section

```
hwlines = addline(dio, 0: 7, 'out');
```

The output argument `hwlines` is a *line object* that reflects the line array contained by `dio`. You can display the class of `hwlines` with the `whos` command.

```
whos hwlines
      Name          Size          Bytes  Class
      hwlines       8x1           536  dioline object
```

```
Grand total is 13 elements using 536 bytes
```

You can use `hwlines` to easily access lines. For example, you can easily configure or return property values for one or more lines. As described in XXX, you can also access lines with the `Line` property.

Once you add lines to a digital I/O object, the properties listed below are automatically assigned values. These properties provide descriptive information about the lines based on their class type and ID.

Table 7-2: Descriptive Digital I/O Line Properties

Property Name	Description
HwLi ne	Specify the hardware line ID.
I ndex	Indicate the MATLAB index of a hardware line.
Parent	Indicate the parent (device object) of a line.
Type	Indicate a line.

You can display the values of these properties for chans with the `get` function.

```
get(hwLines, {'HwLi ne', 'I ndex', 'Parent', 'Type'})
ans =
    [0]    [1]    [1x1 di gi tal i o]    'Li ne'
    [1]    [2]    [1x1 di gi tal i o]    'Li ne'
    [2]    [3]    [1x1 di gi tal i o]    'Li ne'
    [3]    [4]    [1x1 di gi tal i o]    'Li ne'
    [4]    [5]    [1x1 di gi tal i o]    'Li ne'
    [5]    [6]    [1x1 di gi tal i o]    'Li ne'
    [6]    [7]    [1x1 di gi tal i o]    'Li ne'
    [7]    [8]    [1x1 di gi tal i o]    'Li ne'
```

Line and Port Characteristics

As described in the preceding section, when you add lines to a digital I/O object, they must be configured for either input or output. You read values from an input line and write values to an output line. Whether a given hardware line is addressable for input or output depends on the port it resides on. You can classify digital I/O ports into these two groups based on your ability to address lines individually:

- **Port-configurable devices** – You cannot address the lines associated with a port-configurable device individually. This means that you must configure

all the lines for either input or output. If you attempt to mix the two configurations, an error is returned.

You can add any number of available port-configurable lines to a digital I/O object. However, the engine will address all lines behind the scenes. For example, if one line is added to a DIO object, then you automatically get all lines. This means that if a digital I/O object contains lines from a port-configurable device, and you write a value to one or more of those lines, then all the lines are written to even if they are not contained by the device object.

- **Line-configurable devices** – You can address the lines associated with a line-configurable device individually. This means that you can configure individual lines for either input or output. Additionally, you can read and write values on a line-by-line basis. Note that for National Instruments E Series hardware, port 0 is always line-configurable, while all other ports are port-configurable.

You can return the line and port characteristics with the `daqhwinfo` function. For example, National Instruments AT-MIO-16DE-10 board has four ports and eight digital I/O lines per port. To return the digital I/O characteristics for this board

```
hwinfo = daqhwinfo(dio);
```

To display the line and port characteristics for the first port

```
hwinfo.Port(1)
ans =
    ID: 0
    LineIDs: [0 1 2 3 4 5 6 7]
    Direction: 'in/out'
    Config: 'line'
```

To display the line and port characteristics for the second port

```
hwinfo.Port(2)
ans =
    ID: 2
    LineIDs: [0 1 2 3 4 5 6 7]
    Direction: 'in/out'
    Config: 'port'
```

To display the line and port characteristics for the third port

```
hwinfo. Port (3)
ans =
    ID: 3
    LineIDs: [0 1 2 3 4 5 6 7]
    Direction: 'in/out'
    Config: 'port'
```

To display the line and port characteristics for the fourth port

```
hwinfo. Port (4)
ans =
    ID: 4
    LineIDs: [0 1 2 3 4 5 6 7]
    Direction: 'in/out'
    Config: 'port'
```

This information tells you that you can configure all 32 lines for either input or output, and that the first port is line-configurable while the remaining ports are port-configurable.

Referencing Individual Hardware Lines

As described in the preceding section, you can access lines with the `Line` property or with a line object. To reference individual lines, you must specify either MATLAB indices or descriptive line names.

MATLAB Indices

Every hardware line contained by a digital I/O object has an associated MATLAB index that is used to reference the line. When adding lines with the `addline` function, index assignments are made automatically. The line indices start at 1 and increase monotonically up to the number of line group members. The first line indexed in the line group represents the least significant bit (LSB), and the highest indexed line represents the most significant bit (MSB). Unlike adding channels with the `addchannel` function, you cannot manually assign line indices with `addline`.

For example, the digital I/O object `di o` created in the preceding section had the MATLAB indices 1 through 8 automatically assigned to the hardware lines 0 through 7, respectively. To manually swap the hardware line order so that line ID 1 is the LSB, you supply the appropriate index to `hwlines` and use the `HwLine` property.

```
hwlines(1).HwLine = 1;  
hwlines(2).HwLine = 0;
```

Alternatively, you can use the `Line` property.

```
di o.Line(1).HwLine = 1;  
di o.Line(2).HwLine = 0;
```

Descriptive Line Names

Choosing a unique, descriptive name can be a useful way to identify and reference lines – particularly for large line groups. You can associate descriptive names with hardware lines using the `addline` function. For example, suppose you want to add 8 lines to `di o`, and you want to associate the name `Tri gLine` with the first line added.

```
addline(di o, 0, 'out', 'Tri gLine');  
addline(di o, 1: 7, 'out');
```

Alternatively, you can use the `LineName` property.

```
addLine(dio, 0: 7, 'out');
dio.Line(1).LineName = 'TriGLine';
```

You can now use the line name to reference the line.

```
dio.TriGLine.Direction = 'in';
```

Example: Adding Lines for National Instruments Hardware

This example illustrates various ways you can add lines to a digital I/O object associated with a National Instruments AT-MIO-16DE-10 board. This board is a multiport device whose characteristics are described in “Line and Port Characteristics” on page 7-5.

To add eight input lines to `dio` from port 0

```
addLine(dio, 0: 7, 'in');
```

To add four input lines and four output lines to `dio` from port 0

```
addLine(dio, 0: 7, {'in', 'in', 'in', 'in', 'out', 'out', 'out', 'out'});
```

Suppose you want to add the first line from port 0 and port 2 configured for input, and the second line from the same ports configured for output. There are three ways to do this. The first way requires only one call to `addLine` since it uses the hardware line IDs, and not the port IDs.

```
addLine(dio, [0 1 8 9], {'in', 'out', 'in', 'out'});
```

The second way requires two calls to `addLine`, and specifies one line ID and multiple port IDs for each call.

```
addLine(dio, 0, [0 2], 'in');
addLine(dio, 1, [0 2], 'out');
```

The third way requires two calls to `addLine`, and specifies multiple line IDs and one port ID for each call.

```
addLine(dio, 0: 1, 0, {'in', 'out'});
addLine(dio, 0: 1, 2, {'in', 'out'});
```

Writing and Reading Digital I/O Line Values

After you add lines to a digital I/O object, you can:

- Write values to lines
- Read values from lines

Note Unlike analog input and analog output objects, you do not control the behavior of digital I/O objects by configuring properties. This is because buffered digital I/O is not supported, and data is not stored in the engine. Instead, you either write values directly to, or read value directly from one or more hardware lines.

Writing Digital Values

You write values to one or more digital I/O lines with the `putvalue` function. `putvalue` requires the digital I/O object and the values to be written as input arguments. You can specify the values to be written as a decimal value or as a *binary vector* (`binvec`). A binary vector is a logical array that is constructed with the least significant bit (LSB) in the first column and the most significant bit (MSB) in the last column. For example, the decimal value 23 is written in `binvec` notation as `[1 1 1 0 1]` = $2^0 + 2^1 + 2^2 + 2^4$. You may find that `binvecs` are easier to work with than decimal values since there is a clear association between a given line and the value (1 or 0) that is written to it. You can convert decimal values to `binvec` values with the `dec2binvec` function

For example, suppose you create the digital I/O object `di o` and add eight output lines to it from port 0.

```
di o = digitalio('ni daq', 1);  
addline(di o, 0: 7, 'out');
```

To write a value of 23 to the eight lines contained by `di o`, you can write to the device object.

```
data = 23;  
putvalue(di o, data)
```

Alternatively, you can write to individual lines through the `Line` property.

```
putvalue(dio.Line(1:8), data)
```

To write a binary vector of values using the device object and the `Line` property

```
bvdata = dec2binvec(data, 8);  
putvalue(dio, bvdata)  
putvalue(dio.Line(1:8), bvdata)
```

The second input argument supplied to `dec2binvec` specifies the number of bits used to represent the decimal value. Since the preceding commands write to all eight lines contained by `dio`, an eight element binary vector is required. If you do not specify the number of bits, then the minimum number of bits needed to represent the decimal value is used.

Alternatively, you can create the binary vector without using `dec2binvec`.

```
bvdata = logical([1 1 1 0 1 0 0 0]);  
putvalue(dio, bvdata)
```

Rules for Writing Digital Values

Writing values to digital I/O lines follows these rules:

- If the digital I/O object contains lines from a port-configurable device, then the data acquisition engine writes to all port-configurable lines even if they are not contained by the device object.
- When writing decimal values:
 - If the value is too large to be represented by the lines contained by the device object, then an error is returned.
 - You can write to a maximum of 32 lines. To write to more than 32 lines, you must use a `binvec` value.
- When writing `binvec` values:
 - You can write to any number of lines.
 - There must be an element in the binary vector for each line you write to.
- You can always read from a line configured for output. Reading values is discussed on “Reading Digital Values” on page 7-12.
- An error is returned if you write a negative value, or if you write to a line configured for input.

Reading Digital Values

You can read values from one or more digital I/O lines with the `getvalue` function. `getvalue` requires the digital I/O object as an input argument. You can optionally specify an output argument, which represents the returned values as a binary vector. Binary vectors are described in “Writing Digital Values” on page 7-10.

For example, suppose you create the digital I/O object `di o` and add eight input lines to it from port 0.

```
di o = digitalio('ni daq', 1);
addline(di o, 0: 7, 'in');
```

To read the current value of all the lines contained by `di o` and return the result to `out`

```
portval = getvalue(di o)
portval =
    1     1     1     0     1     0     0     0
```

To read the current values of the first four lines contained by `di o` and return the result to `out`

```
lineval = getvalue(di o, Line(1: 5))
lineval =
    1     1     1     0     1
```

You can convert a `binvec` to a decimal value with the `binvec2dec` function. For example, to convert the binary vector `lineval` to a decimal value

```
out = binvec2dec(lineval)
out =
    23
```


Rules for Reading Digital Values

Reading values from digital I/O lines follows these rules:

- If the digital I/O object contains lines from a port-configurable device, then all lines are read even if they are not contained by the device object. However, only values from the lines contained by the digital I/O object are returned.
- You can always read from a line configured for output.
- For National Instruments hardware, lines configured for input return a value of 1 by default.
- `getValue` always returns a binary vector (`binvec`). To convert the `binvec` to a decimal value, use the `binvec2dec` function.

Example: Writing and Reading Digital Values

This example illustrates how to read and write digital values using a line-configurable subsystem. With line-configurable subsystems, you can transfer values on a line-by-line basis.

You can run this example by typing `daqdoc7_1` at the MATLAB command line.

1. Create a device object – Create the digital I/O object `di o` for a National Instruments board. The installed adaptors and hardware IDs are found with `daqhwinfo`.

```
di o = digitalio('ni daq', 1);
```

2. Add lines – Add eight output lines from port 0 (line-configurable).

```
addline(di o, 0: 7, 'out');
```

3. Read and write values – Write a value of 13 to the first four lines as a decimal number and as a binary vector, and read back the values.

```
data = 13;
putvalue(di o.Line(1: 4), data)
val 1 = getvalue(di o);
bvdata = dec2binvec(data);
putvalue(di o.Line(1: 4), bvdata)
val 2 = getvalue(di o);
```

Write a value of 3 to the last four lines as a decimal number and as a binary vector, and read back the values.

```
data = 3;
putvalue(dio, Line(5:8), data)
val3 = getvalue(dio, Line(5:8));
bvdata = dec2binvec(data, 4);
putvalue(dio, Line(5:8), bvdata)
val4 = getvalue(dio, Line(5:8));
```

Read values from the last four lines but switch the most significant bit (MSB) and the least significant bit (LSB).

```
val5 = getvalue(dio, Line(8:-1:5));
```

4. Clean up – When you no longer need `dio`, you should remove it from memory and from the MATLAB workspace.

```
delete(dio)
clear dio
```

Generating Timer Events

The fact that analog input and analog output objects make use of data stored in the engine and clocked I/O leads to the concept of a “running” device object and the generation of events.

However, since the Data Acquisition Toolbox does not support buffered digital I/O operations, digital I/O objects do not store data in the engine. Additionally, reading and writing line values are not clocked at a specific rate in the way that data is sampled by an analog input or analog output subsystem. Instead, values are either written directly to digital lines with `putval ue`, or read directly from digital lines with `getval ue`.

Therefore, the concept of a running digital I/O object does not make sense in the same way that it does for analog input and analog output objects. However, you can “run” a digital I/O object to perform one task: generate timer events. You can use timer events to update and display the state of the digital I/O object. For an example of this, refer to the `di opanel demo`.

Timer Events

The only event supported by digital I/O objects is a timer event. Timer events occur after a specified period of time has passed. Properties associated with generating timer events are given below.

Table 7-3: Digital I/O Timer Event Properties

Property Name	Description
<code>Runni ng</code>	Indicates if the device object is running.
<code>Ti merActi on</code>	Specifies the M-file action function to execute whenever a predefined period of time passes.
<code>Ti merPeri od</code>	Specifies the period of time between timer events.

A timer event is generated whenever the time specified by `Ti merPeri od` passes. This event executes the action function specified for `Ti merActi on`. Time is measured relative to when the device object starts running (`Runni ng is On`). Starting a digital I/O object is discussed in the next section.

Some timer events may not be processed if your system is significantly slowed or if the `TimerPeriod` value is too small. For example, a common application for timer events is to display data. However, since displaying data can be a CPU-intensive task, some of these events may be dropped. For digital I/O objects, timer events are typically used to display the state of the device object. Refer to the `diopanel_demo` for an example that generates timer events to update the state of a digital I/O object.

To see how to construct an action function, refer to “Creating and Executing Action Functions” on page 5-52 or the example below.

Starting and Stopping a Digital I/O Object

You use the `start` function to start a digital I/O object. For example, to start the digital I/O object `di o`

```
start(di o)
```

After `start` is issued, the `Running` property is automatically set to `On`, and timer events can be generated. If you attempt to start a digital I/O object that does not contain any lines or that is already running, then an error is returned.

A digital I/O object will stop executing under these conditions:

- The `stop` function is issued.
- An error occurred in the system.

When the device object stops, `Running` is automatically set to `Off`.

Example: Generating Timer Events

This example illustrates how to generate timer events for a digital I/O object. The action function `daqacti on` displays the event type and device object name. Note that you must issue a stop command to stop the execution of the digital I/O object.

You can run this example by typing `daqdoc7_2` at the MATLAB command line.

1. Create a device object – Create the digital I/O object `di o` for a National Instruments board. The installed adaptors and hardware IDs are found with `daqhwinfo`.

```
di o = digital i o('ni daq', 1);
```

2. Add lines – Add eight input lines from port 0 (line-configurable).

```
addl i ne(di o, 0: 7, 'i n');
```

3. Configure property values – Configure the timer event to call `daqacti on` every five seconds.

```
set(di o, 'Ti merActi on', 'daqacti on')
set(di o, 'Ti merPeri od', 5. 0)
```

Start the digital I/O object. You must issue a stop command when you no longer want to generate timer events.

```
start(di o)
```

4. Clean up – When you no longer need `di o`, you should remove it from memory and from the MATLAB workspace.

```
del ete(di o)
cl ear di o
```

Evaluating the Digital I/O Object Status

At any time after a digital I/O object is created, you can evaluate its status by:

- Returning the value of the `Running` property
- Invoking the display summary

The Display Summary

You can invoke the display summary by typing the digital I/O object at the command line. The displayed information is designed so you can quickly evaluate the status of your data acquisition session. The display is divided into two main sections: general summary information and line summary information.

General Summary Information

The general display summary includes the device object type and the hardware device name, followed by the port parameters. The port parameters include the port ID, and whether the associated lines are configurable for reading or writing.

Line Summary Information

The line display summary includes property values associated with:

- The hardware line mapping
- The line name
- The port ID
- The line direction

The display summary for the example given in “Example: Writing and Reading Digital Values” on page 7-13 is shown below.

General display summary	Display Summary of DigitalIO (DIO) Object Using 'PCI-6024E'.				
	Port Parameters: Port 0 is line configurable for reading and writing				
	Engine status: Engine not required.				
Line display summary	DIO object contains line(s):				
	Index:	LineName:	HwLine:	Port:	Direction:
	1	' '	0	0	' In'
	2	' '	1	0	' In'
	3	' '	2	0	' In'
	4	' '	3	0	' In'
	5	' '	4	0	' In'
	6	' '	5	0	' In'
	7	' '	6	0	' In'
	8	' '	7	0	' In'

You can use the Line property to display only the line summary information.

DIO.Line

Saving and Loading the Session

Overview	8-2
Saving and Loading Device Objects	8-3
Saving Device Objects to an M-File	8-3
Saving Device Objects to a MAT-File	8-5
Logging Information to Disk	8-6
Specifying a Filename	8-7
Retrieving Logged Information	8-8
Example: Logging and Retrieving Information	8-10

Overview

During a data acquisition session, you may want to save your work to one or more disk files in order to document or debug your application, post-process data, and so on. The components of a data acquisition session that you can save to disk are divided into two separate parts:

- Device objects and their associated property values
- Information associated with acquiring data with an analog input object including:
 - Acquired data
 - Event information
 - Device object and channel information
 - Hardware information

You can save specific components of your session to an M-file using the `obj2mfile` function, to a MAT-file using the `save` command, or to a binary ".daq" file using the `LoggingMode` property.

Saving and Loading Device Objects

You can save a device object to disk using two possible formats:

- As an M-file using the `obj2mfile` function
- As a MAT-file using the `save` command

For analog input objects, you can also save acquired data, hardware information, and so on to a log file. Refer to “Logging Information to Disk” on page 8-6 for more information.

Saving Device Objects to an M-File

You can save a device object to an M-file using the `obj2mfile` function. `obj2mfile` provides you with these options:

- Save all property values, or save only those property values that differ from their default values.

Read-only property values are not saved. Therefore, read-only properties use their default values when you load the device object into the MATLAB workspace. To determine if a property is read-only, use the `propinfo` function or examine the property reference pages.

- Save property values using the `set` syntax, the dot notation, or named referencing (if defined).

If the `UserData` property is not empty, or if action properties are set to a cell array of values, then the data stored in these properties is written to a MAT-file when the device object is saved. The MAT-file has the same name as the M-file containing the device object code.

For example, suppose you create the analog input object `ai` for a sound card, add two channels to it, and configure several property values.

```
ai = analoginput('winsound');
addchannel(ai, 1:2, {'Temp1'; 'Temp2'});
time = now;
set(ai, 'SampleRate', 11025, 'TriggerRepeat', 4)
set(ai, 'TriggerAction', {'myaction', time})
start(ai)
```

The following command saves `ai` and the modified property values to the M-file `myai.m`. Since the `TriggerAction` property is set to a cell array of values, its value is automatically written to the MAT-file `myai.mat`.

```
obj2mfile(ai, 'myai.m');
```

```
Created: d:\v6\myfiles\myai.m
```

```
Created: d:\v6\myfiles\myai.mat
```

Use the `type` command to display `myai.m` at the command line.

Loading the Device Object

To load a device object that was saved as an M-file into the MATLAB workspace, type the name of the M-file at the command line. For example, to load `ai` from the M-file `myai.m`

```
ai = myai
```

Note that the read-only properties such as `SamplesAcquired` and `SamplesAvailable` are restored to their default values.

```
get(ai, {'SamplesAcquired', 'SamplesAvailable'})
ans =
     [0]     [0]
```

When loading `ai` into the workspace, the MAT-file `myai.mat` is automatically loaded and the `TriggerAction` property value is restored.

```
ai.TriggerAction
ans =
    'myaction'    [7.3071e+005]
```

Saving Device Objects to a MAT-File

You can save a device object to a MAT-file just as you would any workspace variable – using the `save` command. For example, to save the analog input object `ai` and the variable `time` defined in the preceding section to the MAT-file `myai1.mat`

```
save myai1 ai time
```

Read-only property values are not saved. Therefore, read-only properties use their default values when you load the device object into the MATLAB workspace. To determine if a property is read-only, use the `propinfo` function or examine the property reference pages.

Loading the Device Object

To load a device object that was saved to a MAT-file into the MATLAB workspace, use the `load` command. For example, to load `ai` and `time` from MAT-file `myai1.mat`

```
load myai1
```


Logging Information to Disk

While an analog input object is running, you can log this information to a disk file:

- Acquired data
- Event information
- Device object and channel information
- Hardware information

Logging information to disk provides a permanent record of your data acquisition session, and is an easy way to debug your application.

As shown below, you can think of the logged information as a stream of data and events.



The properties associated with logging information to a disk file are given below.

Table 8-1: Analog Input Logging Properties

Property Name	Description
LogFileName	Specify the name of the disk file to which information is logged.
Logging	Indicate if data is being logged.

Table 8-1: Analog Input Logging Properties (Continued)

Property Name	Description
Loggi ngMode	Specify the destination for acquired data.
LogToDi skMode	Specify whether data, device object information, and hardware information is saved to one disk file or to multiple disk files.

You can initiate logging by setting Loggi ngMode to Di sk or Di sk&Memory. A new log file is created each time you issue the start function, and each different analog input object must log information to a separate log file. Writing to disk is performed as soon as possible after the current data block is filled.

You can choose whether a log file is overwritten or if multiple log files are created with the LogToDi skMode property. If LogToDi skMode is Overwri te, the log file is overwritten. If LogToDi skMode is Index, new log files are created, each with an indexed name based on the value of LogFi l eName.

Specifying a Filename

You specify the name of the log file with the LogFi l eName property. You can specify any value for LogFi l eName, including a directory path, provided the filename is supported by your operating system. Additionally, if LogToDi skMode is Index, then the log filename also follows these rules:

- Indexed log filenames are identified by a number. This number precedes the filename extension and is incremented by one for successive log files.
- If no number is specified as part of the initial log filename, then the first log file does not have a number associated with it. For example, if LogFi l eName is myfi l e. daq, then myfi l e. daq is the name of the first log file, myfi l e01. daq is the name of the second log file, and so on.
- LogFi l eName is updated after the log file is written (after the stop event occurs).
- If the specified log filename already exists, then the existing file is overwritten.

Retrieving Logged Information

You retrieve logged information with the `daqread` function. You can retrieve any part of the information stored in a log file with one call to `daqread`. However, you will probably use `daqread` in one of these two ways:

- Retrieving data and time information
- Retrieving event, device object, channel, and hardware information

Retrieving Data and Time Information

You can characterize logged data by the sample number or the time the sample was acquired. To retrieve data and time information, you use the syntax shown below

```
[data, time, abstime] = daqread('file', 'P1', V1, 'P2', V2, ...);
```

where:

- `data` is the retrieved data. Data is returned as an `m`-by-`n` matrix where `m` is the number of samples and `n` is the number of channels.
- `time` (optional) is the relative time associated with the retrieved data. Time is returned as an `m`-by-1 matrix where `m` is the number of samples.
- `abstime` (optional) is the absolute time of the first trigger. Absolute time is returned as a clock vector.
- `file` is the name of the log file.
- `'P1', V2, 'P2', V2, ...` (optional) are the property name/property value pairs, which allow you to select the amount of data to retrieve, among other things (see below).

`daqread` returns data and time information in the same format as `getdata`. If data from multiple triggers is retrieved, each trigger is separated by a NaN.

You select the amount of data returned and the format of that data with the properties given below.

Table 8-2: daqread Properties

Property Name	Description
Sampl es	Specify the sample range.
Ti me	Specify the relative time range.
Tri ggers	Specify the trigger range.
Channel s	Specify the channel range. Channel names can be specified as a cell array.
DataFormat	Specify the data format as doubl es or nati ve.
Ti meFormat	Specify the time format as vector or matrix.

The Sampl es, Ti me, and Tri ggers properties are mutually exclusive. If none of these three properties is specified, then all the data is returned.

Retrieving Event, Device Object, Channel, and Hardware Information

You can retrieve event, device object, channel, and hardware information, along with data and time information, using the syntax shown below.

```
[data, time, abstime, events, daqi nfo] =
daqread(' file' , ' P1' , V1, ' P2' , V2, ... );
```

events is a structure containing event information associated with the logged data. The events retrieved depend on the value of the Sampl es, Ti me, or Tri ggers property. At a minimum, the trigger event associated with the selected data is returned. The entire event log is returned to events only if Sampl es, Ti me, or Tri ggers is not specified.

daqi nfo is a structure that stores device object, channel, and hardware information in two fields: Obj I nfo and HwI nfo. Obj I nfo is a structure containing property values for the device object and any channels it contains. The property values are returned in the same format as returned by get. HwI nfo is a structure containing hardware information. The hardware information is identical to the information returned by daqhwi nfo(obj) .

Alternatively, you can return only object, channel, and hardware information with the command

```
daqinfo = daqread('file','info');
```

Note When you retrieve object information, the entire event log is returned to `daqinfo.ObjInfo.EventLog` regardless of the number of samples retrieved.

Example: Logging and Retrieving Information

This example illustrates how to log information to a disk file and then retrieve the logged information to MATLAB using various calls to `daqread`.

A sound card is configured for stereo acquisition, data is logged to memory and to a disk file, four triggers are issued, and 2 seconds of data are collected for each trigger at a sampling rate of 8 kHz. You can run this example by typing `daqdoc8_1` at the MATLAB command line.

1. Create a device object – Create the analog input object `ai` for a sound card. The installed adaptors and hardware IDs are found with `daqhwinfo`.

```
ai = analoginput('winsound');  
%ai = analoginput('ni daq', 1);  
%ai = analoginput('cbi', 1);
```

2. Add channels – Add two hardware channels to `ai`.

```
ch = addchannel(ai, 1:2);  
%ch = addchannel(ai, 0:1); % For NI and CBI
```

3. Configure property values – Define a 2 second acquisition for each trigger, set the trigger to repeat three times, and log information to the file `file00.daq`.

```
duration = 2; % Two seconds of data for each trigger  
set(ai, 'SampleRate', 8000)  
ActualRate = get(ai, 'SampleRate');  
set(ai, 'SamplesPerTrigger', duration*ActualRate)  
set(ai, 'TriggerRepeat', 3)  
set(ai, 'LogFileName', 'file00.daq')  
set(ai, 'LoggingMode', 'Disk&Memory')
```

4. Acquire data – Start `ai`, wait for `ai` to stop running, and extract all the data stored in the log file as sample-time pairs.

```
start(ai)
while strcmp(ai.Running, 'On')
end
[data, time] = daqread('file00.daq');
```

Plot the data and label the figure axes.

```
subplot(211), plot(data)
title('Logging and Retrieving Data')
xlabel('Samples'), ylabel('Signal (Volts)')
subplot(212), plot(time, data)
xlabel('Time (seconds)'), ylabel('Signal (Volts)')
```

5. Clean up – When you no longer need `ai`, you should remove it from memory and from the MATLAB workspace.

```
delete(ai)
clear ai
```

Retrieving Data Based on Samples

You can retrieve data based on samples using the `Samples` property. To retrieve samples 1000 to 2000 for both sound card channels

```
[data, time] = daqread('file00.daq', 'Samples', [1000 2000]);
```

Plot the data and label the figure axes.

```
subplot(211), plot(data);
xlabel('Samples'), ylabel('Signal (Volts)')
subplot(212), plot(time, data);
xlabel('Time (seconds)'), ylabel('Signal (Volts)')
```

Retrieving Data Based on Channels

You can retrieve data based on channels using the `Channels` property. To retrieve samples 1000 to 2000 for the second sound card channel

```
[data, time] = daqread('file00.daq', 'Samples', [1000 2000],
'Channels', 2);
```

Plot the data and label the figure axes.

```
subplot(211), plot(data);  
xlabel('Samples'), ylabel('Signal (Volts)')  
subplot(212), plot(time, data);  
xlabel('Time (seconds)'), ylabel('Signal (Volts)')
```

Alternatively, you can retrieve data for the second sound card channel by specifying the channel name.

```
[data, time] = daqread('file00.daq', 'Samples', [1000 2000],  
    'Channels', {'Right'});
```

Retrieving Data Based on Triggers

You can retrieve data based on triggers using the `Triggers` property. To retrieve all the data associated with the second and third triggers for both sound card channels

```
[data, time] = daqread('file00.daq', 'Triggers', [2 3]);
```

Plot the data and label the figure axes.

```
subplot(211), plot(data);  
xlabel('Samples'), ylabel('Signal (Volts)')  
subplot(212), plot(time, data);  
xlabel('Time (seconds)'), ylabel('Signal (Volts)')
```

Retrieving Data Based on Time

You can retrieve data based on time using the `Time` property. `Time` must be specified in seconds and `Time=0` corresponds to the first logged sample. To retrieve the first 25% of the data acquired for the first trigger

```
[data, time] = daqread('file00.daq', 'Time', [0 0.5]);
```

Plot the data and label the figure axes.

```
subplot(211), plot(data);  
xlabel('Samples'), ylabel('Signal (Volts)')  
subplot(212), plot(time, data);  
xlabel('Time (seconds)'), ylabel('Signal (Volts)')
```

Retrieving Event, Object, Channel, and Hardware Information

You can retrieve event, object, channel, and hardware information by specifying the appropriate arguments to `daqread`. For example, to retrieve all event information, you must return all the logged data.

```
[data, time, abstime, events, info] = daqread('file00.daq');
{events.Type}
ans =
'Start' 'Trigger' 'Trigger' 'Trigger' 'Trigger' 'Stop'
```

If you retrieve part of the data, then only the events associated with the requested data are returned.

```
[data, time, abstime, events, info] = daqread('file00.daq',
'Triigger', [1 3]);
{events.Type}
ans =
'Triigger' 'Triigger' 'Triigger'
```

You can retrieve the entire event log as well as object and hardware information by including `info` as an input argument to `daqread`.

```
daqinfo = daqread('file00.daq', 'info')
daqinfo =
    ObjInfo: [1x1 struct]
    HwInfo: [1x1 struct]
```

To return the event log information

```
eventinfo = daqinfo.ObjInfo.EventLog
eventinfo =
6x1 struct array with fields:
    Type
    Data
```


Function Reference

Overview	9-2
Getting Command Line Function Help	9-2
Functions Grouped by Category	9-4
Functions Listed Alphabetically	9-8

Overview

This section provides descriptions of all toolbox functions that you can use directly. In “Functions Grouped by Category” on page 9-4, the functions are summarized according to usage. In “Functions Listed Alphabetically” on page 9-8, the functions are described in detail.

Getting Command Line Function Help

To get command line function help, you should use the `daqhelp` function. For example, to get help for the `addchannel` function, type

```
daqhelp addchannel
```

Alternatively, you can use the `help` command.

```
help addchannel
```

However, the Data Acquisition Toolbox provides “overloaded” versions of several MATLAB functions. That is, it provides toolbox-specific implementations of these functions using the same function name.

To get command line help for an overloaded toolbox function using the `help` command, you must supply one of two possible class directories to `help`.

```
help daqdevice/function_name  
help daqchild/function_name
```

Note that the same help information is returned regardless of the class directory specified.

For example, the Data Acquisition Toolbox provides an overloaded version of the `delete` function. To obtain help for the MATLAB version of this function, type

```
help delete
```


You can determine if a function is overloaded by examining the last section of the help. For `delete`, the help contains the following overloaded versions (not all are shown).

```
Overloaded methods
help char/delete.m
help scribehandle/delete.m
help scribehobj/delete.m
.
.
.
help daqdevice/delete.m
help daqchild/delete.m
```

So, to obtain help on the toolbox version of this function, type

```
help daqdevice/delete
```

For more information on overloaded functions and class directories, refer to “MATLAB Classes and Objects” in the Help browser.

Functions Grouped by Category

This section contains brief descriptions of all toolbox functions. The functions and the device objects they are associated with are categorized according to usage. The supported device objects include analog input (AI), analog output (AO), and digital I/O (DIO).

A number of other MATLAB M-file helper functions are provided with this toolbox to support the functions listed below. These helper functions are not documented since they are not intended for direct use.

Creating Device Objects		AI	AO	DIO
analoginput	Create an analog input object.	✓		
analogoutput	Create an analog output object.		✓	
digitalio	Create a digital I/O object.			✓

Adding Channels and Lines		AI	AO	DIO
addchannel	Add hardware channels to an analog input or analog output object.	✓	✓	
addline	Add hardware lines to a digital I/O object.			✓
addmuxchannel	Add hardware channels when using a multiplexer board.	✓		

Getting and Setting Properties		AI	AO	DIO
get	Return device object properties.	✓	✓	✓
set	Configure or display device object properties.	✓	✓	✓
setverify	Configure and return the specified property.	✓	✓	✓

Executing the Object		AI	AO	DIO
start	Start a device object.	✓	✓	✓
stop	Stop a device object.	✓	✓	✓
trigger	Manually execute a trigger.	✓	✓	
wait til stop	Wait for the device object to stop running.	✓	✓	

Working with Data		AI	AO	DIO
flushdata	Remove data from the data acquisition engine.	✓		
getdata	Extract data, time, and event information from the data acquisition engine.	✓		
getsampl e	Immediately acquire one sample.	✓		
getval ue	Read values from lines.			✓
peekdata	Preview most recent acquired data.	✓		
putdata	Queue data in the engine for eventual output.		✓	
putsampl e	Immediately output one sample.		✓	
putval ue	Write values to lines.			✓

Getting Information and Help		AI	AO	DIO
daqhel p	Display help for device objects, constructors, adaptors, functions, and properties.	✓	✓	✓
daqhwi nfo	Display data acquisition hardware information.	✓	✓	✓

Getting Information and Help (Continued)		AI	AO	DIO
daqpropedit	Invoke the property editor graphical user interface.	✓	✓	✓
propinfo	Return property characteristics for device objects, channels, or lines.	✓	✓	✓

General Purpose		AI	AO	DIO
binvec2dec	Convert binary vector to decimal value.			✓
clear	Remove device objects from the MATLAB workspace.	✓	✓	✓
daqaction	An action function that displays event information for the specified event.	✓	✓	✓
daqfind	Return device objects, channels, or lines from the data acquisition engine to the MATLAB workspace.	✓	✓	✓
daqmem	Allocate or display memory resources.	✓	✓	
daqread	Read a Data Acquisition Toolbox (.daq) file.	✓		
daqschool	Interface for displaying toolbox tutorials.	✓	✓	✓
daqregister	Register or unregister a hardware driver adaptor.	✓	✓	✓
daqreset	Remove device objects and data acquisition DLLs from memory.	✓	✓	✓
dec2binvec	Convert decimal value to binary vector.			✓
delete	Remove device objects, channels, or lines from the data acquisition engine.	✓	✓	✓
disp	Display summary information for device objects, channels, or lines.	✓	✓	✓
ischannel	Check for channels.	✓	✓	✓
isdioline	Check for lines.	✓	✓	✓

General Purpose (Continued)		AI	AO	DIO
<code>isvalid</code>	Determine if device objects, channels, or lines are valid.	✓	✓	✓
<code>length</code>	Return the length of a device object, channel group, or line group.	✓	✓	✓
<code>load</code>	Load device objects, channels, or lines into the MATLAB workspace.	✓	✓	✓
<code>makenames</code>	Generate a list of descriptive channel or line names.	✓	✓	✓
<code>muxchannel</code>	Return multiplexed scanned channel index.	✓		
<code>obj2mfile</code>	Convert device objects, channels, or lines to MATLAB code.	✓	✓	✓
<code>save</code>	Save device objects to a MAT-file.	✓	✓	✓
<code>showdaqevents</code>	Display event log information.	✓	✓	
<code>size</code>	Return the size of a device object, channel group, or line group.	✓	✓	✓

Functions Listed Alphabetically

This section contains detailed descriptions of all toolbox functions. Each function reference page contains some or all of this information:

- The function name
- The purpose of the function
- The function syntax

All valid input argument and output argument combinations are shown. In some cases, an ellipsis (. . .) is used for the input arguments. This means that all preceding input argument combinations are valid for the specified output argument(s).

- A description of each argument
- A description of the function
- Additional remarks about usage
- An example of usage
- Related functions or properties

Purpose Add hardware channels to an analog input or analog output object

Syntax

```
chans = addchannel (obj , hwch)
chans = addchannel (obj , hwch, i ndex)
chans = addchannel (obj , hwch, ' names' )
chans = addchannel (obj , hwch, i ndex, ' names' )
```

Arguments

obj	An analog input or analog output object.
hwch	Specifies the numeric IDs of the hardware channels added to the device object. Any MATLAB vector syntax can be used.
i ndex	The MATLAB indices to associate with the hardware channels. Any MATLAB vector syntax can be used provided the vector elements are monotonically increasing.
' names'	A descriptive channel name or cell array of descriptive channel names.
chans	A column vector of channels with the same length as hwch.

Description

chans = addchannel (obj , hwch) adds the hardware channels specified by hwch to the device object obj . The MATLAB indices associated with the added channels are assigned automatically. chans is a column vector of channels.

chans = addchannel (obj , hwch, i ndex) adds the hardware channels specified by hwch to the device object obj . i ndex specifies the MATLAB indices to associate with the added channels.

chans = addchannel (obj , hwch, ' names') adds the hardware channels specified by hwch to the device object obj . The MATLAB indices associated with the added channels are assigned automatically. names is a descriptive channel name or cell array of descriptive channel names.

chans = addchannel (obj , hwch, i ndex, ' names') adds the hardware channels specified by hwch to the device object obj . i ndex specifies the MATLAB indices to associate with the added channels. names is a descriptive channel name or cell array of descriptive channel names.

Remarks

Rules for Adding Channels

- The numeric values you supply for `hwch` depend on the hardware you access. For National Instruments and ComputerBoards hardware, channels are “zero-based” (begin at zero). For Agilent Technologies hardware and sound cards, channels are “one-based” (begin at one).
- Hardware channel IDs are stored in the `HwChannel` property and the associated MATLAB indices are stored in the `Index` property.
- You can add individual hardware channels to multiple device objects.
- For for sound cards and Agilent Technologies devices, you cannot add a hardware channel multiple times to the same device object.
- For Agilent Technologies devices, added channels must be in increasing order.
- You can configure sound cards in one of two ways: mono mode or stereo mode. For mono mode, `hwch` must be 1. For stereo mode, the first `hwch` value specified must be 1.

Note If you are using National Instruments AMUX-64T multiplexer boards, you must use the `addmuxchannel` function to add channels.

More About MATLAB Indices

Every hardware channel contained by a device object has an associated MATLAB index that is used to reference the channel. Index assignments are made either automatically by `addchannel` or explicitly with the `index` argument and follow these rules:

- If `index` is not specified and no hardware channels are contained by the device object, then the assigned indices automatically start at one and increase monotonically. If hardware channels have already been added to the device object, then the assigned indices automatically start at the next highest index value and increase monotonically.
- If `index` is specified but the indices are previously assigned, then the requested assignment takes precedence and the previous assignment is reindexed to the next available values. If the lengths of `hwch` and `index` are

not equal, then an error is returned and no channels are added to the device object.

- The resulting indices begin at one and increase monotonically up to the size of the channel group.
- If you are using scanning hardware, then the indices define the scan order.
- Sound cards cannot be reindexed.

Rules for Adding Channels to National Instruments 1200 Series Boards

When using National Instruments 1200 Series hardware, you need to modify the above rules in these ways:

- Channel IDs are given in reverse order with `addchannel`. For example, to add eight single-ended channels to the analog input object `ai`
`addchannel (ai , 7: - 1: 0) ;`
- The scan order is from the highest ID to the lowest ID (which must be 0).
- There cannot be any gaps in the channel group.
- When channels are configured in differential mode, the hardware IDs are 0, 2, 4, and 6.

More About Descriptive Channel Names

You can assign hardware channels descriptive names, which are stored in the Channel Name property. Choosing a unique descriptive name can be a useful way to identify and reference channels. For a single call to `addchannel`, you can:

- Specify one channel name that applies to all channels that are to be added
- Specify a different name for each channel to be added

If the number of names specified in a single `addchannel` call is more than one but not equal to the number of channels to be added, then an error is returned. If a channel is to be referenced by its name, then that name must not contain symbols. If you are naming a large number of channels, then the `makenames` function may be useful. If a channel is not assigned a descriptive name, then it must be referenced by index.

A sound card configured in mono mode is automatically assigned the name `Mono`, while a sound card configured in stereo mode is automatically assigned

addchannel

the names `Left` for the first channel and `Right` for the second channel. You can change these default channel names when the device object is created, or any time after the channel is added.

Example

National Instruments

Suppose you create the analog input object `AI 1` for a National Instruments board, and add the first four hardware channels (channels 0-3) to it.

```
AI 1 = analoginput('ni daq', 1);  
addchannel(AI 1, 0:3);
```

The channels are automatically assigned the indices 1-4. If you want to add the first four hardware channels to `AI 1` and assign descriptive names to the channels

```
addchannel(AI 1, 0:3, {'chan1', 'chan2', 'chan3', 'chan4'});
```

Note that you can use the `makenames` function to create a cell array of channel names. If you add channels 4, 5, and 7 to the existing channel group

```
addchannel(AI 1, [4 5 7]);
```

the new channels are automatically assigned the indices 5-7. Suppose instead you add channels 4, 5, and 7 to the channel group and explicitly assign them indices 1-3.

```
addchannel(AI 1, [4 5 7], 1:3);
```

The new channels are assigned the indices 1-3, and the previously defined channels are reindexed as indices 4-7. However, if you assigned channels 4, 5, and 7 to indices 6-8, an error is returned since there is a gap in the indices (index 5 has no associated hardware channel).

Sound Card

Suppose you create the analog input object AI 1 for a sound card. Most sound cards have only two channels that can be added to a device object. To configure the sound card to operate in mono mode, you must specify hwch as 1.

```
AI 1 = analog input ('winsound');  
addchannel (AI 1, 1);
```

The Channel Name property is automatically assigned the value Mono. You can now configure the sound card to operate in stereo mode by adding the second channel.

```
addchannel (AI 1, 2);
```

The Channel Name property is assigned the values Left and Right for the two hardware channels. Alternatively, you can configure the sound card to operate in stereo mode with one call to addchannel .

```
addchannel (AI 1, 1: 2);
```

See Also

Functions

delete, makenames

Properties

Channel Name, HwChannel, Index

addline

Purpose Add hardware lines to a digital I/O object

Syntax

```
lines = addline(obj, hwline, 'direction')  
lines = addline(obj, hwline, port, 'direction')  
lines = addline(obj, hwline, 'direction', 'names')  
lines = addline(obj, hwline, port, 'direction', 'names')
```

Arguments

obj	A digital I/O object.
hwline	The numeric IDs of the hardware lines added to the device object. Any MATLAB vector syntax can be used.
'direction'	The line directions can be In or Out, and can be specified as a single value or a cell array of values.
port	The numeric IDs of the digital I/O port.
'names'	A descriptive line name or cell array of descriptive line names.
lines	A row vector of lines with the same length as hwline.

Description

`lines = addline(obj, hwline, 'direction')` adds the hardware lines specified by `hwline` to the digital I/O object `obj`. `direction` configures the lines for either input or output. `lines` is a row vector of lines.

`lines = addline(obj, hwline, port, 'direction')` adds the hardware lines specified by `hwline` from the port specified by `port` to the digital I/O object `obj`.

`lines = addline(obj, hwline, 'direction', 'names')` adds the hardware lines specified by `hwline` to the digital I/O object `obj`. `names` is a descriptive line name or cell array of descriptive line names.

`lines = addline(obj, hwline, port, 'direction', 'names')` adds the hardware lines specified by `hwline` from the port specified by `port` to the digital I/O object `obj`. `direction` configures the lines for either input or output. `names` is a descriptive line name or cell array of descriptive line names.

Remarks**Rules for Adding Lines**

- The numeric values you supply for `hwLine` depend on the hardware you access. For National Instruments and ComputerBoards hardware, line IDs are “zero-based” (begin at zero).
- You can add a line only once to a given digital I/O object.
- Hardware line IDs are stored in the `HwLine` property and the associated MATLAB indices are stored in the `Index` property.
- For a single call to `addline`, you can add multiple lines from one port or the same line ID from multiple ports. You cannot add multiple lines from multiple ports.
- If a port ID is not explicitly referenced, lines are added first from port 0, then from port 1, and so on.
- You can specify the line directions as a single value or a cell array of values. If a single direction is specified, then all added lines have that direction. If supported by the hardware, you can configure individual lines by supplying a cell array of directions.

More About MATLAB Indices

Every hardware line contained by a device object has an associated MATLAB index that is used to reference the line. Index assignments are made automatically by `addline` and follow these rules:

- If no hardware lines are contained by the device object, then the assigned indices automatically start at one and increase monotonically. If hardware lines have already been added to the device object, then the assigned indices automatically start at the next highest index value and increase monotonically.
- The resulting indices begin at one and increase monotonically up to the size of the line group.
- The first indexed line represents the least significant bit (LSB) and the highest indexed line represents the most significant bit (MSB).

More About Descriptive Line Names

You can assign hardware lines descriptive names, which are stored in the `LineName` property. Choosing a unique descriptive name can be a useful way to identify and reference lines. For a single call to `addline`, you can:

- Specify one line name that applies to all lines that are to be added
- Specify a different name for each line to be added

If the number of names specified in a single `addline` call is more than one but differs from the number of lines to be added, then an error is returned. If a line is to be referenced by its name, then that name must not contain symbols. If you are naming a large number of lines, then the `makenames` function may be useful. If a line is not assigned a descriptive name, then it must be referenced by index.

Example

Create the digital I/O object `di o` and add the first four hardware lines (line IDs 0-3) from port 0.

```
di o = digitalio('ni daq', 1);  
addline(di o, 0: 3, 'in');
```

These lines are automatically assigned the indices 1-4. If you want to add the first four hardware lines to `di o` and assign descriptive names to the lines

```
addline(di o, 0: 3, 'in', {'line1', 'line2', 'line3', 'line4'});
```

Note that you can use the `makenames` function to create a cell array of line names. You can add the first four hardware lines (line IDs 0-3) from port 1 to the existing line group.

```
addline(di o, 0: 3, 1, 'out');
```

The new lines are automatically assigned the indices 5-8.

See Also

Functions

`delete`, `makenames`

Properties

`Hardware`, `Index`, `LineName`

Purpose	Add hardware channels when using a multiplexer board	
Syntax	<pre>addmuxchannel (obj) addmuxchannel (obj , chani ds) chans = addmuxchannel (. . .)</pre>	
Arguments	obj	An analog input object associated with a National Instruments board.
	chani ds	The hardware channel IDs.
	chans	The channels that are added to obj .
Description	<p><code>addmuxchannel (obj)</code> adds as many channels to <code>obj</code> as is physically possible based on the number of National Instruments AMUX-64T multiplexer (mux) boards specified by the <code>NumMuxBoards</code> property. For one mux board, 64 channels are added. For two mux boards, 128 channels are added. For four mux boards, 256 channels are added.</p> <p><code>addmuxchannel (obj , chani ds)</code> adds the channels specified by <code>chani ds</code> to <code>obj</code> . <code>chani ds</code> refers to the hardware channel IDs of the data acquisition board.</p> <p>The actual number of channels added to <code>obj</code> depends on the number of mux boards used. For example, suppose you are using a data acquisition board with 16 channels connected to one mux board. If <code>chani d</code> is 0, then <code>addmuxchannel</code> adds four channels. Refer to the <i>AMUX-64T User Manual</i> for more information about adding mux channels based on hardware channel IDs and the number of mux boards used.</p> <p><code>chans = addmuxchannel (. . .)</code> returns the channels added to <code>chans</code>.</p>	
Remarks	<p>Before using <code>addmuxchannel</code> , you must set the <code>NumMuxBoards</code> property to the appropriate value. You can use as many as four mux boards with one analog input object. <code>addmuxchannel</code> deletes all channels contained by <code>obj</code> before new channels are added.</p>	
See Also	<p>Functions</p> <p><code>muxchani dx</code></p>	

analoginput

Purpose	Create an analog input object	
Syntax	AI = analoginput('adaptor') AI = analoginput('adaptor', ID)	
Arguments	'adaptor'	Name of the hardware driver adaptor. The supported adaptors are ni daq, cbi, hpe1432, and wi nsound.
	ID	The hardware device identifier. ID is optional if the device object is associated with a sound card having an ID of 0.
	AI	The analog input object.
Description	AI = analoginput('adaptor') creates the analog input object AI for a sound card having an ID of 0 (adaptor must be wi nsound). This is the only case where ID is not required. AI = analoginput('adaptor', ID) creates the analog input object AI for the specified adaptor and for the hardware device with device identifier ID. ID can be specified as an integer or a string.	
Remarks	More About Creating Analog Input Objects <ul style="list-style-type: none">• When an analog input object is created, it does not contain any hardware channels. To execute the device object, hardware channels must be added with the addchannel function.• You can create multiple analog input objects that are associated with a single analog input subsystem. However, you can typically execute only one of them at a time.• The analog input object exists in the data acquisition engine and in the MATLAB workspace. If you create a copy of the device object, it references the original device object in the engine.• If ID is a numeric value, then you can specify it as an integer or a string. If ID contains any non-numeric characters, then you must specify it as a string (see the Agilent Technologies example below).• The Name property is automatically assigned a descriptive name that is produced by concatenating adaptor, ID, and - AI. You can change this name at any time.	

More About the Hardware Device Identifier

When data acquisition devices are installed, they are assigned a unique number which identifies the device in software. The device identifier is typically assigned automatically and can usually be manually changed using a vendor-supplied device configuration utility. National Instruments refers to this number as the device number while Agilent Technologies refers to it as the device ID.

For sound cards, the device identifier is typically not exposed to you through the Microsoft Windows environment. However, the Data Acquisition Toolbox automatically associates each sound card with an integer ID value. There are two cases to consider:

- If you have one sound card installed, then ID is 0. You are not required to specify ID when creating an analog input object associated with this device.
- If you have multiple sound cards installed, the first one installed has an ID of 0, the second one installed has an ID of 1, and so on. You must specify ID when creating analog input objects associated with devices not having an ID of 0.

There are two ways you can determine the ID for a particular device:

- Type `daqhwinfo('adaptor')`.
- Execute the vendor-supplied device configuration utility.

Example

National Instruments

To create an analog input object for a National Instruments board defined as device number 1

```
AI = analoginput('ni daq', 1);
```

Agilent Technologies

To create an analog input object for an Agilent Technologies module with device identifier 1 residing in VXI chassis 0

```
AI = analoginput('hpe1432', 'vxi 0: : 1: : instr');
```

Alternatively, you can use the syntax

```
AI = analoginput('hpe1432', 1, 0);
```

analoginput

The HP driver allows you to span multiple hardware devices. To create an analog input object that spans two HP devices with device identifiers 1 and 2 residing in VXI chassis 0

```
AI = analoginput('hpe1432', 'vxi0:1,2:instr');
```

Alternatively, you can use the syntax

```
AI = analoginput('hpe1432', [1, 2], 0);
```

See Also

Functions

`addchannel`, `daqhwinfo`

Properties

`Name`

Purpose	Create an analog output object						
Syntax	<pre>A0 = analogoutput('adaptor') A0 = analogoutput('adaptor', ID)</pre>						
Arguments	<table> <tr> <td>'adaptor'</td><td>Name of the hardware driver adaptor. The supported adaptors are ni daq, cbi, hpe1432, and wi nsound.</td></tr> <tr> <td>ID</td><td>The hardware device identifier. ID is optional if the device object is associated with a sound card having an ID of 0.</td></tr> <tr> <td>A0</td><td>The analog output object.</td></tr> </table>	'adaptor'	Name of the hardware driver adaptor. The supported adaptors are ni daq, cbi, hpe1432, and wi nsound.	ID	The hardware device identifier. ID is optional if the device object is associated with a sound card having an ID of 0.	A0	The analog output object.
'adaptor'	Name of the hardware driver adaptor. The supported adaptors are ni daq, cbi, hpe1432, and wi nsound.						
ID	The hardware device identifier. ID is optional if the device object is associated with a sound card having an ID of 0.						
A0	The analog output object.						
Description	<p><code>A0 = analogoutput('adaptor')</code> creates the analog output object A0 for a sound card having an ID of 0 (<i>adaptor</i> must be wi nsound). This is the only case where ID is not required.</p> <p><code>A0 = analogoutput('adaptor', ID)</code> creates the analog output object A0 for the specified adaptor and for the hardware device with device identifier ID. ID can be specified as an integer or a string.</p>						
Remarks	<p>More About Creating Analog Output Objects</p> <ul style="list-style-type: none"> • When an analog output object is created, it does not contain any hardware channels. To execute the device object, hardware channels must be added with the <code>addchannel</code> function. • You can create multiple analog output objects that are associated with a single analog output subsystem. However, you can typically execute only one of them at a time. • The analog output object exists in the data acquisition engine and in the MATLAB workspace. If you create a copy of the device object, it references the original device object in the engine. • If ID is a numeric value, then you can specify it as an integer or a string. If ID contains any non-numeric characters, then you must specify it as a string (see the Agilent Technologies example below). • The Name property is automatically assigned a descriptive name that is produced by concatenating <i>adaptor</i>, ID, and -A0. You can change this name at any time. 						

More About the Hardware Device Identifier

When data acquisition devices are installed, they are assigned a unique number which identifies the device in software. The device identifier is typically assigned automatically and can usually be manually changed using a vendor-supplied device configuration utility. National Instruments refers to this number as the device number while Agilent Technologies refers to it as the device ID.

For sound cards, the device identifier is typically not exposed to you through the Microsoft Windows environment. However, the Data Acquisition Toolbox automatically associates each sound card with an integer ID value. There are two cases to consider:

- If you have one sound card installed, then ID is 0. You are not required to specify ID when creating an analog output object associated with this device.
- If you have multiple sound cards installed, the first one installed has an ID of 0, the second one installed has an ID of 1, and so on. You must specify ID when creating analog output objects associated with devices not having an ID of 0.

There are two ways you can determine the ID for a particular device:

- Type `daqhwinfo('adaptor')`.
- Execute the vendor-supplied device configuration utility.

Example

National Instruments

To create an analog output object for a National Instruments board defined as device number 1

```
A0 = analogoutput('ni daq', 1);
```

Agilent Technologies

To create an analog output object for an Agilent Technologies module with device identifier 1 residing in VXI chassis 0

```
A0 = analogoutput('hpe1432', 'vxi 0: : 1: : instr');
```

Alternatively, you can use the syntax

```
A0 = analogoutput('hpe1432', 1, 0);
```

The HP driver allows you to span multiple hardware devices. To create an analog output object that spans two HP devices with device identifiers 1 and 2 residing in VXI chassis 0

```
A0 = analogoutput('hpe1432', 'vxi 0: : 1, 2: : instr');
```

Alternatively, you can use the syntax

```
A0 = analogoutput('hpe1432', [1, 2], 0);
```

See Also

Functions

`addchannel`, `daqhwinfo`

Properties

Name

binvec2dec

Purpose	Convert binary vector to decimal value				
Syntax	<code>out = binvec2dec(bin)</code>				
Arguments	<table><tr><td><code>bin</code></td><td>A binary vector.</td></tr><tr><td><code>out</code></td><td>A double array.</td></tr></table>	<code>bin</code>	A binary vector.	<code>out</code>	A double array.
<code>bin</code>	A binary vector.				
<code>out</code>	A double array.				
Description	<code>out = binvec2dec(bin)</code> converts the binary vector <code>bin</code> to the equivalent decimal number and stores the result in <code>out</code> . All nonzero binary vector elements are interpreted as a 1.				
Remarks	A binary vector (<code>binvec</code>) is constructed with the least significant bit (LSB) in the first column and the most significant bit (MSB) in the last column. For example, the decimal number 23 is written as the <code>binvec</code> value [1 1 1 0 1].				
Example	<p>To convert the <code>binvec</code> value [1 1 1 0 1] to a decimal value</p> <pre>binvec2dec([1 1 1 0 1]) ans = 23</pre>				
See Also	Functions <code>dec2binvec</code>				

Purpose	Remove device objects from the MATLAB workspace	
Syntax	<code>clear obj</code> <code>clear obj.Channel(index)</code> <code>clear obj.Line(index)</code>	
Arguments	<code>obj</code>	A device object or array of device objects.
	<code>obj.Channel(index)</code>	One or more channels contained by <code>obj</code> .
	<code>obj.Line(index)</code>	One or more lines contained by <code>obj</code> .
Description	<code>clear obj</code> removes <code>obj</code> and all associated channels or lines from the MATLAB workspace, but not from the data acquisition engine. <code>clear obj.Channel(index)</code> removes the specified channels contained by <code>obj</code> from the MATLAB workspace, but not from the data acquisition engine. <code>clear obj.Line(index)</code> removes the specified lines contained by <code>obj</code> from the MATLAB workspace, but not from the data acquisition engine.	
Remarks	Clearing device objects, channels, and lines follows these rules: <ul style="list-style-type: none">• <code>clear</code> does not remove device objects, channels, or lines from the data acquisition engine. Use the <code>delete</code> function for this purpose.• If multiple references to a device object exist in the workspace, clearing one reference will not invalidate the remaining references.• You can restore cleared device objects to the MATLAB workspace with the <code>daqind</code> function. If you use the <code>help</code> command to display the M-file help for <code>clear</code> , then you must supply the pathname shown below. <code>help daq/private/clear</code>	

clear

Example

Create the analog input object `ai`, copy `ai` to a new variable `ai copy`, and then clear the original device object from the MATLAB workspace.

```
ai = analoginput('winsound');  
ch = addchannel(ai, 1:2);  
ai copy = ai;  
clear ai
```

Retrieve `ai` from the engine with `daqfind`, and demonstrate that `ai` is identical to `ai copy`.

```
ai new = daqfind;  
isequal(ai copy, ai new)  
ans =  
    1
```

See Also

Functions

`daqfind`, `delete`

Purpose	An action function that displays event information for the specified event	
Syntax	<code>daqaction(obj, event)</code>	
Arguments	<code>obj</code>	A device object.
	<code>event</code>	A variable that captures the event information contained by the <code>EventLog</code> property.
Description	<code>daqaction(obj, event)</code> is an example action function that displays information to the MATLAB command window. For all events, the information includes the event type and the name of the device object that caused the event to occur. For events that record the absolute time in <code>EventLog</code> , the event time is also displayed. For run-time error events, the error message is also displayed.	
Remarks	You specify <code>daqaction</code> as the action function to be executed for any event by specifying it as the value for the associated action property. For analog input objects, <code>daqaction</code> is the default value for the <code>DataMissedAction</code> and <code>RuntimeErrorAction</code> properties. For analog output objects, <code>daqaction</code> is the default value for the <code>RuntimeErrorAction</code> property.	
	You can use the <code>showdaqevents</code> function to easily display event information captured by the <code>EventLog</code> property.	
Example	Create the analog input object <code>ai</code> and call the action function <code>daqaction</code> when a trigger event occurs. <pre>ai = analoginput('winsound'); addchannel(ai, 1); set(ai, 'TriggerRepeat', 3) set(ai, 'TriggerAction', 'daqaction') start(ai)</pre>	
See Also	Functions <code>showdaqevents</code>	
	Properties <code>DataMissedAction</code> , <code>EventLog</code> , <code>RuntimeErrorAction</code>	

daqfind

Purpose Return device objects, channels, or lines from the data acquisition engine to the MATLAB workspace

Syntax

```
out = daqfind
out = daqfind('PropertyName', PropertyValue, ...)
out = daqfind(S)
out = daqfind(obj, 'PropertyName', PropertyValue, ...)
```

Arguments

<i>'PropertyName'</i>	A device object, channel, or line property name.
<i>PropertyVal ue</i>	A device object, channel, or line property value.
<i>obj</i>	A device object, array of device objects, channels, or lines.
<i>S</i>	A structure with field names that are property names and field values that are property values.
<i>out</i>	An array or cell array of device objects, channels, or lines.

Description

`out = daqfind` returns all device objects that exist in the data acquisition engine. The output `out` is an array.

`out = daqfind('PropertyName', PropertyValue, ...)` returns all device objects, channels, or lines that exist in the data acquisition engine and have the specified property names and property values. The property name/property value pairs can be specified as a cell array.

`out = daqfind(S)` returns all device objects, channels, or lines that exist in the data acquisition and have the property names and property values specified by `S`. `S` is a structure with field names that are property names and field values that are property values.

`out = daqfind(obj, 'PropertyName', PropertyValue, ...)` returns all device object, channels, or lines listed by `obj` that have the specified property names and property values.

Remarks**More About Finding Device Objects, Channels, or Lines**

daqfind is particularly useful in these circumstances:

- A device object is cleared from the MATLAB workspace, and it needs to be retrieved from the data acquisition engine.
- You need to locate device objects, channels, or lines that have particular property names and property values.

Rules for Specifying Property Names and Property Values

- You can use property name/property value string pairs, structures, and cell array pairs in the same call to daqfind. However, in a single call to daqfind, you can specify only device object properties or channel/line properties.
- You must use the same format as returned by get. For example, if get returns the Channel Name property value as Left, you must specify Left as the property value in daqfind (case matters). However, case does not matter when you specify enumerated property values. For example, daqfind will find a device object with a Running property value of On or on.

Example

You can use daqfind to return a cleared device object.

```
ai = analoginput('winsound');
ch = addchannel(ai, 1:2);
set(ch, {'Channel Name'}, {'Joe'; 'Jack'})
clear ai
ai_new = daqfind;
```

To return the channel associated with the descriptive name Jack

```
ch2 = daqfind(ai_new, 'Channel Name', 'Jack');
```

To return the device object with a sampling rate of 8000 Hz and the descriptive name winsound0-AI, you can pass a structure to daqfind.

```
S.Name = 'winsound0-AI';
S.SampleRate = 8000;
daqobj = daqfind(S);
```

See Also**Functions**

clear, get, propinfo

daqhelp

Purpose	Display help for device objects, constructors, adaptors, functions, and properties	
Syntax	<pre>daqhelp out = daqhelp('name') out = daqhelp(obj) out = daqhelp(obj, 'name')</pre>	
Arguments	'name'	A device object, constructor, adaptor, function, or property name.
	obj	A device object.
	out	Contains the specified help text.
Description	<p>daqhelp displays a complete listing of Data Acquisition Toolbox constructors and functions along with a brief description of each.</p> <p>out = daqhelp('name') returns help for the device object, constructor, adaptor, function, or property specified by name. The help text is returned to out.</p> <p>out = daqhelp(obj) returns a complete listing of functions and properties for the device object obj to out. Help for obj's constructor is also displayed.</p> <p>out = daqhelp(obj, 'name') returns help for name for the specified device object obj to out. name can be a constructor, adaptor, property, or function name.</p>	
Remarks	More About Displaying Help <ul style="list-style-type: none">• When displaying property help, the names in the See Also section that contain all upper-case letters are function names. The names that contain a mixture of upper and lower-case letters are property names.• When displaying function help, the See Also section contains only function names.	

Rules for Specifying Names

For the `daqhelp('name')` syntax:

- If `name` is the name of a device object constructor, a complete listing of the device object's functions and properties is displayed along with a brief description of each function and property. The help for the device object's constructor is also displayed.
- You can display object-specific function information by specifying `name` as `object/function`. For example, to display the help for an analog input object's `getData` function, `name` is `analoginput/getData`.
- You can display object-specific property information by specifying `name` as `obj.property`. For example, to display the help for an analog input object's `SampleRate` property, `name` is `analoginput.SampleRate`.

For the `daqhelp(obj, 'name')` syntax:

- If `name` is the name of a device object constructor and the `.m` extension is included, the constructor help is displayed.
- If `name` is the name of a function or property, the function or property help is displayed.

Example

The following commands are some of the ways you can use `daqhelp` to obtain help on device objects, constructors, adaptors, functions, and properties.

```
daqhelp('analogoutput');
out = daqhelp('analogoutput.m');
daqhelp set
daqhelp analoginput/peekdata
daqhelp analoginput.TriggerDelayUnits
```

The following commands are some of the ways you can use `daqhelp` to obtain information about functions and properties for an existing device object.

```
ai = analoginput('winsound');
daqhelp(ai, 'InitialTriggerTime')
out = daqhelp(ai, 'getsample');
```

See Also

Functions
`propinfo`

daqhwinfo

Purpose Display data acquisition hardware information

Syntax

```
out = daqhwinfo
out = daqhwinfo('adaptor')
out = daqhwinfo(obj)
out = daqhwinfo(obj, 'FieldName')
```

Arguments

<i>'adaptor'</i>	Name of the hardware driver adaptor. The supported adaptors are ni daq, cbi, hpe1432, and wi nsound.
<i>obj</i>	A device object or array of device objects.
<i>'FieldName'</i>	A single field name or a cell array of field names.
<i>out</i>	A structure containing the requested hardware information.

Description `out = daqhwinfo` returns general hardware-related information as a structure to `out`. The returned information includes installed adaptors, the toolbox and MATLAB version, and the toolbox name.

`out = daqhwinfo('adaptor')` returns hardware-related information for the specified *adaptor*. The returned information includes the adaptor DLL name, the board names and IDs, and the device object constructor syntax.

`out = daqhwinfo(obj)` returns hardware-related information for the device object `obj`. If `obj` is an array of device objects, then `out` is a 1-by-`n` cell array of structures where `n` is the length of `obj`. The returned information depends on the device object type, and may include the maximum and minimum sampling rates, the channel gains, the hardware channel or line IDs, and the vendor driver version.

`out = daqhwinfo(obj, 'FieldName')` returns the hardware-related information specified by *FieldName* for the device object `obj`. *FieldName* can be a single field name or a cell array of field names. `out` is an `m`-by-`n` cell array where `m` is the length of `obj` and `n` is the length of *FieldName*. You can return a list of valid field names with the `daqhwinfo(obj)` syntax.

Example

To display all installed adaptors

```
out = daqhwi nfo;
out.Instal ledAdaptors
ans =
    ' cbi '
    ' hpe1432'
    ' ni daq'
    ' wi nsound'
```

To display the device object constructor names for all installed wi nsound devices

```
out = daqhwi nfo(' wi nsound' );
out. Obj ectConstructorName
ans =
    ' anal ogi nput (' wi nsound' , 0) '
    ' anal ogout put (' wi nsound' , 0) '
```

Create the analog input object ai for a sound card. To display the input ranges for ai

```
ai = anal ogi nput (' wi nsound' );
out = daqhwi nfo(ai );
out. InputRanges
ans =
    - 1      1
```

To display the minimum and maximum sampling rates for ai

```
out = daqhwi nfo(ai , { ' Mi nSampl eRate' , ' MaxSampl eRate' })
out =
    [ 8000]      [ 44100]
```

daqmem

Purpose Allocate or display memory resources

Syntax

```
out = daqmem
out = daqmem(obj)
daqmem(obj, maxmem)
```

Arguments

obj	A device object or array of device objects.
maxmem	The amount of memory to allocate.
out	A structure containing information about memory resources.

Description out = daqmem returns the structure out, which contains several fields describing the memory resources associated with your platform and the Data Acquisition Toolbox. The fields are described below.

Field	Description
MemoryLoad	Specifies a number between 0 and 100 that gives a general idea of current memory utilization. 0 indicates no memory use and 100 indicates full memory use.
Total Phys	Indicates the total number of bytes of physical memory.
Avai l Phys	Indicates the number of bytes of physical memory available.
Total PageFi l e	Indicates the total number of bytes that can be stored in the paging file. Note that this number does not represent the actual physical size of the paging file on disk.
Avai l PageFi l e	Indicates the number of bytes available in the paging file.
Total Vi rtual	Indicates the total number of bytes that can be described in the user mode portion of the virtual address space of the calling process.

Field	Description
AvailVirtual	Indicates the number of bytes of unreserved and uncommitted memory in the user mode portion of the virtual address space of the calling process.
UsedDaq	The total memory used by all device objects.

Note that all the above fields, except for UsedDaq, are identical to the fields returned by Windows' MemoryStatus function.

`out = daqmem(obj)` returns a 1-by-N structure `out` containing two fields: `UsedBytes` and `MaxBytes` for the device object `obj`. `N` is the number of device objects specified by `obj`. `UsedBytes` returns the number of bytes used by `obj`. `MaxBytes` returns the maximum number of bytes that can be used by `obj`.

`daqmem(obj, maxmem)` sets the maximum memory that can be allocated for `obj` to the value specified by `maxmem`.

Remarks

More About Allocating and Displaying Memory Resources

- For analog output objects, `daqmem(obj, maxmem)` controls the value of the `MaxSamplesQueued` property.
- If you manually configure the `BufferingConfig` property, then this value supersedes the values specified by `daqmem(obj, maxmem)` and the `MaxSamplesQueued` property.

Example

Create the analog input object `aiwin` for a sound card and the analog input object `aini` for a National Instruments board, and add two channels to each device object.

```
aiwin = analoginput('winsound');
addchannel(aiwin, 1:2);
aini = analoginput('ni daq', 1);
addchannel(aini, 0:1);
```

daqmem

To display the total memory used by all existing device objects

```
out = daqmem;  
out.UsedDaq  
ans =  
    69120
```

To configure the maximum memory used by aiwi n to 640 KB

```
daqmem(aiwi n, 640000)
```

To configure the maximum memory used by each device object with one call to daqmem

```
daqmem([aiwi n ai ni ], [640000 480000])
```


See Also

Properties

BufferingConfig, MaxSamplesQueued

Purpose	Invoke the property editor graphical user interface
Syntax	<code>daqpropedit</code> <code>daqpropedit (obj)</code>
Arguments	<code>obj</code> A device object.
Description	<code>daqpropedit (obj)</code> invokes the property editor graphical user interface (GUI) for the device object <code>obj</code> . The channels or lines contained by <code>obj</code> are displayed as well.
Remarks	<p>Using <code>daqpropedit</code>, you can perform these tasks for device objects, channels, and lines:</p> <ul style="list-style-type: none">• List all existing device objects as well as the channels or lines they contain. This is identical to using the <code>daqfind</code> function.• Configure property values. This is identical to using the <code>set</code> function or the dot notation.• Display the names and current values for configurable properties. This is identical to using the <code>get</code> function or the dot notation.• Display property characteristics. This is identical to using the <code>propinfo</code> function.• Display property help. This is identical to using the <code>daqhelp</code> function.
Example	<p>Create the analog input object <code>ai</code> for a sound card and add one channel to it.</p> <pre>ai = analoginput('winsound'); addchannel(ai, 1:2);</pre> <p>To configure property values for <code>ai</code> using <code>daqpropedit</code></p> <pre>daqpropedit(ai)</pre>

The Data Acquisition Property Editor is shown below.



See Also

Functions
daqfind, daqhelp, get, propinfo, set

Purpose	Read a Data Acquisition Toolbox (.daq) file	
Syntax	<pre> data = daqread('file') data = daqread('file', 'PropertyName', PropertyValue, ...) [data, time] = daqread(...) [data, time, abstime] = daqread(...) [data, time, abstime, events] = daqread(...) [data, time, abstime, events, daqinfo] = daqread(...) daqinfo = daqread('file', 'info') </pre>	
Arguments	'file'	A Data Acquisition Toolbox (.daq) file.
	'PropertyName'	A daqread property name.
	PropertyValue	A daqread property value.
	'info'	Specifies that device object and hardware information are to be returned.
	data	An m-by-n array where m is the number of samples and n is the number of channels.
	time	An m-by-1 array of relative time values where m is the number of samples.
	abstime	The absolute time of the first trigger.
	events	A structure containing event information.
	daqinfo	A structure containing device object and hardware information.
Description	<p>data = daqread('file') reads all the data from file. data is an m-by-n data matrix where m is the number of samples and n is the number of channels. If data includes data from multiple triggers, then m is increased by the number of triggers due to the addition of NaN's.</p>	

daqread

`data = daqread('file', 'PropertyName', PropertyValue, ...)` reads the specified data from `file`. The amount of data returned and the format of the data is specified with the properties shown below.

Property Name	Description
<code>Samples</code>	Specify the sample range.
<code>Time</code>	Specify the relative time range.
<code>Triggers</code>	Specify the trigger range.
<code>Channels</code>	Specify the channel range. Channel names can be specified as a cell array.
<code>DataFormat</code>	Specify the data format as <code>doubles</code> or <code>native</code> .
<code>TimeFormat</code>	Specify the time format as <code>vector</code> or <code>matrix</code> .

The `Samples`, `Time`, and `Triggers` properties are mutually exclusive.

`[data, time] = daqread(...)` returns sample-time pairs. `time` is a vector with the same length as `data` and contains the relative time for each sample. Relative time is measured with respect to the first trigger that occurs.

`[data, time, abstime] = daqread(...)` returns sample-time pairs and the absolute time of the first trigger. `abstime` is returned as a clock vector.

`[data, time, abstime, events] = daqread(...)` returns sample-time pairs, the absolute time of the first trigger, and a log of events. `events` contains the appropriate events based on the `Samples`, `Time`, or `Triggers` value specified. The entire event log is returned only if `Samples`, `Time`, or `Triggers` is not specified.

`[data, time, abstime, events, daqinfo] = daqread(...)` returns sample-time pairs, the absolute time, the event log, and the structure `daqinfo`, which contains two fields: `ObjInfo` and `HwInfo`. `ObjInfo` is a structure containing property name/property value pairs and `HwInfo` is a structure containing hardware information. The entire event log is returned to `daqinfo.ObjInfo.EventLog`.

`daqinfo = daqread('file', 'info')` returns the structure `daqinfo`, which contains two fields: `ObjInfo` and `HwInfo`. `ObjInfo` is a structure containing property name/property value pairs and `HwInfo` is a structure containing hardware information. The entire event log is returned to `daqinfo.ObjInfo.EventLog`.

Remarks

More About .daq Files

- The format used by `daqread` to return data, relative time, absolute time, and event information is identical to the format used by the `getdata` function.
- If data from multiple triggers is read, then the size of the resulting data array is increased by the number of triggers issued since each trigger is separated by a NaN.
- `ObjInfo.EventLog` always contains the entire event log regardless of the value specified by `Samples`, `Time`, or `Triggers`.
- Data Acquisition Toolbox (.daq) files are created by specifying a value for the `LogFileNames` property (or accepting the default value), and configuring the `LoggingMode` property to `Disk` or `Disk&Memory`.

Example

Suppose you configure the analog input object `ai` for a National Instruments board as shown below. The object acquires one second of data for four channels, and saves the data to the output file `data.daq`.

```
ai = analoginput('ni daq', 1);
chans = addchannel(ai, 0:3);
set(ai, 'SampleRate', 1000)
ActualRate = get(ai, 'SampleRate');
set(ai, 'SamplesPerTrigger', ActualRate)
set(ai, 'LoggingMode', 'Disk&Memory')
set(ai, 'LogFileName', 'data.daq')
start(ai)
```

After the data has been collected and saved to a disk file, you can retrieve the data and other acquisition-related information using `daqread`. To read all the sample-time pairs from `data.daq`

```
[data, time] = daqread('data.daq');
```

daqread

To read samples 500 to 1000 for all channels from data.daq

```
data = daqread('data.daq', 'Samples', [500 1000]);
```

To read the first 0.5 seconds of data for channels 1 and 2 from data.daq

```
data = daqread('data.daq', 'Time', [0 0.5], 'Channels', [1 2]);
```

To obtain the channel property information from data.daq

```
daqinfo = daqread('data.daq', 'info');  
chaninfo = daqinfo.ObjInfo.Channel;
```

To obtain a list of event types and event data contained by data.daq

```
daqinfo = daqread('data.daq', 'info');  
events = daqinfo.ObjInfo.EventLog;  
event_type = {events.Type};  
event_data = {events.Data};
```

See Also

Functions

getdata

Properties


EventLog, LogFileName, LoggingMode, LogToDiskMode

Purpose	Register or unregister a hardware driver adaptor
Syntax	<pre>daqregister('adaptor') daqregister('adaptor', 'unload') out = daqregister(...)</pre>
Arguments	<p>'<i>adaptor</i>' Name of the hardware driver adaptor. The supported adaptors are ni daq, cbi, hpe1432, and wi nsound.</p> <p>'unload' Specifies that the hardware driver adaptor is to be unloaded.</p> <p>out Captures the message returned by daqregister.</p>
Description	<p>daqregister('adaptor') registers the hardware driver adaptor specified by <i>adaptor</i>. For third-part adaptors, <i>adaptor</i> must include the full pathname.</p> <p>daqregister('adaptor', 'unload') unregisters the hardware driver adaptor specified by <i>adaptor</i>. For third-part adaptors, <i>adaptor</i> must include the full pathname.</p> <p>out = daqregister(...) captures the resulting message in out.</p>
Remarks	<p>A hardware driver adaptor must be registered so the data acquisition engine can make use of its services. Unless an adaptor is unloaded, registration is required only once.</p> <p>For adaptors that are included with the toolbox, registration occurs automatically when you first create a device object. However, you may need to register third-party adaptors manually. In either case, you must install the associated hardware driver before registration can occur.</p>
Example	<p>The following command registers the sound card adaptor provided with the toolbox.</p> <pre>daqregister('wi nsound');</pre> <p>The following command registers the third-party adaptor <code>myadaptor.dll</code>. Note that you must supply the full pathname to <code>daqregister</code>.</p> <pre>daqregister('D: /MATLABR12/tool box/daq/myadaptors/ myadaptor.dll');</pre>

daqreset

Purpose	Remove device objects and data acquisition DLLs from memory
Syntax	<code>daqreset</code>
Description	<p><code>daqreset</code> removes all device objects existing in the engine, and unloads all data acquisition DLLs loaded by the engine (including the adaptor and engine DLL files).</p> <p>You should use <code>daqreset</code> to return MATLAB to the known initial state of having no device objects and no data acquisition DLLs loaded in memory. When MATLAB is returned to this state, the data acquisition hardware is reset.</p>
See Also	Functions <code>clear</code> , <code>delete</code>

Purpose	Interface for displaying toolbox tutorials
Syntax	daqschool
Description	daqschool launches the Data Acquisition Toolbox Tutorials interface, which is shown below.



Refer to “Demos” on page 2-13 for a list of demos included with daqschool .

dec2binvec

Purpose	Convert decimal value to binary vector				
Syntax	<pre>out = dec2bi nvec(dec) out = dec2bi nvec(dec, bi ts)</pre>				
Arguments	dec	A decimal value. dec must be nonnegative.			
	bi ts	Number of bits used to represent the decimal number.			
	out	A logical array containing the binary vector.			
Description	out = dec2bi nvec(dec) converts the decimal value dec to an equivalent binary vector and stores the result as a logical array in out.				
	out = dec2bi nvec(dec, bi ts) converts the decimal value dec to an equivalent binary vector consisting of at least the number of bits specified by bi ts.				
Remarks	More About Binary Vectors				
	A binary vector (binvec) is constructed with the least significant bit (LSB) in the first column and the most significant bit (MSB) in the last column. For example, the decimal number 23 is written as the binvec value [1 1 1 0 1].				
	More About Specifying the Number of Bits				
	<ul style="list-style-type: none">• If bi ts is greater than the minimum number of bits required to represent the decimal value, then the result is padded with zeros.• If bi ts is less than the minimum number of bits required to represent the decimal value, then the minimum number of required bits is used.• If bi ts is not specified, then the minimum number of bits required to represent the number is used.				
Example	To convert the decimal value 23 to a binvec value				
	<pre>dec2bi nvec(23) ans = 1 1 1 0 1</pre>				

To convert the decimal value 23 to a binvec value using six bits

```
dec2binvec(23, 6)
```

```
ans =
```

```
1    1    1    0    1    0
```

To convert the decimal value 23 to a binvec value using four bits, then the result uses five bits. This is the minimum number required to represent the number.

```
dec2binvec(23, 4)
```

```
ans =
```

```
1    1    1    0    1
```

See Also

Functions

binvec2dec

delete

Purpose Remove device objects, channels, or lines from the data acquisition engine

Syntax

```
delete(obj)
delete(obj.Channel(index))
delete(obj.Line(index))
```

Arguments

obj	A device object or array of device objects.
obj.Channel(index)	One or more channels contained by obj.
obj.Line(index)	One or more lines contained by obj.

Description

delete(obj) removes the device object specified by obj from the engine. If obj contains channels or lines, they are removed as well. If obj is the last object accessing the driver, then the driver and associated adaptor are unloaded.

delete(obj.Channel(index)) removes the channels specified by index and contained by obj from the engine. As a result, the remaining channels may be reindexed.

delete(obj.Line(index)) removes the lines specified by index and contained by obj from the engine. As a result, the remaining lines may be reindexed.

Remarks

Deleting device objects, channels, and lines follows these rules:

- delete removes device objects, channels, or lines from the data acquisition engine but not from the MATLAB workspace. To remove variables from the workspace, use the clear function.
- If multiple references to a device object exist in the workspace, then removing one device object from the engine invalidates the remaining references. These remaining references should be cleared from the workspace with the clear function.
- If you delete a device object while it is running, then a warning is issued before it is deleted. You cannot delete a device object while it is logging or sending data.

You should use delete at the end of a data acquisition session. You can quickly delete all existing device objects with the command delete(daqlind).

If you use the `help` command to display the M-file help for `delete`, then you must supply the pathname shown below.

```
help daq/daqdevice/delete
```

Example

National Instruments

Create the analog input object `ai` for a National Instruments board, add hardware channels 0-7 to it, and make a copy of hardware channels 0 and 1.

```
ai = analoginput('ni daq', 1);  
addchannel(ai, 0:7);  
ch = ai.Channel(1:2);
```

To delete hardware channels 0 and 1

```
delete(ch)
```

These channels are deleted from the data acquisition engine and are no longer associated with `ai`. The remaining channels are reindexed such that the indices begin at 1 and increase monotonically to 6. To delete `ai`

```
delete(ai)
```

Sound Card

Create the analog input object `AI 1` for a sound card, and configure it to operate in stereo mode.

```
AI 1 = analoginput('winsound');  
addchannel(AI 1, 1:2);
```

You can now configure the sound card for mono mode by deleting hardware channel 2.

```
delete(AI 1.Channel(2))
```

If hardware channel 1 is deleted instead, an error is returned.

See Also

Functions

`clear`, `daqreset`

Purpose	Create a digital I/O object	
Syntax	<code>DI0 = digitalio('adaptor', ID)</code>	
Arguments	'adaptor'	Name of the hardware driver adaptor. The supported adaptors are ni daq and cbi .
	ID	The hardware device identifier.
	DI0	The digital I/O object.
Description	DI0 = digitalio('adaptor', ID) creates the digital I/O object DI0 for the specified adaptor and for the hardware device with device identifier ID. ID can be specified as an integer or a string.	
Remarks	<p>More About Creating Digital I/O Objects</p> <ul style="list-style-type: none">• When a digital I/O object is created, it does not contain any hardware lines. To execute the device object, hardware lines must be added with the addline function.• You can create multiple digital I/O objects that are associated with a single digital I/O subsystem. However, you can typically execute only one of them at a time.• The digital I/O object exists in the data acquisition engine and in the MATLAB workspace. If you create a copy of the device object, it references the original device object in the engine.• The Name property is automatically assigned a descriptive name that is produced by concatenating adaptor, ID, and - DI0. You can change this name at any time. <p>More About the Hardware Device Identifier</p> <p>When data acquisition devices are installed, they are assigned a unique number, which identifies the device in software. The device identifier is typically assigned automatically and can usually be manually changed using a vendor-supplied device configuration utility. National Instruments refers to this number as the device number.</p>	

There are two ways you can determine the ID for a particular device:

- Type `daqhwinfo('adaptor')`.
- Execute the vendor-supplied device configuration utility.

Example

To create a digital I/O object for a National Instruments device defined as device number 1

```
DI0 = digitalio('niDAQ', 1);
```

See Also**Functions**

`addline`, `daqhwinfo`

Properties

Name

disp

Purpose	Display summary information for device objects, channels, or lines	
Syntax	<code>di sp(obj)</code> <code>di sp(obj . Channel (i ndex))</code> <code>di sp(obj . Li ne(i ndex))</code>	
Arguments	<code>obj</code>	A device object.
	<code>obj . Channel (i ndex)</code>	One or more channels contained by <code>obj</code> .
	<code>obj . Li ne(i ndex)</code>	One or more lines contained by <code>obj</code> .
Description	<p><code>di sp(obj)</code> displays summary information for the specified device object <code>obj</code> , and any channels or lines contained by <code>obj</code> . Typing <code>obj</code> at the command line produces the same summary information.</p> <p><code>di sp(obj . Channel (i ndex))</code> displays summary information for the specified channels contained by <code>obj</code> . Typing <code>obj . Channel (i ndex)</code> at the command line produces the same summary information.</p> <p><code>di sp(obj . Li ne(i ndex))</code> displays summary information for the specified lines contained by <code>obj</code> . Typing <code>obj . Li ne(i ndex)</code> at the command line produces the same summary information.</p>	
Remarks	You can invoke <code>di sp</code> by typing the device object at the MATLAB command line or by excluding the semicolon when creating a device object, when adding channel or lines, or when configuring property values using the dot notation.	
Example	All the commands shown below produce summary information for the device object <code>AI</code> or the channels contained by <code>AI</code> . <pre>AI = anal ogi nput (' wi nsound') chans = addchannel (AI , 1: 2) AI . Sampl eRate = 44100 AI . Channel (1) . Channel Name = ' CH1' chans</pre>	

Purpose	Remove data from the data acquisition engine
Syntax	<code>flushdata(obj)</code> <code>flushdata(obj, 'mode')</code>
Arguments	<p><code>obj</code> An analog input object or array of analog input objects.</p> <p><code>'mode'</code> Specifies how much data is removed from the engine.</p>
Description	<p><code>flushdata(obj)</code> removes all data from the data acquisition engine and resets the <code>SamplesAvailable</code> property to zero.</p> <p><code>flushdata(obj, 'mode')</code> removes data from the data acquisition engine depending on the value of <code>mode</code>. If <code>mode</code> is <code>all</code>, all data is removed from the engine and the <code>SamplesAvailable</code> property is set to 0. This is the same as <code>flushdata(obj)</code>. If <code>mode</code> is <code>triggers</code>, then data is removed from the engine only when the data corresponds to an integral number of triggers. <code>triggers</code> is a valid choice only when the <code>TriggerRepeat</code> property is greater than 0 and the <code>SamplesPerTrigger</code> property is not <code>inf</code>.</p>
Example	<p>Create the analog input object <code>ai</code> for a National Instruments board and add hardware channels 0-7 to it.</p> <pre>ai = analoginput('ni daq', 1); addchannel(ai, 0:7);</pre> <p>A two-second acquisition is configured and the device object is executed.</p> <pre>set(ai, 'SampleRate', 2000) duration = 2; ActualRate = get(ai, 'SampleRate'); set(ai, 'SamplesPerTrigger', ActualRate*duration) start(ai)</pre> <p>Four thousand samples will be acquired for each channel group member. To extract 1000 samples from the data acquisition engine for each channel</p> <pre>data = getdata(ai, 1000);</pre>

flushdata

You can use `flushdata` to remove the remaining 3000 samples from the data acquisition engine.

```
flushdata(ai)
ai.SamplesAvailable
ans =
    0
```

See Also

Functions

`getdata`

Properties

`SamplesAvailable`, `SamplesPerTrigger`, `TriggerRepeat`

Purpose	Return device object properties	
Syntax	<pre> out = get(obj) out = get(obj.Channel(index)) out = get(obj.Line(index)) out = get(obj, 'PropertyName') out = get(obj.Channel(index), 'PropertyName') out = get(obj.Line(index), 'PropertyName') get(...) </pre>	
Arguments	obj	A device object or array of device objects.
	obj.Channel(index)	One or more channels contained by obj.
	obj.Line(index)	One or more lines contained by obj.
	'PropertyName'	A property name or a cell array of property names.
Description	<p>out = get(obj) returns the structure out, where each field name is the name of a property of obj and each field contains the value of that property.</p> <p>out = get(obj.Channel(index)) returns the structure out, where each field name is the name of a channel property of obj and each field contains the value of that property.</p> <p>out = get(obj.Line(index)) returns the structure out, where each field name is the name of a line property of obj and each field contains the value of that property.</p> <p>out = get(obj, 'PropertyName') returns the value of the property specified by <i>PropertyName</i> to out. If <i>PropertyName</i> is replaced by a 1-by-n or n-by-1 cell array of strings containing property names, then get returns a 1-by-n cell array of values to out. If obj is an array of data acquisition objects, then out will be an m-by-n cell array of property values where m is equal to the length of obj and n is equal to the number of properties specified.</p> <p>out = get(obj.Channel(index), 'PropertyName') returns the value of <i>PropertyName</i> to out for the specified channels contained by obj. If multiple</p>	

channels and multiple property names are specified, then out is an m-by-n cell array where m is the number of channels and n is the number of properties.

`out = get(obj.Line(index), 'PropertyName')` returns the value of *PropertyName* to out for the specified lines contained by obj. If multiple lines and multiple property names are specified, then out is an m-by-n cell array where m is the number of lines and n is the number of properties.

`get(...)` displays all property names and their current values for the specified device object, channel, or line. Base properties are displayed first followed by device-specific properties.

Remarks

If you use the `help` command to display the M-file help for `get`, then you must supply the pathname shown below.

```
help daq/daqdevice/get
```

Example

Create the analog input object `ai` for a sound card and configure it to operate in stereo mode.

```
ai = analoginput('winsound');  
addchannel(ai, 1:2);
```

The commands shown below are some of the ways you can use `get` to return property values.

```
chan = get(ai, 'Channel');  
out = get(ai, {'SampleRate', 'TriggerDelayUnits'});  
out = get(ai);  
get(chan(1), 'Units')  
get(chan, {'Index', 'HwChannel', 'ChannelName'})
```

See Also

Functions

`set`, `setverify`

Purpose Extract data, time, and event information from the data acquisition engine

Syntax

```
data = getdata(obj)
data = getdata(obj, samples)
data = getdata(obj, 'type')
data = getdata(obj, samples, 'type')
[data, time] = getdata(...)
[data, time, abstime] = getdata(...)
[data, time, abstime, events] = getdata(...)
```

Arguments

<code>obj</code>	An analog input object.
<code>samples</code>	The number of samples to extract. If <code>samples</code> is not specified, the number of samples extracted is given by the <code>SamplesPerTrigger</code> property.
<code>'type'</code>	Specifies the format of the extracted data as <code>double</code> (the default) or as <code>native</code> .
<code>data</code>	An m-by-n array where <code>m</code> is the number of samples extracted and <code>n</code> is the number of channels contained by <code>obj</code> .
<code>time</code>	An m-by-1 array of relative time values where <code>m</code> is the number of samples extracted. Relative time is measured with respect to the first sample logged by the engine.
<code>abstime</code>	The absolute time of the first trigger returned as a clock vector. This value is identical to the value stored by the <code>InitialTriggerTime</code> property.
<code>events</code>	A structure containing a list of events that occurred up to the time of the <code>getdata</code> call.

Description

`data = getdata(obj)` extracts the number of samples specified by the `SamplesPerTrigger` property for each channel contained by `obj`. `data` is an m-by-n array where `m` is the number of samples extracted and `n` is the number of channels.

`data = getdata(obj, samples)` extracts the number of samples specified by `samples` for each channel contained by `obj`.

`data = getdata(obj, 'type')` extracts data in the specified format. If *type* is specified as `native`, the data is returned in the native data format of the device. If *type* is specified as `double` (the default), the data is returned as doubles.

`data = getdata(obj, samples, 'type')` extracts the number of samples specified by `samples` in the format specified by *type* for each channel contained by `obj`.

`[data, time] = getdata(...)` returns data as sample-time pairs. `time` is an `m-by-1` array of relative time values where `m` is the number of samples returned. Relative time is measured with respect to the first sample logged by the engine.

`[data, time, abstime] = getdata(...)` extracts data as sample-time pairs and returns the absolute time of the trigger. The absolute time is returned as a clock vector and is identical to the value stored by the `InitialTriggerTime` property.

`[data, time, abstime, events] = getdata(...)` extracts data as sample-time pairs, returns the absolute time of the trigger, and returns a structure containing a list of events that occurred up to the `getdata` call. The possible events that can be returned are identical to those stored by the `EventLog` property.

Remarks

More About getdata

- In most circumstances, `getdata` returns all requested data and does not miss any samples. In the unlikely event that the engine cannot keep pace with the hardware device, it is possible that data is missed. If data is missed, the `DataMissedAction` property is called and the device object is stopped.
- `getdata` is a *blocking* function since it returns execution control to the MATLAB workspace only when the requested number of samples are extracted from the engine for each channel group member.
- You can issue `^C` (**Control+C**) while `getdata` is blocking. This will not stop the acquisition but will return control to MATLAB.
- The amount of data that you can extract from the engine is given by the `SamplesAvailable` property.

More About Extracting Data From the Engine

- Once the requested data is extracted from the engine, the `SamplesAvailable` property value is automatically reduced by the number of samples returned.
- If the requested number of samples is greater than the samples to be acquired, then an error is returned.
- If the requested data is not returned in the expected amount of time, an error is returned. The expected time to return data is given by the time it takes the engine to fill one data block plus the time specified by the `Timeout` property.
- If multiple triggers are included in a single `getdata` call, a NaN is inserted into the returned data and time arrays and the absolute time returned is given by the first trigger.
- MATLAB supports math operations only for the double data type. Therefore, to use math functions on native data, you must convert it to doubles.

Example

Create the analog input object `ai` for a National Instruments board and add hardware channels 0-3 to it.

```
ai = analoginput('ni daq', 1);
addchannel(ai, 0:3);
```

Configure a one second acquisition with `SampleRate` set to 1000 samples per second and `SamplesPerTrigger` set to 1000 samples per trigger.

```
set(ai, 'SampleRate', 1000)
set(ai, 'SamplesPerTrigger', 1000)
start(ai)
```

The following `getdata` command blocks execution control until all sample-time pairs, the absolute time of the trigger, and any events that occurred during the `getdata` call are returned.

```
[data, time, abstime, events] = getdata(ai);
```

`data` is returned as a 1000-by-4 array of doubles, `time` is returned as a 1000-by-1 vector of relative times, `abstime` is returned as a clock vector, and `events` is returned as a 3-by-1 structure array.

getdata

The three events returned are the start event, the trigger event, and the stop event. To return specific event information about the stop event, you must access the `Type` and `Data` fields.

```
EventType = events(3).Type;  
EventData = events(3).Data;
```

See Also

Functions

`flushdata`, `getsample`, `peekdata`

Properties

`DataMissedAction`, `EventLog`, `SamplesAvailable`, `SamplesPerTrigger`, `Timeout`

Purpose	Immediately acquire one sample
Syntax	<code>sample = getsample(obj)</code>
Arguments	<p><code>obj</code> An analog input object.</p> <p><code>sample</code> A row vector containing one sample for each channel contained by <code>obj</code>.</p>
Description	<code>sample = getsample(obj)</code> immediately returns a row vector containing one sample for each channel contained by <code>obj</code> .
Remarks	<p>Using <code>getsample</code> is a good way to test your analog input configuration. Additionally:</p> <ul style="list-style-type: none"> • <code>getsample</code> does not store samples in, or extract samples from, the data acquisition engine. • You can execute <code>getsample</code> at any time after channels have been added to <code>obj</code>. • Except for sound cards, you can use <code>getsample</code> on an analog input object that is not running (<code>Running</code> is <code>Off</code>). For sound cards, the device object must be running.
Example	<p>Create the analog input object <code>ai</code> and add eight channels to it.</p> <pre>ai = analoginput('ni daq', 1); ch = addchannel(ai, 0:7);</pre> <p>The following command returns one sample for each channel.</p> <pre>sample = getsample(ai);</pre>
See Also	<p>Functions</p> <p><code>getdata</code>, <code>peekdata</code></p>

getvalue

Purpose Read values from lines

Syntax
`out = getvalue(obj)`
`out = getvalue(obj.Line(index))`

Arguments

<code>obj</code>	A digital I/O object.
<code>obj.Line(index)</code>	One or more lines contained by <code>obj</code> .
<code>out</code>	A binary vector.

Description

`out = getvalue(obj)` returns the current value from all lines contained by `obj` as a binary vector to `out`.

`out = getvalue(obj.Line(index))` returns the current value from the lines specified by `obj.Line(index)`.

Remarks **More About Reading Values from Lines**

- By default, `out` is returned as a binary vector (binvec). A binvec value is constructed with the least significant bit (LSB) in the first column and the most significant bit (MSB) in the last column. For example, the decimal number 23 is written as the binvec value [1 1 1 0 1].
- You can convert a binvec value to a decimal value with the `binvec2dec` function.
- If `obj` contains lines from a port-configurable device, the data acquisition engine will automatically read from all the lines even if they are not contained by the device object.

Example Create the digital I/O object `di o` and add eight input lines to it.

```
di o = digitalio('ni daq', 1);  
lines = addline(di o, 0: 7, 'i n');
```

To return the current values from all lines contained by `di o` as a binvec value

```
out = getvalue(di o);
```

See Also **Functions**
`binvec2dec`

Purpose	Check for channels
Syntax	<code>out = ischannel (obj . Channel (i ndex))</code>
Arguments	<p><code>obj . Channel (i ndex)</code> One or more channels contained by <code>obj .</code></p> <p><code>out</code> A logical value.</p>
Description	<code>out = ischannel (obj . Channel (i ndex))</code> returns a logical 1 to <code>out</code> if <code>obj . Channel (i ndex)</code> is a channel. Otherwise, a logical 0 is returned.
Remarks	<p><code>ischannel</code> does not determine if channels are valid (associated with hardware). To check for valid channels, use the <code>isvalid</code> function.</p> <p>Typically, you use <code>ischannel</code> directly only when you are creating your own M-files.</p>
Example	<p>Suppose you create the function <code>myfunc</code> for use with the Data Acquisition Toolbox. If <code>myfunc</code> is passed one or more channels as an input argument, then the first thing you should do in the function is check if the argument is a channel.</p> <pre>function myfunc(chan) % Determine if a channel was passed. if ~ischannel (chan) error('The argument passed is not a channel.');</pre> <p><code>end</code></p> <p>You can examine the Data Acquisition Toolbox M-files for examples that use <code>ischannel</code>.</p>
See Also	<p>Functions</p> <p><code>isvalid</code></p>

isdioline

Purpose	Check for lines	
Syntax	<code>out = isdioline(obj.Line(index))</code>	
Arguments	<code>obj.Line(index)</code>	One or more lines contained by <code>obj</code> .
	<code>out</code>	A logical value.
Description	<code>out = isdioline(obj.Line(index))</code> returns a logical 1 to <code>out</code> if <code>obj.Line(index)</code> is a line. Otherwise, a logical 0 is returned.	
Remarks	<code>isdioline</code> does not determine if lines are valid (associated with hardware). To check for valid lines, use the <code>isvalid</code> function.	
	Typically, you use <code>isdioline</code> directly only when you are creating your own M-files.	
Example	Suppose you create the function <code>myfunc</code> for use with the Data Acquisition Toolbox. If <code>myfunc</code> is passed one or more lines as an input argument, then the first thing you should do in the function is check if the argument is a line.	
	<pre>function myfunc(line) % Determine if a line was passed. if ~isdioline(line) error('The argument passed is not a line.');</pre>	
See Also	You can examine the Data Acquisition Toolbox M-files for examples that use <code>isdioline</code> .	
	Functions	
	<code>isvalid</code>	

Purpose	Determine if device objects, channels, or lines are valid								
Syntax	<pre> out = isvalid(obj) out = isvalid(obj.Channel(index)) out = isvalid(obj.Line(index)) </pre>								
Arguments	<table> <tr> <td>obj</td><td>A device object or array of device objects.</td></tr> <tr> <td>obj.Channel(index)</td><td>One or more channels contained by obj.</td></tr> <tr> <td>obj.Line(index)</td><td>One or more lines contained by obj.</td></tr> <tr> <td>out</td><td>A logical array.</td></tr> </table>	obj	A device object or array of device objects.	obj.Channel(index)	One or more channels contained by obj.	obj.Line(index)	One or more lines contained by obj.	out	A logical array.
obj	A device object or array of device objects.								
obj.Channel(index)	One or more channels contained by obj.								
obj.Line(index)	One or more lines contained by obj.								
out	A logical array.								
Description	<p>out = isvalid(obj) returns a logical 1 to out if obj is a valid device object. Otherwise, a logical 0 is returned.</p> <p>out = isvalid(obj.Channel(index)) returns a logical 1 to out if the channels specified by obj.Channel(index) are valid. Otherwise, a logical 0 is returned.</p> <p>out = isvalid(obj.Line(index)) returns a logical 1 to out if the lines specified by obj.Line(index) are valid. Otherwise, a logical 0 is returned.</p>								
Remarks	<p>Invalid device objects, channels, and lines are no longer associated with any hardware and should be cleared from the workspace with the clear function.</p> <p>Typically, you use isvalid directly only when you are creating your own M-files.</p>								
Example	<p>Create the analog input object ai for a National Instruments board and add eight channels to it.</p> <pre> ai = analoginput('ni daq', 1); ch = addchannel(ai, 0:7); </pre> <p>To verify the device object is valid</p> <pre> isvalid(ai) ans = 1 </pre>								

To verify the channels are valid

```
i s v a l i d ( c h ) '  
ans =  
      1      1      1      1      1      1      1      1
```

If you delete a channel, then `i s v a l i d` returns a logical 0 in the appropriate location

```
d e l e t e ( a i . C h a n n e l ( 3 ) )  
i s v a l i d ( c h ) '  
ans =  
      1      1      0      1      1      1      1      1
```

Typically, you will use `i s v a l i d` directly only when you are creating your own M-files. Suppose you create the function `myfunc` for use with the Data Acquisition Toolbox. If `myfunc` is passed the previously defined device object `ai` as an input argument

```
myfunc(ai)
```

the first thing you should do in the function is check if `ai` is a valid device object.

```
function myfunc(obj)  
% Determine if an invalid handle was passed.  
if ~isvalid(obj)  
    error('Invalid data acquisition object passed. ');  
end
```

You can examine the Data Acquisition Toolbox M-files for examples that use `i s v a l i d`.

See Also

Functions

`clear`, `delete`, `ischannel`, `isdioline`

Purpose Return the length of a device object, channel group, or line group

Syntax

```
out = length(obj)
out = length(obj.Channel)
out = length(obj.Line)
```

Arguments

obj	A device object or array of device objects.
obj.Channel	The channels contained by obj.
obj.Line	The lines contained by obj.
out	A double.

Description

out = length(obj) returns the length of the device object obj to out.

out = length(obj.Channel) returns the length of the channel group contained by obj.

out = length(obj.Line) returns the length of the line group contained by obj.

Example

Create the analog input object ai for a National Instruments board and add eight channels to it.

```
ai = analoginput('ni daq', 1);
ai ch = addchannel(ai, 0:7);
```

Create the analog output object ao for a National Instruments board, add one channel to it, and create the device object array ai ao.

```
ao = analogoutput('ni daq', 1);
ao ch = addchannel(ao, 0);
ai ao = [ai ao]
```

Index:	Subsystem:	Name:
1	Analog Input	ni daq1- AI
2	Analog Output	ni daq1- AO

length

To find the length of ai ao

```
length(ai ao)
ans =
    2
```

To find the length of the analog input channel group

```
length(ai ch)
ans =
    8
```

See Also

Functions

size

Purpose	Load device objects, channels, or lines into the MATLAB workspace
Syntax	<pre>load file load file obj 1 obj 2... out = load('file','obj 1','obj 2',...)</pre>
Arguments	<p><code>file</code> The MAT-file name.</p> <p><code>obj 1 obj 2...</code> Device objects, an array of device objects, channels, or lines.</p> <p><code>out</code> A structure containing the loaded device objects.</p>
Description	<p><code>load file</code> returns all variables from the MAT-file <code>file</code> into the MATLAB workspace.</p> <p><code>load file obj 1 obj 2...</code> returns the specified device objects from the MAT-file <code>file</code> into the MATLAB workspace.</p> <p><code>out = load('file','obj 1','obj 2',...)</code> returns the specified device objects from the MAT-file <code>file</code> as a structure to <code>out</code> instead of directly loading them into the workspace. The field names in <code>out</code> match the names of the loaded device objects. If no device objects are specified, then all variables existing in the MAT-file are loaded.</p>
Remarks	<p>Loading device objects follows these rules:</p> <ul style="list-style-type: none">• Unique device objects are loaded into the MATLAB workspace as well as the engine.• If a loaded device object already exists in the engine but not the MATLAB workspace, the loaded device object automatically reconnects to the engine device object.• If a loaded device object already exists in the workspace or the engine but has different properties than the loaded object, then these rules are followed:<ul style="list-style-type: none">- The read-only properties are automatically reset to their default values.- All other property values are given by the loaded object and a warning is issued stating that property values of the workspace object have been updated.

load

- If the workspace device object is running, then it is stopped before loading occurs.
- If identical device objects are loaded, then they point to the same device object in the engine. For example, if you saved the array
`x = [ai 1 ai 1 ai 2]`
only `ai 1` and `ai 2` are created in the engine, and `x(1)` will equal `x(2)`.
- Values for read-only properties are restored to their default values upon loading. For example, the `EventLog` property is restored to an empty vector. Use the `propinfo` function to determine if a property is read only.
- Values for the `BufferingConfig` property when the `BufferingMode` property is set to `Auto`, and the `MaxSamplesQueued` property may not be restored to the same value since both these property values are based on available memory.

Note `load` is not used to read in acquired data that has been saved to a log file. You should use the `daqread` function for this purpose.

If you use the `help` command to display the M-file help for `load`, then you must supply the pathname shown below.

```
help daq/private/load
```

Example

This example illustrates the behavior of `load` when the loaded device object has properties that differ from the workspace object.

```
ai = analoginput('winsound');  
addchannel(ai, 1:2);  
save ai  
ai.SampleRate = 10000;  
load ai  
Warning: Loaded object has updated property values.
```

See Also

Functions

`daqread`, `propinfo`, `save`

Purpose	Generate a list of descriptive channel or line names	
Syntax	<code>names = makenames(' prefix', index)</code>	
Arguments	<code>' prefix'</code>	A string that constitutes the first part of the name.
	<code>index</code>	Numbers appended to the end of <code>prefix</code> – any MATLAB vector syntax can be used to specify <code>index</code> as long as the numbers are positive.
	<code>names</code>	An m-by-1 cell array of channel names where m is the length of <code>index</code> .
Description	<code>names = makenames(' prefix', index)</code> generates a cell array of descriptive channel or line names by concatenating <code>prefix</code> and <code>index</code> .	
Remarks	You can pass <code>names</code> as an input argument to the <code>addchannel</code> or <code>addline</code> function.	
	If <code>names</code> contains more than one descriptive name, then the size of <code>names</code> must agree with the number of hardware channels specified in <code>addchannel</code> , or the number of hardware lines specified in <code>addline</code> .	
	If the channels or lines are to be referenced by name, then <code>prefix</code> must begin with a letter and contain only letters, numbers, and underscores. Otherwise the names can contain any character.	
Example	<p>Create the analog input object <code>AI</code>. You can use <code>makenames</code> to define descriptive names for each channel that is to be added to <code>AI</code>.</p> <pre>AI = analoginput('ni daq', 1); names = makenames('chan', 1:8);</pre> <p><code>names</code> is an eight-element cell array of channel names <code>chan1</code>, <code>chan2</code>, ..., <code>chan8</code>. You can now pass <code>names</code> as an input argument to the <code>addchannel</code> function.</p> <pre>addchannel(AI, 0:7, names);</pre>	
See Also	Functions <code>addchannel</code> , <code>addline</code>	

muxchanidx

Purpose	Return multiplexed scanned channel index	
Syntax	<code>scani dx = muxchani dx(obj , muxboard, muxi dx)</code> <code>scani dx = muxchani dx(obj , absmuxi dx)</code>	
Arguments	obj	An analog input object associated with a National Instruments board.
	muxboard	The multiplexer board.
	muxi dx	The index number of the multiplexed channel.
	absmuxi dx	The absolute index number of the multiplexed channel.
	scani dx	The scanning index number of the multiplexed channel.
Description	<code>scani dx = muxchani dx(obj , muxboard, muxi dx)</code> returns the scanning index number of the multiplexed channel specified by <code>muxi dx</code> . The multiplexer (mux) board is specified by <code>muxboard</code> . For each mux board, <code>muxi dx</code> can range from 0-31 for differential inputs and 0-63 for single-ended inputs. <code>muxboard</code> and <code>muxi dx</code> are vectors of equal length.	
	<code>scani dx = muxchani dx(obj , absmuxi dx)</code> returns the scanning index number of the multiplexed channel specified by <code>absmuxi dx</code> . <code>absmuxi dx</code> is the absolute index of the channel independent of the mux board.	
	For single-ended inputs, the first mux board has absolute index values that range between 0 and 63, the second mux board has absolute index values that range between 64 and 127, the third mux board has absolute index values that range between 128 and 191, the fourth mux board has absolute index values that range between 192 and 255. For example, the absolute index value of the second single-ended channel on the fourth mux board (<code>muxboard</code> is 4 and <code>muxi dx</code> is 1) is 193.	
Remarks	<code>scani dx</code> identifies the column number of the data returned by <code>getdata</code> and <code>peekdata</code> . Refer to the <i>AMUX-64T User Manual</i> for more information about adding mux channels based on hardware channel IDs and the number of mux boards used.	

Example

Create the analog input object `ai` for a National Instruments board that is connected to four AMUX-64T multiplexers, and add 256 channels to `ai` using `addmuxchannel`.

```
ai = analoginput('ni daq', 1);  
ai.InputType = 'SingleEnded';  
ai.NumMuxBoards = 4;  
addmuxchannel(ai);
```

The following two commands return a scanned index value of 14.

```
scandx = muxchandise(ai, 4, 1);  
scandx = muxchandise(ai, 193);
```

See Also**Functions**

`addmuxchannel`

obj2mfile

Purpose Convert device objects, channels, or lines to MATLAB code

Syntax

```
obj2mfile(obj, 'file')  
obj2mfile(obj, 'file', 'syntax')  
obj2mfile(obj, 'file', 'all')  
obj2mfile(obj, 'file', 'syntax', 'all')
```

Arguments

<code>obj</code>	A device object, array of device objects, channels, or lines.
<code>'file'</code>	The file that the MATLAB code is written to. The full pathname can be specified. If an extension is not specified, the .m extension is used.
<code>'syntax'</code>	Syntax of the converted MATLAB code. By default, the set syntax is used. If dot is specified, then the subscripted referencing syntax is used. If named is specified, then named referencing is used (if defined).
<code>'all'</code>	If all is specified, all properties are written to <code>file</code> . If all is not specified, only properties that are not set to their default values are written to <code>file</code> .

Description `obj2mfile(obj, 'file')` converts `obj` to the equivalent MATLAB code using the set syntax and saves the code to `file`. By default, only those properties that are not set to their default values are written to `file`.

`obj2mfile(obj, 'file', 'syntax')` converts `obj` to the equivalent MATLAB code using *syntax* and saves the code to `file`. The values for *syntax* can be *set*, *dot*, or *named*. *set* uses the set syntax, *dot* uses subscripted assignment (dot notation), and *named* uses named referencing (if defined).

`obj2mfile(obj, 'file', 'all')` converts `obj` to the equivalent MATLAB code using the set syntax and saves the code to `file`. **all** specifies that all properties are written to `file`.

`obj2mfile(obj, 'file', 'syntax', 'all')` converts `obj` including all of `obj`'s properties to the equivalent MATLAB code using *syntax* and saves the code to `file`.

Remarks

If UserData is not empty or if any of the action properties are set to a cell array of values, then the data stored in those properties is written to a MAT-file when the object is converted and saved. The MAT-file has the same name as the M-file containing the object code (see the example below).

You can recreate the saved device objects by typing the name of the M-file at the command line. You can also recreate channels or lines, by typing the name of the M-file with a device object as the only input.

Example

Create the analog input object `ai` for a sound card, add two channels, and set values for several properties.

```
ai = analoginput('winsound');
addchannel(ai, 1:2);
set(ai, 'Tag', 'myai', 'TriggerRepeat', 4)
set(ai, 'StartAction', {'myaction', 2, magic(10)})
```

The following command writes MATLAB code to the files `myai.m` and `myai.mat`.

```
obj2mfile(ai, 'myai.m', 'dot')
```

`myai.m` contains code that recreates the analog input code shown above using the dot notation for all properties that have their default values changed. Since `StartAction` is set to a cell array of values, this property appears in `myai.m` as

```
ai.StartAction = startaction;
```

and is saved in `myai.mat` as

```
startaction = {'myaction', 2, magic(10)};
```

To recreate `ai` and assign the device object to a new variable `ai_new`

```
ai_new = myai;
```

The associated MAT-file, `myai.mat`, is automatically loaded.

peekdata

Purpose Preview most recent acquired data

Syntax `data = peekdata(obj, samples)`

Arguments	<code>obj</code>	An analog input object.
	<code>samples</code>	The number of samples to preview for each channel contained by <code>obj</code> .
	<code>data</code>	An m-by-n matrix where m is the number of samples and n is the number of channels.

Description `data = peekdata(obj, samples)` returns the latest number of samples specified by `samples` to `data`.

Remarks **More About Using peekdata**

- Unlike `getdata`, `peekdata` is a *nonblocking* function that immediately returns control to MATLAB. Since `peekdata` does not block execution control, data may be missed or repeated.
- `peekdata` takes a “snapshot” of the most recent acquired data and does not remove samples from the data acquisition engine. Therefore, the `SamplesAvailable` property value is not affected when `peekdata` is called.

Rules for Using peekdata

- You can call `peekdata` before a trigger executes. Therefore, `peekdata` is useful for previewing data before it is logged to the engine or to a disk file.
- In most cases, you will call `peekdata` while the device object is running. However, you can call `peekdata` once after the device object stops running.
- If `samples` is greater than the number of samples currently acquired, all available samples are returned with a warning message stating that the requested number of samples were not available.

Example

Create the analog input object `ai` for a National Instruments board, add eight input channels, and configure `ai` for a two second acquisition.

```
ai = analoginput('ni daq', 1);  
addchannel(ai, 0:7);  
set(ai, 'SampleRate', 2000)  
set(ai, 'SamplesPerTrigger', 4000)
```

After issuing the `start` function, you can preview the data.

```
start(ai)  
data = peekdata(ai, 100);
```

`peekdata` returns 100 samples to `data` for all eight channel group members. If 100 samples are not available, then whatever samples are available will be returned and a warning message is issued. The data is not removed from the data acquisition engine.

See Also**Functions**

`getdata`, `getsample`

Properties

`SamplesAvailable`

propinfo

Purpose

Return property characteristics for device objects, channels, or lines

Syntax

```
out = propinfo(obj)
out = propinfo(obj, 'PropertyName')
```

Arguments

obj

A device object, channels, or lines.

'PropertyName'

A valid obj property name.

out

A structure whose field names are the property names for obj (if *PropertyName* is not specified).

Description

out = propinfo(obj) returns the structure out whose field names are the property names for obj. Each property name in out contains the fields shown below.

Field Name	Description
Type	The property data type. Possible values are action function, any, double, and string.
Constraint	The type of constraint on the property value. Possible values are action function, bounded, enum, and none.
ConstraintValue	The property value constraint. The constraint can be a range of valid values or a list of valid string values.
DefaultValue	The property default value.
ReadOnly	If the property is read-only, a 1 is returned. Otherwise, a 0 is returned.
ReadOnlyRunning	If the property is read-only while the device object is running, a 1 is returned. Otherwise, a 0 is returned.
DeviceSpecific	If the property is device-specific, a 1 is returned. If a 0 is returned, the property is supported for all device objects of a given type.

`out = propinfo(obj, 'PropertyName')` returns the structure `out` for the property specified by *PropertyName*. If *PropertyName* is a cell array of strings, a cell array of structures is returned for each property.

Example

Create the analog input object `ai` for a sound card and configure it to operate in stereo mode.

```
ai = analoginput('winsound');
addchannel(ai, 1:2);
```

To capture all property information for all common `ai` properties

```
out = propinfo(ai);
```

To display the default value for the `SampleRate` property

```
out.SampleRate.DefaultValue
ans =
    8000
```

To display all the property information for the `InputRange` property

```
propinfo(ai.Channel, 'InputRange')
ans =
    Type: 'double'
    Constraint: 'Bounded'
    ConstraintValue: [-1 1]
    DefaultValue: [-1 1]
    ReadOnly: 0
    ReadOnlyRunning: 1
    DeviceSpecific: 0
```

See Also

Functions

`daqhelp`

putdata

Purpose	Queue data in the engine for eventual output	
Syntax	<code>putdata(obj , data)</code>	
Arguments	<code>obj</code>	An analog output object.
	<code>data</code>	The data to be queued in the engine.
Description	<code>putdata(obj , data)</code> queues the data specified by <code>data</code> in the engine for eventual output to the analog output subsystem. <code>data</code> must consist of a column of data for each channel contained by <code>obj</code> .	
Remarks	More About Queueing Data	
	<ul style="list-style-type: none">• Data must be queued in the engine before <code>obj</code> is executed.• <code>putdata</code> is a <i>blocking</i> function since it returns execution control to the MATLAB workspace only when the requested number of samples are queued in the engine for each channel group member.• If the value of the <code>RepeatOutput</code> property is greater than 0, then all queued data is automatically requeued until the <code>RepeatOutput</code> value is reached. <code>RepeatOutput</code> must be configured before <code>start</code> is issued.• After <code>obj</code> executes, you can continue to queue data unless <code>RepeatOutput</code> is greater than 0.• You can queue data in the engine until the value specified by the <code>MaxSamplesQueued</code> property is reached, or the limitations of your hardware or computer are reached.	
	More About Outputting Data	
	<ul style="list-style-type: none">• Data is output as soon as a trigger occurs.• An error is returned if a NaN is included in the data stream.• You can specify <code>data</code> as the native data type of the hardware. Note that MATLAB supports math operations only for the double data type. Therefore, to use math functions on native data, you must convert it to doubles.• If the output data is not within the range specified by the <code>OutputRange</code> property, then the data is clipped.	

- The `SamplesOutput` property keeps a running count of the total number of samples that have been output per channel.
- The `SamplesAvailable` property tells you how many samples are ready to be output from the engine per channel. After data is output, `SamplesAvailable` is automatically reduced by the number of samples sent to the hardware.

Example

Create the analog output object `ao` for a National Instruments board, add two output channels to it, and generate 10 seconds of data to be output.

```
ao = analogoutput('ni daq', 1);
ch = addchannel(ao, 0:1);
set(ao, 'SampleRate', 1000)
data = linspace(0, 1, 10000)';
```

Before you can output data, it must be queued in the engine using `putdata`.

```
putdata(ao, [data data])
start(ao)
```

See Also

Functions

`putsample`

Properties

`MaxSamplesQueued`, `OutputRange`, `RepeatOutput`, `SamplesAvailable`, `SamplesOutput`, `Timeout`, `UnitsRange`

putsample

Purpose	Immediately output one sample	
Syntax	<code>putsample(obj, data)</code>	
Arguments	<code>obj</code>	An analog output object.
	<code>data</code>	The data to be queued in the engine.
Description	<code>putsample(obj, data)</code> immediately outputs the row vector <code>data</code> , which consists of one sample for each channel contained by <code>obj</code> .	
Remarks	Using <code>putsample</code> is a good way to test your analog output configuration. Additionally:	
	<ul style="list-style-type: none">• <code>putsample</code> does not store samples in the data acquisition engine.• <code>putsample</code> can be executed at any time after channels have been added to <code>obj</code>.• <code>putsample</code> is not supported for sound cards.	
Example	Create the analog output object <code>ao</code> for a National Instruments board and add two hardware channels to it.	
	<pre>ao = analogoutput('ni daq', 1); ch = addchannel(ao, 0:1);</pre> <p>To call <code>putsample</code> for <code>ao</code></p> <pre>putsample(ao, [1 1])</pre>	
See Also	Functions	
	<code>putdata</code>	

Purpose	Write values to lines						
Syntax	<pre>putvalue(obj, data) putvalue(obj.Line(index), data)</pre>						
Arguments	<table> <tr> <td><code>obj</code></td><td>A digital I/O object.</td></tr> <tr> <td><code>obj.Line(index)</code></td><td>One or more lines contained by <code>obj</code>.</td></tr> <tr> <td><code>data</code></td><td>A decimal value or binary vector.</td></tr> </table>	<code>obj</code>	A digital I/O object.	<code>obj.Line(index)</code>	One or more lines contained by <code>obj</code> .	<code>data</code>	A decimal value or binary vector.
<code>obj</code>	A digital I/O object.						
<code>obj.Line(index)</code>	One or more lines contained by <code>obj</code> .						
<code>data</code>	A decimal value or binary vector.						
Description	<p><code>putvalue(obj, data)</code> writes <code>data</code> to the hardware lines contained by the digital I/O object <code>obj</code>.</p> <p><code>putvalue(obj.Line(index), data)</code> writes <code>data</code> to the hardware lines specified by <code>obj.Line(index)</code>.</p>						
Remarks	<p>More About Writing Values to Lines</p> <ul style="list-style-type: none"> • You can specify <code>data</code> as either a decimal value or a binary vector. A binary vector (or <i>binvec</i>) is constructed with the least significant bit (LSB) in the first column and the most significant bit (MSB) in the last column. For example, the decimal number 23 is written as the binary vector [1 1 1 0 1]. • If <code>obj</code> contains lines from a port-configurable device, then all lines will be written to even if they are not contained by the device object. • An error will be returned if <code>data</code> is written to an input line. • An error is returned if you attempt to write a negative value. • If a decimal value is written to a digital I/O object and the value is too large to be represented by the hardware, then an error is returned. 						
Example	<p>Create the digital I/O object <code>di o</code> and add four output lines to it.</p> <pre>di o = digitalio('ni daq', 1); lines = addline(di o, 0:3, 'out');</pre> <p>Write the value 8 as a decimal value and as a binary vector</p> <pre>putvalue(di o, 8) putvalue(di o, [0 0 0 1])</pre>						

save

Purpose	Save device objects to a MAT-file	
Syntax	<code>save file</code> <code>save file obj1 obj2...</code>	
Arguments	<code>file</code>	The MAT-file name.
	<code>obj1 obj2...</code>	One or more device objects or an array of device objects.
Description	<code>save file</code> saves all MATLAB variables to the MAT-file <code>file</code> . If an extension is not specified for <code>file</code> , then a <code>.MAT</code> extension is used. <code>save file obj1 obj2...</code> saves the specified device objects to <code>file</code> .	
Remarks	<p>Saving device objects follows these rules:</p> <ul style="list-style-type: none">• You can use <code>save</code> in the functional form as well as the command form shown above. When using the functional form, you must specify the filename and device objects as strings.• Samples associated with a device object are not stored in the MAT-file. You can bring these samples into the MATLAB workspace with the <code>getdata</code> function, and then save them to the MAT-file using a separate variable name. You can also log samples to disk by configuring the <code>LoggingMode</code> property to <code>Disk</code> or <code>Disk&Memory</code>.• Values for read-only properties are restored to their default values upon loading. For example, the <code>EventLog</code> property is restored to an empty vector. Use the <code>propinfo</code> function to determine if a property is read only.• Values for the <code>BufferingConfig</code> property when the <code>BufferingMode</code> property is set to <code>Auto</code>, and the <code>MaxSamplesQueued</code> property may not be restored to the same value since both these property values are based on available memory. <p>If you use the <code>help</code> command to display the M-file help for <code>save</code>, then you must supply the pathname shown below.</p> <pre>help daq/private/save</pre>	
See Also	Functions <code>getdata</code> , <code>load</code> , <code>propinfo</code>	

Purpose	Configure or display device object properties	
Syntax	<pre> set(obj) props = set(obj) set(obj, 'PropertyName') props = set(obj, 'PropertyName') set(obj, 'PropertyName', PropertyValue, ...) set(obj, PN, PV) set(obj, S) </pre>	
Arguments	obj	A device object, array of device objects, channels, or lines.
	'PropertyName'	A property name.
	PropertyValue	A property value.
	PN	A cell array of property names.
	PV	A cell array of property values.
	S	A structure whose field names are device object, channel, or line properties.
	props	A structure array whose field names are the property names for obj, or a cell array of possible values.
Description	<p>set(obj) displays all configurable properties for obj. If a property has a finite list of possible string values, then these values are also displayed.</p> <p>props = set(obj) returns all configurable properties to props. props is a structure array with fields given by the property names, and possible property values contained in cell arrays. if the property does not have a finite set of possible values, then the cell array is empty.</p> <p>set(obj, 'PropertyName') displays the valid values for the property specified by <i>PropertyName</i>. <i>PropertyName</i> must have a finite set of possible values.</p> <p>props = set(obj, 'PropertyName') returns the valid values for <i>PropertyName</i> to props. props is a cell array of possible values or an empty cell array if the property does not have a finite set of possible values.</p>	

`set(obj, 'PropertyName', PropertyValue, ...)` sets multiple property values with a single statement. Note that you can use structures, property name/property value string pairs, and property name/property value cell array pairs in the same call to `set`.

`set(obj, PN, PV)` sets the properties specified in the cell array of strings `PN` to the corresponding values in the cell array `PV`. `PN` must be a vector. `PV` can be `m-by-n` where `m` is equal to the specified number of device objects, channels, or lines and `n` is equal to the length of `PN`.

`set(obj, S)` where `S` is a structure whose field names are device object properties, sets the properties named in each field name with the values contained in the structure.

Remarks

If you use the `help` command to display the M-file help for `set`, then you must supply the pathname shown below.

```
help daq/daqdevice/set
```

Example

Create the analog input object `ai` for a sound card and configure it to operate in stereo mode.

```
ai = analoginput('winsound');  
addchannel(ai, 1:2);
```

To display all of `ai`'s configurable properties and their valid values

```
set(ai)
```

To set the value for the `SampleRate` property to 10000

```
set(ai, 'SampleRate', 10000)
```

The following two commands set the value for the `SampleRate` and `InputType` properties using one call to `set`.

```
set(ai, 'SampleRate', 10000, 'TriggerType', 'Manual')  
set(ai, {'SampleRate', 'TriggerType'}, {10000, 'Manual'})
```

You can also set different channel property values for multiple channels.

```
ch = ai.Channel(1:2);  
set(ch, {'UnitsRange', 'Channel Name'}, {[ -1 1] 'Name1'; [-2 2]  
      'Name2'})
```

See Also

Functions

get, setverify

setverify

Purpose	Configure and return the specified property	
Syntax	<pre>Actual = setverify(obj, 'PropertyName', PropertyValue) Actual = setverify(obj.Channel(index), 'PropertyName', PropertyValue) Actual = setverify(obj.Line(index), 'PropertyName', PropertyValue)</pre>	
Arguments	obj	A device object or array of device objects.
	'PropertyName'	A property name.
	PropertyValue	A property value.
	obj.Channel(index)	One or more channels contained by obj.
	obj.Line(index)	One or more lines contained by obj.
	Actual	The actual value for the specified property.
Description	<p>Actual = setverify(obj, 'PropertyName', PropertyValue) sets <i>PropertyName</i> to PropertyValue for obj, and returns the actual property value to Actual.</p> <p>Actual = setverify(obj.Channel(index), 'PropertyName', PropertyValue) sets <i>PropertyName</i> to PropertyValue for the channels specified by index, and returns the actual property value to Actual.</p> <p>Actual = setverify(obj.Line(index), 'PropertyName', PropertyValue) sets <i>PropertyName</i> to PropertyValue for the lines specified by index, and returns the actual property value to Actual.</p>	
Remarks	<p>setverify is equivalent to the commands</p> <pre>set(obj, 'PropertyName', PropertyValue) Actual = get(obj, 'PropertyName')</pre> <p>Using setverify is not required for setting property values, but it does provide a convenient way to verify the actual property value set by the data acquisition engine.</p>	

setverify is particularly useful when setting the SampleRate, InputRange, and OutputRange properties since these properties can only be set to specific values accepted by the hardware. You can use the propinfo function to obtain information about the valid values for these properties.

If a property value is specified but does not match a valid value, then:

- If the specified value is within the range of supported values:
 - For the SampleRate and InputRange properties, the value is automatically rounded up to the next highest supported value.
 - For all other properties, the value is automatically selected to be the nearest supported value.
- If the value is not within the range of supported values, an error is returned and the current property value remains unchanged.

Example

Create the analog input object ai for a National Instruments AT-MIO-16DE-10 board, add eight hardware channels to it, and set the sample rate to 10,000 Hz using setverify.

```
ai = analoginput('ni daq', 1);
ch = addchannel(ai, 0:7);
ActualRate = setverify(ai, 'SampleRate', 10000);
```

Suppose you use setverify to set the input range for all channels contained by ai to -8 to 8 volts.

```
ActualInputRange = setverify(ai.Channel, 'InputRange', [-8 8]);
```

The InputRange value was actually rounded up to -10 to 10 volts.

```
ActualInputRange{1}
ans =
    -10     10
```

See Also

Functions

get, propinfo, set

Properties

InputRange, OutputRange, SampleRate

showdaqevents

Purpose Display event log information

Syntax `showdaqevents(obj)`
 `showdaqevents(obj, index)`
 `showdaqevents(struct)`
 `showdaqevents(struct, index)`
 `out = showdaqevents(...)`

Arguments

<code>obj</code>	An analog input or analog output object.
<code>index</code>	The event index.
<code>struct</code>	An event structure.
<code>out</code>	A one column cell array of event information.

Description

`showdaqevents(obj)` displays a summary of the event log for `obj`.

`showdaqevents(obj, index)` displays a summary of the events specified by `index` for `obj`.

`showdaqevents(struct)` displays a summary of the events stored in the structure `struct`.

`showdaqevents(struct, index)` displays a summary of the events specified by `index` stored in the structure `struct`.

`out = showdaqevents(...)` outputs the event information to a one column cell array `out`. Each element of `out` is a string that contains the event information associated with that index value.

Remarks

You can pass a structure of event information to `showdaqevents`. This structure can be obtained from the `getdata` function, the `daqread` function, or the `EventLog` property.

Example Create the analog input object ai for a sound card, add two channels to it, and configure ai to execute three triggers.

```
ai = analoginput('winsound');
ch = addchannel(ai, 1:2);
set(ai, 'TriggerRepeat', 2)
```

Start ai and display the trigger event information with showdaqevents.

```
start(ai)
showdaqevents(ai, 2:4)
2 Trigger#1      ( 17:07:06, 0 )      Channel: N/A
3 Trigger#2      ( 17:07:07, 8000 )   Channel: N/A
4 Trigger#3      ( 17:07:08, 16000 )  Channel: N/A
```

See Also Functions
daqread, getdata

Properties
EventLog

size

Purpose Return the size of a device object, channel group, or line group

Syntax

```
d = size(obj)
[m1, m2, m3, ..., mn] = size(obj)
m = size(obj, dim)
d = size(obj.Channel)
[m1, m2, m3, ..., mn] = size(obj.Channel)
m = size(obj.Channel, dim)
d = size(obj.Line)
[m1, m2, m3, ..., mn] = size(obj.Line)
m = size(obj.Line, dim)
```

Arguments

obj	A device object or array of device objects.
dim	The dimension.
obj.Channel	The channels contained by obj.
obj.Line	The lines contained by obj.
d	A two-element row vector containing the number of rows and columns in obj.
m1, m2, m3, . . . , mn	Each dimension of obj is captured in a separate variable.
m	The length of the dimension specified by dim.

Description

`d = size(obj)` returns the two-element row vector `d = [m, n]` containing the number of rows and columns in obj.

`[m1, m2, m3, . . . , mn] = size(obj)` returns the length of the first n dimensions of obj to separate output variables. For example, `[m, n] = size(obj)` returns the number of rows to m and the number of columns to n.

`m = size(obj, dim)` returns the length of the dimension specified by the scalar dim. For example, `size(obj, 1)` returns the number of rows.

`d = size(obj.Channel)` returns the two-element row vector `d = [m, n]` containing the number of rows and columns in the channel group `obj.Channel`.

`[m1, m2, m3, . . . , mn] = size(obj . Channel)` returns the length of the first `n` dimensions of the channel group `obj . Channel` to separate output variables. For example, `[m, n] = size(obj . Channel)` returns the number of rows to `m` and the number of columns to `n`.

`m = size(obj . Channel , di m)` returns the length of the dimension specified by the scalar `di m`. For example, `size(obj . Channel , 1)` returns the number of rows.

`d = size(obj . Li ne)` returns the two-element row vector `d = [m, n]` containing the number of rows and columns in the line group `obj . Li ne`.

`[m1, m2, m3, . . . , mn] = size(obj . Li ne)` returns the length of the first `n` dimensions of the line group `obj . Li ne` to separate output variables. For example, `[m, n] = size(obj . Li ne)` returns the number of rows to `m` and the number of columns to `n`.

`m = size(obj . Li ne, di m)` returns the length of the dimension specified by the scalar `di m`. For example, `size(obj . Li ne, 1)` returns the number of rows.

Example

Create the analog input object `ai` for a National Instruments board and add eight channels to it.

```
ai = analoginput('ni daq', 1);
ch = addchannel(ai, 0:7);
```

To find the size of the device object

```
size(ai)
ans =
     1     1
```

To find the size of the channel group

```
size(ch)
ans =
     8     1
```

See Also

Functions

`length`

start

Purpose	Start a device object
Syntax	start(obj)
Arguments	obj A device object or an array of device objects.
Description	start(obj) initiates the execution of the device object obj.
Remarks	<p>When start is issued for an analog input or analog output object:</p> <ul style="list-style-type: none">• The M-file action function specified for StartAction is executed.• The Running property is set to On.• The start event is recorded in the EventLog property.• Data existing in the engine is flushed. <p>Although an analog input or analog output object may be executing, data logging or sending is not necessarily initiated. Data logging or sending requires a trigger event to occur, and depends on the TriggerType property value.</p> <p>For any device object, you can specify start as the value for an action property.</p> <pre>ai.StopAction = 'start';</pre>
<hr/>	
Note You typically execute a digital I/O object to periodically update and display its state. Refer to the diopanel demo for an example of this.	
<hr/>	
See Also	<p>Functions</p> <p>stop, trigger</p> <p>Properties</p> <p>EventLog, Running, Sending, TriggerType</p>

Purpose	Stop a device object
Syntax	<code>stop(obj)</code>
Arguments	<p><code>obj</code> A device object or an array of device objects.</p>
Description	<code>stop(obj)</code> terminates the execution of the device object <code>obj</code> .
Remarks	<p>An analog input object automatically stops when the requested samples are acquired or data is missed. An analog output object automatically stops when the queued data is output. These two device objects can also stop executing under one of these conditions:</p> <ul style="list-style-type: none"> • The <code>Timeout</code> property value is reached. • A run-time error occurs. <p>For analog input objects, <code>stop</code> must be used when the <code>TriggerRepeat</code> property or <code>SamplesPerTrigger</code> property is set to <code>inf</code>. For analog output objects, <code>stop</code> must be used when the <code>RepeatOutput</code> property is set to <code>inf</code>. When <code>stop</code> is issued for either of these device objects:</p> <ul style="list-style-type: none"> • The <code>Running</code> property is set to <code>Off</code>. • The <code>Logging</code> property or <code>Sending</code> property is set to <code>Off</code>. • The M-file action function specified for <code>StopAction</code> is executed. • The stop event is recorded in the <code>EventLog</code> property. <p>For any device object, you can specify <code>stop</code> as the value for an action property.</p> <pre>ao.TimerAction = 'stop';</pre> <hr/> <p>Note Issuing <code>stop</code> is the only way to stop an executing digital I/O object. You typically execute a digital I/O object to periodically update and display its state. Refer to the <code>diopanel</code> demo for an example.</p> <hr/>

stop

See Also

Functions

start, trigger

Properties

EventLog, Loggi ng, RepeatOut put, Runni ng, Sampl esPerTri gger, Sendi ng, Timeout, Tri ggerRepeat

Purpose	Manually execute a trigger	
Syntax	<code>trigger(obj)</code>	
Arguments	<code>obj</code>	An analog input or analog output object or an array of these device objects.
Description	<code>trigger(obj)</code> manually executes a trigger.	
Remarks	<p>After <code>trigger</code> is issued:</p> <ul style="list-style-type: none">• The absolute time of the trigger event is recorded by the <code>InitialTriggerTime</code> property.• The <code>Logging</code> property or <code>Sending</code> property is set to <code>On</code>.• The M-file action function specified by <code>TriggerAction</code> is executed.• The trigger event is recorded in the <code>EventLog</code> property. <p>You can issue <code>trigger</code> only if <code>TriggerType</code> is set to <code>Manual</code>, <code>Running</code> is <code>On</code>, and <code>Logging</code> is <code>Off</code>.</p> <p>You can specify <code>trigger</code> as the value for an action property.</p> <pre>ai.StartAction = 'trigger';</pre>	
See Also	Functions <code>start</code> , <code>stop</code> Properties <code>InitialTriggerTime</code> , <code>Logging</code> , <code>Running</code> , <code>Sending</code> , <code>TriggerAction</code> , <code>TriggerType</code>	

waittilstop

Purpose	Wait for the device object to stop running				
Syntax	<code>waittilstop(obj, waittime)</code>				
Arguments	<table><tr><td><code>obj</code></td><td>A device object or an array of device objects.</td></tr><tr><td><code>waittime</code></td><td>The maximum time to wait for <code>obj</code> to stop running.</td></tr></table>	<code>obj</code>	A device object or an array of device objects.	<code>waittime</code>	The maximum time to wait for <code>obj</code> to stop running.
<code>obj</code>	A device object or an array of device objects.				
<code>waittime</code>	The maximum time to wait for <code>obj</code> to stop running.				
Description	<p><code>waittilstop(obj, waittime)</code> blocks the MATLAB command line, and waits for <code>obj</code> to stop running. You specify the maximum waiting time, in seconds, with <code>waittime</code>. <code>waittime</code> overrides the value specified for the <code>Timeout</code> property. If <code>obj</code> is an array of device objects, then <code>waittilstop</code> may wait up to the specified time for each device object in the array.</p> <p><code>waittilstop</code> is particularly useful if you want to guarantee that the specified data is acquired before another task is performed.</p>				
Remarks	<p>If <code>obj</code> is not running when <code>waittilstop</code> is issued, or if an error occurs while <code>obj</code> is running, then <code>waittilstop</code> immediately relinquishes control of the command line.</p> <p>When <code>obj</code> stops running, its <code>Running</code> property is automatically set to <code>Off</code>. <code>obj</code> can stop running under one of the these conditions:</p> <ul style="list-style-type: none">• The requested number of samples is acquired (analog input) or sent out (analog output).• The stop function is issued.• A run-time error occurs.• The <code>Timeout</code> property value is reached (<code>waittime</code> supersedes this value). <p>It is not guaranteed that the <code>StopAction</code> property is called before <code>waittilstop</code> returns. The stop action event is recorded by the <code>EventLog</code> property.</p>				

Example

Create the analog input object `ai` for a National Instruments board, add eight channels to it, and configure a 25 second acquisition.

```
ai = analoginput('ni daq', 1);  
ch = addchannel(ai, 0:7);  
ai.SampleRate = 2000;  
ai.TriggerRepeat = 4;  
ai.SamplesPerTrigger = 10000;
```

You can use `waittilstop` to block the MATLAB command line until all the requested data is acquired. Since the expected acquisition time is 25 seconds, the `waittime` argument is 26. If the acquisition does not complete within this time, then a timeout occurs.

```
start(ai)  
waittilstop(ai, 26)
```

See Also

Properties

EventLog, Running, StopAction, Timeout

Base Property Reference

Overview	10-2
Getting Command Line Property Help	10-2
Properties Grouped by Category	10-3
Analog Input Properties	10-3
Analog Output Properties	10-8
Digital I/O Properties	10-12
Properties Listed Alphabetically	10-14

Overview

This section provides descriptions of all toolbox base properties. Base properties apply to all supported hardware subsystems of a given type (analog input, analog output, etc.). For example, the `SampleRate` property is supported for all analog input subsystems regardless of the vendor.

In “Properties Grouped by Category” on page 10-3, the properties are summarized according to subsystem type and usage. In “Properties Listed Alphabetically” on page 10-14, the properties are described in detail.

Getting Command Line Property Help

To get command line property help, you should use the `daqhelp` function. For example, to get help for the `SampleRate` property

```
daqhelp SampleRate
```

Note You can use `daqhelp` without creating a device object.

You can also get property characteristics, such as the default property value, using the `propinfo` function. For example, suppose you create the analog input object `ai` for a sound card and you want to find the default value for the `SampleRate` property.

```
ai = analoginput('winsound');  
out = propinfo(ai, 'SampleRate');  
out.DefaultValue  
ans =  
      8000
```

Properties Grouped by Category

This section contains brief descriptions of all toolbox base properties. The properties are categorized according to these subsystems:

- Analog input properties
- Analog output properties
- Digital I/O properties

Depending on the hardware device you are using, additional property names or property values may be present. The additional property names are described in Chapter 11, “Device-Specific Property Reference.” For example, only analog input and analog output objects associated with a sound card have a `BitsPerSample` property. The additional property values are also device-specific but are included in this chapter. For example, all supported devices have an `InputType` property, but the value `AC-Coupled` is unique to analog input objects associated with a sound card.

Analog Input Properties

Analog input base properties are divided into two main categories: common properties and channel properties. Common properties apply to every channel contained by the analog input object, while channel properties can be configured for individual channels.

Common Properties

The analog input common properties are grouped into the following categories based on usage.

Analog Input Basic Setup Properties	
<code>SamplesPerTrigger</code>	Specify the number of samples to acquire for each channel group member for each trigger that occurs.
<code>SampleRate</code>	Specify the per-channel rate at which analog data is converted to digital data.
<code>TriggerType</code>	Specify the type of trigger to execute.

Analog Input Logging Properties	
LogFi l eName	Specify the name of the disk file to which information is logged.
Loggi ng	Indicate if data is being logged to memory or to a disk file.
Loggi ngMode	Specify the destination for acquired data.
LogToDi skMode	Specify whether data, events, and hardware information are saved to one disk file or to multiple disk files.

Analog Input Trigger Properties	
Ini ti al Tri gger Ti me	Indicate the absolute time of the first trigger.
Manual Tri gger HwOn	Specify that the hardware device starts when a manual trigger is issued.
Tri ggerActi on	Specify the M-file action function to execute when a trigger occurs.
Tri ggerChannel	Specify the channels serving as trigger sources.
Tri ggerCondi ti on	Specify the condition that must be satisfied before a trigger executes.
Tri ggerCondi ti on Val ue	Specify one or more voltage values that must be satisfied before a trigger executes.
Tri ggerDel ay	Specify the delay value for data logging.
Tri ggerDel ay Uni ts	Specify the units in which trigger delay data is measured.

Analog Input Trigger Properties (Continued)	
Tri ggerRepeat	Specify the number of additional times the trigger executes.
Tri ggersExecuted	Indicate the number of triggers that execute.
Tri ggerType	Specify the type of trigger to execute.
Analog Input Status Properties	
Loggi ng	Indicate if data is being logged to memory or to a disk file.
Runni ng	Indicate if the device object is running.
Sampl esAcqui red	Indicate the number of samples acquired per channel.
Sampl esAvai l abl e	Indicate the number of samples available per channel in the engine.
Analog Input Hardware Configuration Properties	
Channel Skew	Specify the time between consecutive scanned hardware channels.
Channel SkewMode	Specify how the channel skew is determined.
Cl ockSource	Specify the clock used to govern the hardware conversion rate.
I nputType	Specify the analog input hardware channel configuration.
Sampl eRate	Specify the per-channel rate at which analog data is converted to digital data.

Analog Input Action Properties	
DataMi ssedActi on	Specify the M-file action function to execute when data is missed.
InputOverRange Acti on	Specify the M-file action function to execute when acquired data exceeds the valid hardware range.
Runti meError Acti on	Specify the M-file action function to execute when a run-time error occurs.
Sampl esAcqui red Acti on	Specify the M-file action function to execute every time a predefined number of samples is acquired for each channel group member.
Sampl esAcqui red Acti onCount	Specify the number of samples to acquire for each channel group member before a samples acquired event is generated.
StartActi on	Specify the M-file action function to execute just before the device object starts running.
StopActi on	Specify the M-file action function to execute just after the device object stops running.
Ti merActi on	Specify the M-file action function to execute whenever a predefined period of time passes.
Ti merPeri od	Specify the period of time between timer events.
Tri ggerActi on	Specify the M-file action function to execute when a trigger occurs.

Analog Input General Purpose Properties	
BufferingConfig	Specify the per-channel allocated memory.
BufferingMode	Specify how memory is allocated.
Channel	Contain hardware channels added to the device object.
EventLog	Store information for specific events.
Name	Specify a descriptive name for the device object.
Tag	Specify a device object label.
Timeout	Specify an additional waiting time to extract data.
Type	Indicate the device object type.
UserData	Store data that you want to associate with a device object.

Channel Properties

The analog input channel properties are given below.

Analog Input Channel Properties	
Channel Name	Specify a descriptive channel name.
HwChannel	Specify the hardware channel ID.
Index	Indicate the MATLAB index of a hardware channel.
InputRange	Specify the range of the analog input subsystem.
NativeOffset	Indicate the offset to use when converting data from native format to doubles.
NativeScaling	Indicate the scaling to use when converting data from native format to doubles.

Analog Input Channel Properties (Continued)	
Parent	Indicate the parent (device object) of a channel.
SensorRange	Specify the range of data you expect from your sensor.
Type	Indicate a channel.
Units	Specify the engineering units label.
UnitsRange	Specify the range of data as engineering units.

Analog Output Properties

Analog output base properties are divided into two main categories: common properties and channel properties. Common properties apply to every channel contained by the analog output object, while channel properties can be configured for individual channels.

Common Properties

The analog output common properties are grouped into the following categories based on usage.

Analog Output Basic Setup Properties	
SampleRate	Specify the per-channel rate at which digital data is converted to analog data.
TriggerType	Specify the type of trigger to execute.

Analog Output Trigger Properties	
InitialTriggerTime	Indicate the absolute time of the first trigger.
TriggerAction	Specify the M-file action function to execute when a trigger occurs.

Analog Output Trigger Properties (Continued)	
TriggersExecuted	Indicate the number of triggers that execute.
TriggerType	Specify the type of trigger to execute.
Analog Output Status Properties	
Running	Indicate if the device object is running.
SamplesAvailable	Indicate the number of samples available per channel in the engine.
SamplesOutput	Indicate the number of samples output per channel from the engine.
Sending	Indicate if data is being sent to the hardware device.
Analog Output Hardware Configuration Properties	
ClockSource	Specify the clock used to govern the hardware conversion rate.
SampleRate	Specify the per-channel rate at which digital data is converted to analog data.
Analog Output Data Management Properties	
MaxSamplesQueued	Indicate the maximum number of samples that can be queued in the engine.

Analog Output Data Management Properties (Continued)	
RepeatOutput	Specify the number of additional times queued data is output.
Timeout	Specify an additional waiting time to queue data.

Analog Output Action Properties	
RuntimeErrorAction	Specify the M-file action function to execute when a run-time error occurs.
SamplesOutputAction	Specify the M-file action function to execute every time a predefined number of samples is output for each channel group member.
SamplesOutputActionCount	Specify the number of samples to output for each channel group member before a samples output event is generated.
StartAction	Specify the M-file action function to execute just before the device object starts running.
StopAction	Specify the M-file action function to execute just after the device object stops running.
TimerAction	Specify the M-file action function to execute whenever a predefined period of time passes.
TimerPeriod	Specify the period of time between timer events.
TriggerAction	Specify the M-file action function to execute when a trigger occurs.

Analog Output General Purpose Properties	
BufferingConfig	Specify the per-channel allocated memory.
BufferingMode	Specify how memory is allocated.
Channel	Contain hardware channels added to the device object.
EventLog	Store information for specific events.
Name	Specify a descriptive name for the device object.
OutOfDataMode	Specify how the value held by the analog output subsystem is determined.
Tag	Specify a device object label.
Type	Indicate the device object type.
UserData	Store data that you want to associate with a device object.

Channel Properties

The analog output channel properties are given below.

Analog Output Channel Properties	
Channel Name	Specify a descriptive channel name.
DefaultChannelValue	Specify the value held by the analog output subsystem.
HwChannel	Specify the hardware channel ID.
Index	Indicate the MATLAB index of a hardware channel.
NativeOffset	Indicate the offset to use when converting data from native format to doubles.

Analog Output Channel Properties (Continued)	
Nati veScal i ng	Indicate the scaling to use when converting data from native format to doubles.
OutputRange	Specify the range of the analog output hardware subsystem.
Parent	Indicate the parent (device object) of a channel.
Type	Indicate a channel.
Uni ts	Specify the engineering units label.
Uni tsRange	Specify the range of data as engineering units.

Digital I/O Properties

Digital I/O base properties are divided into two main categories: common properties and line properties. Common properties apply to every line contained by the digital I/O object, while line properties can be configured for individual lines.

Common Properties

The digital I/O common properties are given below.

Digital I/O Common Properties	
Li ne	Contain hardware lines added to the device object.
Name	Specify a descriptive name for the device object.
Runni ng	Indicate if the device object is running.
Tag	Specify a device object label.
Ti merActi on	Specify the M-file action function to execute whenever a predefined period of time passes.
Ti merPeri od	Specify the period of time between timer events.

Digital I/O Common Properties (Continued)	
Type	Indicate the device object type.
UserData	Store data that you want to associate with a device object.

Line Properties

The digital I/O line properties are given below.

Digital I/O Line Properties	
Direction	Specify whether a line is used for input or output.
HwLine	Specify the hardware line ID.
Index	Indicate the MATLAB index of a hardware line.
LineName	Specify a descriptive line name.
Parent	Indicate the parent (device object) of a line.
Port	Specify the port ID.
Type	Indicate a line.

Properties Listed Alphabetically

This section contains detailed descriptions of all toolbox device-specific properties. Each property reference page contains some or all of this information:

- The property name
 - A description of the property
 - The property characteristics, including:
 - Usage – the device objects the property can be used with, and whether it is a common property, or a channel or line property

Common properties apply to all channels or lines contained by the device object, while channel (line) properties can be set on a per-channel (per-line) basis. The device objects supported by the Data Acquisition Toolbox include analog input (AI), analog output (AO), and digital I/O (DIO) objects.

 - Access – whether the property is read/write or read-only

Read/write property values can be returned with the `get` command and configured with the `set` command. Read-only property values can be returned with the `get` command but cannot be configured with the `set` command.

 - Data type – the property data type

The supported data types include action function, double, string, Channel, Line, and any.

 - Read-only when running – whether a property value can be configured when the device object is running
 - Valid property values including the default value
- When property values are given by a predefined list, the default value is usually indicated by `{}` (curly braces). Default values for some properties are calculated by the data acquisition engine, while others are determined by the hardware driver. If there are device-specific values, they are listed separately.
- An example using the property
 - Related properties and functions

Purpose	Specify the per-channel allocated memory	
Description	<p>BufferingConfig is a two-element vector that specifies the per-channel allocated memory. The first element of the vector specifies the block size, while the second element of the vector specifies the number of blocks. The total allocated memory (in bytes) is given by</p> $(\text{block size}) \cdot (\text{number of blocks}) \cdot (\text{number of channels}) \cdot (\text{native data type})$ <p>You can determine the native data type with <code>daqhwinfo</code>.</p> <p>You can allocate memory automatically or manually. If BufferingMode is Auto, the BufferingConfig values are automatically set by the engine. If BufferingMode is Manual, then you must manually set the BufferingConfig values. If you change the BufferingConfig values, BufferingMode is automatically set to Manual.</p> <p>When memory is automatically allocated by the engine, the block-size value depends on the sampling rate and is typically a binary number. The number of blocks is initially set to a value of 30 but can dynamically increase to accommodate the memory requirements. In most cases, the number of blocks used results in a per-channel memory that is somewhat greater than the SamplesPerTrigger value. When you manually allocate memory, the number of blocks is not dynamic and care must be taken to ensure there is sufficient memory to store the acquired data. If the number of samples acquired or queued exceeds the allocated memory, then an error is returned.</p> <p>You can easily determine the memory allocated and available memory for each device object with the <code>daqmem</code> function.</p>	
Characteristics	Usage	AI, AO, Common
	Access	Read/write
	Data type	Two-element vector of doubles
	Read-only when running	Yes

BufferingConfig

Values The default value is determined by the engine, and is based on the number of channels contained by the device object and the sampling rate. The `BufferingMode` value determines whether the values are automatically updated as data is acquired. For analog output objects, the default number of blocks is zero.

Example Create the analog input object `ai` for a sound card and add two channels to it.

```
ai = analoginput('winsound');  
addchannel(ai, 1:2);
```

The block size and number of blocks are given by `BufferingConfig`, while the native data type for the sound card is given by `daqhwinfo`.

```
ai.BufferingConfig  
ans =  
    512    30  
out = daqhwinfo(ai);  
out.NativeDataType  
ans =  
    int16
```

With this information, the total allocated memory is calculated to be 61,440 bytes. This number is stored by `daqmem`.

```
out = daqmem(ai);  
out.UsedBytes  
ans =  
    61440
```

The allocated memory is more than sufficient to store 8000 two-byte samples for two channels. If more memory was required, then the number of blocks would dynamically grow since `BufferingMode` is set to `Auto`.

See Also

Functions

`daqhwinfo`, `daqmem`

Properties

`BufferingMode`, `SampleRate`, `SamplesPerTrigger`

Purpose	Specify how memory is allocated	
Description	<p>BufferingMode can be set to Auto or Manual . If BufferingMode is set to Auto, the data acquisition engine automatically allocates the required memory. If BufferingMode is set to Manual , you must manually allocate memory with the BufferingConfig property.</p> <p>If BufferingMode is set to Auto and the SampleRate value is changed, then the BufferingConfig values may be recalculated by the engine. Specifically, you can increase (decrease) the block size if SampleRate is increased (decreased). If BufferingMode is set to Auto and you change the BufferingConfig values, then BufferingMode is automatically set to Manual . If BufferingMode is set to Manual , then you cannot set the number of blocks to a value less than three.</p> <p>For most data acquisition applications, you should set BufferingMode to Auto and have memory allocated by the engine since this minimizes the chance of an out-of-memory condition.</p>	
Characteristics	Usage	AI, AO, Common
	Access	Read/write
	Data type	String
	Read-only when running	Yes
Values	{Auto}	Memory is allocated by the data acquisition engine.
	Manual	Memory is allocated manually.
See Also	Functions daqmem	
	Properties BufferingConfig	

Channel

Purpose	Contain hardware channels added to a device object	
Description	<p>Channel is a vector of all the hardware channels contained by an analog input (AI) or analog output (AO) object. Since a newly created AI or AO object does not contain hardware channels, Channel is initially an empty vector. The size of Channel increases as channels are added with the addchannel function, and decreases as channels are removed using the delete function.</p> <p>Channel is used to reference one or more individual channels. To reference a channel, you must know its MATLAB index, which is given by the Index property. For example, you must use Channel with the appropriate indices when configuring channel property values.</p> <p>For scanning hardware, the scan order follows the MATLAB index. Therefore, the hardware channel associated with index 1 is sampled first, the hardware channel associated with index 2 is sampled second, and so on. To change the scan order, you can specify a permutation of the indices with Channel.</p>	
Characteristics	Usage	AI, AO, Common
	Access	Read/write
	Data type	Vector of channels
	Read-only when running	Yes
Values	Values are automatically defined when channels are added to the device object with the addchannel function. The default value is an empty column vector.	
Example	<p>Create the analog input object ai for a National Instruments card and add three hardware channels to it.</p> <pre>ai = analoginput('ni daq', 1); addchannel(ai, 0:2);</pre> <p>To set a property value for the first channel added (ID = 0), you must reference the channel by its index using the Channel property.</p> <pre>chans = ai.Channel(1); set(chans, 'InputRange', [-10 10])</pre>	

Based on the current configuration, the hardware channels are scanned in order from 0 to 2. To swap the scan order of channels 0 and 1, you can specify the appropriate permutation of the MATLAB indices with `Channel`.

```
ai.Channel([1 2 3]) = ai.Channel([2 1 3]);
```

See Also

Functions

`addchannel`, `delete`

Properties

`hwChannel`, `Index`

ChannelName

Purpose	Specify a descriptive channel name	
Description	<p>Channel Name specifies a descriptive name for a hardware channel. If a channel name is defined, then you can reference that channel by its name. If a channel name is not defined, then the channel must be referenced by its index. Channel names are not required to be unique.</p> <p>You can also define descriptive channel names when channels are added to a device object with the addchannel function.</p>	
Characteristics	Usage	AI, AO, Channel
	Access	Read/write
	Data type	String
	Read-only when running	Yes
Values	The default value is an empty string. To reference a channel by name, it must contain only letters, numbers, and underscores and must begin with a letter.	
Example	<p>Create the analog input object ai for a sound card and add two channels to it.</p> <pre>ai = analoginput('winsound'); addchannel(ai, 1:2);</pre> <p>To assign a descriptive name to the first channel contained by ai</p> <pre>Chan1 = ai.Channel(1) set(Chan1, 'Channel Name', 'Joe')</pre> <p>You can now reference this channel by name instead of by index.</p> <pre>set(ai.Joe, 'Units', 'Decibels')</pre>	
See Also	Functions addchannel	

Purpose	Specify the time between consecutive scanned hardware channels	
Description	Channel Skew applies only to scanning hardware and not to simultaneous sample and hold (SS/H) hardware.	
	For scanning hardware, Channel Skew is configurable only when Channel SkewMode is set to Manual . If Channel SkewMode is set to Minimum or Equi sampl e, then Channel Skew is automatically set to the appropriate device-specific read-only value. For SS/H hardware, the only valid Channel Skew value is zero. SS/H hardware includes Agilent Technologies devices and sound cards.	
	The Channel Skew value is specified in seconds.	
Characteristics	Usage	AI, Common
	Access	Read/write (depends on Channel SkewMode value)
	Data type	Double
	Read-only when running	Yes
Values	For all supported hardware, the default value is zero. For SS/H hardware, the only valid value is zero. For scanning hardware, the value depends on Channel SkewMode.	
See Also	Properties Channel SkewMode	

ChannelSkewMode

Purpose	Specify how the channel skew is determined	
Description	<p>For simultaneous sample and hold (SS/H) hardware, Channel SkewMode is None. For scanning hardware, Channel SkewMode can be Minimum, Equi sample, or Manual (National Instruments only). SS/H hardware includes Agilent Technologies devices and sound cards, while scanning hardware includes most National Instruments and ComputerBoards devices. Note that some supported boards from these vendors are SS/H. Examples include the 6110E board and the DAS4020/12 board, respectively.</p> <p>If Channel SkewMode is Minimum, then the minimum channel skew supported by the hardware is used. If Channel SkewMode is Equi sample, the channel skew is given by $[(\text{sampling rate})(\text{number of channels})]^{-1}$. If Channel SkewMode is Manual, then you must specify the channel skew with the Channel Skew property.</p>	
Characteristics	Usage	AI, Common
	Access	Read/write
	Data type	String
	Read-only when running	Yes
Values	National Instruments	
	{ Minimum }	The channel skew is set to the minimum supported value.
	Equi sample	The channel skew is given by $[(\text{sampling rate})(\text{number of channels})]^{-1}$.
	Manual	The channel skew is given by Channel Skew.

ComputerBoards

- { Minimum } The channel skew is set to the minimum supported value.
- Equi sampl e The channel skew is given by $[(\text{sampling rate})(\text{number of channels})]^{-1}$.

Agilent Technologies and Sound Cards

- { None } This is the only supported value for SS/H hardware.

See Also

Properties
Channel Skew, Sampl eRate

ClockSource

Purpose	Specify the clock used to govern the hardware conversion rate	
Description	<p>For all supported hardware except ComputerBoards analog output subsystems, Cl ockSource can be set to Internal , which specifies that the acquisition rate is governed by the internal hardware clock.</p> <p>For subsystems without a hardware clock, you must use software clocking to govern the sampling rate. Software clocking allows a maximum sampling rate of 500 Hz and a minimum sampling rate of 0.0002 Hz. An error is returned if more than 1 sample of jitter is detected. Note that you may not be able to attain rates over 100 Hz on all systems, especially Windows 9X.</p>	
Characteristics	Usage	AI, AO, Common
	Access	Read/write
	Data type	String
	Read-only when running	Yes
Values	National Instruments	
	{ Internal }	The internal hardware clock is used.
	External	Externally control the channel clock (AO only).
	External Sampl e Ctrl	Externally control the channel clock. This value overrides the Channel Skew property value (AI only).
	External Scan Ctrl	Externally control the scan clock. This value overrides the Sampl eRate property value (AI only).
	External Sampl e AndScanCtrl	Externally control the channel and scan clocks. This value overrides the Channel Skew and Sampl eRate property values (AI only).
	<p>For an analog output object, a Cl ockSource value of Internal is analogous to a value of Update.</p>	

ComputerBoards

{ Internal }	The internal hardware clock is used (AI only).
Software	The computer clock is used.
External	Externally control the channel clock (AI only).

Agilent Technologies

{ Internal }	The internal hardware clock is used.
External	The external sample clock, positive true.
Inverted External	The external sample clock, negative true.
VXI Bus/3	The VXI bus clock, divided by 3, provided by some other clock master
VXI BusSample	The VXI bus sample clock

Sound Cards

{ Internal }	The internal hardware clock is used.
--------------	--------------------------------------

See Also

Properties

Channel Skew, SampleRate

DataMissedAction

Purpose	Specify the M-file action function to execute when data is missed	
Description	<p>A data missed event is generated immediately after acquired data is missed. This event executes the action function specified for DataMissedAction. The default value for DataMissedAction is daqacti on, which displays the event type and the device object name.</p> <p>In most cases, data is missed because:</p> <ul style="list-style-type: none">• The engine cannot keep up with the rate of acquisition.• The driver wrote new data into the hardware's FIFO buffer before the previously acquired data was read. You can usually avoid this problem by increasing the size of the memory block with the Bufferi ngConf i g property. <p>Data missed event information is stored in the Type and Data fields of the EventLog property. The Type field value is DataMi ssed. The Data field value is Rel Sampl e, which is the last sample acquired when the event occurred.</p> <p>When a data missed event occurs, the analog input object is automatically stopped.</p>	
Characteristics	Usage	AI, Common
	Access	Read/write
	Data type	String
	Read-only when running	No
Values	The default value is daqacti on.	
See Also	<p>Functions</p> <p>daqacti on</p> <p>Properties</p> <p>EventLog</p>	

Purpose	Specify the value held by the analog output subsystem	
Description	<p>DefaultChannelValue specifies the value to write to the analog output (AO) subsystem when data is finished being output from the engine.</p> <p>DefaultChannelValue is used only when OutOfDataMode is set to DefaultValue. This property guarantees that a known value is held by the AO subsystem if a run-time error occurs. Note that sound cards do not have an OutOfDataMode property.</p>	
Characteristics	Usage	AO, Channel
	Access	Read/write
	Data type	Double
	Read-only when running	Yes
Values	The default value is zero.	
Example	<p>Create the analog output object ao and add two channels to it.</p> <pre>ao = analogoutput('ni daq', 1); addchannel (ao, 0: 1);</pre> <p>You can configure ao so that when it stops outputting data, a value of 1 volt is held for both channels.</p> <pre>ao.OutOfDataMode = 'DefaultValue'; ao.Channel.DefaultChannelValue = 1.0;</pre>	
See Also	Properties OutOfDataMode	

Direction

Purpose	Specify whether a line is used for input or output	
Description	<p>When adding hardware lines to a digital I/O object with <code>addline</code>, you must configure the line direction. The line direction can be <code>In</code> or <code>Out</code>, and is automatically stored in <code>Direction</code>. If a line direction is <code>In</code>, you can only read a value from that line. If a line direction is <code>Out</code>, you can write or read a line value.</p> <p>For line-configurable devices, you can change individual line directions using <code>Direction</code>. For port-configurable devices, you cannot change individual line directions.</p>	
Characteristics	Usage	DIO, Line
	Access	Read/write
	Data type	String
	Read-only when running	Yes
Values	<code>{In}</code>	The line can be read from.
	<code>Out</code>	The line can be read from or written to.
Example	<p>Create the digital I/O object <code>dio</code> and add two input lines and two output lines to it.</p> <pre>dio = digitalio('ni daq', 1); addline(dio, 0: 3, {'In', 'In', 'Out', 'Out'});</pre> <p>To configure all lines for output.</p> <pre>dio.Line(1: 2).Direction = 'Out';</pre>	
See Also	<p>Functions</p> <p><code>addline</code></p>	

Purpose Store information for specific events

Description EventLog is a structure array that stores information related to specific analog input (AI) or analog output (AO) events. Event information is stored in the Type and Data fields of EventLog. Type stores the event type. The logged event types are shown below.

Event Type	Description	AI	AO
Data missed	Data is missed by the engine.	✓	
Input overrange	A signal exceeds the hardware input range.	✓	
Run-time error	A run-time error is encountered. Run-time errors include timeouts and hardware errors.	✓	✓
Start	The start function is issued.	✓	✓
Stop	The device object stops executing.	✓	✓
Trigger	A trigger executes.	✓	✓

Timer events, samples available events (AI), and samples output events (AO) are not logged.

Data stores event-specific information associated with the event type in several fields. For all stored events, Data contains the Rel Sample field, which returns the input or output sample number at the time the event occurred. For the start, stop, run-time error, and trigger events, Data contains the AbsTime field, which returns the absolute time (as a clock vector) the event occurred. Other event-specific fields are included in Data. For a description of these fields, refer to “Configuring Events and Actions” in Chapter 5 for analog input objects, “Configuring Events and Actions” in Chapter 6 for analog output objects, or the appropriate reference pages in this chapter.

EventLog can store a maximum of 1000 events. If this value is exceeded, then the most recent 1000 events are stored. You can use the showdagevents function to easily display stored event information.

EventLog

Characteristics	Usage	AI, AO, Common
	Access	Read-only
	Data type	Structure array
	Read-only when running	N/A

Values Values are automatically added as events occur. The default value is an empty structure array.

Example Create the analog input object `ai` and add four channels to it.

```
ai = analoginput('ni daq', 1);
chans = addchannel(ai, 0:3);
```

Acquire 1 second of data and display the logged event types.

```
start(ai)
events = ai.EventLog;
{events.Type}
ans =
    'Start'    'Trigger'    'Stop'
```

To examine the data associated with the trigger event

```
events(2).Data
ans =
    AbsTime: [1999 2 12 14 54 52.5456]
    Rel Sample: 0
    Channel: []
    Trigger: 1
```

See Also Functions
`showdaqevents`

Purpose	Specify the hardware channel ID	
Description	<p>All channels contained by a device object have a hardware channel ID and an associated MATLAB index. The channel ID is given by <code>HwChannel</code> and the MATLAB index is given by the <code>Index</code> property. The <code>HwChannel</code> value is defined when hardware channels are added to a device object with the <code>addchannel</code> function.</p> <p>The beginning channel ID value depends on the hardware device. For National Instruments hardware, channel IDs are zero-based (begin at zero). For Agilent Technologies hardware and sound cards, channel IDs are one-based (begin at one).</p> <p>For scanning hardware, the scan order follows the MATLAB index. Therefore, the hardware channel associated with index 1 is sampled first, the hardware channel associated with index 2 is sampled second, and so on. To change the scan order, you can assign the channel IDs to different indices using <code>HwChannel</code>.</p>	
Characteristics	Usage	AI, AO, Channel
	Access	Read/write
	Data type	Double
	Read-only when running	Yes
Values	Values are automatically defined when channels are added to the device object with the <code>addchannel</code> function. The default value is one.	
Example	<p>Create the analog input object <code>ai</code> for a National Instruments board and add the first three hardware channels to it.</p> <pre>ai = analoginput('ni daq', 1); addchannel(ai, 0:2);</pre>	

HwChannel

Based on the current configuration, the hardware channels are scanned in order from 0 to 2. To swap the scan order of channels 0 and 1, you can assign these channels to the appropriate indices using `HwChannel`.

```
ai . Channel ( 1 ) . HwChannel = 1;  
ai . Channel ( 2 ) . HwChannel = 0;
```

See Also

Functions

`addchannel`

Properties

`Channel`, `Index`

Purpose	Specify the hardware line ID	
Description	<p>All lines contained by a digital I/O object have a hardware ID and an associated MATLAB index. The hardware ID is given by <code>HwLine</code> and the MATLAB index is given by the <code>Index</code> property. The <code>HwLine</code> value is defined when hardware lines are added to a digital I/O object with the <code>addline</code> function.</p> <p>The beginning line ID value depends on the hardware device. For National Instruments hardware, line IDs are zero-based (begin at zero).</p>	
Characteristics	Usage	DIO, Line
	Access	Read/write
	Data type	Double
	Read-only when running	Yes
Values	Values are automatically defined when lines are added to the digital I/O object with the <code>addline</code> function. The default value is one.	
Example	<p>Suppose you create the digital I/O object <code>di o</code> and add four hardware lines to it.</p> <pre>di o = digitalio('ni daq', 1); addline(di o, 0: 3, 'out');</pre> <p><code>addline</code> automatically assigns the indices 1-4 to these hardware lines. You can swap the hardware lines associated with index 1 and index 2 with <code>HwLine</code>.</p> <pre>di o.Line(1).HwLine = 1; di o.Line(2).HwLine = 0;</pre>	
See Also	<p>Functions</p> <p><code>addline</code></p> <p>Properties</p> <p><code>Line</code>, <code>Index</code></p>	

Index

Purpose	Indicate the MATLAB index of a hardware channel or line	
Description	<p>Every hardware channel (line) contained by a device object has an associated MATLAB index that is used to reference that channel (line). For example, to configure property values for an individual channel, you must reference the channel through the Channel property using the appropriate Index value. Likewise, to configure property values for an individual line, you must reference the line through the Line property using the appropriate Index value.</p> <p>For channels (lines), you can assign indices automatically with the addchannel (addl i ne) function. Channel (line) indices always begin at 1 and increase monotonically up to the number of channels (lines) contained by the device object. For channels, index assignments can also be made manually with the addchannel function.</p> <p>For scanning hardware, the scan order follows the MATLAB index. Therefore, the hardware channel associated with index 1 is sampled first, the hardware channel associated with index 2 is sampled second, and so on. To change the scan order, you can assign the channel IDs to different indices using the HwChannel or Channel property.</p> <p>Index provides a convenient way to access channels and lines programmatically.</p>	
Characteristics	Usage	AI, AO, Channel; DIO, Line
	Access	Read-only
	Data type	Double
	Read-only when running	N/A
Values	Values are automatically defined when channels (lines) are added to the device object with the addchannel (addl i ne) function. The default value is one.	

Example

Create the analog input object `ai` for a sound card and add two hardware channels to it.

```
ai = analoginput('winsound');  
chans = addchannel(ai, 1:2);
```

You can access the MATLAB indices for these channels with `Index`.

```
Index1 = chans(1).Index;  
Index2 = chans(2).Index;
```

See Also

Functions

`addchannel`, `addline`

Properties

`Channel`, `HwChannel`, `HwLine`, `Line`

InitialTriggerTime

Purpose	Indicate the absolute time of the first trigger					
Description	<p>For all trigger types, InitialTriggerTime records the time when Logging or Sending is set to On. The absolute time is recorded as a clock vector.</p> <p>You can return the InitialTriggerTime value with the getdata function, or with the Data.AbsTime field of the EventLog property.</p>					
Characteristics	Usage	AI, AO, Common				
	Access	Read-only				
	Data type	Six-element vector of doubles				
	Read-only when running	N/A				
Values	The value is automatically updated when the trigger executes. The default value is a vector of zeros.					
Example	Create the analog input object ai for a sound card and add two hardware channels to it.					
	<pre>ai = analoginput('winsound'); chans = addchannel(ai, 1:2);</pre>					
	<p>After starting ai, the trigger immediately executes and the trigger time is recorded.</p> <pre>start(ai) abstime = ai.InitialTriggerTime abstime = 1.0e+003 * 1.9990 0.0020 0.0190 0.0130 0.0260 0.0208</pre>					
	To convert the clock vector to a more convenient form					
	<pre>t = fix(abstime); sprintf('%d: %d: %d', t(4), t(5), t(6)) ans = 13: 26: 20</pre>					

See Also

Functions

`getdata`

Properties

`EventLog`, `Loggi ng`, `Sendi ng`

InputOverRangeAction

Purpose Specify the M-file action function to execute when acquired data exceeds the valid hardware range

Description An input overrange event is generated immediately after an overrange condition is detected for any channel group member. This event executes the action function specified for InputOverRangeAct i on.

An overrange condition occurs when an input signal exceeds the range specified by the InputRange property. Overage detection is enabled only if the analog input object is running and an action function is specified for InputOverRangeAct i on.

Input overrange event information is stored in the Type and Data fields of the EventLog property. The Type field value is OverRange. The Data field values are given below.

Data Field Value	Description
Rel Sampl e	The acquired sample number when the event occurred.
Channel	The index of the channel that experienced an overrange signal.
OverRange	Indicates if the channel went from overrange to in range, or from in range to overrange.

Characteristics

Usage	AI, Common
Access	Read/write
Data type	String
Read-only when running	No

Values The default value is an empty string.

See Also Properties
EventLog, Input Range

Purpose Specify the range of the analog input subsystem

Description `InputRange` is a two-element vector that specifies the range of voltages that can be accepted by the analog input (AI) subsystem. You should configure `InputRange` so that the maximum dynamic range of your hardware is utilized.

If an input signal exceeds the `InputRange` value, then an overrange condition occurs. Overrange detection is enabled only if the analog input object is running and a value is specified for the `InputOverRangeAction` property. For many devices, the input range is expressed in terms of the gain and polarity.

AI subsystems have a finite number of `InputRange` values that you can set. If an input range is specified but does not match a valid range, then the next highest supported range is automatically selected by the engine. If `InputRange` exceeds the range of valid values, then an error is returned. Use the `daqhwinfo` function to return the input ranges supported by your board.

Since the engine can set the input range to a value that differs from the value you specify, you should return the actual input range for each channel using the `get` function or the device object display summary. Alternatively, you can use the `setverify` function, which sets the `InputRange` value and then returns the actual value that is set.

Note If your hardware supports a channel gain list, then you can configure `InputRange` for individual channels. Otherwise, `InputRange` must have the same value for all channels contained by the analog input object.

You should use `InputRange` in conjunction with the `SensorRange` property. These two properties should be configured such that the maximum precision is obtained and the full dynamic range of the sensor signal is covered.

InputRange

Characteristics	Usage	AI, Channel
	Access	Read/write
	Data type	Two-element vector of doubles
	Read-only when running	Yes

Values The default value is supplied by the hardware driver.

Example Create the analog input object `ai` for a National Instruments board, and add two hardware channels to it.

```
ai = analoginput('ni daq', 1);
addchannel(ai, 0:1);
```

You can return the input ranges supported by the board with the `InputRanges` field of the `daqhwinfo` function.

```
out = daqhwinfo(ai);
out.InputRanges
ans =
    -0.0500    0.0500
    -0.5000    0.5000
    -5.0000    5.0000
   -10.0000   10.0000
```

To configure both channels contained by `ai` to accept input signals between -10 volts and 10 volts

```
ai.Channel.InputRange = [-10 10];
```

Alternatively, you can use the `setverify` function.

```
ActualRange = setverify(ai.Channel, 'InputRange', [-10 10]);
```

See Also **Functions**
`daqhwinfo`, `setverify`

Properties
`InputOverRangeAction`, `SensorRange`, `Units`, `UnitsRange`

Purpose	Specify the analog input hardware channel configuration	
Description	<p>For National Instruments devices, InputType can be SingleEnded, Differential, or NonReferencedSingleEnded. For ComputerBoards devices, InputType can be SingleEnded, or Differential, for Agilent Technologies devices, InputType can only be Differential. For sound cards, InputType can only be AC-Coupled.</p> <p>If channels have been added to a National Instruments or ComputerBoards analog input object and you change the InputType value, then the channels are automatically deleted if the hardware reduces the number of available channels.</p>	
Characteristics	Usage	AI, Common
	Access	Read/write
	Data type	String
	Read-only when running	Yes
Values	National Instruments	
	{ Differential }	Channels are configured for differential input.
	SingleEnded	Channels are configured for single-ended input.
	NonReferencedSingleEnded	This channel configuration is used when the input signal has its own ground reference, which is tied to the negative input of the instrumentation amplifier.
	ComputerBoards	
	{ Differential }	Channels are configured for differential input.
	SingleEnded	Channels are configured for single-ended input.

InputType

Agilent Technologies

{Differential} Channels are configured for differential input.

Sound Cards

{AC-Coupled} The input is coupled so that constant (DC) signal levels are suppressed.

Purpose	Contain hardware lines added to the device object	
Description	<p>Line is a vector of all the hardware lines contained by a digital I/O (DIO) object. Since a newly created DIO object does not contain hardware lines, Line is initially an empty vector. The size of Line increases as lines are added with the addline function, and decreases as lines are removed with the delete function.</p> <p>You can use Line to reference one or more individual lines. To reference a line, you must know its MATLAB index and hardware ID. The MATLAB index is given by the Index property, while the hardware ID is given by the HwLine property.</p>	
Characteristics	Usage	DIO, Common
	Access	Read/write
	Data type	Vector of lines
	Read-only when running	Yes
Values	Values are automatically defined when lines are added to the DIO object with the addline function. The default value is an empty column vector.	
Example	<p>Create the digital I/O object dio and add four input lines to it.</p> <pre>dio = digitalio('ni daq', 1); addline(dio, 0: 3, 'In');</pre> <p>To set a property value for the first line added (ID = 0), you can reference the line by its index using the Line property.</p> <pre>line1 = dio.Line(1); set(line1, 'Direction', 'Out')</pre>	
See Also	Functions addline, delete	
	Properties HwLine, Index	

LineName

Purpose	Specify a descriptive line name	
Description	<p>LineName specifies a descriptive name for a hardware line. If a line name is defined, then you can reference that line by its name. If a line name is not defined, then the line must be referenced by its index. Line names are not required to be unique.</p> <p>You can also define descriptive line names when lines are added to a digital I/O object with the addLine function.</p>	
Characteristics	Usage	DIO, Line
	Access	Read/write
	Data type	String
	Read-only when running	Yes
Values	The default value is an empty string. To reference a line by name, it must contain only letters, numbers, and underscores and must begin with a letter.	
Example	<p>Create the digital I/O object di o and add four hardware lines to it.</p> <pre>di o = di gi ta l i o(' ni daq' , 1); add l i ne(di o, 0: 3, ' out');</pre> <p>To assign a descriptive name to the first line contained by di o</p> <pre>l i ne1 = di o. Li ne(1); set(l i ne1, ' Li neName' , ' Joe')</pre> <p>You can now reference this line by name instead of index.</p> <pre>set(di o. Joe, ' Di rect i on' , ' I n')</pre>	
See Also	<p>Functions</p> <p>addLine</p>	

Purpose	Specify the name of the disk file to which information is logged	
Description	You can log acquired data, device object property values and event information, and hardware information to a disk file by setting the LoggingMode property to Disk or Disk&Memory.	
	You may specify any value for LogFileName as long as it conforms to the MATLAB naming conventions: the name cannot start with a number and cannot contain spaces. If no extension is specified as part of LogFileName, then daq is used. The default value for LogFileName is logfile.daq.	
	You can choose whether an output file is overwritten or if multiple log files are created with the LogToDiskMode property. Setting LogToDiskMode to Overwrite causes the output file to be overwritten. Setting LogToDiskMode to Index causes new data files to be created, each with an indexed name based on the value of LogFileName.	
Characteristics	Usage	AI, Common
	Access	Read/write
	Data type	String
	Read-only when running	Yes
Values	The default value is logfile.daq.	
See Also	Properties	
	Logging, LoggingMode, LogToDiskMode	

Logging

Purpose	Indicate if data is being logged to memory or to a disk file	
Description	Along with the Runni ng property, Loggi ng reflects the state of an analog input object. Loggi ng can be On or Off.	
	Loggi ng is automatically set to On when a trigger occurs. When Loggi ng is On, acquired data is being stored in memory or to a disk file.	
	Loggi ng is automatically set to Off when the requested samples are acquired, an error occurs, or a stop function is issued. When Loggi ng is Off, you can still preview data with the peekdata function provided Runni ng is On. However, peekdata does not guarantee that all the requested data is returned.	
	To guarantee that acquired data contains no gaps, is must be logged to memory or to a disk file. Data stored in memory is extracted with the getdata function, while data stored to disk is returned with the daqread function. The destination for logged data is controlled with the Loggi ngMode property.	
Characteristics	Usage	AI, Common
	Access	Read-only
	Data type	String
	Read-only when running	N/A
Values	{ Off }	Data is not logged to memory or a disk file.
	On	Data is logged to memory or a disk file.
See Also	Functions daqread, getdata, peekdata, stop	
	Properties Loggi ngMode, Runni ng	

Purpose	Specify the destination for acquired data	
Description	<p>LoggingMode can be set to Disk, Memory, or Disk&Memory. If LoggingMode is set to Disk, then acquired data (as well as device object and hardware information) is streamed to a disk file. If LoggingMode is set to Memory, then acquired data is stored in the data acquisition engine. If LoggingMode is set to Disk&Memory, then acquired data is stored in the data acquisition engine and is streamed to a disk file.</p> <p>When logging to the engine, you must extract the data with the getData function. If the data is not extracted, it may be overwritten.</p> <p>When logging to disk, you can specify the log filename with the LogFileName property, and you can control the number of log files created with the LogToDiskMode property. You can return data stored in a disk file to the MATLAB workspace with the daqread function.</p>	
Characteristics	Usage	AI, Common
	Access	Read/write
	Data type	String
	Read-only when running	Yes
Values	Disk	Acquired data is logged to a disk file.
	{Memory}	Acquired data is logged to memory.
	Disk&Memory	Acquired data is logged to a disk file and to memory.
See Also	<p>Functions</p> <p>daqread, getData</p> <p>Properties</p> <p>LogFileName, LogToDiskMode</p>	

LogToDiskMode

Purpose Specify whether data, events, and hardware information are saved to one disk file or to multiple disk files

Description LogToDiskMode can be set to Overwrite or Index. If LogToDiskMode is set to Overwrite, then the log file is overwritten each time start is issued. If LogToDiskMode is set to Index, a different disk file is created each time start is issued and these rules are followed:

- The first log filename is specified by the initial value of LogFileName.
- If the specified file already exists, it is overwritten and no warning is issued.
- LogFileName is automatically updated with a numeric identifier after each file is written. For example, if LogFileName is initially specified as data.daq, then data.daq is the first filename, data01.daq is the second filename, and so on.

Separate analog input objects are logged to separate files. You can return data stored in a disk file to the MATLAB workspace with the daqread function. If an error occurs during data logging, an error message is returned and data logging is stopped.

Characteristics	Usage	AI, Common
	Access	Read/write
	Data type	String
	Read-only when running	Yes

Values	Index	Multiple log files are written, each with an indexed filename based on the LogFileName property.
	{Overwrite}	The log file is overwritten.

See Also	Functions daqread
	Properties LogFileName, LoggingMode

Purpose	Specify that the hardware device starts when a manual trigger is issued	
Description	<p>Manual TriggerHwOn can be set to Start or Trigger. If Manual TriggerHwOn is Start, then the hardware device associated with your device object starts running after you issue the start function. If Manual TriggerHwOn is Trigger, then the hardware device associated with your device object starts running after you execute a manual trigger with the trigger function. You can use trigger only when you configure the TriggerType property to Manual.</p> <p>You should configure Manual TriggerHwOn to Trigger when you want to synchronize the input and output of data, or you require more control over when an analog input device starts. Note that you cannot use peekdata or acquire pretrigger data when you use this value. Additionally, you should not use this value with repeated triggers since the subsequent behavior is undefined.</p>	
Characteristics	Usage	AI, Common
	Access	Read/write
	Data type	String
	Read-only when running	Yes
Values	{Start}	Start the hardware after the start function is issued.
	Trigger	Start the hardware after the trigger function is issued.
Example	<p>Create the analog input object ai and the analog output object ao for a sound card and add two channels to each device object.</p> <pre>ai = analoginput('winsound'); addchannel(ai, 1:2); ao = analogoutput('winsound'); addchannel(ao, 1:2);</pre>	

ManualTriggerHwOn

To operate the sound card in full duplex mode, and to minimize the time between when ai starts and ao starts, you configure Manual TriggerHwOn to Trigger for ai and TriggerType to Manual for both ai and ao.

```
set([ai ao], 'TriggerType', 'Manual')  
ai.ManualTriggerHwOn = 'Trigger';
```

The analog input and analog output hardware devices will both start after you issue the trigger function. For a detailed example that uses Manual TriggerHwOn, refer to “Starting Multiple Device Objects” on page 6-37.

See Also

Functions

peekdata, start, trigger

Properties

TriggerType

Purpose	Indicate the maximum number of samples that can be queued in the engine	
Description	<p>MaxSampl esQueued indicates the maximum number of samples allowed in the analog output queue. The default value is calculated by the engine, and is based on the memory resources of your system. You can override the default value of MaxSampl esQueued with the daqmem function.</p> <p>The value of MaxSampl esQueued can affect the behavior of put data. For example, if the queued data exceeds the value of MaxSampl esQueued, then put data becomes a blocking function until there is enough space in the queue to add the additional data.</p>	
Characteristics	Usage	AO, Common
	Access	Read-only
	Data type	Double
	Read-only when running	N/A
Values	The value is calculated by the data acquisition engine.	
See Also	Functions daqmem, put data	

Name

Purpose	Specify a descriptive name for the device object	
Description	When a device object is created, a descriptive name is automatically generated and stored in Name. This name is produced by concatenating the name of the adaptor, the device ID, and the device object type. You can change the value of Name at any time.	
Characteristics	Usage	AI, AO, DIO, Common
	Access	Read/write
	Data type	String
	Read-only when running	No
Values	The value is defined after the device object is created.	
Example	Create the analog input object ai for a sound card.	
	<pre>ai = analoginput('winsound');</pre> <p>The descriptive name for ai is given by</p> <pre>ai.Name ans = winsound0-AI</pre>	

Purpose Indicate the offset to use when converting data from native format to doubles

Description `NativeOffset`, along with `NativeScaling`, is used to convert data from your hardware's native format to doubles using the formula

$$\text{data} = (\text{native data})(\text{native scaling}) + \text{native offset}$$

For analog input objects, you return data from the engine in native format using the `getData` function. Additionally, if you log native data to a `.daq` file, then you can read back that data using the `daqread` function. For analog output objects, you queue native data in the engine using the `putData` function. You may want to return or queue data in native format to conserve memory and to increase speed.

The `NativeOffset` value for a given channel may change if you change the `UnitsRange`, `SensorRange`, or `InputRange` property values for that channel.

You can return the native data type of your hardware device with the `daqhwinfo` function.

Characteristics	Usage	AI, AO, Channel
	Access	Read-only
	Data type	Double
	Read-only when running	N/A

Values The default value is device-specific.

Example Create the analog input object `ai` for a National Instruments board, and add eight channels to it.

```
ai = analoginput('ni daq', 1);  
addchannel(ai, 0:7);
```

Start `ai`, collect one second of data for each channel, and extract the data from the engine using the native format of the device.

```
start(ai)  
nativedata = getdata(ai, 1000, 'native');
```

NativeOffset

You can return the native data type of the board with the `daqhwinfo` function.

```
out = daqhwinfo(ai);  
out.NativeDataType  
ans =  
int16
```

Convert the data to doubles using the `NativeScaling` and `NativeOffset` properties.

```
scaling = get(ai.Channel(1), 'NativeScaling');  
offset = get(ai.Channel(1), 'NativeOffset');  
data = double(nativedata)*scaling + offset;
```

See Also

Functions

`daqhwinfo`, `daqread`, `getdata`, `putdata`

Properties

`InputRange`, `NativeScaling`, `SensorRange`, `UnitsRange`

Purpose Indicate the scaling to use when converting data from native format to doubles

Description `NativeScaling`, along with `NativeOffset`, is used to convert data from your hardware's native format to doubles using the formula

$$\text{data} = (\text{native data})(\text{native scaling}) + \text{native offset}$$

For analog input objects, you return data from the engine in native format using the `getData` function. Additionally, if you log native data to a `.daq` file, then you can read back that data using the `daqread` function. For analog output objects, you queue native data in the engine using the `putData` function. You may want to return or queue data in native format to conserve memory and to increase speed.

The `NativeScaling` value for a given channel may change if you change the `UnitsRange`, `SensorRange`, or `InputRange` property values for that channel.

You can return the native data type of your hardware device with the `daqhwinfo` function.

Characteristics	Usage	AI, AO, Channel
	Access	Read-only
	Data type	Double
	Read-only when running	N/A

Values The default value is device-specific.

See Also **Functions**
`daqhwinfo`, `daqread`, `getData`, `putData`

Properties
`InputRange`, `NativeOffset`, `SensorRange`, `UnitsRange`

OutputRange

Purpose	Specify the range of the analog output hardware subsystem	
Description	<p>OutputRange is a two-element vector that specifies the range of voltages that can be output by the analog output (AO) subsystem. You should configure OutputRange so that the maximum dynamic range of your hardware is utilized. For many devices, the output range is expressed in terms of the gain and polarity.</p> <p>AO subsystems have a finite number of OutputRange values that you can set. If an output range is specified but does not match a valid range, then the next highest supported range is automatically selected by the engine. If OutputRange exceeds the range of valid values, then an error is returned. Use the daqhwi nfo function to return the output ranges supported by your board.</p> <p>Since the engine can set the output range to a value that differs from the value you specify, you should return the actual output range for each channel using the get function or the device object display summary. Alternatively, you can use the setveri fy function, which sets the OutputRange value and then returns the actual value that is set.</p>	
Characteristics	Usage	AO, Channel
	Access	Read/write
	Data type	Two-element vector of doubles
	Read-only when running	Yes
Values	The default value is determined by the hardware driver.	
Example	<p>Create the analog output object ao for a National Instruments board and add two hardware channels to it.</p> <pre>ao = analogoutput('ni daq', 1); addchannel(ao, 0: 1);</pre>	

You can return the output ranges supported by the board with the `OutputRanges` field of the `daqhwinfo` function.

```
out = daqhwinfo(ao);  
out.OutputRanges  
ans =  
    0.0000    10.0000  
   -10.0000    10.0000
```

To configure both channels contained by `ao` to output signals between -10 volts and 10 volts

```
ao.Channel.OutputRange = [-10 10];
```

Alternatively, you can use the `setverify` function to configure and return the `OutputRange` value.

```
ActualRange = setverify(ao.Channel, 'OutputRange', [-10 10]);
```

See Also

Functions

`daqhwinfo`, `setverify`

Properties

`Units`, `UnitsRange`

Parent

Purpose	Indicate the parent (device object) of a channel or line	
Description	<p>The parent of a channel (line) is defined as the device object that contains the channel (line).</p> <p>You can create a copy of the device object containing a particular channel or line by returning the value of Parent. You can treat this copy like any other device object. For example, you can configure property values, add channels or lines to it, and so on.</p>	
Characteristics	Usage	AI, AO, Channel; DIO, Line
	Access	Read-only
	Data type	Device object
	Read-only when running	N/A
Values	The value is defined when channels or lines are added to the device object.	
Example	<p>Create the analog input object ai for a National Instruments board and add three hardware channels to it.</p> <pre>ai = analoginput('ni daq', 1); chans = addchannel(ai, 0:2);</pre> <p>To return the parent for channel 2</p> <pre>parent = ai.Channel(2).Parent;</pre> <p>parent is an exact copy of ai.</p> <pre>i = equal(ai, parent) ans = 1</pre>	

Purpose	Specify the port ID	
Description	Hardware lines are often grouped together as a port. Digital I/O subsystems can consist of multiple ports and typically have eight lines per port. When adding hardware lines to a digital I/O object with <code>addl i ne</code> , you can specify the port ID. The port ID is stored in the <code>Port</code> property. If the port ID is not specified, then the smallest port ID value is automatically used.	
Characteristics	Usage	DIO, Line
	Access	Read-only
	Data type	Double
	Read-only when running	N/A
Values	The port ID is defined when line are added to the digital I/O object with <code>addl i ne</code> .	
Example	Create the digital I/O object <code>di o</code> and add two hardware channels to it.	
	<pre>di o = di gi ta l i o(' ni daq' , 1); addl i ne(di o, 0: 1, ' I n');</pre> <p>You can use <code>Port</code> property to return the port IDs associated with the lines contained by <code>di o</code>.</p> <pre>di o. Li ne. Port ans = [0] [0]</pre>	
See Also	Functions <code>addl i ne</code>	

RepeatOutput

Purpose	Specify the number of additional times queued data is output	
Description	<p>To send data to an analog output subsystem, it must first be queued in the data acquisition engine with the <code>putdata</code> function. If you want to continuously output the same data, you can use multiple calls to <code>putdata</code>. However, since each <code>putdata</code> call consumes memory, a long output sequence can quickly bring your system to halt.</p> <p>As an alternative to <code>putdata</code>, you can continuously output previously queued data using <code>RepeatOutput</code>. Since <code>RepeatOutput</code> requeues the data, additional memory resources are not consumed. While the data is being output, you cannot add additional data to the queue.</p>	
Characteristics	Usage	AO, Common
	Access	Read/write
	Data type	Double
	Read-only when running	Yes
Values	The default value is zero.	
Example	<p>Create the analog output object <code>ao</code> for a sound card and add one channel to it.</p> <pre>ao = analogoutput('winsound'); chans = addchannel(ao, 1);</pre> <p>To queue one second of data</p> <pre>data = sin(linspace(0, 10, 8000))'; putdata(ao, data)</pre> <p>To continuously output data for 10 seconds</p> <pre>set(ao, 'RepeatOutput', 9)</pre>	
See Also	<p>Functions</p> <p><code>putdata</code></p>	

Purpose	Indicate if the device object is running	
Description	<p>Along with the Logging or Sending property, Running reflects the state of an analog input or analog output object. Running can be On or Off.</p> <p>Running is automatically set to On once the start function is issued. When Running is On, you can acquire data from an analog input device or send data to an analog output device after the trigger occurs. For digital I/O objects, Running is typically used to indicate if time-based events are being generated.</p> <p>Running is automatically set to Off once the stop function is issued, the specified data is acquired or sent, or a run-time error occurs. When Running is Off, you cannot acquire or send data. However, you can acquire one sample with the getsample function, or send one sample with the putsample function.</p>	
Characteristics	Usage	AI, AO, DIO, Common
	Access	Read-only
	Data type	String
	Read-only when running	N/A
Values	{Off}	The device object is not running.
	On	The device object is running.
See Also	<p>Functions getsample, putsample, start</p> <p>Properties Logging, Sending</p>	

RuntimeErrorAction

Purpose Specify the M-file action function to execute when a run-time error occurs

Description A run-time error event is generated immediately after a run-time error occurs. This event executes the action function specified for `RuntimeErrorAction`. Additionally, a toolbox error message is automatically displayed to the MATLAB workspace. If an error occurs that is not explicitly handled by the toolbox, then the hardware-specific error message is displayed.

The default value for `RuntimeErrorAction` is `daqaction`, which displays the event type, the time the event occurred, and the device object name along with the error message.

Run-time error event information is stored in the `Type` and `Data` fields of the `EventLog` property. The `Type` field value is `Error`. The `Data` field values are given below.

Data Field Value	Description
AbsTime	The absolute time (as a clock vector) the event occurred.
RelSample	The acquired (AI) or output (AO) sample number when the event occurred.
String	The descriptive error message.

Run-time errors include hardware errors and timeouts. Run-time errors do not include configuration errors such as setting an invalid property value.

Characteristics	Usage	AI, AO, Common
	Access	Read/write
	Data type	String
	Read-only when running	No

Values The default value is `daqaction`.

See Also

Functions

`daqaction`

Properties

`EventLog`, `Timeout`

SampleRate

Purpose Specify the per-channel rate at which analog data is converted to digital data, or digital data is converted to analog data

Description `SampleRate` specifies the per-channel rate (in samples/second) that an analog input (AI) or analog output (AO) subsystem converts data. AI subsystems convert analog data to digital data, while AO subsystems convert digital data to analog data.

AI and AO subsystems have a finite (though often large) number of valid sampling rates. If you specify a sampling rate that does not match one of the valid values, then the data acquisition engine automatically selects the next highest supported sampling rate that is within 1% of the requested rate. If there is no valid sampling rate within this tolerance, the next highest value is selected. If a higher value is not available, then an error is returned.

Since the engine can set the sampling rate to a value that differs from the value you specify, you should return the actual sampling rate using the `get` function or the device object display summary. Alternatively, you can use the `set verify` function, which sets the `SampleRate` value and then returns the actual value that is set. To find out the range of sampling rates supported by your board, use the `propinfo` function. Additionally, since the actual sampling rate depends on the number of channels contained by the device object and the `Channel Skew` property value (AI only), `SampleRate` should be the last property you set before starting the device object.

Characteristics	Usage	AI, AO, Common
	Access	Read/write
	Data type	Double
	Read-only when running	Yes

Values The default value is obtained from the hardware driver.

Example

Create the analog input object `ai` for a sound card and add two channels to it.

```
ai = analoginput('winsound');  
addchannel(ai, 1:2);
```

You can find out the range of valid sampling rates with the `ConstraintValue` field of the `propinfo` function.

```
rates = propinfo(ai, 'SampleRate');  
rates.ConstraintValue  
ans =  
      8000      48000
```

To configure the per-channel sampling rate to 48 kHz

```
set(ai, 'SampleRate', 48000)
```

Alternatively, you can use the `setverify` function to configure and return the `SampleRate` value.

```
ActualRate = setverify(ai, 'SampleRate', 48000);
```

See Also

Functions

`propinfo`, `setverify`

Properties

Channel Skew

SamplesAcquired

Purpose	Indicate the number of samples acquired per channel	
Description	SamplesAcquired is continuously updated to reflect the current number of samples acquired by an analog input object. It is reset to zero after a start function is issued.	
	Use the SamplesAvailable property to find out how many samples are available to be extracted from the engine.	
Characteristics	Usage	AI, Common
	Access	Read-only
	Data type	Double
	Read-only when running	N/A
Values	The value is continuously updated to reflect the current number of samples acquired. The default value is zero.	
See Also	Functions start	
	Properties SamplesAvailable	

Purpose Specify the M-file action function to execute every time a predefined number of samples is acquired for each channel group member

Description A samples acquired event is generated immediately after the number of samples specified by the `SamplesAcquiredActionCount` property is acquired for each channel group member. This event executes the action function specified for `SamplesAcquiredAction`.

You should use `SamplesAcquiredAction` if you must access each sample that is acquired. If you do not have this requirement, you may want to use the `TimerAction` property.

Samples acquired event information is not stored in the `EventLog` property.

Characteristics	Usage	AI, Common
	Access	Read/write
	Data type	String
	Read-only when running	No

Values The default value is an empty string.

See Also **Properties**
`EventLog`, `SamplesAcquiredActionCount`, `TimerAction`

SamplesAcquiredActionCount

Purpose Specify the number of samples to acquire for each channel group member before a samples acquired event is generated

Description A samples acquired event is generated immediately after the number of samples specified by SamplesAcquiredActionCount is acquired for each channel group member. This event executes the action function specified by the SamplesAcquiredAction property.

Characteristics	Usage	AI, Common
	Access	Read/write
	Data type	Double
	Read-only when running	Yes

Values The default value is 1024.

See Also Properties
SamplesAcquiredAction

Purpose	Indicate the number of samples available per channel in the engine	
Description	<p>For analog input (AI) objects, SamplesAvailable indicates the number of samples that can be extracted from the engine for each channel group member with the getData function. For analog output (AO) objects, SamplesAvailable indicates the number of samples that have been queued with the putData function, and can be sent (output) to each channel group member.</p> <p>After data has been extracted (AI) or output (AO), the SamplesAvailable value is reduced by the appropriate number of samples. For AI objects, SamplesAvailable is reset to zero after a start function is issued.</p> <p>For AI objects, use the SamplesAcquired property to find out how many samples have been acquired since the start function was issued. For AO objects, use the SamplesOutput property to find out how many samples have been output since the start function was issued.</p>	
Characteristics	Usage	AI, AO, Common
	Access	Read-only
	Data type	Double
	Read-only when running	N/A
Values	The value is automatically updated based on the number of samples acquired (analog input) or sent (analog output). The default value is zero.	
See Also	Functions start	
	Properties SamplesAcquired, SamplesOutput	

SamplesOutput

Purpose	Indicate the number of samples output per channel from the engine	
Description	<p>Sampl esOutput is continuously updated to reflect the current number of samples output by an analog output object. It is reset to zero after the device objects stops and data has been queued with the put data function.</p> <p>Use the Sampl esAvai l abl e property to find out how many samples are available to be output from the engine.</p>	
Characteristics	Usage	AO, Common
	Access	Read-only
	Data type	Double
	Read-only when running	N/A
Values	The value is continuously updated to reflect the current number of samples output. The default value is zero.	
See Also	<p>Functions</p> <p>put data</p> <p>Properties</p> <p>Sampl esAvai l abl e</p>	

Purpose	Specify the M-file action function to execute every time a predefined number of samples is output for each channel group member	
Description	<p>A samples output event is generated immediately after the number of samples specified by the <code>SamplesOutputActionCount</code> property is output for each channel group member. This event executes the action function specified for <code>SamplesOutputAction</code>.</p> <p>Samples output event information is not stored in the <code>EventLog</code> property.</p>	
Characteristics	Usage	AO, Common
	Access	Read/write
	Data type	String
	Read-only when running	No
Values	The default value is an empty string.	
See Also	<p>Properties</p> <p><code>EventLog</code>, <code>SamplesOutputActionCount</code></p>	

SamplesOutputActionCount

Purpose	Specify the number of samples to output for each channel group member before a samples output event is generated	
Description	A samples output event is generated immediately after the number of samples specified by SamplesOutputActionCount is output for each channel group member. This event executes the action function specified by the SamplesOutputAction property.	
Characteristics	Usage	AO, Common
	Access	Read/write
	Data type	Double
	Read-only when running	Yes
Values	The default value is 1024.	
See Also	Properties SamplesOutputAction	

Purpose	Specify the number of samples to acquire for each channel group member for each trigger that occurs	
Description	<p>Sampl esPerTri gger specifies the number of samples to acquire for each analog input channel group member for each trigger that occurs. If Sampl esPerTri gger is set to Inf, then the analog input object continually acquires data until a stop function is issued or an error occurs.</p> <p>The default value of Sampl esPerTri gger is calculated by the data acquisition engine such that one second of data is acquired. This calculation is based on the value of Sampl eRate.</p>	
Characteristics	Usage	AI, Common
	Access	Read/write
	Data type	Double
	Read-only when running	Yes
Values	The default value is set by the engine such that one second of data is acquired.	
Example	<p>Create the analog input object ai for a sound card and add two channels to it.</p> <pre>ai = anal ogi nput (' wi nsound'); addchannel (ai , 1: 2);</pre> <p>By default, a one second acquisition in which 8000 samples are acquired for each channel is defined. To define a two second acquisition at the same sampling rate</p> <pre>set (ai , ' Sampl esPerTri gger' , 16000)</pre>	
See Also	<p>Functions</p> <p>stop</p> <p>Properties</p> <p>Sampl eRate</p>	

Sending

Purpose	Indicate if data is being sent to the hardware device	
Description	Along with the Runni ng property, Sendi ng reflects the state of an analog output object. Sendi ng can be On or Off.	
	Sendi ng is automatically set to On when a trigger occurs. When Sendi ng is On, queued data is being output to the analog output subsystem.	
	Sendi ng is automatically set to Off when the queued data has been output, an error occurs, or a stop function is issued. When Sendi ng is Off, data is not being output to the analog output subsystem although you can output a single sample with the put sampl e function.	
Characteristics	Usage	AO, Common
	Access	Read-only
	Data type	String
	Read-only when running	N/A
Values	{ Off }	Data is not being sent to the analog output hardware.
	On	Data is being sent to the analog output hardware.
See Also	Functions put sampl e	
	Properties Runni ng	

Purpose	Specify the range of data you expect from your sensor	
Description	<p>You use <code>SensorRange</code> to scale your data to reflect the range you expect from your sensor. You can find the appropriate sensor range from your sensor's specification sheet. For example, an accelerometer might have a sensor range of ± 5 volts, which corresponds to ± 50 g's ($1\text{ g} = 9.80\text{ m/s/s}$).</p> <p>The data is scaled while it is extracted from the engine with the <code>getData</code> function according to the formula</p> $\text{scaled value} = (\text{A/D value})(\text{units range})/(\text{sensor range})$ <p>The A/D value is constrained by the <code>InputRange</code> property, which reflects the gain and polarity of your hardware channels. The units range is given by the <code>UnitsRange</code> property</p>	
Characteristics	Usage	AI, Channel
	Access	Read/write
	Data type	Two-element vector of doubles
	Read-only when running	No
Values	The default value is determined by the default value of the <code>InputRange</code> property.	
See Also	Functions <code>getData</code>	
	Properties <code>InputRange</code> , <code>Units</code> , <code>UnitsRange</code>	

StartAction

Purpose Specify the M-file action function to execute just before the device object starts running

Description A start event is generated immediately after the start function is issued. This event executes the action function specified for StartAction. When the action function has finished executing, Runni ng is automatically set to 0n and the device object and hardware device begin executing.

Start event information is stored in the Type and Data fields of the EventLog property. The Type field value is Start. The Data field values are given below.

Data Field Value	Description
AbsTi me	The absolute time (as a cl ock vector) the event occurred.
Rel Sampl e	The acquired (AI) or output (AO) sample number when the event occurred.

Characteristics

Usage	AI, AO, Common
Access	Read/write
Data type	String
Read-only when running	No

Values The default value is an empty string.

See Also

Functions
start

Properties
EventLog, Runni ng

Purpose Specify the M-file action function to execute just after the device object stops running

Description A stop event is generated immediately after the device object and hardware device stop executing. This occurs when:

- A stop function is issued.
- For analog input (AI) objects, the requested number of samples to acquire was reached or data was missed. For analog output (AO) objects, the requested number of samples to output was reached.
- A run-time error occurred.

A stop event executes the action function specified for StopActi on. Under most circumstances, the action function is not guaranteed to complete execution until sometime after the device object and hardware device stop, and the Runni ng property is set to Off.

Stop event information is stored in the Type and Data fields of the EventLog property. The Type field value is Stop. The Data field values are given below.

Data Field Value	Description
AbsTi me	The absolute time (as a cl ock vector) the event occurred.
Rel Sampl e	The acquired (AI) or output (AO) sample number when the event occurred.

Characteristics	Usage	AI, AO, Common
	Access	Read/write
	Data type	String
	Read-only when running	No

Values The default value is an empty string.

StopAction

See Also

Functions

stop

Properties

EventLog, Runni ng

Purpose	Specify a device object label	
Description	Tag provides a means to identify device objects with a label. Using the <code>daqfind</code> function and the Tag value, you can identify and retrieve a device object that was cleared from the MATLAB workspace.	
Characteristics	Usage	AI, AO, DIO, Common
	Access	Read/write
	Data type	String
	Read-only when running	No
Values	The default value is an empty string.	
Example	Create the analog input object <code>ai</code> for a sound card and add two channels to it.	
	<pre>ai = analoginput('winsound'); addchannel(ai, 1:2);</pre>	
	Assign <code>ai</code> a label using Tag.	
	<pre>set(ai, 'Tag', 'Sound')</pre>	
	If <code>ai</code> is cleared from the workspace, you can use <code>daqfind</code> and the Tag value to identify and retrieve the device object.	
See Also	<pre>clear ai ai_cell = daqfind('Tag', 'Sound'); ai = ai_cell{1};</pre>	
	Functions	
	<code>daqfind</code>	

Timeout

Purpose	Specify an additional waiting time to extract or queue data	
Description	The <code>Timeout</code> value (in seconds) is added to the time required to extract data from the engine or queue data to the engine. Since data is extracted with the <code>getData</code> function, and queued with the <code>putData</code> function, <code>Timeout</code> is associated only with these two “blocking” functions.	
	If the requested data is not extracted or queued after waiting the required time, then a timeout condition occurs and control is immediately returned to MATLAB. A timeout is one of the conditions for stopping an acquisition. When a timeout occurs, the action function specified by <code>RuntimeErrorAction</code> is called.	
	<code>Timeout</code> is not associated with hardware timeout conditions. Possible hardware timeout conditions include:	
	<ul style="list-style-type: none">• Triggering on a voltage level and that level never occurs• Externally clocking an acquisition and the external clock signal never occurs• Losing the hardware connection	
Characteristics	To check for hardware timeouts, you may need to poll the appropriate property value.	
	Usage	AI, AO, Common
	Access	Read/write
	Data type	Double
	Read-only when running	Yes
Values	The default value is one second.	
See Also	Functions <code>getData</code> , <code>putData</code>	
	Properties <code>RuntimeErrorAction</code>	

Purpose Specify the M-file action function to execute whenever a predefined period of time passes

Description A timer event is generated whenever the time specified by the `TimerPeriod` property passes. This event executes the action function specified for `TimerAction`. Time is measured relative to when the device object starts running.

Some timer events may not be processed if your system is significantly slowed or if the `TimerPeriod` value is too small. For example, a common application for timer events is to display data. However, since displaying data is a CPU-intensive task, some of these events may be dropped. To guarantee that events are not dropped, you may want to use the `SamplesAcquiredAction` property (analog input) or the `SamplesOutputAction` property (analog output). For digital I/O objects, timer events are typically used to update and display the state of the device object.

Timer event information is not stored in the `EventLog` property.

Characteristics	Usage	AI, AO, DIO, Common
	Access	Read/write
	Data type	String
	Read-only when running	No

Values The default value is an empty string.

See Also **Properties**
`EventLog`, `SamplesAcquiredAction`, `SamplesOutputAction`, `TimerPeriod`

TimerPeriod

Purpose	Specify the period of time between timer events	
Description	<p>Ti merPeri od specifies the time, in seconds, that must pass before the action function specified for Ti merAct i on is called. Time is measured relative to when the hardware device starts running.</p> <p>Some timer events may not be processed if your system is significantly slowed or if the Ti merPeri od value is too small. For example, a common application for timer events is to display data. However, since displaying data is a CPU-intensive task, some of these events may be dropped.</p>	
Characteristics	Usage	AI, AO, DIO, Common
	Access	Read/write
	Data type	Double
	Read-only when running	No
Values	The default value is 0.1 second.	
See Also	Properties Ti merAct i on	

Purpose Specify the M-file action function to execute when a trigger occurs

Description A trigger event is generated immediately after a trigger occurs. This event executes the action function specified for `TriggerAction`. Under most circumstances, the action function is not guaranteed to complete execution until sometime after `Logging` is set to `On` for analog input (AI) objects, or `Sending` is set to `On` for analog output (AO) objects.

Trigger event information is stored in the `Type` and `Data` fields of the `EventLog` property. The `Type` field value is `Trigger`. The `Data` field values are given below.

Data Field Value	Description
AbsTime	The absolute time (as a clock vector) the event occurred.
RelSample	The acquired (AI) or output (AO) sample number when the event occurred.
Channel	The index number for each input channel serving as a trigger source (AI only).
Trigger	The trigger number (AI only).

Characteristics	Usage	AI, AO, Common
	Access	Read/write
	Data type	String
	Read-only when running	No

Values The default value is an empty string.

See Also Functions
`trigger`

Properties
`EventLog`, `Logging`

TriggerChannel

Purpose	Specify the channels serving as trigger sources	
Description	<p>Tri ggerChannel specifies the channels serving as trigger sources. The trigger channel(s) must be specified before the trigger type. You may need to configure the Tri ggerCondi ti on and Tri ggerCondi ti onVal ue properties in conjunction with Tri ggerChannel .</p> <p>For all supported vendors, if Tri ggerType is Software, then you must acquire data from the channel being used for the trigger source. For National Instruments hardware, if Tri ggerType is HwAnal ogChannel , then Tri ggerChannel must be the first element of the channel group. The exception is if you are using the PCI-6110 or PCI-6111 board. In this case, you can specify any channel for the Tri ggerChannel value.</p>	
Characteristics	Usage	AI, Common
	Access	Read/write
	Data type	Vector of channels
	Read-only when running	Yes
Values	The default value is an empty vector.	
Example	<p>Create the analog input object ai , add two channels, and define the trigger source as channel 2.</p> <pre>ai = anal ogi nput (' wi nsound'); ch = addchannel (ai , 1: 2); set (ai , ' Tri ggerChannel ' , ch(2)) set (ai , ' Tri ggerType' , ' Software')</pre>	
See Also	<p>Properties</p> <p>Tri ggerCondi ti on, Tri ggerCondi ti onVal ue, Tri ggerType</p>	

Purpose	Specify the condition that must be satisfied before a trigger executes											
Description	<p>The available trigger conditions depend on the value of Tri ggerType. If Tri ggerType is Immedi ate or Manual , the only available Tri ggerCondi ti on is None. If Tri ggerType is Software, then Tri ggerCondi ti on can be Ri si ng, Fal li ng, Leavi ng, or Enteri ng. These trigger conditions require one or more voltage values to be specified for the Tri ggerCondi ti onVal ue property.</p> <p>Based on the hardware you are using, additional trigger conditions may be available.</p>											
Characteristics	Usage	AI, Common										
	Access	Read/write										
	Data type	String										
	Read-only when running	Yes										
Values	<p>All Supported Hardware</p> <p>The following trigger condition is used when Tri ggerType is Immedi ate or Manual .</p> <table><tr><td>{ None }</td><td>No trigger condition is required.</td></tr></table> <p>The following trigger conditions are available when Tri ggerType is Software.</p> <table><tr><td>{ Ri si ng }</td><td>The trigger occurs when the signal has a positive slope when passing through the specified value.</td></tr><tr><td>Fal li ng</td><td>The trigger occurs when the signal has a negative slope when passing through the specified value.</td></tr><tr><td>Leavi ng</td><td>The trigger occurs when the signal leaves the specified range of values.</td></tr><tr><td>Enteri ng</td><td>The trigger occurs when the signal enters the specified range of values.</td></tr></table>		{ None }	No trigger condition is required.	{ Ri si ng }	The trigger occurs when the signal has a positive slope when passing through the specified value.	Fal li ng	The trigger occurs when the signal has a negative slope when passing through the specified value.	Leavi ng	The trigger occurs when the signal leaves the specified range of values.	Enteri ng	The trigger occurs when the signal enters the specified range of values.
{ None }	No trigger condition is required.											
{ Ri si ng }	The trigger occurs when the signal has a positive slope when passing through the specified value.											
Fal li ng	The trigger occurs when the signal has a negative slope when passing through the specified value.											
Leavi ng	The trigger occurs when the signal leaves the specified range of values.											
Enteri ng	The trigger occurs when the signal enters the specified range of values.											

TriggerCondition

National Instruments

The following trigger condition is used when `TriggerType` is `HwDigital`.

`{None}` No trigger condition is required.

The following trigger conditions are available when `TriggerType` is `HwAnalogChannel` or `HwAnalogPin`.

- `{AboveHighLevel}` The trigger occurs when the analog signal is above the specified value.
- `BelowLowLevel` The trigger occurs when the analog signal is below the specified value.
- `InsideRegion` The trigger occurs when the analog signal is inside the specified region.
- `LowHysteresis` The trigger occurs when the analog signal is less than the specified low value with hysteresis given by the specified high value.
- `HighHysteresis` The trigger occurs when the analog signal is greater than the specified high value with hysteresis given by the specified low value.

ComputerBoards

The following trigger conditions are available when `TriggerType` is `HwDigital`.

<code>GateHigh</code>	The trigger occurs as long as the digital signal is high.
<code>GateLow</code>	The trigger occurs as long as the digital signal is low.
<code>Tri gHigh</code>	The trigger occurs when the digital signal is high.
<code>Tri gLow</code>	The trigger occurs when the digital signal is low.
<code>Tri gPosEdge</code>	The trigger occurs when the positive (rising) edge of the digital signal is detected.
<code>{Tri gNegEdge}</code>	The trigger occurs when the negative (falling) edge of the digital signal is detected.

The following trigger conditions are available when `TriggerType` is `HwAnalog`.

<code>{TriggerAbove}</code>	The trigger occurs when the analog signal makes a transition from below the specified value to above.
<code>TriggerBelow</code>	The trigger occurs when the analog signal makes a transition from above the specified value to below.
<code>GateNegHys</code>	The trigger occurs when the analog signal is more than the specified high value. The acquisition stops if the analog signal is less than the specified low value.
<code>GatePosHys</code>	The trigger occurs when the analog signal is less than the specified low value. The acquisition stops if the analog signal is more than the specified high value.
<code>GateAbove</code>	The trigger occurs as long as the analog signal is more than the specified value.
<code>GateBelow</code>	The trigger occurs as long as the analog signal is less than the specified value.
<code>GateInWindow</code>	The trigger occurs as long as the analog signal is within the specified range of values.
<code>GateOutWindow</code>	The trigger occurs as long as the analog signal is outside the specified range of values.

Agilent Technologies

The following trigger conditions are available when `TriggerType` is `HwDigital`.

<code>{PositiveEdge}</code>	The trigger occurs when the positive (rising) edge of a digital signal is detected.
<code>NegativeEdge</code>	The trigger occurs when the negative (falling) edge of a digital signal is detected.

TriggerCondition

The following trigger conditions are available when `TriggerType` is `HwAnalog`.

<code>{ Rising }</code>	The trigger occurs when the analog signal has a positive slope when passing through the specified range of values.
<code>Falling</code>	The trigger occurs when the analog signal has a negative slope when passing through the specified range of values.
<code>Leaving</code>	The trigger occurs when the analog signal leaves the specified range of values.
<code>Entering</code>	The trigger occurs when the analog signal enters the specified range of values.

Note that when `TriggerType` is `HwAnalog`, the trigger condition values are all specified as two-element vectors. Setting two trigger levels prevents the module from triggering repeatedly due to a noisy signal.

See Also

Properties

`TriggerChannel` , `TriggerConditionValue`, `TriggerType`

Purpose	Specify one or more voltage values that must be satisfied before a trigger executes	
Description	<p>TriggerConditionValue is used when TriggerType is Software, and is ignored when TriggerCondition is None.</p> <p>To execute a software trigger, the values specified for TriggerCondition and TriggerConditionValue must be satisfied. When TriggerCondition is Rising or Falling, TriggerConditionValue accepts a single value. When TriggerCondition is Entering or Leaving, TriggerConditionValue accepts a two-element vector of values.</p>	
Characteristics	Usage	AI, Common
	Access	Read/write
	Data type	Double (or a two-element vector of doubles)
	Read-only when running	Yes
Values	The default value is zero.	
Example	<p>Create the analog input object ai and add one channel to it.</p> <pre>ai = analoginput('winsound'); ch = addchannel(ai, 1);</pre> <p>The trigger executes when a signal with a negative slope passing through 0.2 volts is detected on channel 1.</p> <pre>set(ai, 'TriggerChannel', ch) set(ai, 'TriggerType', 'Software') set(ai, 'TriggerCondition', 'Falling') set(ai, 'TriggerConditionValue', 0.2)</pre>	
See Also	Properties TriggerCondition, TriggerType	

TriggerDelay

Purpose	Specify the delay value for data logging	
Description	<p>You can define both pretriggers and postriggers. Pretriggers are specified with a negative <code>TriggerDelay</code> value while postriggers are specified with a positive <code>TriggerDelay</code> value. You can delay a trigger in units of time or samples with the <code>TriggerDelayUnits</code> property. Pretriggers are not defined for hardware triggers or when <code>TriggerType</code> is <code>Immediate</code>.</p> <p>Pretrigger samples are included as part of the total samples acquired per trigger as specified by the <code>SamplesPerTrigger</code> property. If sample-time pairs are returned to the workspace with the <code>getdata</code> function, then the pretrigger samples are identified with negative time values.</p>	
Characteristics	Usage	AI, Common
	Access	Read/write
	Data type	Double
	Read-only when running	Yes
Values	The default value is zero.	
Example	<p>Create the analog input object <code>ai</code> and add one channel to it.</p> <pre>ai = analoginput('winsound'); ch = addchannel(ai, 1);</pre> <p>Configure <code>ai</code> to acquire 44,100 samples per trigger with 11,025 samples (0.25 seconds) acquired as pretrigger data.</p> <pre>set(ai, 'SampleRate', 44100) set(ai, 'TriggerType', 'Manual') set(ai, 'SamplesPerTrigger', 44100) set(ai, 'TriggerDelay', -0.25)</pre>	
See Also	<p>Properties</p> <p><code>SamplesPerTrigger</code>, <code>TriggerDelayUnits</code></p>	

Purpose	Specify the units in which trigger delay data is measured	
Description	<p>TriggerDelayUnits can be Seconds or Samples. If TriggerDelayUnits is Seconds, then data logging is delayed by the specified time for each channel group member. If TriggerDelayUnits is Samples, then data logging is delayed by the specified number of samples for each channel group member.</p> <p>The trigger delay value is given by the TriggerDelay property.</p>	
Characteristics	Usage	AI, Common
	Access	Read/write
	Data type	String
	Read-only when running	Yes
Values	{Seconds}	The trigger is delayed by the specified number of seconds.
	Samples	The trigger is delayed by the specified number of samples.
See Also	Properties TriggerDelay	

TriggerRepeat

Purpose	Specify the number of additional times the trigger executes	
Description	<p>You can configure a trigger to occur once (one-shot acquisition) or multiple times. If <code>TriggerRepeat</code> is set to its default value of zero, then the trigger executes once. If <code>TriggerRepeat</code> is set to a positive integer value, then the trigger executes the specified number of times. If <code>TriggerRepeat</code> is set to <code>inf</code> then the trigger executes continuously until a stop function is issued or an error occurs.</p> <p>You can quickly evaluate how many triggers have executed by examining the <code>TriggersExecuted</code> property or by invoking the display summary for the device object. The display summary is invoked by typing the device object name at the MATLAB command line.</p>	
Characteristics	Usage	AI, Common
	Access	Read/write
	Data type	Double
	Read-only when running	Yes
Values	The default value is zero.	
See Also	Functions <code>disp</code> , <code>stop</code>	
	Properties <code>TriggersExecuted</code> , <code>TriggerType</code>	

Purpose	Indicate the number of triggers that execute	
Description	You can find out how many triggers executed by returning the value of TriggersExecuted. The trigger number for each trigger executed is also recorded by the Data.Trigger field of the EventLog property.	
Characteristics	Usage	AI, AO, Common
	Access	Read-only
	Data type	Double
	Read-only when running	N/A
Values	The default value is zero.	
Example	Create the analog input object ai and add one channel to it.	
	<pre>ai = analoginput('winsound'); ch = addchannel(ai, 1);</pre>	
	Configure ai to acquire 40,000 samples with five triggers using the default sampling rate of 8000 Hz.	
	<pre>set(ai, 'TriggerRepeat', 4) start(ai)</pre>	
	TriggersExecuted returns the number of triggers executed.	
See Also	<pre>ai.TriggersExecuted ans = 5</pre>	
	Properties	EventLog

TriggerType

Purpose Specify the type of trigger to execute

Description TriggerType can be Immediate, Manual, or Software. If TriggerType is Immediate, the trigger occurs immediately after the start function is issued. If TriggerType is Manual, the trigger occurs immediately after the trigger function is issued. If TriggerType is Software, the trigger occurs when the associated trigger condition is satisfied (AI only).

For a given hardware device, additional trigger types may be available. Some trigger types require trigger conditions and trigger condition values. Trigger conditions are specified with the TriggerCondition property, while trigger condition values are specified with the TriggerConditionValue property.

When a trigger occurs for an analog input object, data logging is initiated and the Logging property is automatically set to On. When a trigger occurs for an analog output object, data sending is initiated and the Sending property is automatically set to On.

Characteristics	Usage	AI, AO, Common
	Access	Read/write
	Data type	String
	Read-only when running	Yes

Values	All Supported Hardware	
	{Immediate}	The trigger executes immediately after start is issued. Pretrigger data cannot be captured.
	Manual	The trigger executes immediately after the trigger function is issued.
	Software	The trigger executes when the associated trigger condition is satisfied. Trigger conditions are given by the TriggerCondition property. (AI only).

National Instruments

- HwDi gi tal The trigger source is the falling edge of an external digital signal. Pretrigger data cannot be captured.
- HwAnal ogChannel The trigger source is an external analog signal (AI only).
- HwAnal ogPi n The trigger source is a low-range external analog signal (AI only).

For 1200 Series hardware, HwDi gi tal is the only device-specific Tri ggerType value for analog input subsystems. Analog output subsystems do not support any device-specific Tri ggerType values.

ComputerBoards

- HwDi gi tal The trigger source is an external digital signal (AI only). Pretrigger data cannot be captured.
- HwAnal og The trigger source is an external analog signal (AI only).

Agilent Technologies

- HwDi gi tal The trigger source is an external digital signal (AI only). Pretrigger data cannot be captured.
- HwAnal og The trigger source is an external analog signal (AI only).
- HwDi gi tal Pos The trigger source is the positive edge of an external digital signal (AO only).
- HwDi gi tal Neg The trigger source is the negative edge of an external digital signal (AO only).

See Also

Functions

start, trigger

Properties

Loggi ng, Sendi ng, Tri ggerCondi ti on, Tri ggerCondi ti onVal ue

Type

Purpose	Indicate the device object type, a channel, or a line							
Description	<p>Type is associated with device objects, channels, and lines. For device objects, Type can be Analog Input, Analog Output, or Digital I/O. Once a device object is created, the value of Type is automatically defined.</p> <p>For channels, the only value of Type is Channel . For lines, the only value of Type is Line. The value is automatically defined when channels or lines are added to the device object.</p>							
Characteristics	Usage	AI, AO, Common, Channel; DIO, Common, Line						
	Access	Read-only						
	Data type	String						
	Read-only when running	N/A						
Values	<p>Device Objects</p> <p>For device objects, Type has these possible values</p> <table><tr><td>Analog Input</td><td>The device object type is analog input.</td></tr><tr><td>Analog Output</td><td>The device object type is analog output.</td></tr><tr><td>Digital I/O</td><td>The device object type is digital I/O.</td></tr></table> <p>The value is automatically defined after the device object is created.</p> <p>Channels and Lines</p> <p>For channels, the only value of Type is Channel . For lines, the only value of Type is Line. The value is automatically defined when channels or lines are added to the device object.</p>		Analog Input	The device object type is analog input.	Analog Output	The device object type is analog output.	Digital I/O	The device object type is digital I/O.
Analog Input	The device object type is analog input.							
Analog Output	The device object type is analog output.							
Digital I/O	The device object type is digital I/O.							

Purpose	Specify the engineering units label	
Description	Uni ts is a string that specifies the engineering units label to associate with your data. You should use Uni ts in conjunction with the Uni tsRange property.	
Characteristics	Usage	AI, AO, Channel
	Access	Read/write
	Data type	String
	Read-only when running	No
Values	The default value is Vol ts.	
See Also	Properties Uni tsRange	

UnitsRange

Purpose	Specify the range of data as engineering units	
Description	You use <code>UnitsRange</code> to scale your data to reflect particular engineering units.	
	For analog input objects, the data is scaled while it is extracted from the engine with the <code>getData</code> function according to the formula	
	$\text{scaled value} = (\text{A/D value})(\text{units range})/(\text{sensor range})$	
	The A/D value is constrained by the <code>InputRange</code> property, which reflects the gain and polarity of your analog input channels. The sensor range is given by the <code>SensorRange</code> property, which reflects the range of data you expect from your sensor.	
	For analog output objects, the data is scaled when it is queued in the engine with the <code>putData</code> function according to the formula	
Characteristics	$\text{scaled value} = (\text{original value})(\text{output range})/(\text{units range})$	
	The output range is constrained by the <code>OutputRange</code> property, which specifies the gain and polarity of your analog output channels.	
	For both objects, you can also use the <code>Units</code> property to associate a meaningful label with your data.	
	Usage	AI, AO, Channel
	Access	Read/write
Values	Data type	Two-element vector of doubles
	Read-only when running	No
	The default value is determined by the default value of the <code>InputRange</code> or the <code>OutputRange</code> property.	
	See Also	
	Functions <code>getData</code> , <code>putData</code>	
	Properties <code>InputRange</code> , <code>OutputRange</code> , <code>SensorRange</code> , <code>Units</code>	

Purpose	Store data that you want to associate with a device object	
Description	UserData stores data that you want to associate with the device object.	
Characteristics	Usage	AI, AO, DIO, Common
	Access	Read/write
	Data type	Any type
	Read-only when running	No
Values	The default value is an empty vector.	
Example	Create the analog input object ai and add two channels to it.	
	<pre>ai = analoginput('ni daq', 1); addchannel(ai, 0:1);</pre> <p>Suppose you want to access filter coefficients during the acquisition. You can create a structure to store these coefficients, which can then be stored in UserData.</p> <pre>coeff.a = 1.0; coeff.b = -1.25; set(ai, 'UserData', coeff)</pre>	

Device-Specific Property Reference

Overview	11-2
Getting Command Line Property Help	11-2
 Properties Listed by Vendor	11-3
National Instruments Properties	11-3
ComputerBoards Properties	11-4
Agilent Technologies Properties	11-4
Sound Card Properties	11-5
 Properties Listed Alphabetically	11-6

Overview

The following sections provide descriptions of all toolbox device-specific properties. Device-specific properties apply only to hardware devices of a specific type or from a specific vendor. For example, the `BitsPerSample` property is supported only for sound cards, while the `NumMuxBoards` property is supported only for National Instruments devices.

In “Properties Listed by Vendor” on page 11-3, the properties are summarized according to vendor. In “Properties Listed Alphabetically” on page 11-6, the properties are described in detail.

Getting Command Line Property Help

To get command line property help, you should use the `daqhelp` function. For example, to get help for the sound card’s `BitsPerSample` property

```
daqhelp BitsPerSample
```

Note You can use `daqhelp` without creating a device object.

You can also get property characteristics, such as the default property value, using the `propinfo` function. For example, suppose you create the analog input object `ai` for a sound card and want to find the default value for the `BitsPerSample` property.

```
ai = analoginput('winsound');  
out = propinfo(ai, 'BitsPerSample');  
out.DefaultValue  
ans =  
    16
```


Properties Listed by Vendor

This section contains brief descriptions of all toolbox device-specific properties. The properties are grouped according to these supported vendors:

- National Instruments properties
- ComputerBoards properties
- Agilent Technologies properties
- Sound card properties (multi-vendor)

You can display device-specific properties with the `set` function. The device-specific properties are displayed after the base properties.

Note Some device-specific property values are not available for all devices. Refer to your hardware documentation for detailed information about device-specific behavior.

National Instruments Properties

The device-specific National Instruments properties for analog input (AI) and analog output (AO) objects are given below.

Property Name	Description	Device Objects
DriveAISenseToGround	Specify if AISENSE is driven to onboard ground.	AI
NumMuxBoards	Specify the number of external multiplexer devices connected.	AI
OutOfDataMode	Specify how the value held by the analog output subsystem is determined	AO
TransferMode	Specify how data is transferred from the data acquisition device to system memory.	AI, AO

ComputerBoards Properties

The device-specific ComputerBoards (Measurement Computing Corporation) properties for analog input (AI) and analog output (AO) objects are given below.

Property Name	Description	Device Objects
OutOfDataMode	Specify how the value held by the analog output subsystem is determined	AO
TransferMode	Specify how data is transferred from the data acquisition device to system memory.	AI, AO

Agilent Technologies Properties

The device-specific Agilent Technologies properties for analog input (AI) and analog output (AO) objects are given below.

Property Name	Description	Device Objects
COLA	Specify whether the source constant-level output is enabled or disabled.	AO
Coupling	Specify the input coupling mode.	AI
GroundingMode	Specify the input channel grounding mode.	AI
InputMode	Specify the channel input mode.	AI
InputSource	Specify the input to the A/D converter.	AI
RampRate	Specify the source ramp-up and ramp-down rate.	AO
SourceMode	Specify the source mode.	AO
SourceOutput	Specify the source output.	AO

Property Name	Description	Device Objects
Span	Specify the measurement bandwidth in Hz.	AI, AO
Sum	Specify whether the source sum input is enabled or disabled.	AO

Sound Card Properties

The device-specific sound card properties for analog input (AI) and analog output (AO) objects are given below.

Property Name	Description	Device Objects
BitsPerSample	Specify the number of bits the sound card uses to represent each sample.	AI, AO
StandardSample Rates	Specify whether the valid sample rates snap to a small set of standard values, or if you can set the sample rate to any value within the allowed bounds.	AI, AO

Properties Listed Alphabetically

This section contains detailed descriptions of all toolbox device-specific properties. Each property reference page contains some or all of this information:

- The property name
- A description of the property
- The property characteristics, including:
 - Vendor – the vendors that support the property
 - Usage – whether it is a common property, or a channel or line property and which device objects the property is associated with

Common properties apply to all channels or lines contained by the device object. Channel (line) properties can be set on a per-channel (per-line) basis. The device objects supported by the Data Acquisition Toolbox include analog input (AI), analog output (AO), and digital I/O (DIO) objects.
 - Access – whether the property is read/write or read-only

Read/write property values can be returned with the `get` command, and configured with the `set` command.

Read-only property values can be returned with the `get` command, but cannot be configured with the `set` command.
 - Data type – the property data type

The supported data types include action function, double, string, Channel, Line, and any.
 - Read-only when running – whether a property value can be configured while the device object is running
- Valid property values including the default value

When property values are given by an enumerated list, the default value is usually indicated by `{}`. Default values for some properties are calculated by the data acquisition engine, while others are determined by the hardware driver. If there are device-specific enumerated values, they are listed separately.
- An example using the property
- Related properties and functions

Purpose	Specify the number of bits the sound card uses to represent each sample	
Description	<p>BitsPerSample can be 8 or 16. If BitsPerSample is 8, the sound card represents each sample with 8 bits. This means that each sample can be represented by any number between 0 and 255. If BitsPerSample is 16, the sound card represents each sample with 16 bits. This means that each sample can be represented by any number between 0 and 65,535.</p> <p>For older Sound Blaster cards configured for full duplex operation, you may not be able to set BitsPerSample to 16 bits for both the analog input and analog output subsystems. Instead, you need to set one subsystem for 8 bits, and the other subsystem for 16 bits.</p> <hr/> <p>Note Sound cards with greater than 16 bits are not supported by the toolbox.</p> <hr/>	
Characteristics	Vendor	Sound cards
	Usage	AI, AO, Common
	Access	Read/write
	Data type	Double
	Read-only when running	Yes
Values	8	Represent data with 8 bits.
	{ 16 }	Represent data with 16 bits.

COLA

Purpose Specify whether the source constant-level output is enabled or disabled

Description COLA can be Off or On. If COLA is Off, the source constant level output is disabled. If COLA is ON, the source constant level output is enabled.

For the Option 1D4 single-channel source, the source COLA output is shared with the source sum input. Only one of these two may be enabled at any one time. For prototype Option 1D4 sources only, one of the two must be enabled at all times, and the default is for the constant-level output to be enabled.

Characteristics	Vendor	Agilent Technologies
	Usage	AO, Channel
	Access	Read/write
	Data Type	String
	Read-only while Running	Yes

Values	{Off}	The source constant level output is disabled.
	On	The source constant level output is enabled.

Purpose	Specify the input coupling mode	
Description	<p>Coupling can be DC or AC. If Coupling is DC, the input is connected directly to the amplifier. If Coupling is AC, a series capacitor is inserted between the input connector and the amplifier. For source channels, Coupling is generally not used since it is usually not possible to AC couple the output of a source.</p> <p>After a hardware reset, Coupling is automatically set to DC.</p>	
Characteristics	Vendor	Agilent Technologies
	Usage	AI, Channel
	Access	Read/write
	Data type	String
	Read-only while Running	Yes
Values	{DC}	The input is connected directly to the amplifier.
	AC	A series capacitor is inserted between the input connector and the amplifier.

DriveAISenseToGround

Purpose	Specify if AISENSE is driven to onboard ground	
Purpose	<p>DriveAISenseToGround can be On or Off. If DriveAISenseToGround is Off, then AISENSE is not driven to onboard ground. If DriveAISenseToGround is On, then AISENSE is driven to onboard ground. AISENSE serves as the hardware reference node when channels are configured in nonreferenced single-ended (NRSE) mode.</p> <p>Channels are configured with the InputType property. If InputType is NonReferencedSingleEnded, then the hardware device uses AISENSE for the negative input of the amplifier, regardless of the DriveAISenseToGround value. If InputType is Differential or SingleEnded, then the hardware device drives AISENSE to onboard ground only if DriveAISenseToGround is On.</p>	
Characteristics	Vendor	National Instruments
	Usage	AI, Common
	Access	Read/write
	Data type	String
	Read-only when running	Yes
Values	{ Off }	Do not drive AISENSE to onboard ground.
	On	Drive AISENSE to onboard ground.
See Also	Properties InputType	

Purpose	Specify the input channel grounding mode	
Description	<p>GroundingMode can be Grounded or Floating. If GroundingMode is Grounded, the low side of the channel is grounded. If GroundingMode is Floating, the low side of the channel floats thereby making the input a differential input. GroundingMode can be set only for input channels. Source channels are never floating, and are always grounded.</p> <p>If a smart break-out box is attached to the channel, then the grounding mode is automatically set to the appropriate value. If a dumb break-out box (or no break-out box) is attached to the channel, the grounding mode is given by the GroundingMode value. In this case, no hardware settings are changed and no errors are generated.</p>	
Characteristics	Vendor	Agilent Technologies
	Usage	AI, Channel
	Access	Read/write
	Data type	String
	Read-only while Running	Yes
Values	{ Grounded }	The input channel is grounded.
	Floating	The input channel is floating.

InputMode

Purpose	Specify the channel input mode
Description	<p>InputMode can be set to Vol t age, ICP, Charge, Mi c, or 200Vol tMi c. You can set InputMode to Charge only if a charge break-out box is attached to the specified channel. Your can set InputMode to Mi c or 200Vol tMi c only if a microphone break-out box is attached to the specified channel. For each input mode, the full-scale setting is configured with the InputRange property.</p> <p>There is no ICP current source inside the E1432 device. Instead, you can attach a break-out box containing an ICP current source to the input. When this ICP break-out box is attached, setting InputMode to ICP enables the ICP current source in the break-out box. If there is no ICP break-out box attached to the input, then setting InputMode to ICP does nothing.</p>
	<p>Note If a channel is not connected to a smart break-out box, then changing its input mode causes the input mode for all channels within the device to change. If there is a smart break-out box present, then you can set the input mode on a per-channel basis.</p>

Characteristics	Vendor	Agilent Technologies
	Usage	AI, Channel
	Access	Read/write
	Data Type	String
	Read-only while Running	Yes

Values	{Vol t age}	The input mode is set to volts.
	ICP	The input mode is set to ICP.
	Charge	The input mode is set to charge-amp.
	Mic	The input mode is set to microphone.
	200Vol tMi c	The input mode is set to microphone with 200 volt supply turned on.

See Also **Properties**
Coupl i ng, I nputRange, I nputType

InputSource

Purpose	Specify the input to the A/D converter	
Description	<p>For input channels, InputSource can be Swi tchBox, CALIN, Ground, and BOBCALIN. The BOBCALIN value is valid only when a smart break-out box (BoB) is connected to the input. Smart BoB's include a charge break-out box or a microphone break-out box. When a smart break-out box is attached to an E1432 or E1433 module, additional input modes are available. These additional input modes are available through the InputMode property.</p> <p>After a hardware reset, InputSource is automatically set to Swi tchBox.</p>	
Characteristics	Vendor	Agilent Technologies
	Usage	AI, Channel
	Access	Read/write
	Data Type	String
	Read-only while Running	Yes
Values	{ Swi tchBox }	Select the front panel connector.
	CALIN	Select the module's CALIN line.
	Ground	Ground the input.
	BOBCALIN	Select the module's CALIN line via the CAL connection in a break-out box.
See Also	Properties InputMode	

Purpose	Specify the number of external multiplexer devices connected	
Description	NumMuxBoards specifies the number of AMUX-64T multiplexer devices connected to your hardware. NumMuxBoards can be 0, 1, 2, or 4. If you are using a 1200 Series board, then NumMuxBoards can only be 0.	
Characteristics	Vendor	National Instruments
	Usage	AI, Common
	Access	Read/write
	Data type	Double
	Read-only when running	No
Values	{ 0 }, 1, 2, or 4	The number of AMUX-64T multiplexer devices connected.

OutOfDataMode

Purpose	Specify how the value held by the analog output subsystem is determined	
Description	<p>When queued data is output to the analog output (AO) subsystem, the hardware typically holds a value. For National Instruments and ComputerBoards devices, the value held is determined by OutOfDataMode.</p> <p>OutOfDataMode can be Hold or DefaultValue. If OutOfDataMode is Hold, then the last value output is held by the AO subsystem. If OutOfDataMode is DefaultValue, then the value specified by the DefaultChannelValue property is held by the AO subsystem.</p>	
Characteristics	Vendor	National Instruments, ComputerBoards
	Usage	AO, Common
	Access	Read/write
	Data type	String
	Read-only while Running	Yes
Values	{Hold}	Hold the last output value.
	DefaultValue	Hold the value specified by DefaultChannelValue.
Example	<p>Create the analog output object ao and add two channels to it.</p> <pre>ao = analogoutput('ni daq', 1); addchannel(ao, 0:1);</pre> <p>You can configure ao so that when queued data is finished being output, a value of 1 volt is held for both channels.</p> <pre>ao.OutOfDataMode = 'DefaultValue'; ao.Channel.DefaultChannelValue = 1.0;</pre>	
See Also	<p>Properties</p> <p>DefaultChannelValue</p>	

Purpose Specify the source ramp-up and ramp-down rate

Purpose For input channels, RampRate is not generally used. For source channels, RampRate is usually used to ensure that the source signal starts and stops smoothly.

Characteristics	Vendor	Agilent Technologies
	Usage	AO, Channel
	Access	Read/write
	Data Type	Double
	Read-only while Running	Yes

Values You can set RampRate to any value between 0 and 100 seconds.

SourceMode

Purpose	Specify the source mode	
Description	<p>If SourceMode is set to Arbitrary, the host program must provide the data to use for the arbitrary source signal.</p> <p>Note that there is no Off source mode. To turn a source channel off, you must make it inactive. When the source is inactive, it is normally low impedance to ground. To make the source high impedance, set the SourceOutput property to Open.</p>	
Characteristics	Vendor	Agilent Technologies
	Usage	AO, Channel
	Access	Read/write
	Data Type	String
	Read-only while Running	Yes
Values	{Arbitrary}	An arbitrary source signal.
See Also	Properties SourceOutput	

Purpose	Specify the source output	
Description	SourceOutput can be Normal , Grounded, Open, CALOUT, or SRC&CALOUT.	
	If SourceOutput is Normal , the normal source output is used. This output is defined by the source mode and other source parameters.	
	If SourceOutput is Grounded, the source output connector remains grounded while the source D/A converter is internally connected to the CALOUT line in the module.	
	If SourceOutput is Open, the source remains open-circuited even when the source is started. The impedance on the output is only about 1 kilohm because the power-fail decay circuit is still connected to the output.	
	If SourceOutput is CALOUT, the source output is connected to the module's internal CALOUT line. This allows the module's CALOUT line to be driven by an external signal applied at the source output connector.	
Characteristics	If SourceOutput is SRC&CALOUT, the source output is connected to the module's internal CALOUT line, and the source D/A converter is also connected to the CALOUT line. This is a combination of the Grounded and CALOUT values, and is useful for multi-mainframe calibration.	
	Vendor	Agilent Technologies
	Usage	AO, Channel
	Access	Read/write
	Data Type	String
	Read-only while Running	Yes

SourceOutput

Values	{ Normal }	Normal source output.
	Grounded	The source output connector remains grounded while the source D/A converter is internally connected to the CALOUT line in the module.
	Open	The source remains open-circuited even when the source is started.
	CALOUT	The source output is connected to the module's internal CALOUT line.
	SRC&CALOUT	The source output connector remains grounded while the source D/A converter is internally connected to the CALOUT line in the module. The source output is also connected to the module's internal CALOUT line.

Purpose	Specify the measurement bandwidth (in Hz)
Description	<p>For an input channel, span specifies the maximum frequency at which valid alias-protected data is received. Frequencies above this value are filtered out. For a source channel, Span specifies the maximum frequency at which the output signal will correctly track the signal that the source is attempting to generate.</p> <p>The valid values for Span depend of the current clock frequency. You should set the clock frequency before setting Span. Normally, the maximum valid span is the clock frequency divided by 2.56. Valid spans are given by the maximum span divided by powers of two, and the maximum span divided by five and by powers of two. The ratio between the span and the maximum span is called the <i>decimation factor</i>.</p> <p>For the E1432 module, the maximum number of decimate-by-two passes allowed is nine. Therefore, the maximum decimation factor is $5 \cdot 2^9$, and the minimum valid span is $(\text{clock frequency}) / (2.56 \cdot 5 \cdot 2^9)$. If the clock frequency is larger than 51.2 kHz, then the module is unable to do a decimation factor of one. In this case, the minimum decimation factor is two and the maximum valid span is $(\text{clock frequency}) / 5.12$.</p> <p>For the E1433 module, the maximum number of decimate-by-two passes allowed is 12, so the maximum decimation factor is $5 \cdot 2^{12}$. Due to limits in the module's DSP processor, when the clock frequency is set higher than 102,400 Hz, it is unable to do any decimation. In this case, the only valid span is $(\text{clock frequency}) / 2.56$. If you attempt to use decimation when the clock frequency is above 102,400 Hz, then an error may occur when the measurement starts.</p> <p>For the Option 1D4 source board, the maximum number of decimate-by-two passes allowed is 16, and the maximum decimation factor is $5 \cdot 2^{16}$.</p> <p>The effective sample rate is defined as the rate at which data is received from an input or used by a source, and is normally equal to 2.56 times the span. If the data is oversampled, then the effective sample rate is 5.12 times the span.</p> <p>If the digital filters in a module have a cutoff that is sharper than $1/2.56$, then some of the frequencies above the maximum span may contain valid alias-protected data. This is the case with the E1432 and E1433 modules, which have a top span filter cutoff of $(\text{clock frequency}) / 2.226$, which is 23 kHz when the clock frequency is 51.2 kHz, 88.3 kHz when the clock frequency is</p>

Span

196.608 kHz. However, Span ignores the extra bandwidth so that the maximum span is always 1/2.56 times the effective sample rate.

Span applies to an entire E1432 module rather than to one of its channels. After a hardware reset, each module is automatically set to the maximum legal span.

Characteristics	Vendor	Agilent Technologies
	Usage	AI, AO, Common
	Access	Read/write
	Data Type	Double
	Read-only while Running	Yes

Values Normally, the maximum valid span is given by the clock frequency divided by 2.56. Valid spans are given by the maximum span divided by powers of two, and the maximum span divided by five and by powers of two. The value set for Span automatically updates the SampleRate value.

See Also Properties
SampleRate

Purpose Specify whether the valid sample rates snap to a small set of standard values, or if you can set the sample rate to any value within the allowed bounds

Description StandardSampleRates can be On or Off. If StandardSampleRates is Off, then it is possible to set the sample rate to any value within the bounds supported by the hardware. For most sound cards, the lower bound is 8.000 kHz, while the upper bound is 44.0 kHz. For newer sound cards, an upper bound of 96.0 kHz may be supported. The specified sample rate is rounded up to the next integer value.

If StandardSampleRates is On, then the available sample rates snap to a small set of standard values. The standard values are 8.000 kHz, 11.025 kHz, 22.050 kHz, and 44.100 kHz. If you specify a sampling rate that is within one percent of a standard value, then the sampling rate snaps to that standard value. If you specify a sampling rate that is not within one percent of a standard value, then the sampling rate rounds up to the closest standard value.

Regardless of the StandardSampleRates value, if you specify a sampling rate that is outside the allowed limits, then an error is returned.

Characteristics	Vendor	Sound cards
	Usage	AI, AO, Common
	Access	Read/write
	Data type	String
	Read-only when running	Yes

Values	{On}	The sample rate can be set only to a small set of standard values.
	Off	If supported by the hardware, the sample rate can be set to any value within the allowed bounds, up to a maximum of 96.0 kHz.

Sum

Purpose	Specify whether the source sum input is enabled or disabled	
Description	<p>Sum can be Off or On. If Sum is Off, the sum input is disabled. If Sum is On, the sum input is enabled. The signal on the sum input is internally added to the output that the source would otherwise produce.</p> <p>For the Option 1D4 single-channel source, the source sum input is shared with the source COLA output. Only one of these two may be enabled at any one time. For prototype Option 1D4 sources, one of the two must be enabled at all times. By default, the constant-level output is enabled and the sum input is disabled.</p>	
Characteristics	Vendor	Agilent Technologies
	Usage	AO, Channel
	Access	Read/write
	Data Type	String
	Read-only while Running	No
Values	{Off}	Disable the source sum input.
	On	Enable the source sum input.

Purpose	Specify how data is transferred from the data acquisition device to system memory
Description	<p>For National Instruments hardware, TransferMode can be Interrupts or SingleDMA for both analog input and analog output subsystems. If TransferMode is Interrupts, then data is transferred from the hardware first-in, first-out memory buffer (FIFO) to system memory using interrupts. If TransferMode is SingleDMA, then data is transferred from the hardware FIFO to system memory using a single direct memory access (DMA) channel. Some boards also support a TransferMode of Dual DMA for analog input subsystems. For example, the AT-MIO-16E-1 board supports this transfer mode. If TransferMode is Dual DMA, then data is transferred from the hardware FIFO to system memory using two DMA channels. Depending on your system resources, data transfer via interrupts can significantly degrade system performance.</p> <p>For ComputerBoards hardware, TransferMode can be Default, InterruptPerPoint, DMA, or InterruptPerBlock. If TransferMode is Default, the transfer mode is automatically selected by the driver based on the board type and the sampling rate. If TransferMode is InterruptPerPoint, a single conversion is transferred for each interrupt. You should use this property value if your sampling rate is less than 5 kHz or you specify a small block size for memory buffering (as defined by the BufferingConfig property). If TransferMode is DMA, data is transferred using a single DMA channel. If TransferMode is InterruptPerBlock, a block of data is transferred for each interrupt. You should use this property value if your sampling rate is greater than 5 kHz and you are using a board that has a fast maximum sampling rate. Note that a data block is defined by the board, and usually corresponds to half the FIFO size.</p>

Note If your sampling rate is greater than around 5 kHz, you should avoid using interrupts if possible. The recommended TransferMode setting for your application will be described in your hardware documentation, and depends on the specific board you are using and your platform configuration.

TransferMode

Characteristics	Vendor	National Instruments, ComputerBoards
	Usage	AI, AO, Common
	Access	Read/write
	Data type	String
	Read-only when running	Yes

Values	National Instruments	
	Interrupts	Transfer data using interrupts.
	SingleDMA	Transfer data using a single DMA channel.
	Dual DMA	Transfer data using two DMA channels.

This default property value is supplied by the driver. For most devices that support data transfer via interrupts and DMA, SingleDMA is the default value.

ComputerBoards		
{Default}	The transfer mode is automatically selected by the driver based on the board type and the sampling rate.	
InterruptPerPoint	Transfer single data points using interrupts.	
DMA	Transfer data using a single DMA channel (AI only).	
InterruptPerBlock	Transfer a block of data using interrupts (AI only).	

Troubleshooting Your Hardware

Overview	A-2
Troubleshooting National Instruments Hardware	A-3
What Driver Are You Using?	A-3
Is Your Hardware Functioning Properly?	A-5
Troubleshooting ComputerBoards Hardware	A-6
What Driver Are You Using?	A-6
Is Your Hardware Functioning Properly?	A-8
Troubleshooting Agilent Technologies Hardware	A-9
What Driver Are You Using?	A-9
Is Your Hardware Functioning Properly?	A-10
Troubleshooting Sound Cards	A-12
Microphone and Sound Card Types	A-15
Testing with a Microphone	A-16
Testing with a CD Player	A-16
Running in Full Duplex Mode	A-18
Other Things to Try	A-19
Registering the Hardware Driver Adaptor	A-19
Contacting The MathWorks	A-20

Overview

This appendix describes simple tests you can perform to troubleshoot your National Instruments, ComputerBoards, Agilent Technologies, and Windows sound card hardware. The tests involve using software provided by the vendor or the operating system (sound cards), and does not involve using the Data Acquisition Toolbox.

To accurately test your hardware, you should use these tools to match the requirements of your data acquisition session. For example, you should select the appropriate sampling rate, number of channels, acquisition mode (continuous or single-point), and input range.

Note If you cannot access your board using the vendor's software, then you will not be able to do so with the Data Acquisition Toolbox.

If these tests do not help you, then you may need to register the hardware driver adaptor or contact The MathWorks for support. Contact information is provided in “Contacting The MathWorks” on page A-20 as well as in the beginning of this guide. If the problem is with your hardware, then you should contact the hardware vendor.

Troubleshooting National Instruments Hardware

If you are having trouble using the Data Acquisition Toolbox with a supported National Instruments board, the reason may be that:

- You are using a hardware driver that is incompatible with the toolbox.
- Your hardware is not functioning properly.

What Driver Are You Using?


The Data Acquisition Toolbox is compatible only with specific versions of the NI-DAQ driver, and is not guaranteed to work with any other versions. For a list of the NI-DAQ driver versions that are compatible with the Data Acquisition Toolbox, refer to the product page on the MathWorks Web site at <http://www.mathworks.com/products/daq/>.

If you think your driver is incompatible with the Data Acquisition Toolbox, then you should verify that your hardware is functioning properly before updating drivers. If your hardware is functioning properly, then you are probably using unsupported drivers. Visit the National Instruments Web site at <http://www.natinst.com/> for the latest NI-DAQ drivers.

You can find out which version of NI-DAQ you are using with National Instruments' Measurement & Automation Explorer. You should be able to access this program through the Windows Desktop. The driver version is available through the **Help** menu.

Help->About Measurement & Automation Explorer->System Info->Software


For example, the version of NI-DAQ used by a PCI-6024E board is shown below.



Is Your Hardware Functioning Properly?

To troubleshoot your National Instruments hardware, you should use the **Test Panel**. The **Test Panel** allows you to test each subsystem supported by your board, and is installed as part of the NI-DAQ driver software. You can access the **Test Panel** by right-clicking on the appropriate device in the Measurement & Automation Explorer and choosing **Test Panel**.

For example, suppose you want to verify that the analog input subsystem on your PCI-6024E board is operating correctly. To do this, you should connect a known signal – such as that produced by a function generator – to one or more channels using a screw terminal panel. The result of such a test is shown below for channel 1.



If the **Test Panel** does not provide you with the expected results for the subsystem under test, and you are sure that your test setup is configured correctly, then the problem is probably with the hardware.

To get support for your National Instruments hardware, visit their Web site at <http://www.natinst.com/>.

Troubleshooting ComputerBoards Hardware

If you are having trouble using the Data Acquisition Toolbox with a supported ComputerBoards (Measurement Computing Corporation) board, the reason may be that:

- You are using a hardware driver that is incompatible with the toolbox.
- Your hardware is not functioning properly.

What Driver Are You Using?

The Data Acquisition Toolbox is compatible only with specific versions of the the Universal Library drivers or the associated release of the InstaCal software, and is not guaranteed to work with any other versions. For a list of the driver versions that are compatible with the Data Acquisition Toolbox, refer to the product page on the MathWorks Web site at <http://www.mathworks.com/products/daq/>.

If you think your driver is incompatible with the Data Acquisition Toolbox, then you should verify that your hardware is functioning properly before updating drivers. If your hardware is functioning properly, then you are probably using unsupported drivers. Visit the ComputerBoards Web site at <http://www.computerboards.com/> for the latest drivers.


Suppose you are using InstaCal software with your hardware. You can access this software through the Windows **Start** button.

Start->Programs->ComputerBoards->InstaCal

The driver version is available through the **Help** menu.

Help->About InstaCal


For example, the version of InstaCal used by a PCI-DAS4020/12 board is shown below.



Is Your Hardware Functioning Properly?

To troubleshoot your ComputerBoards hardware, you should use the test feature provided by InstaCal. To access this feature, select the board you want to test from the PC Board List, and select **Analog** from the **Test** menu.

For example, suppose you want to verify that the analog input subsystem on your PCI-DAS4020/12 board is operating correctly. To do this, you should connect a known signal – such as that produced by a function generator – to one of the channels using a BNC cable. The result of such a test is shown below for channel 0.



If InstaCal does not provide you with the expected results for the subsystem under test, and you are sure that your test setup is configured correctly, then the problem is probably with the hardware.

To get support for your ComputerBoards hardware, visit their Web site at <http://www.computerboards.com/>.

Troubleshooting Agilent Technologies Hardware

If you are having trouble using the Data Acquisition Toolbox with a supported Agilent Technologies device, the reason may be that:

- You are using a hardware driver that is incompatible with the toolbox.
- Your hardware is not functioning properly.

What Driver Are You Using?

The Data Acquisition Toolbox is compatible only with specific versions of the HP E1432 driver, and is not guaranteed to work with any other versions. You can find out which driver version you are using with the Soft Front Panel, which is described in the next section.

If you think your driver is incompatible with the Data Acquisition Toolbox, then you should verify that your hardware is functioning properly before updating drivers. If your hardware is functioning properly, then you are probably using unsupported drivers. Visit the Agilent Web site at <http://agilent.com/> for the latest drivers.


For a list of the HP E1432 driver versions that are compatible with the Data Acquisition Toolbox, refer to the product page on the MathWorks Web site at <http://www.mathworks.com/products/daq/>.

Is Your Hardware Functioning Properly?


To troubleshoot your Agilent hardware, you should use the HP E1432 Soft Front Panel. The Soft Front Panel allows you to test each module supported by the HP E1432 driver software, and is installed as part of this software. You can access the Soft Front Panel through the Windows **Start** button.

Start->Programs->hpe1432->HP E1432 Soft Front Panel

For example, suppose you want to verify that the HP E1432 module is operating correctly. To do this, you should connect a known signal – such as that produced by a function generator – to the module. You then configure the input parameters as shown below.



The result of such a test is shown below for channel 1.



If the Soft Front Panel does not provide you with the expected results for the module under test, and you are sure that your test setup is configured correctly, then the problem is probably with the hardware.

To get support for your Agilent Technologies hardware, visit their Web site at <http://www.agilent.com/>.

Troubleshooting Sound Cards


You can verify that your sound card is functioning properly by recording data and then playing back the recorded data. Recording data uses the sound card's analog input subsystem, while playing back data uses the sound card's analog output subsystem. Successful completion of these two tasks indicates your sound card works properly. The data to be recorded can come from two sources:

- A microphone
- A CD player

The first thing you should do is enable your sound card's ability to record and play data. This is done using the Microsoft Windows **Multimedia Properties** dialog box. You can access this dialog box using the Windows **Start** button.

Start->Settings->Control Panel->Multimedia


The **Multimedia Properties** dialog box for a Windows NT 4.0 platform is shown below, and is configured for both playback and recording.



You can record data and then play it back using the Windows **Sound Recorder** panel. To access this application

Start->Programs->Accessories->Multimedia->Sound Recorder


The figure below shows how to record and play data.




You must also make sure that your microphone or CD player is enabled for recording and playback using the Windows **Volume Control** panel. To access this application

Start->Programs->Accessories->Multimedia->Volume Control


The **Volume Control** panel is shown below. The CD, microphone, and line devices are enabled for playback when the **Mute** check box is unchecked for the **CD Balance**, **Microphone Balance**, and **Line Balance** volume controls, respectively. You can play . WAV files by leaving the **Mute** check box unchecked for the **Wave Balance** volume control.



If the CD, microphone, or Wave Output controls do not appear in the **Volume Control** panel, you must modify the playback properties by selecting **Properties** from the **Options** menu. The **Properties** dialog box is shown below for playback devices. Select the appropriate device check box to enable playback.



To check if the CD and microphone are enabled for recording, select the **Recording** radio button in the **Properties** dialog box, and then select the appropriate device check box to enable recording. The **Properties** dialog box is shown below for recording devices.




The **Recording Control** panel is shown below. You enable the CD or microphone for recording when the **Select** check box is checked for the **CD Balance** or **Microphone Balance** controls, respectively.



Microphone and Sound Card Types

Your microphone will be one of two possible types: powered or unpowered. You can use powered microphones only with Sound Blaster or Sound Blaster-compatible microphone inputs. You can use unpowered microphones with any sound card microphone input. Many laptops must use unpowered microphones since they do not have Sound Blaster compatible sound cards.

As shown below, you can easily identify these two microphone types by their jacks.



You can find out which sound card brand you have installed by selecting the **Devices** tab on the **Multimedia Properties** dialog box. Refer to “Troubleshooting Sound Cards” on page A-12 for a picture of this dialog box.

Testing with a Microphone

To test your sound card with a microphone, follow these steps:

- 1 Plug the microphone into the appropriate sound card jack. For a Sound Blaster sound card, this jack is labeled MIC IN.
- 2 Record audio data by selecting the **Record** button on the **Sound Recorder** and then speak into the microphone. While recording, the green line in the **Sound Recorder** should indicate that data is being captured. If this is the case, then the analog input subsystem on your sound card is functioning properly.
- 3 After recording the audio data, save it to disk. The data is automatically saved as a . WAV file.
- 4 Play the saved . WAV file. While playing, the green line in the **Sound Recorder** should indicate that data is being captured. If this is the case, then the analog output subsystem on your sound card is functioning properly.

If you are not able to record or play data, make sure that the sound card and input devices are enabled for recording and playback as described in the beginning of this section.

Testing with a CD Player


To test your sound card with a CD player, follow these steps:

- 1 Check that your CD is physically connected to your sound card.
 - Open your computer and locate the back of the CD player.
 - If there is a wire connecting the Audio Out CD port with the sound card, you can record audio data from your CD. If there is no wire connecting your CD and sound card, you must either make this connection or use the microphone to record data.

- 2 Put an audio CD into your CD player. The Windows **CD Player** application should be automatically invoked and begin playing the CD. If this doesn't occur, then you must access the application manually.

Start->Programs->Accessories->Multimedia->CD Player

The figure below shows how to play a CD with the **CD Player** application.



- 3 While the CD is playing, record audio data by selecting the **Record** button on the **Sound Recorder**. While recording, the green line in the **Sound Recorder** should indicate that data is being captured. If this is the case, then the analog input subsystem on your sound card is functioning properly. Note that the CD player converts digital audio data to analog audio data. Therefore, the CD sends analog data to the sound card.
- 4 After recording the audio data, save it to disk. The data is automatically saved as a . WAV file.
- 5 Play the saved . WAV file. While playing, the green line in the **Sound Recorder** should indicate that data is being captured. If this is the case, then the analog output subsystem on your sound card is functioning properly.

If you are not able to record or play data, make sure that the sound card and input devices are enabled for recording and playback as described in the beginning of this section.

Running in Full Duplex Mode

The term *full duplex* refers to a system that can send and receive information simultaneously. For sound cards, full duplex means that the device can acquire input data via an analog input subsystem while outputting data via an analog output subsystem at the same time.

Note that *full* tells you nothing about the bit resolution or the number of channels used in each direction. Therefore, sound cards can simultaneously receive and send data using 8 or 16 bits while in mono or stereo mode. A common restriction of full duplex mode is that both subsystems must be configured for the same sampling rate.

If you try to run your card in full duplex mode and the following error is returned

```
?? Error using ==> daqdevice/start
Device 'Winsound' already in use.
```

then your sound card is not configured properly, it does not support this mode, or you don't have the correct driver installed.

If your card supports full duplex mode, then you may need to enable this feature through the **Multimedia Properties** dialog box. Refer to “Troubleshooting Sound Cards” on page A-12 for a picture of this dialog box. If you are unsure about the full duplex capabilities of your sound card, refer to its specification sheet or user manual. It is usually very easy to update your hardware drivers to the latest version by visiting the vendor's Web site.

Other Things to Try

If troubleshooting your hardware does not help you, then you may need to register the hardware driver adaptor or contact The MathWorks for support.

Registering the Hardware Driver Adaptor

When you first create a device object, the associated hardware driver adaptor is automatically registered so that the data acquisition engine can make use of its services.

The hardware driver adaptors included with the toolbox are all located in the `daq/private` directory. The full name for each adaptor is shown below.

Table A-1: Supported Vendors and Full Adaptor Names

Vendor	Full Adaptor Name
National Instruments	<code>mwni daq. dll</code>
ComputerBoards	<code>mwcbi . dll</code>
Agilent Technologies	<code>mwhe1432. dll</code>
Windows sound cards	<code>mwwi nsound. dll</code>

If for some reason, a toolbox adaptor is not automatically registered, then you need to register it manually using the `daqregister` function. For example, to manually register the sound card adaptor

```
daqregister('wi nsound');
```

If you are using a third-party adaptor, then you may need to register it manually. If so, you must supply the full pathname to `daqregister`. For example, to register the third-party adaptor `myadaptor. dll`

```
daqregister('D: /MATLABR12/tool box/daq/myadaptors/myadaptor. dll')
```

Note You must install the associated hardware driver before adaptor registration can occur.

Contacting The MathWorks

If you need support from The MathWorks, visit our Web site at <http://www.mathworks.com/support/> or e-mail us at support@mathworks.com.

Before contacting The MathWorks, you should run the `daqsupport` function. This function returns diagnostic information such as:

- The versions of the MathWorks products you are using
- Your MATLAB path
- The characteristics of your hardware

The output from `daqsupport` is automatically saved to a text file, which you can use to help troubleshoot your problem. For example, if you are having trouble with your sound card, type

```
daqsupport('winound')
```

Managing Your Memory Resources

Overview	B-2
How Much Memory Do You Need?	B-3
Example: Managing Memory Resources	B-4

Overview

When data is acquired from an analog input subsystem or output to an analog output subsystem, it must be temporarily stored in computer memory.

The Data Acquisition Toolbox allocates memory in terms of *data blocks*. A data block is defined as the smallest “slice” of memory that the data acquisition engine can usefully manipulate. For example, acquired data is logged to a disk file using an integral number of data blocks. A representation of allocated memory using *n* data blocks is shown below.



The Data Acquisition Toolbox strives to make memory allocation as simple as possible. For this reason, the data block size and number of blocks are automatically calculated by the engine. This calculation is based on the parameters of your acquisition such as the sampling rate, and is meant to apply to most common data acquisition applications. Additionally, as data is acquired, the number of blocks dynamically increases up to a predetermined limit. However, the engine cannot guarantee that the appropriate block size, number of blocks, or total memory is allocated under these conditions:

- You select certain property values. For example, if the samples to acquire per trigger is significantly less than the FIFO buffer of your hardware.
- You acquire data at the limits of your hardware, your computer, or the toolbox. In particular, if you are acquiring data at very high sampling rates, then the allocated memory must be carefully evaluated to guarantee that samples are not lost.

You are free to override the memory allocation rules used by the engine and manually change the block size and number of blocks provided the device object is not running. However, you should do so only after careful consideration since system performance may be adversely affected, which may result in lost data.

You can manage memory resources using the `BufferingConfig` property and the `daqmem` function. With `BufferingConfig`, you can configure and return the block size and number of blocks used by a device object. With `daqmem`, you can return the current state of the memory resources used by a device object, and configure the maximum memory that one or more device objects can use.

How Much Memory Do You Need?

The memory (in bytes) required for data storage depends on these factors:

- The number of hardware channels you use
- The number of samples you need to store in the engine
- The data type size of each sample

The memory required for data storage is given by the formula

$$\text{memory required} = \text{samples stored} \times \text{channel number} \times \text{data type}$$

Of course, the number of samples you need to store in the engine at any time depends on your particular needs. The memory used by a device object is given by the formula

$$\text{memory used} = \text{block size} \times \text{block number} \times \text{channel number} \times \text{data type}$$

The block size and block number is given by the `BufferingConfig` property. The data type is given by the `NativeDataType` field of the `daqhwinfo` function.

You can display the memory resources used by (and available to) a device object with the `daqmem` function. For analog input objects, memory is used when channels are added. For analog output objects, memory is used when data is queued in the engine. For both device objects, the memory used can dynamically change based on the number of samples acquired or queued.

Example: Managing Memory Resources

Suppose you create the analog input object `ai` for a sound card, add two channels to it, and configure a four second acquisition using a sampling rate of 11.025 kHz.

```
ai = analoginput('winsound');  
addchannel(ai, 1:2);  
set(ai, 'SampleRate', 11025);  
set(ai, 'SamplesPerTrigger', 44100);
```

You return the default block size and number of blocks with the `BufferingConfig` property.

```
get(ai, 'BufferingConfig')  
ans =  
      1024      30
```

You return the memory resources with the `daqmem` function.

```
daqmem(ai)  
ans =  
      UsedBytes: 122880  
      MaxBytes: 18011136
```

The `UsedBytes` field tells you how much memory is currently used by `ai`, while the `MaxBytes` field tells you the maximum memory that `ai` can use to store acquired data. Note that the value returned for `MaxBytes` depends on the total available computer memory, and may be different for your platform.

You can verify the `UsedBytes` value with the formula given in the previous section. However, you must first find the size (in bytes) of each sample using the `daqhwinfo` function.

```
hwinfo = daqhwinfo(ai);  
hwinfo.NativeDataType  
ans =  
      int16
```

The value of the `NativeDataType` field tells you that each sample requires two bytes. Therefore, the initial allocated memory is 122,880 bytes. However, if you want to keep all the acquired data in memory, then 176,400 bytes are required.

The Data Acquisition Toolbox will accommodate this memory requirement by dynamically increasing the number of data blocks after you start ai .

```
start(ai)
```

After all the data is acquired, you can examine the final number of data blocks used by ai .

```
ai.BufferingConfig
ans =
    1024    44
```

The final total memory used is

```
daqmem(ai)
ans =
    UsedBytes: 180224
    MaxBytes: 1801136
```

Note that this was more than enough memory to store all the acquired data

Glossary

Accuracy	A determination of how close a measurement comes to the true value.
Acquiring data	The process of inputting an analog signal from a sensor into an analog input subsystem, and then converting the signal into bits that the computer can read.
Action function	An M-file function that you construct to suit your specific data acquisition needs. If you supply the action function as the value for an action property, then the function is executed when the event associated with the action property occurs.
Action property	A property associated with a specific event type. When an event occurs, the engine examines the associated action property. If an action function is given as the value for the action property, then that function is executed. All event types have an action property.
Actuator	A device that converts data output from your computer into a physical variable.
Adaptor	The interface between the data acquisition engine and the hardware driver. The adaptor's main purpose is to update the engine with properties that are unique to the hardware device.
A/D converter	An analog input subsystem.
Analog input subsystem	Hardware that converts real-world analog input signals into bits that a computer can read. This is also referred to as an AI subsystem, an A/D converter, or an ADC.
Analog output subsystem	Hardware that converts digital data to a real-world analog signal. This is also referred to as an AO subsystem, a D/A converter, or a DAC.
Bandwidth	The range of frequencies present in the signal being measured. You can also think of bandwidth as being related to the rate of change of the signal. A slowly varying signal has a low bandwidth, while a rapidly varying signal has a high bandwidth.
Base property	A property that applies to all supported hardware subsystems of a given type (analog input, analog output, etc.). For example, the <code>SampleRate</code> property is supported for all analog input subsystems regardless of the vendor.
Channel	A component of an analog input subsystem or an analog output subsystem that you read data from, or write data to.

Channel group	The collection of channels contained by an analog input object or an analog output object. For scanning hardware, the channel group defines the scan order.
Channel property	A property that applies to individual channels.
Channel skew	The time gap between consecutively sampled channels. Channel skew exists only for scanning hardware.
Common property	A property that applies to every channel or line contained by a device object.
Configuration	The process of supplying the device object with the resources and information necessary to carry out the desired tasks. Configuration consists of two steps: adding channels or lines, and setting property values to establish the desired behavior
Counter/timer subsystem	Hardware that is used for event counting, frequency and period measurement, and pulse train generation. This subsystem is not supported by the Data Acquisition Toolbox.
D/A converter	A digital to analog subsystem.
Data acquisition session	A process that encompasses all the steps you must take to acquire data using an analog input object, output data using an analog output object, or read values from or write values to digital I/O lines. These steps are broken down into initialization, configuration, execution, and termination.
Data block	The smallest “slice” of memory that the data acquisition engine can usefully manipulate.
Device object	A MATLAB object that allows you to access a hardware device.
Device-specific property	A property that applies only for specific hardware devices. For example, the <code>BitsPerSample</code> property is supported only for sound cards.
Differential input	Input channel configuration where there are two signal wires associated with each input signal — one for the input signal and one for the reference (return) signal. The measurement is the difference in voltage between the two wires, which helps reduce noise and any voltage common to both wires.

Digital I/O subsystem	Hardware that sends or receives digital values (logic levels). This is also referred to as a DIO subsystem.
DMA	Direct memory access (DMA) is a system of transferring data whereby samples are automatically stored in system memory while the processor does something else.
Engine	A MEX-file dynamic link library (DLL) file that stores the device objects and associated property values that control your data acquisition application, controls the synchronization of events, and controls the storage of acquired or queued data.
Engineering units properties	Channel properties that allow you to linearly scale input or output data.
Event	An event occurs at a particular time after a condition is met and may result in one or more actions. Many event types are automatically generated by the toolbox, while others are generated only after you configure specific properties.
Execution	The process of starting the device object and hardware device. While an analog input object is executing, you can acquire data. While an analog output object is executing, you can output data.
FIFO buffer	The first-in first-out (FIFO) memory buffer, which is used by data acquisition hardware to temporarily store data.
Full duplex	A system that can send and receive information simultaneously. For sound cards, full duplex means that the device can acquire input data via an analog input subsystem while outputting data via an analog output subsystem at the same time.
Input range	The span of input values for which an A/D conversion is valid.
Interrupts	The slowest but most common method to move acquired data from the hardware to system memory. Interrupt signals can be generated when one sample is acquired or when multiple samples are acquired.
Line	A component of a digital I/O subsystem that you can read digital values from, or write digital value to.
Line group	The collection of lines contained by a digital I/O object.

Line properties	Properties that are configured for individual lines.
Logging	A state of the Data Acquisition Toolbox where an analog input object stores acquired data to memory or a log file.
Noise	Any measurement that is not part of the phenomena of interest.
Onboard clock	A timer chip on the hardware board which is programmed to generate a pulse train at the desired rate. In most cases, the onboard clock controls the sampling rate of the board.
Output range	The span of output values for which a D/A conversion is valid.
Postrigger data	Data that is acquired and stored in the engine after the trigger event occurs.
Precision	A determination of how exactly a result is determined without reference to what the result means.
Pretrigger data	Data that is acquired and stored in the engine before the trigger event occurs.
Properties	A characteristic of the toolbox or the hardware driver that you can configure to suit your needs. The property types supported by the toolbox include base properties, device-specific properties, common properties, and channel or line properties.
Quantization	The process of converting an infinitely precise analog signal to a binary number. This process is performed by an A/D converter.
Queuing data	The process of storing data in the engine for eventual output to an analog output subsystem.
Running	A state of the Data Acquisition Toolbox where a device object is executing.
Sample rate	The per-channel rate (in samples/second) that an analog input or analog output subsystem converts data.
Sampling	The process whereby an A/D converter or a D/A converter takes a “snapshot” of the data at discrete times. For most applications, the time interval between samples is kept constant (e.g., sample every millisecond) unless externally clocked.

Scanning hardware	Data acquisition hardware that samples a single input signal, converts that signal to a digital value, and then repeats the process for every input channel used.
Sending	A state of the Data Acquisition Toolbox where an analog output object is outputting (sending) data from the engine to the hardware.
Sensor	A device that converts a physical variable into a signal that you can input into your data acquisition hardware.
Signal conditioning	The process of making a sensor signal compatible with the data acquisition hardware. Signal conditioning includes amplification, filtering, electrical isolation, and multiplexing.
Single-ended input	Input channel configuration where there is one signal wire associated with each input signal, and all input signals are connected to the same ground. Single-ended measurements are more susceptible to noise than differential measurements due to differences in the signal paths.
SS/H hardware	Data acquisition hardware that simultaneously samples all input signals, and then holds the values until the A/D converter digitizes all the signals.
Subsystem	A data acquisition hardware component that performs a specific task. The Data Acquisition Toolbox supports analog input, analog output, and digital I/O subsystems.
Trigger event	An analog input trigger event initiates data logging to memory or a disk file. An analog output trigger event initiates the output of data from the engine to the hardware.

The following abbreviations are used in this index:

Agilent – Agilent Technologies

AI – analog input

AO – analog output

CBI – ComputerBoards

DIO – digital I/O

HP – Hewlett-Packard

NI – National Instruments

A

A/D converter 1-7

input range 5-56

sampling rate 4-9

absolute time 5-19

accuracy 1-30

acquiring data 3-22

continuous

samples per trigger 5-23

simultaneous input and output 6-37

trigger repeats 5-30

single point 9-61

action function 5-52

action properties

AI object 5-46

AO object 6-27

saving property values to a MAT-file 8-3

actuator 1-4

adaptor kit 2-7

adaptors

registering A-19

supported hardware 2-7

third-party A-19

addchannel

AI object 4-4

AO object 6-4

addline 7-4

addmuxchannel 9-17

Agilent Technologies hardware

channel configuration 5-5

decimation factor 11-21

driver A-9

properties 11-4

trigger types

AI object 5-44

AO object 6-26

troubleshooting A-9

AISENSE 11-10

alias 1-37

AMUX-64T

adding channels 9-17

channel indices 9-72

analog input object

acquisition

continuous 5-23

single point 9-61

adding channels 4-4

creating 4-3

display summary 4-23

engineering units 5-56

events and actions 5-46

extracting data 5-13

logging

data 4-14

information to disk 8-6

previewing data 5-9

properties

basic setup 4-9

channel 10-7

common 10-3

configuring 3-18

sampling rate 4-9

starting 4-13

- status evaluation 4-22
- stopping 4-14
- triggers
 - configuring 5-20
 - types 4-11
- analog output object
 - adding channels 6-4
 - creating 6-3
 - display summary 6-15
 - engineering units 6-34
 - events and actions 6-27
 - output
 - continuous 6-22
 - single point 9-82
 - properties
 - basic setup 6-6
 - channel 10-11
 - common 10-8
 - configuring 3-18
 - queueing data for output 6-9
 - sampling rate 6-6
 - starting 6-9
 - status evaluation 6-14
 - stopping 6-10
 - triggers
 - configuring 6-21
 - types 6-8
- analog triggers
 - Agilent hardware 5-45
 - CBI hardware 5-43
 - NI hardware 5-39
- analog input 4-3
 - Agilent hardware 9-19
 - NI hardware 4-19
 - sound card 4-16
- analog output 6-3
 - Agilent hardware 9-22

- NI hardware 6-12
- sound card 6-11
- antialiasing filter 1-34
- array
 - data returned by `getdata` 5-13
 - device object 3-5

B

- bandwidth 1-10
- base properties 3-13
- binary vector 7-10
- `binvec2dec` 7-12
- block. *See* data block
- blocking function
 - `getdata` 5-13
 - `putdata` 6-9
- board ID 4-3
- buffer
 - configuration 10-15
 - extracting data 5-14
 - previewing data 5-10
 - queueing data 6-17

C

- calibration 1-3
- Channel 4-6
- channel gain list 4-5
- channel group
 - AI object 4-4
 - AO object 6-4
- channel names 4-7
- channel properties 3-12
 - AI object 10-7
 - AO object 10-11
- channel skew 5-7

- Channel Name 4-7
- channels 3-8
 - adding
 - AI object 4-4
 - AO object 6-4
 - descriptive names 4-6
 - input configuration 5-4
 - mapping to hardware IDs 3-9
 - referencing 4-6
- Channel Skew 5-8
- Channel SkewMode 5-8
- cleaning up the MATLAB environment
 - clear 3-25
 - delete 3-25
 - with daqfind 9-48
- clear 3-25
- clipping 6-35
- clock 5-37
- clocked acquisition 1-23
- common properties 3-12
 - AI object 10-3
 - AO object 10-8
 - DIO object 10-12
- ComputerBoards hardware
 - channel configuration 5-5
 - driver A-6
 - properties 11-4
 - trigger types (AI) 5-41
 - troubleshooting A-6
- configuring property values
 - dot notation 3-18
 - set 3-18
- constructor 3-4
- Contents 2-12
- continuous acquisition
 - example using AI and AO 6-37
 - samples per trigger 5-23

- trigger repeats 5-30
- continuous output 6-22
- creation function 3-4
- custom adaptors 2-7

D

- D/A converter 1-7
 - output range 6-34
 - sampling rate 6-6
- daqacton
 - AI example 5-53
 - default property value
 - data missed event (AI) 5-47
 - run-time error event
 - AI object 5-47
 - AO object 6-28
- daqfind 3-25
- daqhelp 2-20
- daqhwinfo 2-17
- daqmemB-4
- daqpropedit 3-20
- daqread 8-8
- daqregister A-19
- daqreset 9-44
- daqschool 2-13
- daqsupport A-20
- data
 - extracting from engine 5-13
 - previewing 5-9
 - queuing for output 6-17
- data acquisition session 3-2
 - acquiring data (AI) 4-13
 - adding channels
 - AI object 4-4
 - AO object 6-4
 - adding lines 7-4

- cleaning up 3-25
- configuring properties
 - AI object 4-9
 - AO object 6-6
- creating a device object
 - AI object 4-3
 - AO object 6-3
 - DIO object 7-3
- loading 8-3
- outputting data (AO) 6-9
- saving 8-3
- data block B-2
 - polling 5-11
- data flow
 - acquired data 2-4
 - output data 2-6
- data missed event 5-47
- DataMissedAction 5-47
- debugging your hardware A-20
 - daqsupport A-20
- dec2binvec 7-10
- delete 3-25
- demos 2-13
- descriptive names
 - channels 4-6
 - lines 7-8
- device ID 4-3
- device object 3-4
 - array 3-5
 - simultaneous input and output 6-37
 - configuring property values 3-18
 - copying 3-6
 - creating
 - AI object 4-3
 - AO object 6-3
 - DIO object 7-3
 - invalid 3-7
 - loading 8-3
 - saving 8-3
 - specifying property names 3-19
 - starting 3-23
 - stopping 3-24
- device-specific properties 3-13
 - Agilent hardware 11-4
 - CBI hardware 11-4
 - NI hardware 11-3
 - sound cards 11-5
- differential inputs 1-24
- digital I/O object
 - adding lines 7-4
 - creating 7-3
 - display summary 7-18
 - port types 7-5
 - properties
 - common 10-12
 - line 10-13
 - reading values 7-12
 - starting 7-16
 - status evaluation 7-18
 - stopping 7-16
 - writing values 7-10
- digital triggers
 - Agilent hardware
 - AI object 5-45
 - AO object 6-26
 - CBI hardware (AI) 5-42
 - NI hardware
 - AI object 5-39
 - AO object 6-25
- digital values
 - reading 7-12
 - writing 7-10
- digital i o 7-3
- disk logging 10-47

- display summary
 - AI object 4-23
 - AO object 6-15
 - DIO object 7-18
- DMA 1-28
 - NI hardware 11-25
- documentation examples 2-12
- dot notation
 - configuring property values 3-18
 - returning property values 3-17
 - saving property values to an M-file 8-3
- driver
 - Agilent hardware A-9
 - CBI hardware A-6
 - NI hardware A-3
- E**
- E1432 driver A-9
- engine 2-4
 - extracting data from 5-13
 - queuing data to 6-17
- engineering units
 - AI object 5-56
 - AO object 6-34
- event types
 - data missed (AI) 5-47
 - input overrange (AI) 5-47
 - run-time error
 - AI object 5-47
 - AO object 6-28
 - samples acquired (AI) 5-47
 - samples output (AO) 6-28
 - start
 - AI object 5-48
 - AO object 6-28
 - stop
 - AI object 5-48
 - AO object 6-28
 - timer
 - AI object 5-48
 - AO object 6-28
 - trigger
 - AI object 5-48
 - AO object 6-29
- EventLog
 - AI object 5-49
 - AO object 6-29
- events 2-3
 - AI object 5-46
 - AO object 6-27
 - displaying with showdaqevents
 - AI object 5-36
 - AO object 6-24
- example index 2-12
- examples
 - acquiring data
 - NI hardware 4-19
 - sound card 4-15
 - adding lines 7-9
 - generating timer events (DIO) 7-17
 - logging and retrieving information (AI) 8-10
 - outputting data with a sound card 6-11, 6-12
 - performing a linear conversion
 - AI object 5-57
 - AO object 6-35
 - polling the data block (AI) 5-11
 - previewing and extracting data 5-15
 - reading and writing DIO values 7-13
 - retrieving event information
 - AI object 5-51
 - AO object 6-30
 - using action properties
 - AI object 5-53

- AO object 6-31
 - using putdata 6-19
 - voice activation (AI)
 - pretriggers 5-28
 - repeating triggers 5-31
 - software trigger 5-23
- execution
 - AI object 4-13
 - AO object 6-9
 - DIO object 7-15
- external clock 1-23
 - clock sources 10-24
- extracting data 5-13

F

- fft 4-17
- FIFO 1-27
 - TransferMode 11-25
- filtering 1-34
- finding device objects 9-28
- floating signal 1-24
- flow of data
 - acquired 2-4
 - output 2-6
- flushdata 9-53
- full duplex A-18
 - BitsPerSample 11-7
- functions
 - addchannel 9-9
 - addline 9-14
 - addmuxchannel 9-17
 - analoginput 9-18
 - analogoutput 9-21
 - binvec2dec 9-24
 - clear 9-25
 - daqaction 9-27
 - daqfind 9-28
 - daqhelp 9-30
 - daqhwnfo 9-32
 - daqmem 9-34
 - daqpropedit 9-37
 - daqread 9-39
 - daqregister 9-43
 - daqreset 9-44
 - daqschool 9-45
 - dec2binvec 9-46
 - delete 9-48
 - digitalio 9-50
 - disp 9-52
 - flushdata 9-53
 - get 9-55
 - getdata 9-57
 - getsampl e 9-61
 - getvalue 9-62
 - ischannel 9-63
 - isdioline 9-64
 - isvalid 9-65
 - length 9-67
 - load 9-69
 - makenames 9-71
 - muxchannel 9-72
 - obj2mfile 9-74
 - peekdata 9-76
 - propinfo 9-78
 - putdata 9-80
 - putsampl e 9-82
 - putvalue 9-83
 - save 9-84
 - set 9-85
 - setverify 9-88
 - showdaqevents 9-90
 - size 9-92
 - start 9-94

stop 9-95
trigger 9-97
waittilstop 9-98

G

gain 1-22
 engineering units (AI) 5-56
gain list 4-5
get 3-15
getdata 5-13
 event information 9-59
 native data 10-53
 time information 5-18
getsampl e 9-61
getval ue 7-12
graphical property editor 3-20
grounded signal 1-24

H

hardware
 initializing 9-44
 resources 2-17
 scanning 1-17
 setting up 1-4
 simultaneous sample and hold 1-19
 supported vendors 2-7
hardware ID
 channel 4-4
 device (board) 4-3
 line 7-4
 mapping to channels 3-9
 port 7-4
hardware triggers
 AI object 5-37
 AO object 6-25

help 2-20
holding the last output value 11-16
HP E1432 driver A-9
HwChannel
 AI object 4-5
 AO object 6-5
HwLi ne 7-5

I

ID
 channel 4-4
 HwChannel 4-6
 device (board) 4-3
 line 7-4
 HwLi ne 7-8
 mapping to channels 3-9
 port 7-4
immediate trigger
 AI object 5-23
 AO object 6-22
Index
 AI object 4-5
 AO object 6-5
 DIO object 7-5
indexing
 channel array 4-6
 line array 7-8
initializing the hardware 9-44
Ini ti al Tri gger Ti me 5-19
input overrange event 5-47
input range 1-22
 engineering units 5-57
InputOverRangeActi on 5-47
InputType 5-4
InstaCal A-6
 hardware configuration 2-17

- internal clock 1-23
- interrupts 1-27
 - NI hardware 11-25
- invalid device object 3-7
- is channel 9-63
- is d i o l i n e 9-64
- is s n a n 5-35
- is v a l i d 9-65

J

- jitter 1-23

L

- least significant bit (DIO) 7-8
- L i n e 7-8
- line group 7-4
- line names 7-9
- line object 7-4
- line properties 3-12
- line-configurable device 7-6
- L i n e N a m e 7-9
- lines 3-8
 - adding 7-4
 - descriptive names 7-8
 - referencing 7-8
- l o a d 8-5
- loading device objects
 - MAT-file 8-5
 - M-file 8-4
- LogF i l e N a m e 8-7
- Logg i n g 5-20
- logging
 - data to memory 3-22
 - information to disk (AI) 8-6
 - file name specification 8-7

- multiple files 8-7
 - retrieving data with daqread 8-8

- Logg i n g M o d e 8-7

- LogT o D i s k M o d e 8-7

M

- makeNames 9-71
- managing data
 - acquired 5-9
 - output 6-17
- manual trigger
 - AI object 5-23
 - AO object 6-22
- mapping channels to hardware IDs 3-9
- MAT-file
 - device objects, saving to 8-5
 - properties, saving to 8-3
- maximum samples queued 10-51
- MaxS a m p l e s Q u e u e d 6-18
- Measurement and Automation Explorer A-3
 - hardware configuration 2-17
- memory resources B-2
- mono mode 4-7
- most significant bit (DIO) 7-8
- multifunction boards 1-6
- multiple device objects
 - array 3-5
 - starting 6-37
 - stopping 6-38
- multiplexing 1-12
- mux board
 - adding channels 9-17
 - channel indices 9-72
- m u x c h a n i d x 9-72

N

Name

- AI object 4-3
- AO object 6-4
- DIO object 7-3

National Instruments hardware

- AISENSE 11-10
- channel configuration 5-5
- data transfer mechanisms 11-25
- driver A-3
- properties 11-3
- trigger types
 - AI object 5-38
 - AO object 6-25
- troubleshooting A-3

native data

- getdata 9-57
- offset 10-53
- putdata 9-80
- scaling 10-55

NI-DAQ driver A-3

noise 1-33

Nyquist frequency 4-16

Nyquist theorem 1-35

O

obj2mfile 8-3

object constructor 3-4

onboard clock 1-23

one-shot acquisition 5-30

online help 2-20

output range 6-34

outputting data 3-22

- continuous 6-22
- holding the last value 11-16
- single point 9-82

overloaded functions 9-2

overrange condition 1-22

P

PC clock 1-23

peekdata 5-9

polarity 1-22

engineering units (AI) 5-56

polling the data block 5-11

port characteristics 7-5

port-configurable device 7-5

postriggers 5-28

precision 1-31

pretriggers 5-27

previewing data 5-9

properties

- BitsPerSample 11-7
- BufferingConfig 10-15
- BufferingMode 10-17
- Channel 10-18
- Channel Name 10-20
- Channel Skew 10-21
- Channel SkewMode 10-22
- ClockSource 10-24
- COLA 11-8
- Coupling 11-9
- DataMissedAction 10-26
- DefaultChannelValue 10-27
- Direction 10-28
- DriveAISenseToGround 11-10
- EventLog 10-29
- GroundingMode 11-11
- HwChannel 10-31
- HwLine 10-33
- Index 10-34
- InitialTriggerTime 10-36

InputMode 11-12
InputOverRangeAction 10-38
InputRange 10-39
InputSource 11-14
InputType 10-41
Line 10-43
LineName 10-44
Logging 10-46
LoggingMode 10-47
LogToDiskMode 10-48
ManualTriggerHwOn 10-49
MaxSamplesQueued 10-51
Name 10-52
NativeOffset 10-53
NativeScaling 10-55
NumMuxBoards 11-15
OutOfDataMode 11-16
OutputRange 10-56
Parent 10-58
Port 10-59
RampRate 11-17
RepeatOutput 10-60
Running 10-61
SampleRate 10-64
SamplesAcquired 10-66
SamplesAcquiredAction 10-67
SamplesAcquiredActionCount 10-68
SamplesAvailable 10-69
SamplesOutput 10-70
SamplesOutputAction 10-71
SamplesOutputActionCount 10-72
SamplesPerTrigger 10-73
Sending 10-74
SensorRange 10-75
SourceMode 11-18
SourceOutput 11-19
Span 11-21

StandardSampleRates 11-23
StartAction 10-76
StopAction 10-77
Sum 11-24
Tag 10-79
Timeout 10-80
TimerAction 10-81
TimerPeriod 10-82
TransferMode 11-25
TriggerAction 10-83
TriggerChannel 10-84
TriggerCondition 10-85
TriggerConditionValue 10-89
TriggerDelay 10-90
TriggerDelayUnits 10-91
TriggerRepeat 10-92
TriggersExecuted 10-93
TriggerType 10-94
Type 10-96
Units 10-97
UnitsRange 10-98
UserData 10-99
property characteristics 2-21
property types
 base 3-13
 channel 3-12
 AI object 10-7
 AO object 10-11
 common 3-12
 AI object 10-3
 AO object 10-8
 DIO object 10-12
 device-specific 3-13
 Agilent hardware 11-4
 CBI hardware 11-4
 NI hardware 11-3
 sound cards 11-5

- line 3-12
 - DIO object 10-13
- property values
 - configuring 3-18
 - default 3-19
 - graphical property editor 3-20
 - saving 8-3
 - specifying names 3-19
- propinfo 2-21
- putdata 6-17
- putsample 9-82
- putvalue 7-10

Q

- quantization 1-20
- queuing data for output 6-17
 - maximum number of samples 10-51
- Quick Reference Guide 2-13

R

- reading digital values 7-12
- Readme 2-12
- read-only properties 2-21
- registering your adaptor A-19
- relative time 5-18
- repeating triggers 5-30
- RepeatOutput 6-18
- resetting the hardware 9-44
- retrieving data from a log file 8-8
- returning property values
 - dot notation 3-17
 - get 3-15
 - set 3-14
- Runni ng 3-23
- running device objects 3-22

- run-time error event
 - AI object 5-47
 - AO object 6-28
- RuntimeErrorActi on 10-62
 - AI object 5-47
 - AO object 6-28

S

- SampleRate 5-6
- samples acquired event 5-47
- samples output event 6-28
- SamplesAcqui red 5-15
- SamplesAcqui redActi on 5-47
- SamplesAcqui redActi onCount 5-47
- SamplesAvai lable 5-15
- SamplesOutput 6-14
- SamplesOutputActi on 6-28
- SamplesOutputActi onCount 6-28
- SamplesPerTri gger 4-12
 - postrigger data 5-28
 - pretrigger data 5-27
- sampling 1-16
- sampling rate
 - AI subsystem 5-6
 - AO subsystem 6-6
- saturation 6-35
- save 8-5
- saving
 - device objects
 - MAT-file 8-5
 - M-file 8-3
 - information to disk (AI) 8-6
 - property values to a MAT-file 8-3
- scaling the data
 - AI object 5-56
 - AO object 6-34

- scanning hardware 1-17
 - channel order
 - channels
 - scan order 4-5
 - Sendi ng 6-21
 - sending data 3-22
 - sensors 1-8
 - range 10-75
 - session 3-2
 - loading 8-3
 - saving 8-3
 - set
 - configuring property values 3-18
 - returning property values 3-14
 - saving property values to an M-file 8-3
 - settling time 1-31
 - setverify 4-10
 - showdaqevents 5-36
 - AI object 5-51
 - AO object 6-31
 - signal conditioning 1-11
 - simultaneous input and output 6-37
 - simultaneous sample and hold hardware 1-19
 - single-ended inputs 1-25
 - single-point
 - acquisition 9-61
 - output 9-82
 - skew 5-7
 - software clock 1-23
 - CBI hardware 10-24
 - software trigger 5-23
 - sound cards
 - channel configuration 5-5
 - device-specific properties 11-5
 - mono mode 4-7
 - standard sample rates 11-23
 - stereo mode 4-8
 - troubleshooting A-12
 - start 3-23
 - start event
 - AI object 5-48
 - AO object 6-28
 - StartActi on
 - AI object 5-48
 - AO object 6-28
 - starting multiple device objects 6-37
 - state
 - logging 3-22
 - running 3-22
 - sending 3-22
 - status evaluation
 - AI object 4-22
 - AO object 6-14
 - DIO object 7-18
 - stereo mode 4-8
 - stop 3-24
 - stop event
 - AI object 5-48
 - AO object 6-28
 - StopActi on
 - AI object 5-48
 - AO object 6-28
 - synchronizing triggers 10-49
- T**
- third-party adaptors A-19
 - time
 - absolute 5-19
 - initial trigger 5-19
 - relative 5-18
 - Timeout 5-31
 - timer event
 - AI object 5-48

- AO object 6-28
 - DIO object 7-15
 - TimerAction
 - AI object 5-48
 - AO object 6-28
 - DIO object 7-15
 - TimerPeriod
 - AI object 5-48
 - AO object 6-28
 - DIO object 7-15
 - toolbox components
 - data acquisition engine 2-4
 - hardware driver adaptor 2-7
 - M-files 2-3
 - transducer 1-4
 - trigger
 - AI object 5-23
 - AO object 6-22
 - trigger event
 - AI object 5-48
 - AO object 6-29
 - TriggerAction
 - AI object 5-48
 - AO object 6-29
 - TriggerCondition 5-21
 - TriggerConditionValue 5-21
 - TriggerDelay 5-26
 - TriggerDelayUnits 5-26
 - triggered acquisition 5-20
 - triggered output 6-21
 - TriggerRepeat 5-30
 - triggers
 - delays 5-26
 - postriggers 5-28
 - pretriggers 5-27
 - repeating 5-30
 - samples acquired for each trigger 4-12
 - synchronizing for AI and AO 10-49
 - times
 - AI object 5-37
 - AO object 6-24
 - initial trigger 5-19
 - trigger conditions (AI) 5-21
 - trigger types
 - AI object 5-21
 - AO object 6-22
 - TriggersExecuted
 - AI object 5-36
 - AO object 6-24
 - TriggerType
 - AI object 5-21
 - AO object 6-22
 - Type
 - AI object 4-3, 4-5
 - AO object 6-4, 6-5
 - DIO object 7-3, 7-5
- ## U
- undersampling 1-35
 - Universal Library driver A-6
 - UserData
 - saving values to a MAT-file 8-3
- ## V
- verifying property values 4-10
 - voice activation example 5-23
- ## W
- waittillstop 9-98
 - writing digital values 7-10

